

## Time to fly

Hear what you want to hear by amplifying the sounds and voices most important to you.

- 1 ) Plug your earphone or headphone into the headphone jack.
- 2 ) Position the internal microphone near the preferred sound.
- 3 ) Adjust the volume to your comfort.



## Take Flight

Transform your senses. Listen to life, to the sounds you choose, with the new Pocketalker® 2.0 personal amplifier from Williams Sound.

Designed with you in mind, the 2.0 exudes beauty and grace both inside and out. Clean lines designed for comfort and style. Cutting-edge technology clearly in tune with your life.

So, spread those wings. Personally tailor your listening experiences with the 2.0's superior sound and unparalleled frequency range.

**Amplify your life with Pocketalker 2.0**

## Wherever you land



One-on-one discussions



Television listening



Conversing in the car



Restaurants

PERSONAL AMPLIFIER

1.800.843.3544 (USA)  
+1.952.943.2252 (World)

info@williamssound.com (USA)

info-intl@williamssound.com (World)

www.williamssound.com

 WILLIAMS SOUND®

©Copyright 2015 - Williams Sound, LLC BROM 052A



[Click here to view detailed product descriptions and prices.](#)  
**Call 800-397-5899 to order today!**



Beauty is **more** than skin deep.



Ergonomic, lightweight design



Extended high and low frequency range



Superior sound quality



Internal microphone



Built-in T-coil loop compatible



Rechargeable ready



### Hearing aids vs. **Pocketalker**

A hearing aid is a customized medical device that amplifies **all** sounds.

The Pocketalker personal amplifier is used to amplify **specific** sounds or voices. It helps you hear what you want to hear by reducing distracting background noise.



HELPING PEOPLE HEAR