



- Designed to resist accidental activation with you in control
- Safely encases the needle as the catheter is advanced
- Patented needle gasket designed to minimize blood exposure
- Choice of catheter biomaterials — OCRILON® polyurethane and FEP polymer
- Choice of needlepoint and hub design eases your transition to safety



Encased safety that sets the standard

Catheter Gauge and Length	PROTECTIV Safety I.V. Catheter (FEP Polymer)	PROTECTIV-W Safety I.V. Catheter (FEP Polymer)	PROTECTIV PLUS Safety I.V. Catheter (OCRILON Polyurethane)	PROTECTIV PLUS-W Safety I.V. Catheter (OCRILON Polyurethane)	Flow Rate ML/Min
14G x 1 1/4"	3048	3078	3068	3088	325
16G x 1 1/4"	3042	3072	3062	3082	215
18G x 1 1/4"	3055	3075	3065	3085	110
20G x 1"	3057	3077	3067	3087	65
20G x 1 1/4"	3056	3076	3066	3086	63
22G x 1"	3050	3070	3060	3080	38
24G x 5/8"		3073		3083	24
24G x 3/4"	3053		3063		24
50 units per box, 4 boxes per case					

www.vascularaccess.com

THE DETAILS GIVEN IN THIS LEAFLET ARE CORRECT AT THE TIME OF GOING TO PRESS. THE COMPANY RESERVES THE RIGHT TO IMPROVE THE EQUIPMENT SHOWN.

Smiths Medical ASD, Inc.

160 Weymouth Street, Rockland, MA 02370 USA

T: 800 848 1757 F: 800 621 2590 www.smiths-medical.com

Smiths Medical, part of the global technology business Smiths Group

This device is designed to reduce the risk of accidental needlesticks. However, care must still be taken to avoid needlesticks. Please refer to the Instructions for Use supplied with the product for detailed instructions, warnings and cautions. Universal Precautions must be adhered to, in accordance with CDC/OSHA standards for bloodborne pathogens, when starting or maintaining any I.V. catheter, to avoid the risk of exposure to contaminated blood. This device is PVC, DEHP and latex free.

Caution: Federal (USA) law restricts these devices to sale by or on the order of a physician.

Smiths Medical design mark, JELCO, PROTECTIV and OCRILON are trademarks of the Smiths Medical family of companies. The symbol ® indicates it is registered in the U.S. Patent and Trademark Office and certain other countries. ©2006 Smiths Medical family of companies. All rights reserved. VA505 2006.01

smiths medical