



2 YEAR WARRANTY
REGISTER ONLINE AT:
RESPONDO2.COM/WARRANTY



Ideal For Home & Travel

Lightweight & Portable!

Includes dual purpose cleaning and travel case for at home or on the go sanitizing!

Easy To Use

Long lasting rechargeable lithium-ion battery with USB charging adapter.

Maintenance Free

No additional cleaning solutions or filters required.



261 Wolfner Drive. | St. Louis, MO 63026
636-600-4030 | Toll Free: 866-333-4030
Fax: 636-600-4035

www.purifyo3.com
www.respondo2.com

Purify O₃™

Portable Ozone CPAP/BiPAP Sanitizer



*Don't take a chance
with your health!*



Sanitize with *Purify O₃*
The simple, effective way to clean
your CPAP/BiPAP mask & equipment.

www.purifyo3.com

REST EASY WITH

Purify O₃



Disinfect With Just The Touch Of A Button

Simple, Easy To Use Controls: One button push to start with automatic shut-off.

Quiet, Worry-Free Cleaning: Whisper quiet design, audible alert when cycle is complete.

Universal Design: Fits all 22mm CPAP/ BiPAP machines and masks.

Safe & Effective: Provides hospital grade sanitation in just 35 minutes.

All-In-One Solution

Patient friendly sanitizing kit includes storage/sanitizing case.



Simplifies Daily Mask Cleaning

Wipe residue from mask, push button to activate cleaning cycle, and zip mask and purify O3 into provided sanitizing bag.

Why Sanitize With Ozone Technology?

Natural: Activated oxygen, known as ozone (O₃), is a naturally formed gas and is widely used as a disinfecting agent in a number of industries world-wide.

Effective: Purify O₃ provides hospital grade sanitation and quickly, permanently kills 99.99% of viruses and bacteria.

Affordable: No pre-wash, replacement bags or extra filters required. Purify O₃ includes all you need to sanitize your device in one affordable kit!