



ROHO® QUADTRO SELECT® HIGH PROFILE® Cushion

The QUADTRO SELECT® HIGH PROFILE® features four distinct compartments, each controlled by our exclusive ISO FLO® Memory Control. This gives you ultimate control over the inflation levels across your ROHO cushion, making it easy to ensure you are positioned properly. The interconnected 4.25 in. (10.5 cm) tall air cells envelop your body, delivering maximum immersion and positioning, while protecting you from skin and soft tissue breakdown.



Pressure injury prevention and treatment



Stability



In-place customizability with ISO FLO® Memory Control



Maximum positioning control



Four independent compartments for individualized fit

Dimensions

Width: 12.75 in. to 26.25 in. / 32.5 cm to 67 cm

Depth: 11.75 in. to 67 in. / 30 cm to 51.5 cm

Height: 4.25 in. / 10.5 cm

Weight: Approx. 3.5 lbs. / 1.6 kg


Weight limit

No weight limit.
However, the cushion must be properly sized to the individual.

Limited warranty

Cushion: 36 months

Cover: 6 months



Our evidence shows the QUADTRO SELECT HIGH PROFILE is ideal for people who:

- Currently have, have a history of or are at risk for skin/soft tissue breakdown and/or any category/stage pressure injury (including deep tissue injury)
- Require accommodation or correction of pelvic asymmetry, up to approx. 3 in. (7.5 cm)
- Require pelvic and lower extremity positioning
- Have asymmetrical posture and need front-to-back and side-to-side positioning
- Need additional stability while seated
- Tend to lean or slide in their wheelchair



ROHO® QUADTRO SELECT® MID PROFILE™ Cushion

The QUADTRO SELECT® MID PROFILE™ features four distinct compartments, each controlled by our exclusive ISOFLOR® Memory Control. This gives you ultimate control over the inflation levels across your ROHO cushion, making it easy to ensure you are positioned properly. The interconnected 3.25 in. (8.5 cm) air cells deliver the ideal combination of customizable positioning and stability, all while protecting you from skin and soft tissue breakdown.

- ✓ Pressure injury prevention and treatment
- ✓ Four independent compartments for individualized fit
- ✓ In-place individualization with ISOFLOR® Memory Control
- ✓ Blend of positioning and stability

Dimensions

Width: 14.5 in. to 19.5 in. /
36.5 cm to 49.5 cm

Depth: 13.5 in. to 18.5 in. /
34.5 cm to 47 cm

Height: 3.25 in. / 8.5 cm

Weight: Approx. 3 lbs. /
1.3 kg

Weight limit

No weight limit.
However, the cushion
must be properly sized to
the individual.

Limited warranty

Cushion: 36 months
Cover: 6 months



Our evidence shows the QUADTRO SELECT MID PROFILE is ideal for people who:

- Currently have, have a history of, or are at risk for skin/soft tissue breakdown and/or any category/stage pressure injury (including deep tissue injury)
- Require pelvic and lower extremity positioning
- Have asymmetrical posture and need front-to-back and side-to-side positioning



ROHO® QUADTRO SELECT® LOW PROFILE® Cushion

Designed for those with a high activity level, the QUADTRO SELECT® LOW PROFILE® features four distinct compartments, each controlled by our exclusive ISO FLO® Memory Control. This gives you ultimate control over the inflation levels across your ROHO cushion, making it easy to ensure you are positioned properly. The interconnected 2.25 in. (5.5 cm) air cells deliver the ideal level of sitting stability, while offering maximum skin and soft tissue protection.

The QUADTRO SELECT LOW PROFILE Cushion is ideal for use as a secondary cushion, in an office chair, a sport wheelchair, in the car or anywhere you need extra cushioning and protection.



Four independent compartments for individualized fit



In-place individualization with ISO FLO® Memory Control

Dimensions

Width: 11 in. to 24.25 in. /
28 cm to 62.5 cm

Depth: 11.75 in. to 20.25 in. /
30 cm to 51.5 cm

Height: 2.25 in. / 5.5 cm

Weight: Approx. 2.25 lbs. /
1 kg

Weight limit

No weight limit.
However, the cushion
must be properly sized to
the individual.

Limited warranty

Cushion: 36 months
Cover: 6 months



Our evidence shows the QUADTRO SELECT LOW PROFILE is ideal for people who:

- Currently have, have a history of, or are at risk for skin/soft tissue breakdown and/or any category/stage pressure injury (including deep tissue injury)
- Have normal or impaired sensation
- Have asymmetrical posture and need front-to-back and side-to-side positioning