


# ROHO® ADAPTOR® PAD

The ADAPTOR Pad turns any flat or contoured surface into a safe, cushioned, protective environment. Since the air cells are not linked to one another, the ADAPTOR can be custom-cut and fit to address a particular need.

It's not a seat; it's designed to protect other pressure points. It can be cut to fit arm rests, trays, foot rests or anywhere protection is needed.

- ✓ Individually sealed air cells
- ✓ Skin/soft tissue protection
- ✓ Can be cut to size
- ✓ Absorbs shock



#### Dimensions

Height: 1 in. / 2 cm  
Width: 9.75 in. / 25 cm  
Depth: 13.5 in. / 34 cm  
Weight: 0.5 lbs. / 0.2 kg

#### Weight limit

None, if properly sized  
Limited warranty  
24 months