

TLCannula™ is now available in the following popular models of Salter-Style® Cannulas

1600TLC...A unique, patented design which sets the standard for comfort. All Salter-Style cannulas are anatomically shaped with soft, curved nasal prongs.

1606TLC...All the comfort features of the 1600 model, connecting tubing is not included.

1600QTLC...All the comfort features of the 1600 model. However, the curved face piece and nasal prongs are larger to accommodate higher flows (up to 8 LPM). It is considerably quieter and provides a more soothing environment for the patient.


1616TLC...Micro...all the comfort features of the 1600 model. This version uses a smaller molded face piece, nasal prongs and cannula tubes. It is lighter in weight, less conspicuous and is intended for low flow (up to 3 LPM) applications.

4804TLC and 4807TLC...Demand TLCannulas for use with dual port oxygen conservation delivery systems. A permanent barrier in the cannula's face piece allows one prong to sense inhalation/exhalation while oxygen is delivered through the other prong.

All TLCannulas incorporate the patented features common to all Salter-Style cannulas. Over the years our cannulas have earned customer acclaim as an industry standard for long term "wearability comfort". The gray colored, closed-cell foam cushions on TLC models are permanently installed on the cannula tubing. The foam cushions slide easily on the head set tubing to achieve the most comfortable fit and to diminish chafing and tissue injury over the ear.

TLCannula

Gray colored, soft foam ear cushions are pre-assembled over the cannula headset tubing for comfort and convenience.



Ordering Information

Salter-Style® TLCannula™

	Units (per case)	Part#
Adult, TLCannula , clear with foam cushions and 7' supply tube	25	1600TLC*
Adult, TLCannula , clear with foam cushions and 1' supply tube	25	1600TLC-1*
Adult, TLCannula , clear with foam cushions and 4' supply tube	25	1600TLC-4*
Adult, TLCannula , clear with foam cushions and 14' supply tube	25	1600TLC-14*
Adult, TLCannula , clear with foam cushions and 15' supply tube	25	1600TLC-15*
Adult, TLCannula , clear with foam cushions and 25' supply tube	25	1600TLC-25*
Adult, Quiet TLCannula , clear large face piece, foam cushions and 7' supply tube (for oxygen flows up to 8 LPM)	25	1600QTLC*
Small Adult/Pediatric, TLCannula , clear with foam cushions and 7' supply tube	25	1602TLC
Adult, TLCannula , clear with foam cushions, no supply tube	25	1606TLC
Adult, Micro TLCannula , clear, foam cushions and 7' supply tube (for oxygen flows up to 3 LPM)	25	1616TLC
Adult, Micro TLCannula , clear, foam cushions and 4' supply tube (for oxygen flows up to 3 LPM)	25	1616TLC-4
Adult, TLCannula , clear with female thread grip connector and 7' supply tube	25	1699TLC
NEW Adult, Demand TLCannula , clear with foam cushions and 4' smooth bore supply tube	10 or 25	4804TLC**
NEW Adult, Demand TLCannula , clear with foam cushions and 7' smooth bore supply tube	10 or 25	4807TLC**
Adult, Conventional Style TLCannula , clear with foam cushions and 7' supply tube	25	1053TLC*
Adult, Conventional Style TLCannula , clear with foam cushions and no supply tube	25	1057TLC

The closed cell foam cushions are available as separate items for self-installation on nasal cannulas.
(see SLC-71 Salter Accessories E-Z Wrap™ P/N 1016).



100 W. Sycamore Road - Arvin, CA 93203
Phone: 661-854-3166
1-800-421-0024
1-800-235-4203
Fax: 661-854-3850
Toll Free Fax: 1-800-628-4690
E-mail: salterlabs@us.salterlabs.com
www.salterlabs.com

ISO 9001 & ISO 13485 **CE** **0482**