

LEGACY TRACHEOSTOMY CARE

Shiley™ Disposable Inner Cannula Tracheostomy Tubes



This tube of choice offers:

- Four configurations for ventilatory support and weaning: cuffed, cuffless, fenestrated, or nonfenestrated
- Disposable inner cannula with snap-lock 15 mm connector[†]
- Soft-swivel flange that moves with the patient
- High-volume, low-pressure cuff that provides a gentle and effective seal
- Jackson sizes 4, 6, 8, and 10

Latex free. Packaged individually, sterile.

[†] Disposable inner cannulae sold separately, in boxes of 10.

