

Item# F-005SLU

The Spa Lift Ultra is an ADA compliant lift designed to provide access to both residential and commercial elevated pools & spas. The Ultra is constructed of quality stainless steel with a durable powder coat finish.

Easy push button hand control operation allows for 360 degree powered rotation. Access your personal backyard or commercial spa environment while enjoying the benefits of relaxation along with increased independence.



FEATURES:

- ADA compliant
- 400 lb (181 kg) weight capacity
- Clearances vary by wall height & width
- Wall height range 21" - 51" (54 - 129 cm)
- 360 degree powered rotation
- Adjustable pull-out leg rest
- Flip-up armrests
- 5 year structural warranty
- 5 year pro-rated electronics warranty



Available in a variety of accent colors
(blue plastics standard)



OPTIONS/ACCESSORIES:

ITEM

DESCRIPTION


F-SPAAK8	Anchor, standard tri-point anchor system
F-107UEWDA	Anchor, wood deck applications
F-450SLUC	Cover, standard lockable, blue
F-450SLUCS	Cover, for use with solar charger
F-422PLH-2	Headrest
F-045SCH	Solar charging station

PROUDLY SOLD BY:

ADA COMPLIANT


The Spa Lift Ultra

PRODUCT SPECIFICATIONS



APPLICABLE POOL PROFILES

RAISED POOL / SPA



To confirm specific application go to www.aquacreek.com for our **Pool Lift Selector**



Become a member of our new **Dealer and Architect** page. Access to CAD blocks, CSI specs and sales materials.