



ABLE LIFE™
Comfortably Independent

NEW & IMPROVED

SPACE SAVER WALKER™

**FOLDS
QUICK AND
EASY FOR
STRESS FREE
TRAVEL**



**WEIGHS ONLY 7 LBS
SUPPORTS UP TO 400 LBS**



GO ANYWHERE WITH THE MOST PORTABLE WALKER

**PRODUCT # 4200-BW
4200-RR
4200-CB**

EXPERIENCE **ABLE** LIFE



**Adjusts to fit users
4'10" - 6'8"**



**Weighs only 7 lbs
supports up to 400 lbs**



**Easily folds 4 times
smaller than the
average walker**



6" wheels

FEATURES & BENEFITS


- Portable - Quickly and easily folds up for storage in small spaces
- 6" wheels and rear Easy-Glide feet allows for trouble-free maneuvering over all surfaces
- Height adjustable, and supports up to 400 lbs

SPECIFICATIONS

- **Weight:** 7 lbs
- **Height Adjustment:** 32" - 38.5" for users 4' 10" to 6'8"
- **Wheels:** 6" diameter wheels
- **Open Handle Width:** 18"
- **Open Base Depth:** 21.5"
- **Open Base Width:** 25.5"
- **Collapsed Product Dimension:** 36" H X 5.25" W x 8" D"
- **Package Dimensions:** 36.25"L X 8.5"W X 6"H, 8 lbs
- **Case of 3 Dimensions:** 37.25"L X 9.5"W X 19.25"H, 28 lbs
- **Made Of:** Anodized Aluminum
- **Colors:** Black Walnut, Cobalt Blue, Regal Rose
- **Weight Capacity:** 400 lbs



4200-BW-Black Walnut



4200-RR-Regal Rose



4200-CB-Cobalt Blue