

## Splash Defense Transfer Bench with Curtain Guard Protection

HCPCS:




FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G



FIGURE H

- New "U" shape cut-out design enables the shower curtain liner to stay within the bath tub and avoiding water spillage outside onto the bathroom floor
- Reversible to accommodate any bathroom
- Tool-free assembly of back, legs and arm
- Height adjusts in  $\frac{1}{2}$ " increments with unique "Dual Column" extension legs
- 400 lb weight capacity
- A-frame construction provides stability
- Durable blow-molded plastic bench and backrest
- Pinch-free lever allows for push pins to be depressed without pinching fingers
- Made using recycled materials

PRODUCTS	Description
RTL12032KDR	Splash Defense Transfer Bench
SPECIFICATIONS	Description
Armrest Length	11.75"
Height Adjustment	1/2" increments
Height Setting (armrest to floor)	27.25"
Outside Width Between Legs	33.5"
Seat Depth	16"
Seat To Floor Height	16"
Seat Width	28"
Warranty	Limited Lifetime