

# Splash!™ Series



## Lift Color



GRAY MIST

Due to printing technology  
actual color may differ.

## Splash! variations include:

- Hi/Lo – Useful for facilities that have both in-ground pools and raised spas
- Spa – Designed for raised spas
- Extended Reach – For use on pools with unusual gutter configurations requiring a longer reach
- ER Hi/Lo – For facilities that have both in-ground pools and raised spas, requiring a longer reach than the Splash Hi/Lo

This fixed lifting system is designed for swimming facilities seeking to provide user-friendly access to their swimming pools.

- Third-party tested & verified ADA compliant
- Integrated armrests
- State of California compliant
- 400 lb/181.5 kg lifting capacity\*
- LiftOperator® Intelligent Controller
- 344° rotation helps ensure a safe transfer area
- Powder-coated stainless steel and aluminum construction

\* Lifting capacity of Extended Reach models are 300 lb/136 kg



Model No.	Description Pool Lift with Armrests	Shipping – (Class 85)			
		Weight	Length	Width	Height
300-0000	Splash!	280 lbs/127kg	40"/102cm	48"/122cm	60"/152cm
350-0000	Splash! Hi/Lo	280 lbs/127kg	40"/102cm	48"/122cm	60"/152cm
370-0000	Splash! Extended Reach	335 lbs/152kg	40"/102cm	48"/122cm	78"/198cm
375-0000	Splash! Spa	280 lbs/127kg	40"/102cm	48"/122cm	60"/152cm
385-0000	Splash! ER Hi/Lo	305 lbs/138kg	40"/102cm	48"/122cm	88"/224cm
390-0000	Splash! Package*	275 lbs/125kg	40"/102cm	48"/122cm	78"/198cm
390-2000	Splash! ER Package*	290 lbs/177kg	40"/102cm	48"/122cm	78"/198cm

### Notes:

- Includes battery, charger, battery console cover, water-resistant hand control, stainless steel anchor socket with cover, footrest, seat belt assembly and armrests.
- To order a lift with an Activation Key (for LiftOperator), add a K to the end of the part number
- To order without an anchor jig add an N to end of part number

\* Package includes mast cover and caddy

## Splash! 300

The Splash! 300 has the same features as the regular Splash!, but with a 300lb/159kg lifting capacity. In addition, it offers easier core drilling installation, with no additional footings required on a standard 4" reinforced deck.

Model No.	Description Pool Lift with Armrests	Shipping – (Class 85)			
		Weight	Length	Width	Height
300-3000	Splash! 300	280 lbs/127kg	40"/102cm	48"/122cm	60"/152cm
350-3000	Splash! 300 Hi/Lo	280 lbs/127kg	40"/102cm	48"/122cm	60"/152cm
375-3000	Splash! 300 Spa	280 lbs/127kg	40"/102cm	48"/122cm	60"/152cm

### Notes:

- Includes battery, charger, battery console cover, water resistant hand control, stainless steel anchor socket with cover, footrest, seat belt assembly and armrests
- To order a lift with an Activation Key (for LiftOperator), add a K to the end of the part number
- To order without an anchor jig add an N to end of part number

## Parts & Accessories:

- 1001495 Battery
- 500-1000 Spineboard Attachment
- 920-5000 Cover
- 970-5000 Seat Saver Cover
- 170-3000A Armrest Replacement (pair)
- 400-0000 Caddy
- 900-6000 Stability Strap
- 300-6200L Locking Anchor & Cap
- 900-1000 Seat Belt

A completed Deck Profile form is required with your pool lift order.