

# The Portable Outlet CPAP Battery

## Rechargeable Battery with AC and USB Outlets

Item#  
**CAP1008**

### Includes

- Main unit
- Wall charger cord  
(100-240V AC adapter/power supply brick)
- USB cable
- Travel bag



### Features

- Three-prong AC outlet
- Two USB ports
- LED power indication meter
- Travel bag for convenience
- Meets TSA/FAA battery guidelines
- 1-year warranty
- Rated for over 500 charge cycles
- Recharges in less than 4 hours
- May be used for other devices powered by AC or USB
- AC and USB ports may be used simultaneously

### Specs

- 2 lbs
- 5in wide x 7in long x 1.5in tall
- Rated Capacity: 155Wh/11Ah
- Rated Voltage: 14.5V
- Charge Input: 16.8V/2A
- USB Output: 5V/4.2A
- AC Output: 110v AC 60hz/100w
- AC Adapter Input: 100V-240V/50-60Hz



**NO ADDITIONAL ADAPTERS NEEDED**  
for use in the U.S.A