

The SureTemp® Plus 690—
 setting the standard
 of care even higher

SureTemp® Plus 690 specifications

| | |
|--|--|
| Calibration Accuracy..... | ±0.2 °F (±0.1 °C) (Monitor Mode) |
| Temperature Predict Time: | |
| Oral..... | Approx. 4-6 secs |
| Adult Axillary (18 years and older)..... | Approx. 15 secs |
| Pediatric Axillary (17 years and younger)..... | Approx. 10 secs |
| Rectal..... | Approx. 10 secs |
| Unit Dimensions..... | 8.46" x 3.18" x 2.43" (215 mm x 81 mm x 62 mm) |
| Weight..... | 12.6 oz (357 g) |
| Ambient Temperature Range..... | 50.0 °F (10 °C) to 104 °F (40.0 °C) |
| Patient Temperature Range..... | 80.0 °F (26.7 °C) to 110 °F (43.3 °C) |
| Power Source..... | (3) 1.5 Vdc AA batteries |
| Battery Life..... | Approx. 6,000 temperatures |
| Warranty..... | 2-year limited instrument, 1-year limited probe, and 90-day limited probe well |



ordering information

thermometers

- 01690-200** SureTemp Plus Electronic Thermometer (Model 690), Oral Probe with Oral Probe Well, 2-Year Limited Warranty
- 01690-201** SureTemp Plus Electronic Thermometer (Model 690), Rectal Probe with Rectal Probe Well, 2-Year Limited Warranty
- 01690-300** SureTemp Plus Electronic Thermometer (Model 690), One Per Room, Wall Mount, 9ft Oral Probe with Oral Probe Well, 2-Year Limited Warranty (*must be locked to the wall*)
- 01690-301** SureTemp Plus Electronic Thermometer (Model 690), One Per Room, Wall Mount, 9ft Rectal Probe with Rectal Probe Well, 2-Year Limited Warranty (*must be locked to the wall*)
- 01690-400** SureTemp Plus Electronic Thermometer (Model 690), Wall Mount, Oral Probe with Oral Probe Well, 2-Year Limited Warranty
- 01690-700** SureTemp Plus Electronic Thermometer (Model 690), Wall Mount, 9ft Oral Probe with Oral Probe Well, 2-Year Limited Warranty, with Rolling Stand and Holder for Spare Probe & Well

accessories

- 21326-0000** Wall Holder
- 21330-0000** Holder for Spare Probe and Well
- 06138-000** Calibration Key Assembly, 690/692
- 02892-000** Probe and Well Kit, 4 ft, Rectal
- 02892-100** Probe and Well Kit, 9 ft, Rectal
- 02893-000** Probe and Well Kit, 4 ft, Oral
- 02893-100** Probe and Well Kit, 9 ft, Oral
- 02891-0000** Probe Well, Blue
- 02891-1000** Probe Well, Red
- 08273-000** Thermometer Stand with Holder for Spare Probe and Well
- 21333-0000** Bracket, M600 Stand Upgrade, with Holder for Spare Probe and Well
- 98024-0000** Extended Warranty M690

replacement probe covers

- 05031-101** Case of 1,000 Probe Covers

Call Welch Allyn to order the SureTemp® Plus 690
 4341 State Street Road, PO Box 220
 Skaneateles Falls, NY 13153-0220 USA
 (p) 800.535.6663 • (f) 315.685.3361
www.welchallyn.com

WelchAllyn®

Advancing Frontline Care™

SureTemp is a registered trademark of Welch Allyn © 2006

SM2546 Rev D



Introducing the SureTemp® Plus 690 Electronic Thermometer

Now taking temperatures is easier, safer, and more accurate.
 And taking care of all your patients is more convenient.

WelchAllyn®

Advancing Frontline Care™



The SureTemp® Plus 690—
with the speed and accuracy
you require, and the versatility
you need for use with all your
patients, all the time.

- **Easy to use, intuitive operation**
- **Removable probe well reduces the risk of cross-contamination**
- **Improved accuracy for even more reliable readings**
- **Greater versatility for use with all your patients, all the time**

When we set out to improve the Welch Allyn SureTemp, we knew we had our work cut out for us. After all, when you start with the thermometer that helped set the standard for speed and accuracy, the expectations run pretty high.

So how did we do? See for yourself by taking a look at the innovative SureTemp Plus 690—easier, safer and more accurate than the original, and still built to withstand constant use in demanding environments.

Easy: An easy-to-read LCD display with flashing icons that clearly communicate in any language, and an intuitive user interface and ergonomic shape make using the SureTemp Plus 690 a breeze.

Safe: Patient safety is improved thanks to an interchangeable probe well that reduces the risk of cross-contamination.

Fast and Accurate: In addition to improved accuracy, the SureTemp Plus 690 is one of the fastest thermometers on the market, providing approximately 4-6 second oral, 10-second pediatric axillary, 15-second adult axillary, and 10-second rectal measurements.

Versatile: To top it all off, you can now use the SureTemp Plus 690 with all your patients, from infants to adults—whether you need to take oral, axillary or rectal temperatures.



Optional wall mount has storage space for 50 additional probe covers

Interchangeable rectal probe and well reduce risk of cross-contamination

The SureTemp® Plus 690

WelchAllyn®

Advancing Frontline Care™