

Ventyv[®]

Nitrile PF Exam
PLUS 3.5 Violet Blue



Ventyv[®]
Outsmart Infection[®]



sritrangusa.com | ventyv.com | Hello@sritrangusa.com | 1.888.483.6898 |
+1 (813) 606-4301 | 5820 West Cypress Street, Suite H | Tampa, FL 33607

Ventyv®

Nitrile PF Exam
PLUS 3.5 Violet Blue



Ventyv® Nitrile Powder-Free Plus 3.5 Violet Blue (Elephant) Exam Glove is ideal for use in a large variety of clinical settings. This high end glove is a favorite amongst healthcare workers for its strength and polymer coating* to make it easy to put on and take off.



Elephant



Reorder #	Size	Gloves/ Box ¹	Boxes/ Case
10333101	XS	200	10
10334101	S	200	10
10335101	M	200	10
10336101	L	200	10
10337101	XL	180	10

PRODUCT ATTRIBUTES

- Violet Blue
- Nitrile²
- Powder-free
- Polymer Coated*
- Tested for Use with Chemotherapy Drugs⁴
- Fingertip Textured

THICKNESS (Single Wall)

- Finger: 0.05 mm / 1.96 mils
- Palm: 0.05 mm / 1.96 mils

PHYSICAL ATTRIBUTES³

- Length: 220 mm (XS, S)
230 mm (M, L, XL)
- Width: 70 ± 10 mm (XS)
80 ± 10 mm (S)
95 ± 10 mm (M)
110 ± 10 mm (L)
120 ± 10 mm (XL)

PHYSICAL PROPERTIES

- Tensile strength before aging: 14 MPa
- Tensile strength after aging: 14 MPa
- Elongation before aging: 500%
- Elongation after aging: 400%



80%

Inner boxes made with 80% recycled material

90%

Outer cartons made with 90% recycled material

1. By weight
2. This product is not made with Natural Rubber Latex
3. Physical Attributes meet or exceed ASTM standards for exam gloves
4. CAUTION: Gloves used for protection against chemotherapy drug exposure must be selected specifically for the type of chemicals used. Review material safety data sheets for the drug being used to determine the required level of protection. Contact Ventyv for specific testing information.
* Polymer coating optional

Ventyv®
Outsmart Infection®



sritrangusa.com | ventyv.com | Hello@sritrangusa.com | 1.888.483.6898 |
+1 (813) 606-4301 | 5820 West Cypress Street, Suite H | Tampa, FL 33607