

## Calcium Alginate Dressing with Antimicrobial Silver

**Calcium Alginate Dressings with Antimicrobial Silver:** Assist with infection reduction. Primary dressing for wounds with moderate to heavy exudate. Silver ions are released in the presence of wound fluid and protect the dressing from bacterial colonization and provide an effective barrier to bacterial penetration. Requires a secondary dressing.

**Indications:** Effective barrier to bacterial penetration in moderate to heavily exuding wounds such as diabetic foot ulcers, leg ulcers, pressure ulcers, partial and full thickness wounds, donor sites and traumatic and surgical wounds.

**Contraindications:** Third degree burns, patients with a known sensitivity to alginates or silver and controlling of heavy bleeding. Should not be used with products that create an ionic exchange in the wound site, i.e., do not use with enzymatic debriding agents.

**Change Frequency:** May be left in place for up to seven days. Change dressing as directed by a healthcare professional.

**Recommended Secondary Dressings:** Cover with an appropriate secondary dressing depending on the amount of exudate and frequency of dressing change.



### ORDERING INFORMATION

DESCRIPTION	MFR #	QTY	HCPCS	EXUDATE LEVEL*
Calcium Alginate Dressing, Sheet, with Antimicrobial Silver, Sterile, 2" x 2"	3557	10ea/bx	A6196	●●○
Calcium Alginate Dressing, Sheet, with Antimicrobial Silver, Sterile, 4" x 4¾"	3558	10ea/bx	A6196	●●●
Calcium Alginate Dressing, Sheet, with Antimicrobial Silver, Sterile, 4" x 8"	3559	5ea/bx	A6197	●●●
Calcium Alginate Dressing, Rope, with Antimicrobial Silver, Sterile, ¾" x 12"	3560	5ea/bx	A6199	●●●

EXUDATE LEVEL KEY: ●○● = Light ●●○ = Moderate ●●● = Heavy

\*Exudate level indicates highest level of exudate that the dressing can accommodate

**Satisfaction Guaranteed**

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

### Never Compromise on Care<sup>SM</sup>

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them – making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.