

Cardinal Health™ Sterilization Wrap

Your trusted source for ease and assurance



Everything you need, wrapped in a single solution

Every day, you're tasked to do more with less. To improve productivity with less time, budget and staff. And do it all while you strive to keep your patients safe and improve the quality of care.

That's why we created a new choice in sterilization wrap: a full selection of products that are durable, reliable and cost effective.

All the selection you want

Whatever your clinicians need, now you have a single product line that can take care of it all. Choose the wrap that best meets your goals of maximizing efficiency, safety and quality of care.


Two Color:

The blue and green layers help make it fast and easy to differentiate, identify and inspect trays

Single Layer:

Supports sequential wrapping methodologies

- Made from SMS polypropylene, non-woven material for strength and sterility maintenance.
- Designed with a simple, "cross-stitch" bond pattern for clear, visual feedback of breaches that may occur during handling, storage and/or sterilization.
- Available in six basis weights and a wide range of sizes to meet the needs of today's Central Sterile Department.



Cardinal Health™ Sterilization Wraps

Two Color

- Made from two layers of SMS — one layer blue, the other green — bonded together along opposing sides
- Delivers the protection of double-wrapping in a single sheet (simultaneous wrapping) — resulting in time and productivity efficiencies when wrapping and opening^{1,2}
- Select the color you want — blue or green, to help differentiate and identify instrument trays by type or use

Single Layer

- Made from one layer of blue SMS
- Supports simultaneous double or sequential wrapping

Two Color

	Catalog no.	Size/Dimensions	Case qty.		Catalog no.	Case qty.		
CH100	CH1G0012	12 x 12 in. (30 x 30cm)	480	CH100	CH110012	1000		
	CH1G0015	15 x 15 in. (38 x 38cm)	480		CH110015	1000		
	CH1G0018	18 x 18 in. (45 x 45cm)	480		CH110018	1000		
	CH1G0020	20 x 20 in. (50 x 50cm)	480		CH110020	1000		
	CH1G0024	24 x 24 in. (60 x 60cm)	240		CH110024	500		
	CH1G0030	30 x 30 in. (76 x 76cm)	144		CH110030	300		
	CH1G0036	36 x 36 in. (91 x 91cm)	144		CH110036	300		
	CH1G0040	40 x 40 in. (101 x 101cm)	120		CH110040	250		
	CH1G0045	45 x 45 in. (114 x 114cm)	120		CH110045	250		
	CH1G0048	48 x 48 in. (121 x 121cm)	120		CH110048	250		
	CH1G0054	54 x 54 in. (137 x 137cm)	96		CH110054	100		
	n/a	54 x 72 in. (137 x 182cm)	n/a		CH110072	100		
CH200	CH2G0012	12 x 12 in. (30 x 30cm)	480	CH200	CH210012	1000		
	CH2G0018	18 x 18 in. (45 x 45cm)	480		CH210018	1000		
	CH2G0024	24 x 24 in. (60 x 60cm)	240		CH210024	500		
	CH2G0030	30 x 30 in. (76 x 76cm)	144		CH210030	300		
	CH2G0036	36 x 36 in. (91 x 91cm)	144		CH210036	300		
	CH2G0040	40 x 40 in. (101 x 101cm)	120		CH210040	250		
	CH2G0048	48 x 48 in. (121 x 121cm)	48		CH210048	100		
	CH2G0054	54 x 54 in. (137 x 137cm)	48		CH210054	100		
	CH2G0072	54 x 72 in. (137 x 182cm)	48		CH210072	100		
	CH300	CH3G0024	24 x 24 in. (60 x 60cm)		240	CH300	CH310024	500
CH3G0030		30 x 30 in. (76 x 76cm)	144	CH310030	300			
CH3G0036		36 x 36 in. (91 x 91cm)	144	CH310036	300			
CH3G0040		40 x 40 in. (101 x 101cm)	96	CH310040	250			
CH3G0045		45 x 45 in. (114 x 114cm)	48	CH310045	100			
CH3G0048		48 x 48 in. (121 x 121cm)	96	CH310048	100			
CH3G0054		54 x 54 in. (137 x 137cm)	48	CH310054	100			
CH3G0072		54 x 72 in. (137 x 182cm)	48	CH310072	100			
n/a	18 x 18 in. (45 x 45cm)	n/a	CH410018	600				
CH400	CH4G0024	24 x 24 in. (60 x 60cm)	120	CH400	CH410024	250		
	n/a	30 x 30 in. (76 x 76cm)	n/a		CH410030	200		
	CH4G0036	36 x 36 in. (91 x 91cm)	72		CH410036	150		
	CH4G0040	40 x 40 in. (101 x 101cm)	72		CH410040	150		
	CH4G0045	45 x 45 in. (114 x 114cm)	48		CH410045	100		
	CH4G0048	48 x 48 in. (121 x 121cm)	48		CH410048	50		
	CH4G0054	54 x 54 in. (137 x 137cm)	48		CH410054	50		
	CH4G0072	54 x 72 in. (137 x 182cm)	24		CH410072	50		
	CH500	CH5G0018	18 x 18 in. (45 x 45cm)		288	CH500	CH510018	600
		CH5G0024	24 x 24 in. (60 x 60cm)		120		CH510024	250
CH5G0030		30 x 30 in. (76 x 76cm)	120	CH510030	250			
CH5G0036		36 x 36 in. (91 x 91cm)	72	CH510036	150			
CH5G0045		45 x 45 in. (114 x 114cm)	48	CH510045	100			
CH5G0048		48 x 48 in. (121 x 121cm)	48	CH510048	50			
CH5G0054		54 x 54 in. (137 x 137cm)	48	CH510054	50			
CH5G0060		60 x 60 in. (152 x 152cm)	24	CH510060	50			
CH5G0072		54 x 72 in. (137 x 182cm)	24	CH510072	50			
CH5G0090		54 x 90 in. (137 x 228cm)	24	CH510090	50			
CH600	CH6G0036	36 x 36 in. (91 x 91cm)	72	CH600	CH610036	144		
	CH6G0040	40 x 40 in. (101 x 101cm)	48		CH610040	96		
	CH6G0045	45 x 45 in. (114 x 114cm)	48		CH610045	96		
	CH6G0048	48 x 48 in. (121 x 121cm)	48		CH610048	48		
	CH6G0054	54 x 54 in. (137 x 137cm)	48		CH610054	48		
	CH6G0072	54 x 72 in. (137 x 182cm)	24		CH610072	48		

1. Webster, Joan et al. Barrier properties and cost implications of a single versus a double wrap for storing sterile instrument packs. American Journal of Infection Control, Volume 33, Issue 6, 348 - 352.

2. Rutala WA, Weber DJ. Choosing a sterilization wrap. Infection Control Today 2000;4:64,70.

Cardinal Health™ Sterilization Wraps (continued)

Paper

- Made from wetlaid cellulose — blue in color
- Single layer format to support sequential or simultaneous double wrapping techniques
- Available in one basis weight and 16 sizes

Paper

Catalog no.	Size/Dimensions	Case qty.
4008	8 x 8 in. (20 x 20cm)	1000
4009	9 x 9 in. (23 x 23cm)	1000
4012	12 x 12 in. (30 x 30cm)	1000
4015	15 x 15 in. (38 x 38cm)	500
4018	18 x 18 in. (45 x 45cm)	500
4020	20 x 20 in. (50 x 50cm)	500
4024	24 x 24 in. (60 x 60cm)	500
4030	30 x 30 in. (76 x 76cm)	250
4036	36 x 36 in. (91 x 91cm)	250
4040	40 x 40 in. (101 x 101cm)	250
4045	45 x 45 in. (114 x 114cm)	250
4048	48 x 48 in. (121 x 121cm)	100
4054	54 x 54 in. (137 x 137cm)	100
4060	40 x 60 in. (101 x 152cm)	100
4065	60 x 60 in. (152 x 152cm)	100
4072	54 x 72 in. (137 x 182cm)	100

**We are your
trusted source for
everything you need**

to be efficient,
productive —
and save money
at the same time.

Ask your Cardinal Health representative
for a trial today or call customer service
at 877.CARDINAL.



cardinalhealth.com

Spunbond/Meltblown/Spunbond (SMS): Consists of meltblown polypropylene between layers of spunbonded polypropylene. It is strong, fluid-resistant, breathable and provides durable protection.

© 2016 Cardinal Health. All Rights Reserved. CARDINAL HEALTH, the Cardinal Health LOGO and ESSENTIAL TO CARE are trademarks or registered trademarks of Cardinal Health.
Lit. No. 2PERI16-529512 (05/2016) PRC (03/2016)

Cardinal Health
Waukegan, Illinois