

carefix[®]
Creating functional comfort



CarePocket™
Leg bag holder



TYTEX[®]

CarePocket™ is specially designed to fixate a urine leg bag, without impeding appearance. The weight of a full leg bag is supported by the whole product. This way, the tension required to hold a bag in place is distributed evenly around the leg, rather than being concentrated in two narrow rings at the top and the bottom of the bag. This design not only reduces the risk of any constriction or pressure sores, but also prevents the bag from moving around and helps to dampen any undesired sounds, thus providing discretion.

Features:

- Easy to use
- Safe retention of leg bag for optimum support
- Tension is distributed evenly around the leg
- Helps to reduce unpleasant noises
- Reduces the risk of pressure sores
- Provides extra security during physical activity
- Very soft, stretchy and breathable material, that is comfortable to wear
- Fits left and right leg equally well
- The inner pocket has leg bag opening as well as bottom tap opening
- Can hold up to 750 ml/25 oz
- Practical colour code for size indication
- Washable at 60°C/140°F
- Latex-free
- Oeko-Tex certified

Indications:

For fixation of urinary leg bags.

Composition:

80% Polyester, 20% Elastane.

Washing instructions:



Recommended sizes:

Based on leg circumference.

SIZE RANGE	S		M		L		XL	
CarePocket™	CM	INCHES	CM	INCHES	CM	INCHES	CM	INCHES
	30 - 45	12 - 17½	45 - 65	17½ - 25½	55 - 80	21½ - 31½	65 - 90	25½ - 35½

Packaging:



1 pc. or 5 pcs. per bag.

Instruction for use:



Pull on the leg bag holder (ensure the color size code is at the top).



Place the leg bag in the pocket and pull the tap through the lower whole.



The CarePocket™ is in place.

For further information please contact:

Dealer: