

## Bariatric Aluminum Rehab Shower Commode Chair with Four Locking Casters



- Can be positioned over a standard toilet or used as a portable self-contained commode
- Floor to seat height can be adjusted 4" without tools
- White aluminum frame is attractive and rust resistant
- Padded seat back and arms
- Height adjustable
- Comes with swing-away footrests

PRODUCTS	Description
NRS185008	White and Blue, Back

SPECIFICATIONS	Description
Overall Length (with Riggings)	35.04"
Overall Width	26"
Overall Height	40" - 42"
Seat to Floor Height	20.5" - 22.4"
Seat Dimensions	22" (W) x 17.1" (H)
Backrest Height	20.1"
Net Weight	36.8 lbs.
Gross Weight	43.3 lbs.
Weight Capacity	500 lbs.
Warranty	Limited Lifetime

ACCESSORIES	Description
RTL12085	Biodegradable Sanitary Commode Liner, 12/bx, 6 bx/cs

PARTS	Description
NRS185008	Back
NRS185006-02	Seat Clip, 1/ea
NRS185006-04	Arm Pads, 1/pr
NRS185006-05	Left Armrest Bracket, 1/ea
NRS185006-06	Right Armrest Bracket, 1/ea
NRS185006-07	Armrest Tube, 1/ea
NRS185006-08	Caster with Leg, 1/ea
NRS185006-09	Frame Clips, 4/pkg
NRS185006-11	Swing Away Footrest, 1/pr
NRS185006-12	Footplates, 1/pr
NRS185006-13	Angle Adj. Footplate Mnt. & Hdwr., 1/ea
NRS185006-14	Down Tube & Plastic Insert, 1/ea
NRS185006-15	Footrest Height Adj. Knob & Bolt, 1/ea

PARTS	Description
NRS185006-16	Footrest Pivot Saddle & Bolts, 1/ea
NRS0185006-17	Footrest Frame Latch
NRS185006-21	Back Locking Knob, 1/pr
NRS185006-22	Anti-Rattle Insert, 1/ea
NRS185006-23	Frame Plugs, 6/pkg
11106-1	Pail, 1/ea
NRS185008-SEAT	Seat, 1/ea
NRS185008-Back	Backrest Only, 1/ea
STDS851-SA	Velcro Seat Belt, 1/ea
NRS185006-17	FootrestFrameLatchForNRS185006/7, 1/ea
NRS185006	Self Propel Shower Commode Chair, 1/cs