



**Urology Products**  
**1-800-962-9888**



## Urology Products

### Table of Contents

#### Dover™ Silver Products

Foley Catheters ..... 1  
 Urine Meter Foley Trays ..... 2  
 Foley Trays ..... 3

#### Dover™ 100% Silicone Products

Foley Catheters ..... 4  
 Urine Meter Foley Trays ..... 5  
 Foley Trays ..... 6

#### Dover™ Latex Products

Foley Catheters .....7-8  
 Urine Meter Foley Trays ..... 9  
 Foley Trays ..... 10

#### Dover™ Insertion Trays

Silver Insertion Trays ..... 11  
 Silicone Insertion Trays ..... 11  
 Latex Insertion Trays..... 11  
 Universal Trays..... 11

#### Dover™ Collection Products

Urine Meters ..... 12  
 Add-A-Foley Urine Meter Trays ..... 12  
 Drainage Bags ..... 13  
 Add-A-Foley Drainage Bag Trays..... 14  
 Miscellaneous Accessories ..... 14

#### Dover™ Urethral Products

Latex Free Catheters ..... 15  
 Latex Catheters ..... 16  
 Urethral Trays..... 17

#### Patient Care Products

Male External Catheters ..... 18  
 Leg Bags & Accessories ..... 18  
 Irrigation Trays..... 19  
 Midstream Catch Kits ..... 19  
 Specimen Containers..... 20  
 Enema and Ancillary Products ..... 20

#### Index

Index ..... 21–25

## Dover™ Silver Foley Catheters



- 100% silicone with silver hydrogel coating
- 16" length
- Rounded tip, 2 opposing eyes
- 100% silicone Foley catheter eliminates latex allergy concerns
- Outer and inner catheter and balloon surfaces contain a patented silver and hydrogel coating which provides a consistent release of ionic silver particles
- Hydrogel coating eases insertion to prevent urethra trauma

#### 2-Way

Code	Size (Fr)	Ship Carton
------	-----------	-------------

##### 5 cc, 2-way

605122IC	12 Fr	10
605148IC	14 Fr	10
605163IC	16 Fr	10
605189IC	18 Fr	10
605205IC	20 Fr	10
605221IC	22 Fr	10
605247IC	24 Fr	10
605262IC	26 Fr	10

##### 30 cc, 2-way

630187IC	18 Fr	10
630203IC	20 Fr	10
630229IC	22 Fr	10
630245IC	24 Fr	10
630260IC	26 Fr	10

#### 3-Way

Code	Size (Fr)	Ship Carton
------	-----------	-------------

##### 30 cc, 3-way

665183IC	18 Fr	10
665209IC	20 Fr	10
665225IC	22 Fr	10
665241IC	24 Fr	10

#### Specialty

Code	Size (Fr)	Ship Carton
------	-----------	-------------

##### 5 cc, 2-way, Coudé Tip

20512ICC	12 Fr	10
20514ICC	14 Fr	10
20516ICC	16 Fr	10
20518ICC	18 Fr	10
20520ICC	20 Fr	10
20522ICC	22 Fr	10

##### 30 cc, 2-way, Coudé Tip

23018ICC	18 Fr	10
23020ICC	20 Fr	10

##### 5 cc, 2-way, 400 Series Temperature Sensing

50516ICT	16 Fr	12
50518ICT	18 Fr	12



## Dover™ Silicone Foley Catheters



- 100% silicone Foley catheter eliminates latex allergy concern
- 16" length
- Biocompatible, inorganic, non-porous, non-toxic, non-adhering material reduces incidence of encrustation, urethritis, and stricture vs. latex catheters
- Packaged in an outer peel-back pouch and an easy-to-open inner sheath

### 2-Way

Code	Size (Fr)	Ship Carton
------	-----------	-------------

#### 3 cc, 2-way (Pediatric) 12"

8887603085	8 Fr	10
8887603101	10 Fr	10
6108	8 Fr (Luer Lock)	12
6010	10 Fr (Luer Lock)	12

#### 5 cc, 2-way

8887605122	12 Fr	10
8887605148	14 Fr	10
8887605163	16 Fr	10
8887605189	18 Fr	10
8887605205	20 Fr	10
8887605221	22 Fr	10
8887605247	24 Fr	10
8887605262	26 Fr	10
8887605282	28 Fr	10
8887605303	30 Fr	10

#### 30 cc, 2-way

8887630142	14 Fr	10
8887630163	16 Fr	10
8887630187	18 Fr	10
8887630203	20 Fr	10
8887630229	22 Fr	10
8887630245	24 Fr	10
8887630260	26 Fr	10
8887630284	28 Fr	10
8887630305	30 Fr	10

### 3-Way

Code	Size (Fr)	Ship Carton
------	-----------	-------------

#### 5 cc, 3-way

8887664161	16 Fr	10
8887664182	18 Fr	10
8887664203	20 Fr	10
8887664224	22 Fr	10
8887664245	24 Fr	10
8887664266	26 Fr	10
8887664287	28 Fr	10

#### 30 cc, 3-way

8887665161	16 Fr	10
8887665183	18 Fr	10
8887665209	20 Fr	10
8887665225	22 Fr	10
8887665241	24 Fr	10
8887665262	26 Fr	10
8887665283	28 Fr	10

**Catheter Kits, 5 cc, 2-way Foley**  
Includes 100% Silicone Foley catheter with lubricating jelly (5 g) and 10 cc prefilled water syringe

Code	Size (Fr)	Ship Case
8887606146	14 Fr	20
8887606161	16 Fr	20
8887606187	18 Fr	20
8887606203	20 Fr	20

### Specialty

Code	Size (Fr)	Ship Carton
------	-----------	-------------

#### 5 cc, 2-way, 400 Series Temperature Sensing

50514T	14 Fr	12
50516T	16 Fr	12
50518T	18 Fr	12

#### 5 cc, 2-way, Coudé tip

20512C	12 Fr	10
20514C	14 Fr	10
20516C	16 Fr	10
20518C	18 Fr	10
20520C	20 Fr	10
20522C	22 Fr	10
20524C	24 Fr	10

#### 30 cc, 2-way, Coudé tip

23018C	18 Fr	10
23020C	20 Fr	10
23022C	22 Fr	10
23024C	24 Fr	10
23026C	26 Fr	10

#### 5 cc, 2-way, Council tip

40516L	16 Fr	10
40518L	18 Fr	10
40520L	20 Fr	10
40522L	22 Fr	10
40524L	24 Fr	10

## Dover™ Silicone Urine Meter Trays



- Standard Foley prep tray (100% latex-free)  
Sterile tray contains: CSR wrap, fenestrated drape, underpad, blue nitrile exam gloves, forceps, prepping balls, 10 cc prefilled sterile water inflation syringe, unique bellows lubricant syringe (10 g), povidone-iodine solution, specimen container with label, Foley reminder label and patient education pamphlet
- 100% silicone Foley catheter
- Extra long 60" tubing



Code	Size (Fr)	Pre-attached Securement Device	Clampable Tubing with Top Vent	Sample Port			Temperature Sensing 400 Series	Tube Management with Sheet Clip	Spring Meter Port	Drain		Self-leveling Bed Hook	Ship Case
				Luer Lock	Needle-less	Needle				Enclosed Spout	Drain Tube		

#### 400 mL Meter

P4P14SD	14 Fr	•	•	•				•	•	•		•	10
P4P16SD	16 Fr	•	•	•				•	•	•		•	10
P4P18SD	18 Fr	•	•	•				•	•	•		•	10
P4P14TSD	14 Fr	•	•	•		•		•	•	•		•	10
P4P16TSD	16 Fr	•	•	•		•		•	•	•		•	10
P4P18TSD	18 Fr	•	•	•		•		•	•	•		•	10
2004LL	14 Fr		•	•				•	•	•		•	10
7006LL	16 Fr		•	•				•	•	•		•	10
7008LL	18 Fr		•	•				•	•	•		•	10

2004	14 Fr			•				•	•	•		•	10
2006	16 Fr			•				•	•	•		•	10
2008	18 Fr			•				•	•	•		•	10

P4P16TS	16 Fr		•	•		•	•	•	•	•		•	10
P4P18TS	18 Fr		•	•		•	•	•	•	•		•	10

#### 200 mL Meter

8122	16 Fr					•			•			•	10
8223	18 Fr					•			•			•	10

#### 400 mL Coudé Tip

P4P16CSD	16 Fr	•	•	•				•	•	•		•	10
P4P18CSD	18 Fr	•	•	•				•	•	•		•	10

#### 400 mL Meter Coudé Tip with Hand Sanitizer and Castile Soap Wipe

P4P16CSDK	16 Fr	•		•				•	•	•		•	10
P4P18CSDK	18 Fr	•		•				•	•	•		•	10

## Dover™ Silicone Foley Trays



- Standard Foley prep tray (100% latex-free)  
Sterile tray contains: CSR wrap, fenestrated drape, underpad, blue nitrile exam gloves, forceps, prepping balls, 10 cc prefilled sterile water inflation syringe (10 g), unique bellows lubricant syringe, povidone-iodine solution, specimen container with label, Foley reminder label and patient education pamphlet
- 100% silicone Foley catheters
- 2000 mL capacity
- 48" tubing



Code	Size (Fr)	Pre-attached Securement Device	Top Vent	Sample Port			Tube Management with Sheet Clip	Sheet Clip	Drain		Anti-Reflex Device	Self-leveling Bed Hook	Ship Case
				Luer Lock	Needle-less	Needle			Enclosed Spout	Tube			
PP14SD	14 Fr	•	•	•			•		•		•	•	10
PP16SD	16 Fr	•	•	•			•		•		•	•	10
PP18SD	18 Fr	•	•	•			•		•		•	•	10
6144LL	14 Fr		•	•			•		•		•	•	10
6146LL	16 Fr		•	•			•		•		•	•	10
6148LL	18 Fr		•	•			•		•		•	•	10
6947LF	14 Fr				•		•	•		•	•	•	10
8949	16 Fr				•		•	•		•	•	•	10
8950	18 Fr				•		•	•		•	•	•	10
6950	18 Fr				•		•	•		•	•	•	10
6070	14 Fr					•	•	•		•	•	•	10
6153	16 Fr					•	•	•		•	•	•	10
6154	18 Fr					•	•	•		•	•	•	10
6140	16 Fr					•	•	•	•	•	•	•	10
6210	18 Fr					•	•	•	•	•	•	•	10
<b>Coudé Tip</b>													
PP16CSD	16 Fr	•	•	•			•		•		•	•	10
PP18CSD	18 Fr	•	•	•			•		•		•	•	10

## Dover™ Latex Foley Catheters



- 16" length
- Coatings ease insertion to help prevent urethra trauma

### Hydrogel Coating

Code	Size (Fr)	Coating	Ship Carton
<b>5 cc, 2-way</b>			
1612	12 Fr	Hydrogel	12
1614	14 Fr	Hydrogel	12
1616	16 Fr	Hydrogel	12
1619	18 Fr	Hydrogel	12
1621	20 Fr	Hydrogel	12
1622	22 Fr	Hydrogel	12
1624	24 Fr	Hydrogel	12
1626	26 Fr	Hydrogel	12
1628	28 Fr	Hydrogel	12
1631	30 Fr	Hydrogel	12

### 30 cc, 2-way

1416	16 Fr	Hydrogel	12
1419	18 Fr	Hydrogel	12
1420	20 Fr	Hydrogel	12
1422	22 Fr	Hydrogel	12
1424	24 Fr	Hydrogel	12
1426	26 Fr	Hydrogel	12
1427	28 Fr	Hydrogel	12
1431	30 Fr	Hydrogel	12

### 5 cc, 3-way

2718	18 Fr	Hydrogel	12
2720	20 Fr	Hydrogel	12
2722	22 Fr	Hydrogel	12
2723	24 Fr	Hydrogel	12

### 30 cc, 3-way

2819	18 Fr	Hydrogel	12
2821	20 Fr	Hydrogel	12
2822	22 Fr	Hydrogel	12
2824	24 Fr	Hydrogel	12
2826	26 Fr	Hydrogel	12

### Catheter Kit

Includes: latex Foley catheter with lubricating jelly (5 g) and 10 cc prefilled water syringe

Code	Size (Fr)	Coating	Ship Case
1814	14 Fr	Hydrogel	12
1817	16 Fr	Hydrogel	12
1818	18 Fr	Hydrogel	12

### Silicone Elastomer Coating

Code	Size (Fr)	Coating	Ship Carton
<b>3 cc, 2-way 12" (Pediatric)</b>			
8887654088	8 Fr	Silicone Elastomer	10
8887654104	10 Fr	Silicone Elastomer	10
<b>5 cc, 2-way</b>			
402712	12 Fr	Silicone Elastomer	10
402714	14 Fr	Silicone Elastomer	10
402716	16 Fr	Silicone Elastomer	10
402718	18 Fr	Silicone Elastomer	10
402720	20 Fr	Silicone Elastomer	10
402722	22 Fr	Silicone Elastomer	10
402724	24 Fr	Silicone Elastomer	10
402726	26 Fr	Silicone Elastomer	10
402728	28 Fr	Silicone Elastomer	10
402730	30 Fr	Silicone Elastomer	10

### 30 cc, 2-way

403712	12 Fr	Silicone Elastomer	10
403714	14 Fr	Silicone Elastomer	10
403716	16 Fr	Silicone Elastomer	10
403718	18 Fr	Silicone Elastomer	10
403720	20 Fr	Silicone Elastomer	10
403722	22 Fr	Silicone Elastomer	10
403724	24 Fr	Silicone Elastomer	10
403726	26 Fr	Silicone Elastomer	10
403728	28 Fr	Silicone Elastomer	10
403730	30 Fr	Silicone Elastomer	10

### 5 cc, 3-way

8887688169	16 Fr	Silicone Elastomer	10
8887688185	18 Fr	Silicone Elastomer	10
8887688201	20 Fr	Silicone Elastomer	10
8887688227	22 Fr	Silicone Elastomer	10
8887688243	24 Fr	Silicone Elastomer	10
8887688268	26 Fr	Silicone Elastomer	10

### 30 cc, 3-way

8887689167	16 Fr	Silicone Elastomer	10
8887689183	18 Fr	Silicone Elastomer	10
8887689209	20 Fr	Silicone Elastomer	10
8887689225	22 Fr	Silicone Elastomer	10
8887689241	24 Fr	Silicone Elastomer	10
8887689266	26 Fr	Silicone Elastomer	10

## Dover™ Latex Foley Catheters (continued)

### Silicone Oil Coating

Code	Size (Fr)	Coating	Ship Carton
<b>5 cc, 2-way</b>			
8887623125	12 Fr	Silicone oil	10
3558	14 Fr	Silicone oil	10
3560	16 Fr	Silicone oil	10
3563	18 Fr	Silicone oil	10
3565	20 Fr	Silicone oil	10
3567	22 Fr	Silicone oil	10
3570	24 Fr	Silicone oil	10
3573	26 Fr	Silicone oil	10
3575	28 Fr	Silicone oil	10
3577	30 Fr	Silicone oil	10

### 30 cc, 2-way,

8887624149	14 Fr	Silicone oil	10
3601	16 Fr	Silicone oil	10
3607	18 Fr	Silicone oil	10
3611	20 Fr	Silicone oil	10
3614	22 Fr	Silicone oil	10
3618	24 Fr	Silicone oil	10
3619	26 Fr	Silicone oil	10
3623	28 Fr	Silicone oil	10
3631	30 Fr	Silicone oil	10

### Specialty Catheters

Code	Size (Fr)	Coating	Ship Carton
<b>5 cc, 2-way, Coudé tip</b>			
1612C	12 Fr	Hydrogel	12
1614C	14 Fr	Hydrogel	12
1616C	16 Fr	Hydrogel	12
1618C	18 Fr	Hydrogel	12
1620C	20 Fr	Hydrogel	12
1622C	22 Fr	Hydrogel	12
1624C	24 Fr	Hydrogel	12



### 5 cc, 2-way, Red Latex, Radiopaque, Coudé tip

1512C	12 Fr	Hydrogel	12
1514C	14 Fr	Hydrogel	12
1516C	16 Fr	Hydrogel	12
1518C	18 Fr	Hydrogel	12
1520C	20 Fr	Hydrogel	12
1522C	22 Fr	Hydrogel	12
1524C	24 Fr	Hydrogel	12

### 30 cc, 2-way, Coudé tip

1416C	16 Fr	Hydrogel	12
1418C	18 Fr	Hydrogel	12
1420C	20 Fr	Hydrogel	12
1422C	22 Fr	Hydrogel	12
1424C	24 Fr	Hydrogel	12

### T.U.R.P. 30 – 50 cc, up to 90 cc, 2-way, 6 Staggered Eyes, Nylon Reinforced Catheter Wall

8887648221	22 Fr	Silicone Elastomer	10
8887648247	24 Fr	Silicone Elastomer	10



### T.U.R.P. 30 – 50 cc, up to 90 cc, 3-way, 4 opposing eyes, Nylon Reinforced Catheter Wall

8887649229	22 Fr	Silicone Elastomer	10
8887649245	24 Fr	Silicone Elastomer	10

## Dover™ Latex Urine Meter Foley Trays



- Standard Foley prep tray. Sterile tray contains: CSR wrap, fenestrated drape, underpad, blue nitrile exam gloves, forceps, prepping balls, 10 cc prefilled sterile water inflation syringe, unique bellows lubricant syringe (10 g), povidone-iodine solution, specimen container with label, Foley reminder label and patient education pamphlet
- Latex Foley catheter with hydrogel coating
- 400 mL urine meter capacity
- Extra long 60" tubing



Code	Size (Fr)	Pre-attached Securement Device	Clampable Tubing with Top Vent	Sample Port		Tube Management with Sheet Clip	Sheet Clip	Spring Meter Port	Enclosed Drain Spout	Self-leveling Bed Hook	Ship Case
				Luer Lock	Needleless						
P4P14ULD	14 Fr	•	•	•		•		•	•	•	10
P4P16ULD	16 Fr	•	•	•		•		•	•	•	10
P4P18ULD	18 Fr	•	•	•		•		•	•	•	10
7014LL	14 Fr		•	•		•		•	•	•	10
7016LL	16 Fr		•	•		•		•	•	•	10
7018LL	18 Fr		•	•		•		•	•	•	10
2014	14 Fr				•		•	•	•	•	10
2016	16 Fr				•		•	•	•	•	10
2018	18 Fr				•		•	•	•	•	10
<b>Metal-free Push Pull Valve</b>											
701672	16 Fr				•		•		•	•	10
701872	18 Fr				•		•		•	•	10
<b>Coudé Tip</b>											
7016C	16 Fr		•	•		•		•	•	•	10

### Dover™ Latex Foley Trays



PP16ULD

- Standard Foley prep tray. Sterile tray contains: CSR wrap, fenestrated drape, underpad, blue nitrile exam gloves, forceps, prepping balls, 10 cc prefilled sterile water inflation syringe, unique bellows lubricant syringe (10 g), povidone-iodine solution, specimen container with label, Foley reminder label and patient education pamphlet

- Latex Foley catheter with Ultramer™ hydrogel coating
- 2000 mL capacity
- 48" tubing



Code	Size (Fr)	Pre-attached Securement Device	Top Vent	Sample Port			Tube Management with Sheet Clip	Sheet Clip	Drain		Anti-Reflex Device	Self-leveling Bed Hook	Ship Case
				Luer Lock	Needle-less	Needle			Enclosed Spout	Tube			
PP14ULD	14 Fr	•	•	•		•		•		•	•	10	
PP16ULD	16 Fr	•	•	•		•		•		•	•	10	
PP18ULD	18 Fr	•	•	•		•		•		•	•	10	
6164LL	14 Fr		•	•		•		•		•	•	10	
6165LL	16 Fr		•	•		•		•		•	•	10	
6166LL	18 Fr		•	•		•		•		•	•	10	
8944	14 Fr		•		•		•	•		•	•	10	
8946	16 Fr		•		•		•	•		•	•	10	
8948	18 Fr		•		•		•	•		•	•	10	
6162	14 Fr		•			•		•	•	•	•	10	
6155	16 Fr		•			•		•	•	•	•	10	
6156	18 Fr		•			•		•	•	•	•	10	
6014	16 Fr		•			•		•	•	•	•	10	
6080	18 Fr		•			•		•	•	•	•	10	
6946	16 Fr		•			•		•	•	•	•	10	
6948	18 Fr		•			•		•	•	•	•	10	
<b>Silicone Oil Coating, 40" Tubing</b>													
3716	16 Fr		•			•		•	•	•	•	10	
3718	18 Fr		•			•		•	•	•	•	10	
<b>Silicone Elastomer with Teardrop Shape and Hook &amp; Loop</b>													
407716	16 Fr					•		•	•	•		10	
407718	18 Fr					•		•	•	•		10	
<b>Coudé Tip</b>													
6165C	16 Fr		•		•		•	•		•	•	10	

### Insertion Trays

#### Dover™ Silver Insertion Trays

##### Dover™ 100% Silicone Foley Catheter with Silver Hydrogel Coating

Sterile tray contains: CSR wrap, underpad, fenestrated drape, PVP solution, prepping balls, forceps, lubricating jelly (5 g), 10cc prefilled water syringe, specimen container with label, graduated basin (1000 mL), blue nitrile exam gloves

Code	Size (Fr)	Ship Case
2116IC	16 Fr	20
2118IC	18 Fr	20



#### Dover™ Silicone Insertion Trays

##### Dover™ 100% Silicone Foley Catheter

Sterile tray contains: CSR wrap, underpad, fenestrated drape, PVP solution, prepping balls, forceps, lubricating jelly (5 g), 10cc prefilled water syringe, specimen container with label, graduated basin (1000 mL), blue nitrile exam gloves

Code	Size (Fr)	Ship Case
8887600776	16 Fr	20
8887600784	18 Fr	20



#### Dover™ Latex Insertion Trays

##### Dover™ Latex Foley Catheter

Sterile tray contains: CSR wrap, underpad, fenestrated drape, PVP solution, prepping balls, forceps, lubricating jelly (5 g), 10cc prefilled water syringe, specimen container with label, graduated basin (1000 mL), blue nitrile exam gloves

Code	Size (Fr)	Ship Case
2101	16 Fr	20
3044	18 Fr	20



#### Universal Trays

All trays feature (unless indicated otherwise): underpad, blue nitrile exam gloves, fenestrated drape, lubricating jelly (5 g) and graduated basin

Code	Antiseptic	Syringe	Cotton balls	Forceps	Specimen Container	Lid	Graduated Basin	Ship Case
5029	PVP Solution	10 cc	•	•	•		1000 mL	20
7100	BZK Swabs	10 cc			•		1000 mL	20
5027	PVP Solution	30 cc	•	•	•		1000 mL	20
76000	PVP Swabs	30 cc				•	1000 mL	20
76010	PVP Swabs	10 cc				•	1000 mL	20
76012	PVP Swabs	10 cc			•	•	1200 mL	20
76020	BZK Swabs	30 cc				•	1000 mL	20
76030	BZK Swabs	10 cc				•	1000 mL	20
8887601220	PVP Solution	10 cc	•	•	•		1000 mL	20



5029

### Collection



7000LL

#### Urine Meter

Designed for precise measurement when output is critically important. Capacity ranges from 150 mL to 400 mL.

#### Next Generation:

- Lowered & centered entry point
- Horizontal face plate with self leveling chamber
- Top vent promotes optimal flow within the collection system
- Luer Lock sample port and meter port helps prevent needle stick injury
- Tube management & sheet clip improve flow dynamics

### Urine Meters



Code	Meter Capacity	Top Vent	Sample Port			Tube Management with Sheet Clip	Spring Meter Port	Drain		Self-leveling Bed Hook	Ship Case
			Luer Lock	Needleless	Needle			Enclosed Spout	Drain Tube		
<b>Extra-long 60" tubing</b>											
6600	150 mL		•						•	•	10
6770	200 mL		•						•	•	10
8000	200 mL				•				•	•	10
7000	400 mL	•		•		•	•	•		•	10
7000LL	400 mL	•	•			•	•	•		•	10
2000	400 mL			•		•	•	•		•	10

### Add-A-Foley Urine Meter Trays

Includes Standard Foley Prep Tray



Code	Meter Capacity	Clampable Tubing with Top Vent	Sample Port			Drain		Ship Case
			Luer Lock	Needleless	Needle	Enclosed Spout	Drain Tube	
7020LL	400 mL	•	•			•		10
2020	400 mL			•		•		10
8324	200 mL				•		•	10

With catheter securement device

P4PAD	400 mL	•	•			•		10
-------	--------	---	---	--	--	---	--	----

### Collection



6308LL

#### Drainage Bag

Capacity offerings in 1000 mL, 2000 mL and 4000 mL.

#### Next Generation:

- Updated enclosed drain spout on a new 45° pedestal
- Pivoting hanger hooks
- Top vent promotes optimal flow within the collection system
- Luer Lock sample port helps in preventing needle stick injury
- Tube management & sheet clip improve flow dynamics

### Drainage Bags



Code	Size	Sample Port			Tube Management with Sheet Clip	Sheet Clip	Drain		Anti-reflux Device	Hanger			Ship Case
		Luer Lock	Needleless	Needle			Enclosed Spout	Drain Tube		Bed Hook	String	Hook & Loop	
<b>48" tubing</b>													
6251	4000 mL			•		•		2	2		•		10
6261	4000 mL			•		•		•	•		•		15
6308LL	2000 mL	•			•		•		•	•			20
6308	2000 mL		•			•	•		•	•			20
6208	2000 mL			•		•	•		•	•			20
6209	2000 mL			•		•		•	•	•			20
6300 *	2000 mL			•		•		•	•	•			20
8206	2000 mL		•			•	•		•	•			20
6206	2000 mL		•			•	•		•	•			20
6310	2000 mL			•		•		•	•			•	20
3057	2000 mL			•		•		•	•	•			20
<b>45" tubing</b>													
6510	1000 mL	•				•		•	•	•			20
<b>40" tubing</b>													
3512V *	2000 mL			•				•	•			•	20
3512 *	2000 mL			•				•	•	•			20
3502 *	2000 mL			•				•		•			20

\*Indicates poly bag packaging – interior of drainage bag and fluid pathway are guaranteed sterile



Collection



**Add-A-Foley 2000 mL Drainage Bag Trays**

Standard Foley prep tray (100% latex-free) Sterile tray contains: CSR wrap, fenestrated drape, underpad, blue nitrile exam gloves, forceps, prepping balls, unique bellows lubricant syringe (10 g), povidone-iodine solution, specimen container and Foley reminder label



Code	Syringe	Sample Port			Drain		Ship Case
		Luer Lock	Needleless	Needle	Enclosed Spout	Drain Tube	
<b>48" Tubing</b>							
6176LL	10 cc	•			•		10
8256	10 cc		•		•		10
6175	10 cc			•		•	10
6255	10 cc			•		•	10
6256	10 cc			•	•		10
<b>With catheter securement device</b>							
PPAD	10 cc	•			•		10
<b>With push-pull valve</b>							
407020	10 cc		•			•	10
<b>40" Tubing</b>							
3515	10 cc			•		•	10
3531	30 cc			•		•	10
3532	30 cc			•		•	10


**Miscellaneous Accessories**



Code	Description	Ship Quantity
54456P	Dover™ Catheter Securement	40/Case
73342	Dover™ Pediatric Catheter Securement	50/Case
1600	Catheter Plug	50/Case
95516	Sterile Catheter Care Kit	50/Case
8887600149	Foley Leg Strap	10/Case
1106	Lubricating Jelly Bellows Syringe (10 g)	280/Case

**Dover™ Urethral Catheters**

**100% Latex Free Urethral Catheters (PVC)**



Code	Size (Fr)	Ship Case
------	-----------	-----------

**Rob-Nel catheters, 16" length, smooth rounded tip**

8888492017	8 Fr	100
8888492025	10 Fr	100
8888492033	12 Fr	100
8888492041	14 Fr	100
8888492058	16 Fr	100
8888492066	18 Fr	100
8888492120	20 Fr	100
8888492122	22 Fr	100
8888492124	24 Fr	100
8888492126	26 Fr	100
8888492128	28 Fr	100
8888492130	30 Fr	100

**Robinson catheters, 16" length, smooth rounded tip, clear vinyl**

400608	8 Fr	100
400610	10 Fr	100
400612	12 Fr	100
400614	14 Fr	100
400616	16 Fr	100
400618	18 Fr	100

**Robinson catheter with 5 gram lubricating jelly packet**

2540	14 Fr	100
------	-------	-----

**Female catheter, 6.5" length**

8887660903	14 Fr	50
------------	-------	----

**Pediatric urethral catheters, 11" length, hydrogel coating and Luer lock adapter**

Code	Size (Fr)	Ship Carton
6006	6 Fr (Ped.)	12

## Dover™ Urethral Catheters (continued)

### Latex Urethral Catheters

Code	Size (Fr)	Ship Carton
<b>Red Rubber Coudé, 16" length, hydrophilic coating</b>		
8403	14 Fr	12
8404	16 Fr	12



Code	Size (Fr)	Ship Carton
<b>Red Rubber, 12" length, hydrophilic coating, smooth rounded tip</b>		
8408	8 Fr	12
8410	10 Fr	12
8412	12 Fr	12



Code	Size (Fr)	Ship Carton
<b>Red Rubber, 16" length, hydrophilic coating, smooth rounded tip</b>		
8414	14 Fr	12
8416	16 Fr	12
8418	18 Fr	12
8420	20 Fr	12
8422	22 Fr	12
8424	24 Fr	12

Code	Size (Fr)	Ship Case
<b>Red Rubber Robinson, 12" length, radiopaque, smooth rounded tip</b>		
8887660085	8 Fr	100
8887660101	10 Fr	100

Code	Size (Fr)	Ship Case
<b>Red Rubber Robinson, 16" length, radiopaque, smooth rounded tip</b>		
8887660127	12 Fr	100
8887660143	14 Fr	100
8887660168	16 Fr	100
8887660184	18 Fr	100
8887660218	20 Fr	100
8887660221	22 Fr	100
8887660242	24 Fr	100
8887660263	26 Fr	100
8887660284	28 Fr	100
8887660305	30 Fr	100



## Dover™ Urethral Trays

### Open Trays

All trays feature: underpad, blue nitrile exam gloves, lubricating jelly, specimen container with label, graduated basin



Code	Size (Fr)	Graduated Basin	Material	Antiseptic	CSR Wrap	Cotton Balls	Ship Case
------	-----------	-----------------	----------	------------	----------	--------------	-----------

#### Hydrogel coated catheters

3217	14 Fr	1000 mL	Red Rubber Latex	PVP Solution	•	•	20
------	-------	---------	------------------	--------------	---	---	----

#### Urethral Catheter

75039	14 Fr.	1000 mL	Vinyl, Latex-free		•	•	20
3141	14 Fr	1000 mL	Vinyl, Latex-free	PVP Solution	•	•	20
3215	16 Fr	1000 mL	Red Rubber Latex	PVP Swab	•		20
406622	14 Fr	1200 mL	Vinyl	BZK Swabs	•		20
8887600057	14 Fr	1000 mL	Red Rubber Latex	PVP Solution	•	•	20
75000	14 Fr	1200 mL	Vinyl	BZK Swabs			20
75005	14 Fr	1200 mL	Red Rubber Latex	BZK Swabs			20
75030	14 Fr	1200 mL	Vinyl	PVP Swabs			20
75035	14 Fr	1200 mL	Red Rubber Latex	PVP Swabs			20
3143	14 Fr	1000 mL	Vinyl	PVP Solution	•	•	20
8887600073	14 Fr	1000 mL	Vinyl	PVP Solution	•	•	20

### Add-A-Cath Tray

Code	Graduated Basin	Antiseptic	Exam Gloves	Lubricating Jelly	Specimen Container	Underpad/ Drape	CSR Wrap	Ship Case
------	-----------------	------------	-------------	-------------------	--------------------	-----------------	----------	-----------

#### Open Urethral Tray

3305	1000 mL	PVP Solution	•	•	•	•	•	20
3307	1000 mL	PVP Solution	•	•	•	•	•	20
75010	1200 mL	PVP Swab Sticks	•	•	•	•		20
75020	1200 mL	BZK Swab Sticks	•	•	•	•		20

### Closed Trays

Trays include (unless indicated otherwise): CSR wrap, underpad, blue nitrile exam gloves, antiseptic, lubricating jelly, specimen container, 1500 mL collection bag



Code	Size (Fr)	Material	Antiseptic	Ship Case
------	-----------	----------	------------	-----------

#### Hydrogel Coated Catheters

3410	14 Fr	Red Rubber Latex	PVP Solution	20
3408	16 Fr	Red Rubber Latex	PVP Solution	20

#### Clear Vinyl Catheters

3450	14 Fr	Vinyl	PVP Solution	20
3450LF	14 Fr	Vinyl	PVP Solution	20
2480	14 Fr	Vinyl	BZK Swabs	20

#### Urethral Catheter with 1000 mL bag only

3170	14 Fr	Vinyl	None	50
------	-------	-------	------	----

Code	Size (Fr)	Material	Antiseptic	Ship Carton
------	-----------	----------	------------	-------------

#### Female Catheterization Kit with 15 mL Specimen Tube and Flip Cap

3411	8 Fr	Vinyl	PVP Swabs	25
------	------	-------	-----------	----

## Dover™ Patient Care Products

### Male External Catheters (Latex)



Code	Size	Attachment	Ship Case
<b>Uri-Drain Male External Catheter</b> Latex rolled sheath, molded one piece, individually packaged			
8884732300	1.3"	Foam Strap	144
8884734600	1.1"	Foam Strap	144
8884732500	1.3"	Two Sided Adhesive	144
8884732000	1.1"	Two Sided Adhesive	144
<b>Texas Catheter</b> Soft, lightweight latex with hard plastic tip			
8884730300	1.3"	Foam Strap, Bulk	144
8884731300	1.3"	Foam Strap, 12/box	144
<b>Adhesive Strap</b>			
8884730600	5.25" x 1"	Foam Strap	500

### Leg Bags



Leg bag features: latex free leg strap, easy drain valve or tubing and anti reflux

Code	Size			Latex	Extension Tubing	Ship Case
	Small 9 oz.	Medium 17 oz. / 500 mL	Large 25 oz. / 730 mL			
<b>Disposable</b>						
3433		•			•	50
145514		•				50
145516			•			50
8887601121		•				20
8887601139			•			20
<b>Reusable</b>						
8884733100	•			•		24
8884733200		•		•		24
8884733900			•	•		24

Code	Description	Ship Case
8884731900	18" Urinary Extension Tube	24

## Dover™ Patient Care Products

### Irrigation Trays



Code	1200 cc Tray	500 cc Container	Antiseptic	CSR Wrap	Sterile Saline	Lid	Tip Protector	Ship Case
<b>60 cc Bulb Syringe Tray</b>								
67800	•	•	•			•	•	20
67803	•	•	•	•			•	20
<b>60 cc Piston Syringe Tray</b>								
68800	•	•	•			•	•	20
68803	•	•	•	•			•	20
68820	•		•		•	•	•	20
3685	•	•				•	•	20

### Syringes Only, 60 cc

Code	Description	Ship Case
67000	Sterile Bulb Syringe	50
68000	Sterile Piston Syringe	50

### Midstream Catch Kits



Code	Capacity	Description	Patient Label	BZK Wipes	Castile Soap	Pneumatic Capable	Funnel	Screw Cap	Sterile	Ship Case
22000	4.5 oz	Funnel Top	•	•			•		•	24
22015	4.5 oz	Funnel Top	•		•		•		•	24
24000	4.5 oz	Wide mouth container	•	•				•	•	50
24015	4.5 oz	Wide mouth container	•		•			•	•	50
24035	7 oz	Wide mouth container	•		•			•	•	50
25000	4.5 oz	Protective cover	•	•				•	•	100
25015	4.5 oz	Protective cover	•		•			•	•	100
5205	4.5 oz	Specimen container	•		•				•	48
5204	4.5 oz	Specimen container	•	•					•	48
2090SA	4 oz	Collection kit	•		•	•	•		•	100
4090SA	4 oz	Collection kit			•	•			•	100
2040SA	2 oz	Collection kit	•		•	•	•			100

Code	Description	Ship Case
2110SA	Calculi Strainer only	100

## Dover™ Ancillary Products

### Specimen Containers



Code	Capacity	Lid	Ship Case
<b>Measuring System</b>			
4014	26 oz		100
2400SA	700 cc	•	60
2500SA	800 cc		100



Code	Capacity	Sterile	Ship Case
<b>Pediatric Urine Collector</b>			
145500	100 mL		500
145501	100 mL	•	300

### Patient Care Products (Latex Free)



Code	Size (Fr)	Description	Tube length	Sterile	Pre-lubricated	Castile Soap	Underpad	Ship Case
145524	24 Fr	Flatus Bag	19"	•	•			50
145617		Vaginal Irrigation (Douche) Bag	60"	•				50
145540	24 Fr	Enema Bag	60"	•		•	•	50
145541	24 Fr	Enema Bag	60"		•			50
145546	24 Fr	Enema Bucket 1400 cc	60"		•	•	•	50
155730	18 Fr	Rectal Tube	20"					50
155731	24 Fr	Rectal Tube	20"	•	•			50

## Index

Code	HCPCS	Latex	Page	Code	HCPCS	Latex	Page
1106	*	No	14	1628	A4338	Yes	7
1416	A4338	Yes	7	1631	A4338	Yes	7
1416C	A4340	Yes	8	1814	A4338	Yes	7
1418C	A4340	Yes	8	1817	A4338	Yes	7
1419	A4338	Yes	7	1818	A4338	Yes	7
1420	A4338	Yes	7	2000	A4315	No	12
1420C	A4340	Yes	8	2004	*	No	5
1422	A4338	Yes	7	2004LL	A4315	No	5
1422C	A4340	Yes	8	2006	A4315	No	5
1424	A4338	Yes	7	2008	A4315	No	5
1424C	A4340	Yes	8	2014	A4314	Yes	9
1426	A43386	Yes	7	2016	A4314	Yes	9
1427	A4338	Yes	7	2018	A4314	Yes	9
1431	A4338	Yes	7	2020	A4354	No	12
145500	A9270	No	20	2040SA	A9270	No	19
145501	A9270	No	20	20512C	A4340	No	4
145514	A4458	No	18	20512ICC	*	No	1
145516	A4358	No	18	20514C	A4340	No	4
145524	A9270	No	20	20514ICC	*	No	1
145540	A4458	No	20	20516C	A4340	No	4
145541	A4358	No	20	20516ICC	*	No	1
145546	A4358	No	20	20518C	A4340	No	4
145617	A9270	Yes	20	20518ICC	*	No	1
1512C	A4340	Yes	8	20520C	A4340	No	4
1514C	A4340	Yes	8	20520ICC	*	No	1
1516C	A4340	Yes	8	20522C	A4340	No	4
1518C	A4340	Yes	8	20522ICC	*	No	1
1520C	A4340	Yes	8	20524C	A4340	No	4
1522C	A4340	Yes	8	2090SA	A9270	No	19
1524C	A4340	Yes	8	2101	A4312	Yes	11
155730	A9270	No	20	2110SA	A9270	No	19
155731	A9270	No	20	2116IC	A4312	No	11
1600	A9270	No	14	2118IC	A4312	No	11
1612	A4338	Yes	7	22000	A9270	No	19
1612C	A4340	Yes	8	22015	A9270	No	19
1614	A4338	Yes	7	23018C	A4340	No	4
1614C	A4340	Yes	8	23018ICC	*	No	1
1616	A4338	Yes	7	23020C	A4340	No	4
1616C	A4340	Yes	8	23020ICC	*	No	1
1618C	A4340	Yes	8	23022C	A4340	No	4
1619	A4338	Yes	7	23024C	A4340	No	4
1620C	A4340	Yes	8	23026C	A4340	No	4
1621	A4338	Yes	7	24000	A9270	No	19
1622	A4338	Yes	7	2400SA	A9270	No	20
1622C	A4340	Yes	8	24015	A9270	No	19
1624	A4338	Yes	7	24035	A9270	No	19
1624C	A4340	Yes	8	2480	A4353	No	17
1626	A4338	Yes	7	25000	A9270	No	19

\* HCPCS code pending verification by PDAC. Please call the Reimbursement Hot Line (800-346-7197 x6076) for update.

**Index**

Code	HCPCS	Latex	Page
2500SA	A9270	No	20
25015	A9270	No	19
2540	A4351; A4402	No	15
2718	A4346	Yes	7
2720	A4346	Yes	7
2722	A4346	Yes	7
2723	A4346	Yes	7
2819	A4346	Yes	7
2821	A4346	Yes	7
2822	A4346	Yes	7
2824	A4346	Yes	7
2826	A4346	Yes	7
3044	A4312	Yes	11
3057	A4357	No	13
3141	A4353	No	17
3143	A4353	No	17
3170	A4353	No	17
3215	A4353	Yes	17
3217	A4353	Yes	17
3305	A4310	No	17
3307	A4310	No	17
3408	A4353	Yes	17
3410	A4353	Yes	17
3411	A9270	No	17
3433	A4358	No	18
3450	A4353	No	17
3450LF	A4353	No	17
3502	A4357	No	13
3512	A4357	No	13
3512V	A4357	No	13
3515	A4354	Yes	14
3531	A4354	Yes	14
3532	A4354	No	14
3558	A4338	Yes	8
3560	A4338	Yes	8
3563	A4338	Yes	8
3565	A4338	Yes	8
3567	A4338	Yes	8
3570	A4338	Yes	8
3573	A4338	Yes	8
3575	A4338	Yes	8
3577	A4338	Yes	8
3601	A4338	Yes	8
3607	A4338	Yes	8
3611	A4338	Yes	8
3614	A4338	Yes	8
3618	A4338	Yes	8
3619	A4338	Yes	8

Code	HCPCS	Latex	Page
3623	A4338	Yes	8
3631	A4338	Yes	8
3685	A4320	No	19
3716	A4314	Yes	10
3718	A4314	Yes	10
400608	A4351	No	15
400610	A4351	No	15
400612	A4351	No	15
400614	A4351	No	15
400616	A4351	No	15
400618	A4351	No	15
4014	A9270	No	20
402712	A4338	Yes	7
402714	A4338	Yes	7
402716	A4338	Yes	7
402718	A4338	Yes	7
402720	A4338	Yes	7
402722	A4338	Yes	7
402724	A4338	Yes	7
402726	A4338	Yes	7
402728	A4338	Yes	7
402730	A4338	Yes	7
403712	A4346	Yes	7
403714	A4346	Yes	7
403716	A4346	Yes	7
403718	A4346	Yes	7
403720	A4346	Yes	7
403722	A4346	Yes	7
403724	A4346	Yes	7
403726	A4346	Yes	7
403728	A4346	Yes	7
403730	A4346	Yes	7
40516L	A4344	No	4
40518L	A4344	No	4
40520L	A4344	No	4
40522L	A4344	No	4
40524L	A4344	No	4
406622	A4353	No	17
407020	*	No	14
407716	*	Yes	10
407718	*	Yes	10
4090SA	A9270	No	19
5027	A4310	No	11
5029	A4310	No	11
50514T	A4340	No	4
50516ICT	A4344 & A9279	No	1
50516T	A4340	No	4
50518ICT	A4344 & A9279	No	1

**Index**

Code	HCPCS	Latex	Page
50518T	A4340	No	4
5204	A9270	No	19
5205	A9270	No	19
54456P	*	No	14
6006	A4351	No	15
6010	A4344	No	4
6014	A4314	Yes	10
605122IC	A4344	No	1
605148IC	A4344	No	1
605163IC	A4344	No	1
605189IC	A4344	No	1
605205IC	A4344	No	1
605221IC	A4344	No	1
605247IC	A4344	No	1
605262IC	A4344	No	1
6070	A4315	No	6
6080	A4314	Yes	10
6108	A4344	No	4
6140	A4354	No	6
6144ICLL	A4314	No	3
6144LL	A4315	No	6
6146ICLL	A4314	No	3
6146LL	A4315	No	6
6148ICLL	A4314	No	3
6148LL	A4315	No	6
6153	A4315	No	6
6154	A4315	No	6
6155	A4314	Yes	10
6156	A4314	Yes	10
6162	A4314	Yes	10
6164LL	A4314	Yes	10
6165C	A4314	Yes	10
6165LL	A4314	Yes	10
6166LL	A4314	Yes	10
6175	A4354	No	14
6176LL	A4354	No	14
6206	A4357	No	13
6208	A4357	No	13
6209	A4357	No	13
6210	A4354	No	6
6251	A4357	No	13
6255	A4354	No	14
6256	A4354	No	14
6261	A4357	No	13
6300	A4357	No	13
630187IC	A4344	No	1
630203IC	A4344	No	1
630229IC	A4344	No	1

Code	HCPCS	Latex	Page
630245IC	A4344	No	1
630260IC	A4344	No	1
6308	*	No	13
6308LL	A4357	No	13
6310	A4357	No	13
6510	A4357	No	13
6600	A4357	No	12
665183IC	A4346	No	1
665209IC	A4346	No	1
665225IC	A4346	No	1
665241IC	A4346	No	1
67000	A4322	No	19
6770	A4357	Yes	12
67800	A4320	No	19
67803	A4320	No	19
68000	A4322	No	19
68800	A4320	No	19
68803	A4320	No	19
68820	A4320	No	19
6946	A4314	Yes	10
6947LF	A4315	No	6
6948	A4314	Yes	10
6950	A4315	No	6
7000	*	No	12
7000LL	A4315	No	12
7004ICLL	A4314	No	2
7006ICLL	A4314	No	2
7006LL	A4315	No	5
7008ICLL	A4314	No	2
7008LL	A4315	No	5
7014LL	A4314	Yes	9
701672	*	Yes	9
7016C	A4314	Yes	9
7016LL	A4314	Yes	9
701872	*	Yes	9
7018LL	A4314	Yes	9
7020LL	A4354	Yes	12
7100	A4310	No	11
73342	*	No	14
75000	A4353	No	17
75005	A4353	Yes	17
75010	A4310	No	17
75020	A4310	No	17
75030	A4353	No	17
75035	A4353	Yes	17
75039	*	No	17
76000	A4310	No	11
76010	A4310	No	11

\* HCPCS code pending verification by PDAC. Please call the Reimbursement Hot Line (800-346-7197 x6076) for update.

**Index**

Code	HCPCS	Latex	Page
76012	A4310	No	11
76020	A4310	No	11
76030	A4310	No	11
8000	A4357	No	12
8122	A4315	No	5
8206	A4357	No	13
8223	A4315	No	5
8256	A4354	No	14
8324	A4354	No	12
8403	A4352	Yes	16
8404	A4352	Yes	16
8408	A4351	Yes	16
8410	A4351	Yes	16
8412	A4351	Yes	16
8414	A4351	Yes	16
8416	A4351	Yes	16
8418	A4351	Yes	16
8420	A4351	Yes	16
8422	A4351	Yes	16
8424	A4351	Yes	16
8884730300	A4349	Yes	18
8884730600	A4333	No	18
8884731300	A4349	Yes	18
8884731900	A4331	No	18
8884732000	A4349	Yes	18
8884732300	A4349	Yes	18
8884732500	A4349	Yes	18
8884733100	A4358	Yes	18
8884733200	A4358	Yes	18
8884733900	A4358	Yes	18
8884734600	A4349	Yes	18
8887600057	A4353	Yes	17
8887600073	A4353	No	17
8887600149	A4334	No	14
8887600776	A4312	No	11
8887600784	A4312	No	11
8887601121	A4358	No	18
8887601139	A4358	No	18
8887601220	A4310	No	11
8887603085	A4344	No	4
8887603101	A4344	No	4
8887605122	A4344	No	4
8887605148	A4344	No	4
8887605163	A4344	No	4
8887605189	A4344	No	4
8887605205	A4344	No	4
8887605221	A4344	No	4
8887605247	A4344	No	4

Code	HCPCS	Latex	Page
8887605262	A4344	No	4
8887605282	A4344	No	4
8887605303	A4344	No	4
8887606146	A4344	No	4
8887606161	A4344	No	4
8887606187	A4344	No	4
8887606203	A4344	No	4
8887623125	A4338	Yes	8
8887624149	A4338	Yes	8
8887630142	A4344	No	4
8887630163	A4344	No	4
8887630187	A4344	No	4
8887630203	A4344	No	4
8887630229	A4344	No	4
8887630245	A4344	No	4
8887630260	A4344	No	4
8887630284	A4344	No	4
8887630305	A4344	No	4
8887648221	A4338	Yes	8
8887648247	A4338	Yes	8
8887649229	A4346	Yes	8
8887649245	A4346	Yes	8
8887654088	A4338	Yes	7
8887654104	A4338	Yes	7
8887660085	A4351	Yes	16
8887660101	A4351	Yes	16
8887660127	A4351	Yes	16
8887660143	A4351	Yes	16
8887660168	A4351	Yes	16
8887660184	A4351	Yes	16
8887660218	A4351	Yes	16
8887660221	A4351	Yes	16
8887660242	A4351	Yes	16
8887660263	A4351	Yes	16
8887660284	A4351	Yes	16
8887660305	A4351	Yes	16
8887660903	A4351	No	15
8887664161	A4346	No	4
8887664182	A4346	No	4
8887664203	A4346	No	4
8887664224	A4346	No	4
8887664245	A4346	No	4
8887664266	A4346	No	4
8887664287	A4346	No	4
8887665161	*	No	4
8887665183	A4346	No	4
8887665209	A4346	No	4
8887665225	A4346	No	4

**Index**

Code	HCPCS	Latex	Page
8887665241	A4346	No	4
8887665262	*	No	4
8887665283	*	No	4
8887688169	A4346	Yes	7
8887688185	A4346	Yes	7
8887688201	A4346	Yes	7
8887688227	A4346	Yes	7
8887688243	A4346	Yes	7
8887688268	A4346	Yes	7
8887689167	A4346	Yes	7
8887689183	A4346	Yes	7
8887689209	A4346	Yes	7
8887689225	A4346	Yes	7
8887689241	A4346	Yes	7
8887689266	A4346	Yes	7
8888492017	A4351	No	15
8888492025	A4351	No	15
8888492033	A4351	No	15
8888492041	A4351	No	15
8888492058	A4351	No	15
8888492066	A4351	No	15
8888492120	A4351	No	15
8888492122	A4351	No	15
8888492124	A4351	No	15
8888492126	A4351	No	15
8888492128	A4351	No	15
8888492130	A4351	No	15
8944	A4314	Yes	10
8946	A4314	Yes	10
8948	A4314	Yes	10
8949	A4315	No	6
8950	A4315	No	6
95516	A4335	No	14
P4P14SD	A4315 & A4333	No	5
P4P14TSD	*	No	5
P4P14ULD	*	Yes	9

Code	HCPCS	Latex	Page
P4P14XSD	A4315 & A4333	No	2
P4P16CSD	*	No	5
P4P16CSDK	*	No	5
P4P16SD	A4315 & A4333	No	5
P4P16TS	A4315	No	5
P4P16TSD	*	No	5
P4P16ULD	*	Yes	9
P4P16XCSDK	*	No	2
P4P16XSD	A4315 & A4333	No	2
P4P16XTS	A4315	No	2
P4P16XTSD	*	No	2
P4P18CSD	*	No	5
P4P18CSDK	*	No	5
P4P18SD	A4315 & A4333	No	5
P4P18TS	A4315	No	5
P4P18TSD	*	No	5
P4P18ULD	*	Yes	9
P4P18XCSDK	*	No	2
P4P18XSD	A4315 & A4333	No	2
P4P18XTS	A4315	No	2
P4P18XTSD	*	No	2
P4PAD	A4354 & A4333	No	12
PP14SD	A4315 & A4333	No	6
PP14ULD	*	Yes	10
PP14XSD	A4315 & A4333	No	3
PP16CSD	*	No	6
PP16SD	A4315 & A4333	No	6
PP16ULD	*	Yes	10
PP16XCSD	*	No	3
PP16XSD	A4315 & A4333	No	3
PP18CSD	*	No	6
PP18SD	A4315 & A4333	No	6
PP18ULD	*	Yes	10
PP18XCSD	*	No	3
PP18XSD	A4315 & A4333	No	3
PPAD	A4354 & A4333	No	14

\* HCPCS code pending verification by PDAC. Please call the Reimbursement Hot Line (800-346-7197 x6076) for update.



Visit us at [www.covidien.com/NursingCare](http://www.covidien.com/NursingCare)

\* HCPCS code pending verification by PDAC. Please call the Reimbursement Hot Line (800-346-7197 x6076) for update.

The Centers for Medicare and Medicaid Services was the primary source of any HCPCS information contained in this document. Covidien does not represent or guarantee that any HCPCS information contained in this document is complete, accurate, or applicable to any particular patient or third-party payer and disclaims all liability for any consequence resulting from any reliance on such information. The final billing-related decisions for any product must be made by the health care provider considering the medical necessity of the product furnished as well as the requirements of third-party payers and any local, state, or federal laws and regulations that apply to the product furnished. Covidien is providing this information in an educational capacity with the understanding that Covidien is not engaged in rendering legal, accounting, or other professional services. Covidien encourages all health care providers to consult with their own advisors regarding coding and payment.

The products shown in this catalog may not be licensed or registered in all parts of the world.  
Please contact your sales representative for availability.

Cover photo credit: Thinkstock

COVIDIEN, COVIDIEN with logo and Covidien logo are U.S. and internationally registered trademarks of Covidien AG.

Other brands are trademarks of a Covidien company. ©2013 Covidien.

H6200URO –10M – 0813

