



DISPOSABLE EXAM PRODUCTS

Coverings for Patients, Equipment and Staff

Patient Comfort. Staff Convenience.



Selection. Convenience. Value.


Whether you're looking to guard against stains and spills, reduce cross-contamination or to provide extra comfort for patients, Medline has you covered. In our full selection, you'll find all the sizes, features and materials you need. We even have environmentally friendly options – look for our greensmart™ logo for recycled exam table paper.

With Medline disposable linens, you'll streamline processes and improve efficiency. There's no reprocessing, restocking or maintenance – just use them once and throw them away. They help eliminate any concerns about inventory control or waste from throwing away soiled reusable linens.

Plus, you can count on the highest-quality American craftsmanship with our tissue and tissue/poly products. Look throughout this catalog for products that are made in the USA. We manufacture and distribute all our disposable linens, so you'll cut out the middleman and save.

Table of Contents

Exam Table Paper	3
Apparel	5
Sheets and Blankets	9
Cubicle Curtains	13
Professional Towels	14
Glossary of Terms	15



Exam Table Paper

Greensmart products go through a rigorous review and certification for sustainability throughout the entire product lifecycle, from production and transport through use and disposal.

Offer a comfortable surface for patients while protecting exam tables from dirt and moisture. Choose from either a traditional smooth finish with machine glazing or a soft crepe finish for more absorbency, less noise and less slipping. All our exam table paper is mechanically metered to exact lengths.

Smooth Deluxe Exam Table Paper (Made from 50% Recycled Pulp)

Individually Poly Bagged

ITEM NUMBER	SIZE	COLOR	PACKAGING
NON24322	18" x 225' (45.7 cm x 68.58 m)	White	12 Rls/cs
NON24323	20" x 225' (50.8 cm x 68.58 m)	White	12 Rls/cs
NON24326	21" x 225' (53.3 cm x 68.58 m)	White	12 Rls/cs
NON24319	14.5" x 225' (36.8 cm x 68.58 m)	White	12 Rls/cs



Crepe Deluxe Exam Table Paper (Made from 100% Recycled Pulp)

Individually Poly Bagged

ITEM NUMBER	SIZE	COLOR	PACKAGING
NON24329	14.5" x 125' (36.8 cm x 38.10 m)	White	12 Rls/cs
NON24324	18" x 125' (45.7 cm x 38.10 m)	White	12 Rls/cs
NON24325	21" x 125' (53.3 cm x 38.10 m)	White	12 Rls/cs

Smooth Standard Exam Table Paper (Made from 50% Recycled Pulp)

Bulk Packaging

ITEM NUMBER	SIZE	COLOR	PACKAGING
NON23319	14.5" x 225' (36.8 cm x 68.58 m)	White	12 Rls/cs
NON23322	18" x 225' (45.7 cm x 68.58 m)	White	12 Rls/cs
NON23323	20" x 225' (50.8 cm x 68.58 m)	White	12 Rls/cs
NON23326	21" x 225' (53.3 cm x 68.58 m)	White	12 Rls/cs
NON18550	18" x 550' (45.7 cm x 167.64 m)	White	12 Rls/cs
NON21260	21" x 260' (53.3 cm x 79.24 m)	White	12 Rls/cs
NON25323	20" x 225' (50.8 cm x 68.58 m)	White	12 Rls/cs

Crepe Standard Exam Table Paper (Made from 100% Recycled Pulp)

Bulk Packaging

ITEM NUMBER	SIZE	COLOR	PACKAGING
NON23324	18" x 125' (45.7 cm x 38.10 m)	White	12 Rls/cs
NON23321	20" x 125' (50.8 cm x 38.10 m)	White	12 Rls/cs
NON23325	21" x 125' (53.3 cm x 38.10 m)	White	12 Rls/cs



Exam Table Paper

Printed Standard Exam Table Paper – Bulk Packaging



ITEM NUMBER	SIZE	PRINT	PACKAGING
NON23324BC	Crepe, 18" x 125' (45.7 cm x 38.10 m)	Beach Comber	6 Rls/cs
NON25202	Smooth, 18" x 225' (45.7 cm x 68.58 m)	Beach Comber	6 Rls/cs
NON25206	Smooth, 21" x 225' (53.3 cm x 68.58 m)	Beach Comber	6 Rls/cs
NON23324FT	Crepe, 18" x 125' (45.7 cm x 38.10 m)	Fish Tales	6 Rls/cs
NON23325FT	Crepe, 21" x 125' (53.3 cm x 38.10 m)	Fish Tales	6 Rls/cs
NON23322FT	Smooth, 18" x 225' (45.7 cm x 68.58 m)	Fish Tales	6 Rls/cs
NON23326FT	Smooth, 21" x 225' (53.3 cm x 68.58 m)	Fish Tales	6 Rls/cs

Crepe Infant Scale Paper – Bulk Packaging

ITEM NUMBER	SIZE	COLOR	PACKAGING
NON253200CS	18" x 24" (45.7 cm x 61 cm)	White	1,000/cs
NON243200	20" x 30" (50.8 cm x 76.2 cm)	White	1,000/cs
NON25331	12" x 20" (30.5 cm x 50.8 cm)	White	600/cs
NON25332	18" x 20" (45.7 cm x 50.8 cm)	White	400/cs



Apparel

Provide patients with extra comfort and the privacy they deserve with Medline disposable apparel. Choose from tissue/poly/tissue and multi-layer materials in a variety of colors.

Tissue / Poly / Tissue Material

- Outer softness for superb patient comfort
- Poly barrier to help protect against fluids

Premium Tissue / Poly / Tissue Patient Capes

ITEM NUMBER	SIZE	COLOR	PACKAGING
NON24259	30" x 21" (76.2 cm x 53.3 cm)	Mauve	100/cs
NON24258	30" x 21" (76.2 cm x 53.3 cm)	White	100/cs
NON24256	30" x 21" (76.2 cm x 53.3 cm)	Blue	100/cs
NON24257	30" x 21" (76.2 cm x 53.3 cm)	Lavender	100/cs



Premium material offers exceptional durability with a heavier basis weight.

Deluxe Tissue/ Poly / Tissue Patient Capes

ITEM NUMBER	SIZE	COLOR	PACKAGING
NON25249	30" x 21" (76.2 cm x 53.3 cm)	Blue	100/cs
NON24248	30" x 21" (76.2 cm x 53.3 cm)	White	100/cs



Deluxe material offers excellent durability with an industry-standard basis weight.

Apparel

Tissue / Poly / Tissue Material

Premium Tissue / Poly / Tissue Patient Gowns

ITEM NUMBER	SIZE	COLOR	PACKAGING
NON24244	Waist Ties, 30" x 42" (76.2 cm x 106.7 cm)	Blue	50/cs
NON24245	Waist Ties, 30" x 42" (76.2 cm x 106.7 cm)	White	50/cs
NON24252	Waist Ties, 30" x 42" (76.2 cm x 106.7 cm)	Mauve	50/cs
NON24243	Waist Ties, 30" x 42" (76.2 cm x 106.7 cm)	Purple	50/cs



Premium material offers exceptional durability with a heavier basis weight.

Deluxe Tissue / Poly / Tissue Patient Gowns

ITEM NUMBER	SIZE	COLOR	PACKAGING
NON24355	Waist Ties, 30" x 42" (76.2 cm x 106.7 cm)	White	50/cs
NON24356	Waist Ties, 30" x 42" (76.2 cm x 106.7 cm)	Blue	50/cs



Deluxe material offers excellent durability with an industry-standard basis weight.

Deluxe Tissue / Poly / X-Ray/Patient Gowns

ITEM NUMBER	SIZE	COLOR	PACKAGING
NON24354	Waist Ties, 30" x 42" (76.2 cm x 106.7 cm)	Blue	50/cs



Deluxe material offers excellent durability with an industry-standard basis weight.

Apparel

Polypropylene Material

- Multi-layer consists of inner layers of meltblown polypropylene sandwiched between outer layers of spunbond polypropylene
- Strong, fluid-resistant and breathable
- As comfortable as fabric
- Single Ply Polypropylene consists of bonding fibers together to form a single layer of breathable, woven-like material
- Choose from sleeveless and short-sleeve gown styles



Sleeveless Multi-Layer Patient Gowns

ITEM NUMBER	SIZE	COLOR	PACKAGING
NON27046	Regular/Large	Blue	50/cs

Short-Sleeved Multi-Layer Patient Gowns

ITEM NUMBER	SIZE	COLOR	PACKAGING
NON27146SL	Regular/Large	Blue	50/cs
NON27146SLXL	X-Large	Blue	50/cs
NON27346SL	Regular/Large, 3 Arm Hole	Blue	50/cs



Multi-Layer Patient Capes

ITEM NUMBER	SIZE	COLOR	PACKAGING
NON27036	Regular/Large	Dark Blue	100/cs



Multi-Layer Robes

ITEM NUMBER	SIZE	COLOR	PACKAGING
NON27148	Regular/Large	Blue	25/cs
NON27148XL	X-Large	Blue	25/cs
NON27148XXL	XX-Large	Blue	25/cs



Polypropylene Telemetry Pouches with Apron Ties and Snap Closure

ITEM NUMBER	SIZE	COLOR	PACKAGING
NON249005	4.5" x 7.5" (11.4 x 19 cm)	Blue	100/cs



Bariatric Patient Gowns with Soft Scrim Material

ITEM NUMBER	SIZE	COLOR	PACKAGING
TID918156	XX-Large	Blue	25/cs
TID918158	XXX-Large	Blue	25/cs

Apparel

Disposable Scrubs

- An alternative to reusable and laundered scrubs
- Made from breathable, fluid-resistant multi-layer material
- Available in blue



* When ordering, add the size code to the end of the item number.

S = Small

M = Medium

L = Large

XL = Extra Large

XXL = XX Large

XXXL = XXX Large

(Example: NON27203XL = Extra-Large Scrub Pants)

Multi-Layer Disposable Scrub Shirts

ITEM NUMBER	DESCRIPTION	COLOR	PACKAGING
NON27202*	V-Neck, Available Small-XXL	Blue	30/cs
NON27212*	Round Neck, Available Small-XXXL	Blue	30/cs

Generously cut shirts have one chest pocket and are available in either round neck or v-neck.

Multi-Layer Disposable Scrub Pants and Shorts

ITEM NUMBER	DESCRIPTION	COLOR	PACKAGING
NON27203*	Scrub Pants, Drawstring Waist, Available Small-XXL	Blue	30/cs
NON27213*	Scrub Pants, Elastic Waist, Available Small-XXXL	Blue	30/cs
NON27209*	Scrub Shorts, Elastic Waist, Available Small-XXL	Blue	30/cs



NON27209

Pants have one rear pocket, reinforced crotch area, and are available with drawstring or elastic waist.

Sheets and Blankets

Full sheet construction offers complete patient privacy. Choose from flat sheets, fitted sheets, impervious sheets and our newest product, super absorbent gurney sheets. (Not recommended for patient transfer.)

Drape Sheets

- Full-sheet construction to aid in patient modesty
- Variety of materials, including 2-ply, 3-ply, Tissue/Poly, Poly/Cellulose and Non-Woven
- Soft and comfortable, yet durable
- Use for draping patients or equipment
- Polyester/cellulose drape sheets provide extra protection from fluids

2-Ply, All-Tissue Drape Sheets

ITEM NUMBER	SIZE	COLOR	PACKAGING
TID918328	13" x 24" (33 cm x 61 cm)	White	500/cs
TID980823	24" x 40" (61 cm x 101.6 cm)	White	200/cs
NON24339	40" x 48" (101.6 cm x 121.9 cm)	White	100/cs
NON23339	40" x 48" (101.6 cm x 121.9 cm), Economical Thinner 2-Ply	White	100/cs
TID918322	40" x 48" (101.6 cm x 121.9 cm)	Peach	100/cs
NON24339M	40" x 48" (101.6 cm x 121.9 cm)	Mauve	100/cs
NON24339A	40" x 60" (101.6 cm x 152.4 cm)	White	100/cs
NON24339B	40" x 72" (101.6 cm x 182.9 cm)	White	50/cs



3-Ply, All Tissue Drape Sheets

ITEM NUMBER	SIZE	COLOR	PACKAGING
NON24336	40" x 48" (101.6 cm x 121.9 cm)	Color	100/cs
NON24336A	40" x 60" (101.6 cm x 152.4 cm)	White	100/cs
NON24339C	40" x 72" (101.6 cm x 182.9 cm)	White	50/cs
NON24336C	40" x 90" (101.6 cm x 228.6 cm)	White	50/cs



Sheets and Blankets

Full sheet construction offers complete patient privacy. Choose from flat sheets, fitted sheets, impervious sheets and our newest product, super absorbent gurney sheets. (Not recommended for patient transfer.)

Drape Sheets

Tissue/Poly Exam Sheets (Made in the USA)

ITEM NUMBER	DESCRIPTION	COLOR	PACKAGING
NON24340	Standard, 40" x 48" (101.6 cm x 121.9 cm)	Color	100/cs
NON24338	Standard, 40" x 60" (101.6 cm x 152.4 cm)	Blue	100/cs
NON24333	Standard, 40" x 72" (101.6 cm x 182.9 cm)	Blue	50/cs
NON24335	Standard, 40" x 90" (101.6 cm x 228.6 cm)	Blue	50/cs
NON24337	Standard, 46" x 96" (116.8 cm x 243.8 cm)	Blue	50/cs



Sheets Polyester/Cellulose Drape Sheets (Made in the USA)

ITEM NUMBER	DESCRIPTION	COLOR	PACKAGING
NON38700	Moisture-Proof, 40" x 84" (101.6 cm x 213.4 cm)	White	100/cs



Extra protection from fluids.

Tissue/Poly Bed Flat Sheet (Made in the USA)

ITEM NUMBER	DESCRIPTION	COLOR	PACKAGING
NON24330A	Bed Sheet, Standard, Tissue/Poly, 58" x 96" (147.3 cm x 243.8 cm)	Blue	25/cs
NON24330B	Bed Sheet, Standard, Tissue/Poly, 58" x 102" (147.3 cm x 259 cm)	Blue	25/cs



Non-Woven Flat Sheet

ITEM NUMBER	DESCRIPTION	COLOR	PACKAGING
NON33025	Sheet, SMMS, 36" x 60" (91.4 cm x 152.4 cm)	Dark Blue	50/cs
NON33100	Sheet, SMMS, (1.4 oz.), 40" x 84" (101.6 cm x 213.4 cm)	Blue	50/cs



Sheets and Blankets

Stretcher / Gurney Sheets and Blankets

- Choose from a variety of sheets, sheet kits and pillow cases
- Draw sheets are soft and comfortable, yet durable
- Fitted stretcher sheet kits provide a low-cost alternative to traditional linens

Draw Sheets

ITEM NO.	DESCRIPTION	COLOR	PACKAGING
NON32100	Rayon-Polyester, 32" x 48" (81.3 cm x 121.9 cm)	White	50/cs
NON32200	Rayon-Polyester, 40" x 60" (101.6 cm x 152.4 cm)	White	50/cs



Soft and comfortable, yet durable.



A low-cost alternative to traditional linens.

Fitted Stretcher Sheet Kit

ITEM NO.	DESCRIPTION	COLOR	PACKAGING
NON37250	Spunbond, Top, Bottom, and Pillowcase	Blue	24/cs
NON37500	Spunbond, Top, Bottom, Pillowcase and Pad	Blue	24/cs
NON47200	Top Sheet (White) and Fitted Heavy Weight (85 Gram) Bottom Sheet (Blue)	Blue/White	50/cs

Fitted Stretcher Sheet

ITEM NO.	DESCRIPTION	COLOR	PACKAGING
NON34000	Spunbond, 65 g (2.29 oz.), 32" x 72" (81.3 cm x 182.9 cm)	Blue	50/cs



Sheets and Blankets

Stretcher / Gurney Sheets and Blankets

Fitted Gurney Sheets

ITEM NO.	DESCRIPTION	COLOR	PACKAGING
NON34600	Spunbond, 65 g (2.29 oz.), 40" x 80" (101.6 cm x 203.2 cm)	Blue	50/cs



Super-Absorbent Fitted Gurney Sheets

ITEM NO.	DESCRIPTION	COLOR	PACKAGING
NON34900	Spunbond, 65 g (2.29 oz.), 41" x 79" (104.1 cm x 200.1 cm), absorbent area measures 20" x 72" (50.8 cm x 182.9 cm)	Blue	24/cs



Emergency Blankets

ITEM NO.	DESCRIPTION	COLOR	PACKAGING
NONDB4080	Poly/Cellulose, 40" x 80" (101.6 cm x 203.2 cm)	Gray	10/cs
LPC980043	Emergency, Poly/foam, 58" x 90" (147.3 cm x 228.6 cm)	Yellow	18/cs
AHS1646015	Mylar®, 52" x 84", (132.1 cm x 213.4 cm)	Silver	5/cs



An excellent warming wrap for patient transport – keep on hand for emergency or disaster situations.

Disposable Pillowcases and Sheets (Made in the USA)

ITEM NO.	DESCRIPTION	COLOR	PACKAGING
NON24345	Tissue/Poly, Standard 21" x 30" (53.3 cm x 76.2 cm)	White	100/cs
NON24346	Tissue/Poly, Standard 21" x 30" (53.3 cm x 76.2 cm)	Blue	100/cs
NON38500	Poly/Cellulose, 20" x 30" (50.8 cm x 76.2 cm)	White	100/cs
TID919350	Polyethylene, 21" x 30" (53.3 cm x 76.2 cm), Impervious	White	100/cs
NON25300	Ultracel®, 21" x 30" (53.3 cm x 76.2 cm), Fluid-Resistant	White	100/cs



SMS (spunbond/melt blown/spunbond) provides excellent fluid repellency as well as breathability.

Disposable SMS Pillowcases and Sheets

ITEM NO.	DESCRIPTION	COLOR	PACKAGING
NON32300	SMS, 20" x 29" (50.8 cm x 73.7 cm)	White	100/cs
NON32500	SMS, 20" x 29" (50.8 cm x 73.7 cm)	Blue	100/cs

Cubicle Curtains

Exclusive Medline designs are great for emergency rooms and isolation rooms, as well as for temporary use while servicing permanent drapes. The complete snap-together system easily attaches to the mesh header. Each panel or the entire system can be quickly changed by one person. The curtains comply with standard 701 of the National Fire Protection Association (NFPA) for flame propagation of textiles and films.

Disposable Cubicle Curtain Panels

ITEM NUMBER	SIZE	COLOR	PACKAGING
NON528472	60" x 72" (152.4 cm x 182.9 cm)	Light Blue	24/cs
NON528772	60" x 72" (152.4 cm x 182.9 cm)	Dark Blue	24/cs
NON528764	60" x 64" (152.4 cm x 162.6 cm)	Dark Blue	24/cs


Disposable Cubicle Mesh Header

ITEM NUMBER	SIZE	COLOR	PACKAGING
NON52060	60" x 22" (152.4 cm x 55.9 cm)	White	24/cs

Mesh-Free Cubicle Curtains


ITEM NUMBER	SIZE	COLOR	PACKAGING
NON527480	60" x 80" (152.4 cm x 203.2 cm)	Light Blue	24/cs

No mesh needed – the eyelets are in the SMS panels, which snap together.



Professional Towels

Professional towels are used as drapes, covers, surface protection, dental bibs and more. Choose from a wide selection of colors in 2-or 3-ply all-tissue or tissue/poly materials.



2-Ply Tissue/Poly Professional Towels

ITEM NUMBER	SIZE	COLOR	PACKAGING
NON24356W	13" x 18" (33 cm x 45.7 cm)	White	500/cs
NON24356B	13" x 18" (33 cm x 45.7 cm)	Blue	500/cs
NON24361	17" x 19" (43.2 cm x 48.3 cm)	White	500/cs

3-Ply Tissue Professional Towels

ITEM NUMBER	SIZE	COLOR	PACKAGING
NON24357W	13" x 18" (33 cm x 45.7 cm)	White	500/cs
NON24357B	13" x 18" (33 cm x 45.7 cm)	Blue	500/cs
NON24359	17" x 19" (43.2 cm x 48.3 cm)	White	500/cs
NON24357M	13" x 18" (33 cm x 45.7 cm)	Mauve	500/cs

3-Ply Tissue/Poly Professional Towels

ITEM NUMBER	SIZE	COLOR	PACKAGING
NON24358W	13" x 18" (33 cm x 45.7 cm)	White	500/cs
NON24358B	13" x 18" (33 cm x 45.7 cm)	Blue	500/cs

Glossary of Terms

SMS

3-ply material consisting of spunbond/melt blown/ spunbond material, for fluid repellency while still being breathable.

SMMS

4-ply material consisting of spunbond/melt blown/ spunbond material, for superb fluid repellency while still being breathable.

PPSB

Polypropylene, spundbond.

2-Ply Tissue

Two sheets of tissue paper, embossed.

3-Ply Tissue

Three sheets of tissue paper, embossed.

Tissue/Poly

One layer of tissue and one layer of poly film, embossed.

Tissue/Poly/Tissue

One layer of poly film between two layers of tissue, embossed.

Smooth

Durable machine-glazed paper with 50% recycled content.

Crepe

Soft, quieter, more absorbent paper with 100% recycled content.



Environmental Services Brochure



TWL Housekeeping, Infection Control & Floor Care

VISIT WWW.MEDLINE.COM

- 1 Download a QR Code reader app
- 2 Launch app
- 3 Hover the phone over this QR Code



Medline Industries, Inc.
 One Medline Place
 Mundelein, IL 60060

Medline United States
1-800-MEDLINE (633-5463)
www.medline.com
info@medline.com

Medline Canada
1-800-396-6996
www.medline.ca
canada@medline.com

Medline México
01-800-831-0898
www.medlinemexico.com
mexico@medline.com

Some products may not be available for sale in Mexico or Canada. We reserve the right to correct any errors that may occur within this brochure.
 ©2011 Medline Industries, Inc. Medline is a registered trademark and greensmart is a trademark of Medline Industries, Inc
 Mylar is a registered trademark of Dupont Teijin Films US. Ultracel is a registered trademark of Bayer Aktiengesellschaft