




UROLOGY SOLUTIONS



Products and resources that support urological care across all settings

You're here for them. We're here for you.

You take responsibility for your patients' care as they move through the continuum of care – from the hospital to skilled nursing to home. So do we. You'll find Medline at every step of your patients' healthcare journey. Listening. Responding. And helping you achieve optimal outcomes.



Highest Quality Products at the Best Value

We offer you a comprehensive range of clinically proven products with manufacturer-direct savings.


- » Patented, nurse-designed urinary catheter trays that align clinical behavior with evidence-based practices
- » Bladder scanners and Foley alternatives to help prevent unnecessary catheterizations
- » Urological care solutions for chronic conditions

Proactive Support

Having a Medline Representative with the clinical resources to create a product portfolio tailored to your needs is invaluable.

- » Experienced clinicians prepared to help you resolve urological challenges
- » Product specialists ready to problem solve with you

E



Staff and Patient Education

Knowledge is power. Medline University, our online learning resource, gives your caregivers and patients immediate access to educational materials to help them take charge.

- » Staff can earn free continuing education credits on a wide variety of subjects to help them reach their full potential
- » Easy-to-understand urological guides offer patients clear descriptions and visual information on their condition and care

Table of Contents

■ Foley Catheter Trays and Bladder Scanner	4
■ Foley Catheters	26
■ Urine Collection	34
■ Catheter Securement Devices	42
■ Alternatives to Foley Catheters	44
■ Pediatrics	54
■ Irrigation and Other Products	62



FOLEY CATHETER TRAYS AND BLADDER SCANNERS

ERASE CAUTI® Solution 5
 One Layer Foley Catheter Trays 8
 Accuryn® Monitoring System 17
 2-Layer Foley Catheter Trays 18
 Foley Catheter Insertion Trays 23
 BioCon™ Bladder Scanners 24



ERASECAUTI®
Comprehensive Care Solution

Removing Barriers to CAUTI Prevention.

Catheter-associated urinary tract infections (CAUTIs) are one of the most common healthcare-associated infections (HAIs)¹. CAUTI protocols exist, but with so many competing priorities, compliance to best practices is often inconsistent.

Medline worked with clinical leaders in hospitals nationwide to develop an effective and easy-to-adopt process that helps drive successful practice change. It combines clinical thinking with product and educational resources to create a CAUTI-reduction solution that can be customized to your facility's specific needs.

ERASE CAUTI® Comprehensive Care Solution



Learn how our ERASE CAUTI Comprehensive Care Solution can make a difference in your facility. Contact your Medline Representative or visit medline.com/erase-cauti.

1. Meddings J, Kreiin SL, Fakhri MG, et al. Reducing Unnecessary Urinary Catheter Use and Other Strategies to Prevent Catheter-Associated Urinary Tract Infections: Brief Update Review. In: Making Health Care Safer II: An Updated Critical Analysis of the Evidence for Patient Safety Practices. Rockville (MD): Agency for Healthcare Research and Quality (US); 2013 Mar. (Evidence Reports/Technology Assessments, No. 211.) Chapter 9. Available from: <http://www.ncbi.nlm.nih.gov/books/NBK133354/> Accessed November 24, 2015.

Part of the
ERASECAUTI®
 Comprehensive Care Solution






Products to Support Practice

Throughout this catalog, you'll find innovatively designed products to help reduce unnecessary indwelling urinary catheter (IUC) use, improve IUC insertion technique and enhance the care and maintenance of these devices. All these products support our ERASE CAUTI solution.

This section features our extensive line of One Layer IUC trays designed by nurses to drive compliance to best insertion practices.

- » Closed-system and Add-a-Cath configurations
- » Regular or Temperature Sensing Foleys
- » Silicone, silver silicone, or latex Foley catheters
- » Foley sizes ranging from 6-18 Fr
- » Drain bags or urine meters

Color-coded packaging identifies the catheter inside the tray, helping to prevent the inconvenience and expense of pulling the wrong tray.


- | | | | |
|---|---|---|--|
|  | Add-A-Cath Trays |  | Select Silicone® Foley Catheter Trays |
|  | Silvertouch® 100% Silicone Foley Catheter Trays |  | Select Silicone Temperature Sensing Foley Catheter Trays |
|  | Simply Latex® Foley Catheter Trays | | |



Look for new innovations like our vented drain bags and meters to improve flow, our elastic-blend silicone catheter, the new SafeSecure securement device, and more.

Changing Behavior by Design.

Created by nurses for nurses, our patented one layer Foley catheter tray helps staff consistently follow aseptic technique. Through its design and the acronym ERASE, the tray intuitively guides clinicians through the insertion process and CAUTI prevention.



One Layer Foley Catheter Trays

Medline's patented* One Layer Foley catheter tray design intuitively leads healthcare workers through established best practices. This unique design also includes labeling and instructions, making it easier to do the right thing.

The Trays Include:

- » Patient education card
- » Foley securement device
- » PVP swabsticks, 3 pack
- » 10 mL prefilled inflation syringe
- » Lubricating jelly syringe
- » Specimen container
- » Wrapped Aloetouch® 3G gloves
- » Fenestrated drape
- » Moisture-proof underpad
- » Foley Insetag®
- » Hand sanitizer packet
- » Sureprep® skin prep pad for securement device placement
- » Alcohol prep pad for securement device placement
- » All Foley Catheter Trays are sterile

*For a current list of patent numbers, visit www.medline.com/patents.



URO170100

One Layer Add-a-Cath Trays

The Add-a-Cath tray allows clinicians to select the Foley catheter based on the patient's needs.

- » Features Medline's One Layer Foley catheter tray to help ensure that best practices are followed
- » Includes 10 mL inflation syringe and all necessary Foley catheter insertion components
- » Not made with natural rubber latex
- » Foley catheter not included

Item No.	Description	HCPCS	Pkg.
One Layer Trays – Add-a-Cath Tray			
DYND160100*	2,000 mL Drain Bag	A4354 + A4333	10/cs
DYND160200	400 mL Urine Meter with 2,500 mL Drain Bag	A4354 + A4333	10/cs
Total One Layer Trays – Include Peri-Care Wipe			
URO170100*	2,000 mL Drain Bag	No Code	10/cs
URO170200	400 mL Urine Meter with 2,500 mL Drain Bag	No Code	10/cs
URO170200S*	400 mL Urine Meter with 2,500 mL Drain Bag	No Code	10/cs
Total One Layer Trays – Includes Peri-Care Wipe, Vented Drainage and SafeSecure			
URO180100*	2,000 mL Drain Bag	No Code	10/cs
URO180200	400 mL Urine Meter with 2,500 mL Drain Bag	No Code	10/cs
URO180200S*	400 mL Urine Meter with 2,500 mL Drain Bag	No Code	10/cs

*Metal-Free product is safe for MRI use



URO170416

One Layer Silvertouch® 100% Silicone Foley Catheter Trays

Our patented process binds silver ions to the 100% silicone Foley catheter. The hydrophilic coating is then covalently bonded to the catheter surface.

- » Features Medline's One Layer Foley Catheter Tray to help ensure best practices are followed
- » Foley catheter hydrates quickly and maintains lubricity to help reduce friction and irritation
- » Not made with natural rubber latex

Item No.	Description	HCPCS	Pkg.
----------	-------------	-------	------

One Layer Trays

DYND160414*	14 Fr 10 mL, with 2,000 mL Drain Bag	A4315 + A4333	10/cs
DYND160416*	16 Fr 10 mL, with 2,000 mL Drain Bag	A4315 + A4333	10/cs
DYND160418*	18 Fr 10 mL, with 2,000 mL Drain Bag	A4315 + A4333	10/cs
DYND160514	14 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	A4315 + A4333	10/cs
DYND160516	16 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	A4315 + A4333	10/cs
DYND160518	18 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	A4315 + A4333	10/cs
DYND160516S*	16 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	No Code	10/cs
DYND160518S*	18 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	No Code	10/cs

Total One Layer Trays – Include Peri-Care Wipe

URO170414*	14 Fr 10 mL, with 2,000 mL Drain Bag	No Code	10/cs
URO170416*	16 Fr 10 mL, with 2,000 mL Drain Bag	No Code	10/cs
URO170418*	18 Fr 10 mL, with 2,000 mL Drain Bag	No Code	10/cs
URO170514	14 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	No Code	10/cs
URO170516	16 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	No Code	10/cs
URO170518	18 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	No Code	10/cs

*Metal-Free product is safe for MRI use

One Layer Foley Catheters

One Layer SimplyLatex® Foley Catheter Trays

Innovative manufacturing process produces the maximum amount of silicone coating available.

- » Features Medline's One Layer tray to help ensure best practices are followed
- » Silicone-elastomer coating on the latex Foley catheter creates a slippery surface to reduce urethral trauma and helps resist encrustation



DYND160116



URO170116

Item No.	Description	HCPCS	Pkg.
----------	-------------	-------	------

One Layer Trays

DYND160114*	14 Fr 10 mL, with 2,000 mL Drain Bag	A4314 + A4333	10/cs
DYND160116*	16 Fr 10 mL, with 2,000 mL Drain Bag	A4314 + A4333	10/cs
DYND160118*	18 Fr 10 mL, with 2,000 mL Drain Bag	A4314 + A4333	10/cs
DYND160214	14 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	A4314 + A4333	10/cs
DYND160216	16 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	A4314 + A4333	10/cs
DYND160218	18 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	A4314 + A4333	10/cs
DYND160216S*	16 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag and Metal-Free Slide-Tap Drainage Port	A4314 + A4333	10/cs
DYND160218S*	18 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag and Metal-Free Slide-Tap Drainage Port	A4314 + A4333	10/cs

Total One Layer Trays – Include Peri-Care Wipe

URO170114*	14 Fr 10 mL, with 2,000 mL Drain Bag	No Code	10/cs
URO170116*	16 Fr 10 mL, with 2,000 mL Drain Bag	No Code	10/cs
URO170118*	18 Fr 10 mL, with 2,000 mL Drain Bag	No Code	10/cs
URO170214	14 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	No Code	10/cs
URO170216	16 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	No Code	10/cs
URO170218	18 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	No Code	10/cs

Total One Layer Trays – Includes Peri-Care Wipe, Vented Drainage and SafeSecure

URO180114*	14 Fr 10 mL, with 2,000 mL Drain Bag	No Code	10/cs
URO180116*	16 Fr 10 mL, with 2,000 mL Drain Bag	No Code	10/cs
URO180118*	18 Fr 10 mL, with 2,000 mL Drain Bag	No Code	10/cs
URO180214	14 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	No Code	10/cs
URO180216	16 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	No Code	10/cs
URO180218	18 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	No Code	10/cs

*Metal-Free product is safe for MRI use



DYND160716C

One Layer Coude Foley Catheter Trays

For difficult insertions, coude-tip Foley catheters feature a bent tip that maneuvers around obstructions leading to the bladder.

- » Featuring Medline's One Layer Foley Catheter Tray to help ensure that best practices are followed
- » Large oval-shaped eye for maximum drainage
- » Smooth reinforced tip resists buckling
- » Available with silicone-elastomer coated latex or 100% silicone coude catheter

Item No.	Description	HCPCS	Pkg.
DYND160116C*	16 Fr 10 mL Latex Coude Foley Catheter Tray with 2,000 mL Drain Bag	A4352 + A4333	10/cs
DYND160716C*	16 Fr 10 mL 100% Silicone Coude Foley Catheter Tray with 2,000 mL Drain Bag	A4352 + A4333	10/cs
URO170816C	16 Fr 10 mL 100% Silicone Coude Foley Catheter Tray with 400 mL Urine Meter with 2,500 mL Drain Bag, Includes Peri-Care Wipe	A4352 + A4333	10/cs

*Metal-Free product is safe for MRI use

One Layer Foley Catheters



DYND160814

One Layer SelectSilicone® Foley Catheter Trays

Maximize safety and comfort with our 100% silicone Foley catheters.

- » Features Medline's One Layer Foley Catheter Tray to help ensure best practices are followed
- » Completely inert for less tissue irritation and encrustation during extended periods of indwelling use
- » Foley catheter has bullet-shaped tip for easier insertion and less buckling
- » Large, smooth eyes for maximum comfort
- » Not made with natural rubber latex

Item No.	Description	HCPCS	Pkg.
One Layer Trays			
DYND160714*	14 Fr 10 mL, with 2,000 mL Drain Bag	A4315 + A4333	10/cs
DYND160716*	16 Fr 10 mL, with 2,000 mL Drain Bag	A4315 + A4333	10/cs
DYND160718*	18 Fr 10 mL, with 2,000 mL Drain Bag	A4315 + A4333	10/cs
DYND160812*	12 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	A4315 + A4333	10/cs
DYND160814	14 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	A4315 + A4333	10/cs
DYND160816	16 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	A4315 + A4333	10/cs
DYND160818	18 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	A4315 + A4333	10/cs
DYND160814S*	14 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag and Metal-Free Slide-Tap Drainage Port	A4315 + A4333	10/cs
DYND160816S*	16 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag and Metal-Free Slide-Tap Drainage Port	A4315 + A4333	10/cs
DYND160818S*	18 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag and Metal-Free Slide-Tap Drainage Port	A4315 + A4333	10/cs

*Metal-Free product is safe for MRI use



URO170716

Item No.	Description	HCPCS	Pkg.
Total One Layer Trays – Include Peri-Care Wipe			
URO170714*	14 Fr 10 mL, with 2,000 mL Drain Bag	No Code	10/cs
URO170716*	16 Fr 10 mL, with 2,000 mL Drain Bag	No Code	10/cs
URO170718*	18 Fr 10 mL, with 2,000 mL Drain Bag	No Code	10/cs
URO170812*	12 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	No Code	10/cs
URO170814	14 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	No Code	10/cs
URO170816	16 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	No Code	10/cs
URO170818	18 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	No Code	10/cs
URO170814S*	14 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag and Metal-Free Slide-Tap Drainage Port	No Code	10/cs
URO170816S*	16 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag and Metal-Free Slide-Tap Drainage Port	No Code	10/cs
URO170818S*	18 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag and Metal-Free Slide-Tap Drainage Port	No Code	10/cs
Total One Layer Trays – Includes Peri-Care Wipe, Vented Drainage and SafeSecure			
URO180714*	14 Fr 10 mL, with 2,000 mL Drain Bag	No Code	10/cs
URO180716*	16 Fr 10 mL, with 2,000 mL Drain Bag	No Code	10/cs
URO180718*	18 Fr 10 mL, with 2,000 mL Drain Bag	No Code	10/cs
URO180714CN*†	14 Fr 10 mL, with 2,000 mL Drain Bag	No Code	10/cs
URO180716CN*†	16 Fr 10 mL, with 2,000 mL Drain Bag	No Code	10/cs
URO180718CN*†	18 Fr 10 mL, with 2,000 mL Drain Bag	No Code	10/cs
URO180814	14 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	No Code	10/cs
URO180816	16 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	No Code	10/cs
URO180816S*	16 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag and Metal-Free Slide-Tap Drainage Port, SafeSecure	No Code	10/cs
URO180818	18 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	No Code	10/cs
URO180818S*	18 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag and Metal-Free Slide-Tap Drainage Port, SafeSecure	No Code	10/cs
URO180814CN†	14 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	No Code	10/cs
URO180816CN†	16 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	No Code	10/cs
URO180818CN†	18 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	No Code	10/cs
URO180814SCN*†	14 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	No Code	10/cs
URO180816SCN*†	16 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	No Code	10/cs
URO180818SCN*†	18 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	No Code	10/cs

* Metal-Free product is safe for MRI use

† Includes Centurion Securement Device

One Layer Foley Catheters



DYND160814T

One Layer SelectSilicone Temperature-Sensing Foley Catheter Trays

Temperature-sensing catheter helps clinicians continuously monitor a patient's core body temperature via the bladder.

- » Features Medline's One Layer Foley Catheter Tray to help ensure best practices are followed
- » 400 Series thermistor
- » Catheter has $\pm 0.2^{\circ}\text{C}$ accuracy
- » Available in pediatric and adult sizes
- » Not made with natural rubber latex

Item No.	Description	HCPCS	Pkg.
----------	-------------	-------	------

One Layer Trays

DYND160814T	14 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	A4315 + A4333	10/cs
DYND160816T	16 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	A4315 + A4333	10/cs
DYND160818T	18 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	A4315 + A4333	10/cs
DYND160814TS*	14 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag and Metal-Free Slide-Tap Drainage Port	No Code	10/cs
DYND160816TS*	16 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag and Metal-Free Slide-Tap Drainage Port	A4315 + A4333	10/cs
DYND160818TS*	18 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag and Metal-Free Slide-Tap Drainage Port	No Code	10/cs

Total One Layer Trays – Includes Peri-Care Wipe

URO170814T	14 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	No Code	10/cs
URO170816T	16 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	No Code	10/cs
URO170818T	18 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	No Code	10/cs
URO170814TS*	14 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag and Metal-Free Slide-Tap Drainage Port	No Code	10/cs
URO170816TS*	16 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag and Metal-Free Slide-Tap Drainage Port	No Code	10/cs
URO170818TS*	18 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag and Metal-Free Slide-Tap Drainage Port	No Code	10/cs

*Conditionally MRI compatible



Item No.	Description	HCPCS	Pkg.
Total One Layer Trays – Includes Peri-Care Wipe, Vented Drainage and SafeSecure			
URO180814T	14 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	No Code	10/cs
URO180816T	16 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	No Code	10/cs
URO180818T	18 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	No Code	10/cs
URO180814TS*	14 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag and Metal-Free Slide-Tap for Drainage Port	No Code	10/cs
URO180816TS*	16 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag and Metal-Free Slide-Tap for Drainage Port	No Code	10/cs
URO180818TS*	18 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag and Metal-Free Slide-Tap for Drainage Port	No Code	10/cs
URO180816TCN†	16 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	No Code	10/cs
URO180816TSCN*†	16 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag and Metal-Free Slide-Tap for Drainage Port	No Code	10/cs

One Layer Pediatric Trays

DYND160806T*	6 Fr 3 mL, with 200 mL Urine Meter with 2,000 mL Drain Bag	No Code	10/cs
DYND160808T*	8 Fr 5 mL, with 200 mL Urine Meter with 2,000 mL Drain Bag	A4315 + A4333	10/cs
DYND160810T*	10 Fr 5 mL, with 200 mL Urine Meter with 2,000 mL Drain Bag	A4315 + A4333	10/cs

Total One Layer Pediatric Trays – Includes Peri-Care Wipe

URO170806T*	6 Fr 3 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	No Code	10/cs
URO170808T*	8 Fr 5 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	A4315 + A4333	10/cs
URO170810T*	10 Fr 5 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	A4315 + A4333	10/cs

Total One Layer Pediatric Trays – Includes Peri-Care Wipe, Vented Drainage and SafeSecure

URO180806T*	6 Fr 3 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	No Code	10/cs
-------------	--	---------	-------

*Conditionally MRI compatible

† Includes Centurion Securement Device

Pediatric One Layer SelectSilicone Foley Catheter Trays

Buddy the Brave® is here to help ease the catheterization process on your pediatric patients. Featured on Medline's new Foley pediatric catheter trays, he provides valuable education for parents and advocates, as well as a fun activity card and sticker for kids.



Kit Includes:

- » Aloetouch 3G gloves
- » Foley Insertag
- » Hand antiseptic packet
- » Under-buttocks drape
- » Fenestrated drape
- » Lubricating jelly syringe
- » Sterile water syringe
- » Specimen container
- » PVP swabsticks, 3 pack
- » Pediatric securement device
- » Exact measure urine meter
- » SurePrep pad for securement device placement
- » Alcohol prep pad for securement device placement
- » Pediatric activity education card



One Layer 100% Silicone Pediatric Foley Catheter Trays

- » Features Medline's One Layer Foley Catheter Tray to help ensure best practices are followed
- » Available with either a 200 mL or 400 mL urine meter with connected drain bag
- » Metal-free and not made with natural rubber latex

Item No.	Description	HCPCS	Pkg.
One Layer Trays			
DYND160806	6 Fr 3 mL, with 200 mL Urine Meter with 2,000 mL Drain Bag	A4315 + A4333	10/cs
DYND160808	8 Fr 5 mL, with 200 mL Urine Meter with 2,000 mL Drain Bag	A4315 + A4333	10/cs
DYND160810	10 Fr 5 mL, with 200 mL Urine Meter with 2,000 mL Drain Bag	A4315 + A4333	10/cs
Total One Layer Trays – Includes Peri-Care Wipe			
URO170806	6 Fr 3 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	A4315 + A4333	10/cs
URO170808	8 Fr 5 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	A4315 + A4333	10/cs
URO170810	10 Fr 5 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	A4315 + A4333	10/cs
Total One Layer Trays – Includes Peri-Care Wipe, Vented Drainage and SafeSecure			
URO180806S	6 Fr 3 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	No Code	10/cs
URO180808S	8 Fr 5 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	No Code	10/cs
URO180810S	10 Fr 5 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	No Code	10/cs

Accuryn™ Monitoring System



Accuryn Monitoring System

Accuryn turns the Foley catheter into a next generation diagnostic tool that provides accurate, real-time measurement of urine output, intra-abdominal pressure and core body temperature—all at the push of a button.

This pioneering technology:

- » Requires virtually no setup
- » Delivers continuous data
- » Automatically clears the drain line
- » Aligns with best practices to reduce CAUTI risk

Item No.	Description	Pkg.
UROACCRN200	Accuryn Monitor with power cable	1/cs
UROACCRNMTSD	Accuryn Floor Stand, with IV pole clamp and stand mount	5/cs
UROACCRNSTND	Accuryn Floor Stand, without mount	1/cs
UROACCRNMNT	Accuryn Mount, IV pole clamp and stand mount	1/cs
URACCPWR	Accuryn Monitor power cable	1/cs
URO190814*	Accuryn Foley catheter tray, 14 Fr, U/O and Temp	10/cs
URO190816	Accuryn Foley catheter tray, 16 Fr, U/O and Temp	10/cs
URO190816P	Accuryn Foley catheter tray, 16 Fr, U/O, Temp and IAP	10/cs
URO190818*	Accuryn Foley catheter tray, 18 Fr, U/O and Temp	10/cs

*Coming in 2019



2-Layer Foley Catheter Trays

A closed system reduces the chance of bacterial introduction and secures the Foley catheter to the drain bag tubing. Medline's 2-Layer Foley Catheter Trays contain all the necessary supplies for catheter insertion.

All Sterile Foley Catheter Trays Include:

- » Compartmentalized prep tray
- » Aloetouch 3G wrapped gloves not made with natural rubber latex
- » Pre-saturated PVP swabsticks, 3 pack
- » Moisture-proof underpad
- » Fenestrated drape
- » Sterile water syringe
- » Lubricating jelly
- » Specimen container with label
- » 2,000 mL drain bag



2-Layer Add-a-Cath Foley Catheter Trays

Add-a-Cath trays contain all of the necessary supplies to prepare for catheter insertion, but allow clinicians to separately select the appropriate Foley catheter according to patient needs. Add-a-Cath Trays contain a blunt cannula sample port and do not contain natural rubber latex.

- » All kits contain 10 mL inflation syringe

Item No.	Description	Inflation Syringe Size	HCPCS	Pkg.
DYND11090	2,000 mL Drain Bag with Anti-Reflux Tower	10 mL	A4354	10/cs
DYND11095	2,000 mL Drain Bag with Anti-Reflux Tower	30 mL	A4354	10/cs
DYND11020	2,000 mL Drain Bag with Anti-Reflux Device	10 mL	A4354	10/cs
DYND11050	2,000 mL Drain Bag with Anti-Reflux Device	30 mL	A4354	10/cs



DYND140218

2-Layer Silvertouch 100% Silicone Foley Catheter Trays

Patented Silvertouch 100% Silicone Foley Trays feature a Silvertouch Foley catheter pre-connected to a Medline drain bag. The brilliant color on the Silvertouch catheter is proof of the silver hydrophilic coating—ensuring a full week's worth of comfort. All the necessary components for proper Foley catheter insertion are included.

» Hydrophilic coating on Silvertouch catheter helps to reduce friction and irritation

Item No.	Description	Sample Port Type	HCPCS	Pkg.
DYND140214	14 Fr 10 mL, with 2,000 mL Drain Bag with Anti-Reflux Device	Luer Lock	A4315	10/cs
DYND140216	16 Fr 10 mL, with 2,000 mL Drain Bag with Anti-Reflux Device	Luer Lock	A4315	10/cs
DYND140218	18 Fr 10 mL, with 2,000 mL Drain Bag with Anti-Reflux Device	Luer Lock	A4315	10/cs
DYND140416	16 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	Luer Lock	A4315	10/cs

2-Layer Foley Catheter Trays



DYND11003

2-Layer Silicone-Elastomer Coated Latex Foley Catheter Trays

Silicone-elastomer coated latex trays feature a latex Foley catheter pre-connected to a drainage bag or urine meter. All necessary components for proper Foley catheter insertion are included.

- » Innovative manufacturing process provides the maximum amount of silicone coating available
- » Internally and externally coated latex catheter in order to resist encrustation and bacterial colonization

Item No.	Description	Sample Port Type	HCPCS	Pkg.
DYND11855	14 Fr 10 mL, with 2,000 mL Drain Bag with Anti-Reflux Device	Blunt Cannula	A4314	10/cs
DYND11860	16 Fr 10 mL, with 2,000 mL Drain Bag with Anti-Reflux Device	Blunt Cannula	A4314	10/cs
DYND11865	18 Fr 10 mL, with 2,000 mL Drain Bag with Anti-Reflux Device	Blunt Cannula	A4314	10/cs
DYND11003	16 Fr 10 mL, with 2,000 mL Drain Bag with Anti-Reflux Tower	Blunt Cannula	A4314	10/cs
DYND11008	18 Fr 10 mL, with 2,000 mL Drain Bag with Anti-Reflux Tower	Blunt Cannula	A4314	10/cs
DYND1516SEC	16 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	Blunt Cannula	A4314	10/cs
DYND1518SEC	18 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	Blunt Cannula	A4314	10/cs
URO1SPLDB14	14 Fr 10 mL, with 2,000 mL Drain Bag with Anti-Reflux Tower	Luer Lock	No Code	10/cs
URO1SPLDB16	16 Fr 10 mL, with 2,000 mL Drain Bag with Anti-Reflux Tower	Luer Lock	No Code	10/cs
URO1SPLDB18	18 Fr 10 mL, with 2,000 mL Drain Bag with Anti-Reflux Tower	Luer Lock	No Code	10/cs
URO1SPLUM14	14 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	Luer Lock	No Code	10/cs
URO1SPLUM16	16 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	Luer Lock	No Code	10/cs
URO1SPLUM18	18 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	Luer Lock	No Code	10/cs



2-Layer 100% Silicone Foley Catheter Trays

100% silicone trays feature a Foley catheter pre-connected to a Medline drain bag or urine meter. All the necessary components for proper Foley catheter insertion are included.

- » 100% silicone Foley catheter designed for maximum drainage
- » Foley catheter has a symmetrical balloon for reliability
- » Not made with natural rubber latex

Item No.	Description	Sample Port Type	HCPCS	Pkg.
DYND11519	16 Fr 10 mL, with 2,000 mL Drain Bag with Anti-Reflux Device	Blunt Cannula	A4315	10/cs
DYND11520	18 Fr 10 mL, with 2,000 mL Drain Bag with Anti-Reflux Device	Blunt Cannula	A4315	10/cs
URO1SPSDB14	14 Fr 10 mL, with 2,000 mL Drain Bag with Anti-Reflux Tower	Luer Lock	No Code	10/cs
DYND11523	16 Fr 10 mL, with 2,000 mL Drain Bag with Anti-Reflux Tower	Luer Lock	A4315	10/cs
DYND11528	18 Fr 10 mL, with 2,000 mL Drain Bag with Anti-Reflux Tower	Luer Lock	A4315	10/cs
URO1SPSUM14	14 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	Luer Lock	No Code	10/cs
DYND11524	16 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	Luer Lock	A4315	10/cs
DYND11529	18 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	Luer Lock	No Code	10/cs
DYND1516SIL	16 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	Blunt Cannula	A4315	10/cs
DYND1518SIL	18 Fr 10 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	Blunt Cannula	A4315	10/cs



DYND11NS14S

Non-Sterile Foley Training Kits

- » Pre-connected Foley tray specifically for simulation use
- » Staining PVP swab sticks replaced by non-staining alcohol swab sticks to mimic aseptic technique
- » 100% silicone catheter for ease of use with manikins
- » Not for patient use

Item No.	Description	HCPCS	Pkg.
DYND11NS12S	12 Fr Silicone Foley Tray	No Code	10/cs
DYND11NS14S	14 Fr Silicone Foley Tray	No Code	10/cs

Foley Catheter Insertion Trays

Medline Foley Catheter Insertion Trays provide the freedom and flexibility to choose the drain bag and catheter to fit the needs of each patient.

All Kits Include:

- » Graduated basin
- » Gloves not made with natural rubber latex
- » Underpad
- » Fenestrated drape
- » Sterile water inflation syringe
- » Lubricating jelly packet
- » All insertion trays are sterile



DYND10160

Item No.	Peel-Lid Tray	CSR Wrapped Tray	10 mL Syringe	30 mL Syringe	PVP Swabstick 3/pk	BZK Swabstick 3/pk	Aloetouch 3G Gloves	Vinyl Gloves	Specimen Container	HCPCS	Pkg.
DYNC1810	•		•		•			•		A4310	20/cs
DYNC1814	•			•		•		•		A4310	20/cs
DYNC1815	•			•	•			•		A4310	20/cs
DYND10160		•	•		•		•		•	A4310	20/cs
DYND10200		•		•	•		•			A4310	20/cs
DYND18100		•	•		•		•			A4310	20/cs
DYND10150*		•	•		•		•		•	A4311	20/cs

*DYND10150 also contains a 16 FR 10 mL Silicone-Elastomer Coated Latex Foley Catheter

BioCon™ 750 Bladder Scanner

Real-Time Bladder Targeting

Quickly and confidently locate the bladder for accurate measurement.

You shouldn't have to do multiple scans to get an acceptable result. BioCon 750 makes it easier than ever to find the bladder with two industry-firsts:

- » Proprietary BladderPoint™ technology puts a bull's-eye in the middle of the bladder
- » Pre-scan shows a live view of the ultrasound image right from the start

Dependable Performance

Reliability is designed into every BioCon 750.

- » Standard five-year warranty includes drop protection
- » Built-in self-diagnostics alert you of any hardware or maintenance issues
- » Tough polycarbonate casing promotes robust and uninterrupted service

Versatile Care

Accommodates various patient populations.

- » One-click, 360-degree measurement requires no probe rocking
- » Probe has a rear grip for easier use on bariatric patients
- » Can be used on pediatric patients

Staff Convenience

With BioCon 750, bladder scanning is quick, simple and efficient.

- » Data and images can be easily uploaded to any EMR system
- » Integrated thermal printer provides documentation at the push of a button
- » Onboard training guide walks staff through scanning technique



New aluminum shell for added probe protection

Item No.	Description	Pkg.
MDSB10750	BioCon 750 Bladder Scanner w/Printer	1/ea
MDSB10700STND	Mobile Stand w/Basket	1/ea
MDSB10750CK	Calibration Kit	1/ea
MDSBSPHANTOM	Phantom Bladder	1/ea
MDSB10500PAPE	Thermal Paper	23/cs
MDS092005	Ultrasound Gel, 8.5-oz. Squeeze Bottle	12/cs

BioCon™ 900 Bladder Scanner

Real-Time Bladder Viewing


Easily locate the bladder position for accurate measurement. You shouldn't have to do multiple scans to get an acceptable result. BioCon 900 makes it easier to find the bladder with two industry-firsts:

- » BladderPoint™ technology puts a bull's-eye in the middle of the bladder
- » Pre-scan shows a live view of the ultrasound image

Built-In Versatility

Bladder scanning is quick, simple and efficient.

- » Measure bladder volume in 3 seconds
- » Increased scan depth to industry best **230 mm** for individuals with higher BMIs
- » Data and images can be easily uploaded to your EMR system or printed from integrated printer
- » Less cleaning necessary for contact isolation patients



Dependable Performance

BioCon 900 features self-diagnosis and five-year warranty.

- » Durable probe construction helps minimize the chance for damage
- » **Five-year** warranty with drop protection coverage comes standard
- » Self-diagnosis quickly tells you if the probe is working correctly
- » No calibration required

Industry-Leading Support

We provide ongoing education, training and service throughout the device's lifespan.

- » Staff is trained on how to use the equipment through online videos, competency testing and onsite in-servicing
- » Toll-free technical service hotline: 1-866-274-2131
- » Service includes free equipment diagnosis and loaner units



Locking mechanism for safe storage and transport.



Item No.	Description	Pkg.
MDSB10900	BioCon 900 Bladder Scanner w/ Printer	1/ea
MDSMRC900	BioCon 900 Stand	1/ea
MDSBSPHANTOM	Phantom Bladder	1/ea
MDSB10500PAPE	Thermal Paper	12/ea, 12/cs
MDS092005	Ultrasound Gel, 8.5-oz. Squeeze Bottle	12/cs



FOLEY CATHETERS

- SelectSilicone Foley Catheters **27**
- Silvertouch® Foley Catheters **29**
- Simply Latex Foley Catheters **30**
- SelectSilicone Temperature-Sensing
Foley Catheters **32**

SelectSilicone Foley Catheters

SelectSilicone Foley Catheters

Maximize safety and comfort with our 100% Silicone Foley Catheters.

- » Completely inert for less tissue irritation and encrustation during extended periods of indwelling use
- » Easy insertion and less buckling with the bullet-shaped tip
- » Larger inner lumen for better drainage when compared to industry-leading Foley
- » Clear for easy visualization of clots, mucous and urine flow
- » Radiopaque for easy visualization during scans
- » Symmetrical balloon shape for improved reliability
- » Not made with natural rubber latex

Item No.	Size	Balloon Size	HCPCS	Pkg.
----------	------	--------------	-------	------

2-Way 100% Silicone Foley Catheters

DYND11552	6 Fr	3 mL	A4344	10/cs
DYND11553	8 Fr	5 mL	A4344	10/cs
DYND11554	10 Fr	5 mL	A4344	10/cs
DYND11500	12 Fr	10 mL	A4344	10/cs
DYND11501	14 Fr	10 mL	A4344	10/cs
DYND11502	16 Fr	10 mL	A4344	10/cs
DYND11503	18 Fr	10 mL	A4344	10/cs
DYND11504	20 Fr	10 mL	A4344	10/cs
DYND11505	22 Fr	10 mL	A4344	10/cs
DYND11506	24 Fr	10 mL	A4344	10/cs
DYND11532	16 Fr	30 mL	A4344	10/cs
DYND11533	18 Fr	30 mL	A4344	10/cs
DYND11534	20 Fr	30 mL	A4344	10/cs
DYND11535	22 Fr	30 mL	A4344	10/cs
DYND11536	24 Fr	30 mL	A4344	10/cs



Symmetrical balloon shape helps with maintaining proper retention.


SelectSilicone Foley Catheters



Pediatric SelectSilicone Foley Catheter + Syringe Kits

Pediatric-sized Foley catheters packaged with the proper amount of sterile water for proper balloon inflation and to prevent over-inflation.

Item No.	Size	Pkg.
2-Way 100% Silicone Foley + Inflation Syringes		
URO2S0615K	6 Fr 1.5 mL Foley, 3 mL Sterile Water Syringe	10/cs
URO2S0803K	8 Fr 3 mL Foley, 5 mL Sterile Water Syringe	10/cs
URO2S1003K	10 Fr 3 mL Foley, 5 mL Sterile Water Syringe	10/cs



Specialty Tiemann/Coude-Tip SelectSilicone Foley Catheters

For difficult insertions, coude-tip Foley catheters feature a bent tip that maneuvers around obstructions leading to the bladder.

- » Large oval-shaped eye for maximum drainage
- » Smooth reinforced tip resists buckling
- » Completely inert material causes less tissue irritation and encrustation during extended periods of indwelling use
- » Not made with natural rubber latex

Item No.	Size	Balloon Size	HCPCS	Pkg.
2-Way Coude 100% Silicone Foley Catheters				
DYND11591	14 Fr	10 mL	A4340	10/cs
DYND11592	16 Fr	10 mL	A4340	10/cs
DYND11593	18 Fr	10 mL	A4340	10/cs
DYND11594	20 Fr	10 mL	A4340	10/cs



3-Way SelectSilicone Foley Catheters

The 3-way Foley catheter adds continuous irrigation capabilities for use after bladder or other urinary surgeries.

- » Medications or other solutions can be continuously introduced into the bladder
- » Large inner drainage lumen – 30% larger than those on latex catheters
- » Radiopaque tip for ease of visualization during scans
- » If irrigation is discontinued, a catheter plug may be used to close the lumen (see DYND12200 on page 36)
- » Not made with natural rubber latex

Item No.	Size	Balloon Size	HCPCS	Pkg.
3-Way 100% Silicone Foley Catheters				
DYND11572	16 Fr	30 mL	A4346	10/cs
DYND11573	18 Fr	30 mL	A4346	10/cs
DYND11574	20 Fr	30 mL	A4346	10/cs
DYND11575	22 Fr	30 mL	A4346	10/cs
DYND11576	24 Fr	30 mL	A4346	10/cs


Silvertouch Foley Catheters

Silver Ions Provide Antimicrobial Action

Silver ions have long been recognized for their broad-spectrum antimicrobial action. Ionic silver is effective against gram-positive and gram-negative bacteria, as well as fungi. Ionic silver hydrophilic catheter coatings can be very effective against the most common and serious UTI pathogens, including *MRSA*, *VRE*, *E. coli*, *Candida albicans*, *Staph aureus*, *Pseudomonas*, *Klebsiella*, and others (see chart below). Silvertouch catheters incorporate the power of silver through a patented UV bonding process that binds silver ions to the catheter's lubricious coating.



The Zone of Inhibition (ZOI) diameter is a measure of the effectiveness of a compound. The larger the clear area, the more effective the compound.



Zone of Inhibition Study Results

	<i>E. COLI</i>	<i>P. AERUGINOSA</i>	<i>K. PNEUMONIAE</i>	<i>C. ALBICANS</i>	<i>MRSA</i>	<i>E. CLOACAE</i>	<i>C. KOSERI</i>	<i>VRE</i>	<i>E. FAECALIS</i>	<i>S. SAPROPHYTICUS</i>
Silvertouch avg. ZOI (mm)	5.13	5.80	4.93	4.93	25.13	5.07	5.00	2.00	2.00	7.73

Silvertouch Foley Catheters

A patented process permanently binds silver ions to the catheter's inner and outer surfaces. The catheter is completely inert and made of 100% silicone, so it poses no risk of allergic reaction to latex.

- » Pure ionic silver provides powerful antimicrobial protection against gram-positive and gram-negative bacteria, as well as fungi
- » Bonded lubricious hydrophilic coating hydrates rapidly and improves patient comfort over time
- » Not made with natural rubber latex


Item No.	Size	Balloon Size	HCPCS	Pkg.
----------	------	--------------	-------	------

2-Way Silvertouch 100% Silicone Foley Catheters

DYND141014	14 Fr	10 mL	A4344	10/cs
DYND141016	16 Fr	10 mL	A4344	10/cs
DYND141018	18 Fr	10 mL	A4344	10/cs
DYND141020	20 Fr	10 mL	A4344	10/cs
DYND141216	16 Fr	30 mL	A4344	10/cs
DYND141218	18 Fr	30 mL	A4344	10/cs
DYND141220	20 Fr	30 mL	A4344	10/cs



Latex Foley Catheters



Symmetrical balloon for added reliability and proper placement in the bladder.

Silicone-Elastomer Coated Latex Foley Catheters

Our innovative manufacturing process ensures the maximum amount of silicone coating available.

- » Silicone-elastomer coating resists encrustation
- » Large, oval-shaped eyes provide maximum drainage
- » Catheter feel is soft, flexible and conformable

Item No.	Size	Balloon Size	HCPCS	Pkg.
----------	------	--------------	-------	------

2-Way Silicone-Elastomer Coated Latex Foley Catheters

DYND11706*	6 Fr	3 mL	A4338	12/cs
DYND11708*	8 Fr	5 mL	A4338	12/cs
DYND11710*	10 Fr	5 mL	A4338	12/cs
DYND11752	12 Fr	10 mL	A4338	12/cs
DYND11754	14 Fr	10 mL	A4338	12/cs
DYND11756	16 Fr	10 mL	A4338	12/cs
DYND11758	18 Fr	10 mL	A4338	12/cs
DYND11760	20 Fr	10 mL	A4338	12/cs
DYND11762	22 Fr	10 mL	A4338	12/cs
DYND11764	24 Fr	10 mL	A4338	12/cs
DYND11766	26 Fr	10 mL	A4338	12/cs
DYND11768	28 Fr	10 mL	A4338	12/cs
DYND11770	30 Fr	10 mL	A4338	12/cs
DYND11772	12 Fr	30 mL	A4338	12/cs
DYND11774	14 Fr	30 mL	A4338	12/cs
DYND11776	16 Fr	30 mL	A4338	12/cs
DYND11778	18 Fr	30 mL	A4338	12/cs
DYND11780	20 Fr	30 mL	A4338	12/cs
DYND11782	22 Fr	30 mL	A4338	12/cs
DYND11784	24 Fr	30 mL	A4338	12/cs
DYND11786	26 Fr	30 mL	A4338	12/cs
DYND11788	28 Fr	30 mL	A4338	12/cs
DYND11790	30 Fr	30 mL	A4338	12/cs

*Includes Stylet

SimplyLatex Foley Catheters

- » 17-inch catheter for added maneuverability
- » Silicone-elastomer coating resists encrustation
- » Symmetrical balloon provides added reliability
- » Large, oval-shaped eyes provide maximum drainage

Item No.	Size	Balloon Size	HCPCS	Pkg.
2-Way Silicone-Elastomer Coated Latex Foley Catheters, 17" (43.2 cm)				
URO2L143	14 Fr	10 mL	A4338	12/cs
URO2L163	16 Fr	10 mL	A4338	12/cs
URO2L183	18 Fr	10 mL	A4338	12/cs
URO2L203	20 Fr	10 mL	A4338	12/cs
URO2L223	22 Fr	10 mL	A4338	12/cs
URO2L243	24 Fr	10 mL	A4338	12/cs



Specialty Coude-Tip Silicone-Elastomer Coated Latex Foley Catheters

For difficult insertions, coude-tip Foley catheters feature a bent tip that maneuvers around obstructions leading to the bladder.

- » Silicone-elastomer coating resists encrustation
- » Symmetrical balloon provides added reliability
- » Large, oval-shaped eyes provide maximum drainage

Item No.	Size	Balloon Size	HCPCS	Pkg.
2-Way Coude-Tip Silicone-Elastomer Coated Latex Foley Catheters				
DYND11214	14 Fr	10 mL	A4340	12/cs
DYND11216	16 Fr	10 mL	A4340	12/cs
DYND11218	18 Fr	10 mL	A4340	12/cs
DYND11220	20 Fr	10 mL	A4340	12/cs



3-Way Silicone-Elastomer Coated Latex Foley Catheters

The 3-way Foley catheter adds continuous irrigation capabilities for use after bladder or other urinary surgeries.

- » Medications or other solutions can be continuously introduced into the bladder
- » If irrigation is discontinued, a catheter plug may be used to close the lumen (see DYND12200 on page 36)

Item No.	Size	Balloon Size	HCPCS	Pkg.
3-Way Silicone-Elastomer Coated Latex Foley Catheters				
DYND11800	16 Fr	30 mL	A4346	12/cs
DYND11801	18 Fr	30 mL	A4346	12/cs
DYND11802	20 Fr	30 mL	A4346	12/cs
DYND11803	22 Fr	30 mL	A4346	12/cs
DYND11804	24 Fr	30 mL	A4346	12/cs
DYND11805	26 Fr	30 mL	A4346	12/cs



SelectSilicone Temperature-Sensing Foley Catheters

ACCURATE

- » Continuous measurements with +/- 0.2°C accuracy

CONVENIENT

- » Saves nursing time monitoring temperature where Foley catheterization is already indicated

COMPATIBLE

- » Works with most hypothermia units, anesthesia machines and patient monitors from GE®, Philips™, Siemens Climate Solutions® and others



Improving Patient Outcomes

1

In the OR

Medline's Temperature-Sensing Foley Catheters can help you comply with SCIP-Inf-10 by monitoring the patient's core temperature with an integrated thermistor.

2

In Therapeutic Hypothermia

There are several types of therapeutic hypothermia, each of which occurs within a narrow temperature range¹; therefore, accurate and continuous monitoring of core temperature is critical. Bladder temperature probes enable continuous monitoring and feedback in conjunction with hypothermia devices. Other techniques, including intermittent tympanic temperature measurements, are less reliable.²

3

In the ICU

The ability to use a temperature-sensing catheter to measure both urine output and continuous body temperature is a distinct advantage, which may eliminate the need to measure temperature at other sites.³

References: 1. Polderman K and Herold I. Therapeutic hypothermia and controlled normothermia in the intensive care unit: practical considerations, side effects, and cooling methods. *Crit Care Med* 2009; 37 (3); 1101-1119. 2. Nolan JP, Morley PT, Vanden Hoek TL, Hickey RW et al. Therapeutic hypothermia after cardiac arrest: an advisory statement by the advanced life support task force of the International Liaison Committee on Resuscitation. *Circulation*. 2003 Jul 8;108(1):118-21. 3. Fallis WM. Monitoring urinary bladder temperature in the intensive care unit: state of the science. *Am J Crit Care*. 2002 Jan;11(1):38-45.

Temperature-Sensing Foley Catheters

400-Series temperature-sensing Foley catheters help with compliance with SCIP-Inf-10. The 100% silicone catheter has an integrated thermistor that aids in monitoring the body's core temperature via the bladder. The connection port has a molded connector for a sealed, moisture-resistant, one-way fit.

- » 8-inch (20 cm) side arm – a sealed electrical connection away from the patient
- » Smooth tip – easy insertion without buckling
- » Accuracy within +/- 0.2°C
- » Fully compatible with most patient monitors
- » Conditionally MRI compatible
- » Not made with natural rubber latex

Item No.	Size	Balloon Size	HCPCS	Pkg.
----------	------	--------------	-------	------

400-Series Temperature-Sensing Foley Catheters

DYND110508	8 Fr	3 mL	A4344	10/cs
DYND110510	10 Fr	3 mL	A4344	10/cs
DYND110512	12 Fr	10 mL	A4344	10/cs
DYND110514	14 Fr	10 mL	A4344	10/cs
DYND110516	16 Fr	10 mL	A4344	10/cs
DYND110518	18 Fr	10 mL	A4344	10/cs



Pediatric SelectSilicone Temperature-Sensing Foley Catheters + Syringe Kits

Pediatric-sized Foley catheters packaged with the proper amount of sterile water for proper balloon inflation and to prevent overinflation.

Item No.	Size	Pkg.
----------	------	------

2-Way 100% Silicone Foley Catheter + Inflation Syringes

URO2S0615TK	6 Fr 1.5 mL Temperature-Sensing Foley, 3 mL Sterile Water Syringe	10/cs
URO2S0803TK	8 Fr 3 mL Temperature-Sensing Foley, 5 mL Sterile Water Syringe	10/cs
URO2S1003TK	10 Fr 3 mL Temperature-Sensing Foley, 5 mL Sterile Water Syringe	10/cs



Temperature-Sensing Cables

Interconnecting cables are made of a durable barrier to minimize interference.

- » Designed for compatibility with most monitors and anesthesia/hypothermia machines
- » All cables are 12-feet (3.7 m) long
- » Not made with natural rubber latex

Item No.	Description	Pkg.
DYND110101	¼" (0.6 cm) Adapter	10/cs
DYND110102	Philips Adapter	10/cs
DYND110103	Siemens Adapter	10/cs



¼" (0.6 cm) Adapter



Philips Adapter



Siemens Adapter



URINE COLLECTION

- Urine Meters **35**
- Drain Bags **36**
- Leg Bags **39**
- Accessories **40**

Urine Meters



DYND17403

DYND17407



DYND27403S

Urine Meters

Medline urine meters offer the lowest cost method of accurately measuring urine output. Each is packaged OR sterile in a Tyvek® peel pouch.

- » Adjustable bed-sheet clamp for securing drain bag tubing
- » Includes 0.31-inch (0.79 cm) inner diameter tubing with either 60-inch or 48-inch (152.4 or 122 cm) length
- » Drain bag capacity: 2,500 mL or 2,000 mL
- » Typically used on critically ill and surgical patients
- » Needle-less sample port, either luer lock or blunt cannula access
- » Not made with natural rubber latex

Item No.	Description	Tubing	Sample Port	HCPCS	Pkg.
DYND17403	400 mL Urine Meter with 2,500 mL Drain Bag, Metal Clamp	60" (152.4 cm)	Blunt Cannula Tip Access	A4357	10/cs
DYND27403	400 mL Urine Meter with 2,500 mL Drain Bag, Metal Clamp	60" (152.4 cm)	Luer Lock	A4357	10/cs
DYND27403S*	400 mL Urine Meter with 2,500 mL Drain Bag, Slide Tap	60" (152.4 cm)	Luer Lock	A4357	10/cs
DYND17407*	200 mL Urine Meter with 2,000 mL Drain Bag, Slide Tap	48" (122 cm)	Luer Lock	A4357	10/cs
DYND17407L*	200 mL Urine Meter with 2,000 mL Drain Bag, Slide Tap	60" (152.4 cm)	Luer Lock	A4357	10/cs
URO27403**	400 mL Urine Meter with 2,500 mL Drain Bag, Vented Tubing	60" (152.4 cm)	Luer Lock	A4357	10/cs
URO27403S**	400 mL Urine Meter with 2,500 mL Drain Bag, Slide Tap, Vented Tubing	60" (152.4 cm)	Luer Lock	A4357	10/cs

*Metal-free product is safe for MRI use. **Includes vented tubing

Drain Bags

Only Medline drain bags have a slide tap for easy one-handed drainage. The needleless sample port was first designed by Medline and provides a safe way to take urine samples. Star tubing (50-inch (127 cm) length, 0.281-inch (0.714 cm) inner diameter) facilitates drainage of thick, fibrous urine and resists kinking.

- » Reinforced hanger with built-in anti-kink drain tube guide
- » Adjustable bed sheet clamp for securing drain bag tubing
- » Capacity: 2,000 mL or 4,000 mL
- » Includes needle free sample port, either luer lock or blunt cannula access
- » Not made with natural rubber latex



DYND25205



Anti-reflux tower minimizes backflow and increases patient safety.



DYND15405

Drain Bags with Anti-Reflux Tower

The urinary drain bag with anti-reflux tower provides reflux protection without mechanical parts to hinder flow.

- » Non-wettable air filter to help prevent vacuum action and facilitate drainage
- » 2,000 mL capacity for standard use or 4,000 mL capacity for use with continuous bladder irrigation

Item No.	Description	Sterility	Sample Port	HCPCS	Pkg.
DYND15205	2,000 mL with Slide Tap Drainage Port	Sterile Fluid Pathway	Blunt Cannula Tip	A4357	20/cs
DYND25205 [†]	2,000 mL with Slide Tap Drainage Port	OR Sterile	Luer Lock	A4357	20/cs
DYND15203	2,000 mL with Metal Clamp Drainage Port	Sterile Fluid Pathway	Blunt Cannula Tip	A4357	20/cs
DYND15405	4,000 mL with Metal Clamp Drainage Port	OR Sterile	Blunt Cannula Tip	A4357	20/cs
DYND25405S [†]	4,000 mL with Slide Tap Drainage Port	OR Sterile	Luer Lock	A4357	20/cs
URO25205V ^{**}	2,000 mL Vented Urinary Drainage Bag	OR Sterile	Luer Lock	A4357	20/cs

*Includes Foley Insetag †Metal free product is safe for MRI use



Anti-reflux control with a 1-way mylar valve.



Drain Bags with Anti-Reflux Device

The urinary drain bag with anti-reflux device has a one-way clear Mylar® valve for visualization of urine flow.

- » Non-wettable air filter to help prevent vacuum action and facilitate drainage
- » Foley InserTag® included for easy tracking of catheter days

Item No.	Description	Sterility	Sample Port	HCPCS	Pkg.
DYND15207	2,000 mL with Slide Tap Drainage Port	Sterile Fluid Pathway	Blunt Cannula Tip	A4357	20/cs
DYND15210	2,000 mL with Metal Clamp Drainage Port	OR Sterile	Blunt Cannula Tip	A4357	20/cs



URO51211C

Pre-Covered Drain Bag

This drain bag with anti-reflux tower comes pre-covered, with patient dignity in mind.

- » All-in-one system makes compliance easy
- » Cover conceals both front and back of drainage bag
- » String tie allows for hanging on a variety of beds
- » Cover is open at the bottom for easy drainage
- » Button closure allows for cover replacement if soiled

Item No.	Description	Sterility	Sample Port	HCPCS	Pkg.
URO51211C	2,000 mL with Slide Tap Drainage Port	Sterile Fluid Pathway	Blunt Cannula Tip	A4357	20/cs

Drain Bags



Slide-tap mechanism allows for easy, one-handed bag drainage.



DYNC1674



UR053211

Drain Bag with Anti-Reflux Valve

The drain bag includes a vinyl flap system that minimizes back flow into the drainage tubing.

- » Allows free flow of fluid into the drain bag
- » Exclusive slide-tap drainage

Item No.	Description	Sterility	Sample Port	HCPCS	Pkg.
DYNC1674	2,000 mL with Slide Tap Drainage Port	Sterile Fluid Pathway	Blunt Cannula Tip	A4357	20/cs
UR053211	2,000 mL with Slide Tap Drainage Port	Sterile Fluid Pathway	Blunt Cannula Tip	A4357	50/cs

See pg. 41 for Medline's exclusive urinary drain bag covers, designed with patient dignity in mind.

Low-Bed Drain Bag

The patented (US 7,846,142 and 8,430,855) drain bag hangs only 6-inches (18 cm) to allow low placement on an adjustable bed without touching the floor.

- » Anti-reflux valve for back-flow prevention
- » 50-inches (127 cm) of tubing to facilitate proper drainage
- » Exclusive slide-tap drainage port allows for minimal splashing during drainage
- » Reversible design for versatile placement on either side of the bed



DYND15215

Item No.	Description	Sterility	Sample Port	HCPCS	Pkg.
DYND15215	1,200 mL with Slide Tap Drainage Port	Sterile Fluid Pathway	Blunt Cannula Tip	A4357	20/cs

Leg Bags



URO5454



DYND12574



DYND12588

Leg Bags with Pre-Attached Straps

Leg bag can be used for those patients that require a catheter, but are able to ambulate on their own. Leg bags available with twist valve, slide-tap or flip drain.

- » Twist valve for easy emptying and minimal spillage
- » Slide tap allows for exceptional ease of use
- » Flip drain is easy to maneuver for limited dexterity users
- » Not made with natural rubber latex
- » Elastic straps come preconnected to each bag
- » Built in anti-reflux valve in each bag helps prevent backflow
- » Sterile fluid pathway

Item No.	Description	Size	HCPCS	Pkg.
DYND12574	Twist-Valve Drainage Port	Medium, 20-oz. (600 mL)	A4358	48/cs
DYND12575	Twist-Valve Drainage Port with Extension Tubing, 18" (46 cm)	Medium, 20-oz. (600 mL)	A4358	48/cs
DYND12578	Twist Valve Drainage Port	Large, 32-oz. (900 mL)	A4358	48/cs
DYND12584	Fabric Back, Exclusive Slide-Tap Drainage Port	Medium, 20-oz. (600 mL)	A4358	48/cs
DYND12588	Fabric Back, Exclusive Slide-Tap Drainage Port	Large, 32-oz. (900 mL)	A4358	48/cs
URO5454	Fabric Back, Flip Drain Spout	Medium, 20-oz. (600 mL)	A4358	48/cs
URO5458	Fabric Back, Flip Drain Spout	Large, 32-oz. (900 mL)	A4358	48/cs

Accessories



DYND12200



DYND12582

DYND12582
(Shown here on DYND12588)

Accessories

- » The plug for catheter drainage protector tube fits securely into the funnel end of most Foley catheters and the included cover protects the drainage bag connector
- » Leg-bag extension tubing can be cut to length to meet most patient needs
- » Leg bag straps with hook and loop fasteners offer increased comfort and durability and are compatible with all Medline leg bags
- » All items listed are not made with natural rubber latex

Item No.	Description	HCPCS	Pkg.
DYND12550	Leg-Bag Extension Tubing, 18" (46 cm), Sterile	A4331	50/cs
DYND12582	Leg Bag Straps	A5114	20 pairs/cs
DYND12200	Plug for Catheter Drainage Protector Tube, Sterile	A9270	100/cs

Drain Bag Covers

Improve patient dignity with Medline's patented (US 8,348,914) urinary drain bag covers.

- » Blue nonwoven material conceals fluid in the drain bag
- » Two closure devices make it easy to affix to a urinary drain bag
- » Design allows for easy-access drainage
- » Design fits most urinary drain bags
- » Low-bed design fits Medline's patented low bed drain bag (DYND15215)
- » Not made with natural rubber latex

Item No.	Description	Pkg.
DYND15200	Urinary Drain Bag Cover	20/cs
DYND15200LOW	Low Bed Drain Bag Cover	20/cs



DYND15200



DYND15200LOW

Drain Bag Holder

- » Made from durable poly/cotton
- » Attaches to a wheel chair frame with 100% polyester straps
- » Three separate snap positions for easy adjustment

Item No.	Description	Pkg.
MDT825150	Drain Bag Holder	1/ea



MDT825150

IV Pole with Drain Bag Hook

- » Tool-Free Urine Bag Hook
- » For use with Medline's IV Pole, item MDS80600

Item No.	Description	HCPCS	Pkg.
MDS80600UBH	Hook for Urine Bag on IV Pole	No Code	10/cs
MDS80600	IV Pole with Quick Release Casters	E0776	2/cs



Tool-free Urine Bag Hook (MDS80600UBH), attaches easily to IV Pole (MDS80600) with a quick release lever.





CATHETER SECUREMENT DEVICES

Foley Catheter Securement 43

Foley Catheter Securement

Adhesive Securement Devices

Keep the catheter securely in place and improve patient safety with our easy-to-use Foley catheter tubing securement devices.

- » Adhesive patch with adjustable hook-and-loop securement prevents unnecessary pulls
- » UROSECURE is packaged with alcohol prep pad and Sureprep® skin prep pad for proper skin preparation before application
- » Pediatric securement device was designed with the patient in mind, incorporating Buddy the Brave, the size of the device is just right for your pediatric patients
- » Packaged individually

Item No.	Description	Size	HCPCS	Pkg.
SAFESECURE	SafeSecure Foley Catheter Securement Device	3 x 4" (7.6 x 10.2 cm)	No Code	50/cs
UROSECURE	Foley Catheter Securement Device	6 x 1.75" (15.2 x 4.5 cm)	A4333	100/cs
DYND7600P	Pediatric Adhesive Securement Device	1.375 x 3.5" (3.5 x 8.9 cm)	A4333	100/cs



SAFESECURE



UROSECURE



DYND7600P

Centurion Adhesive Securement Devices

Item No.	Description	Size	HCPCS	Pkg.
FCS200XT	Centurion Adhesive Securement Device	3 x 4.5" (7.6 x 11.4 cm)	No Code	50/cs
FSC300XT*	Centurion Adhesive Securement Device	3 x 4.5" (7.6 x 11.4 cm)	A4333	100/cs

*Maximum Securement

Adhesive Patch

Medical grade adhesive patch with hook-and-loop holder for properly securing the Foley catheter.

- » Safe and secure method of keeping the catheter in place
- » Non-sterile

Item No.	Description	Size	HCPCS	Pkg.
DYND16700	Adhesive Patch	2 x 3" (5.1 x 7.6 cm)	A4333	50/cs



FSC300XT



DYND16700

Leg Straps

Keep the catheter securely in place by utilizing a leg strap. Available in multiple sizes depending on patient needs.

- » Elastic comfort strap stabilizes the catheter without restricting patient movement
- » Foam catheter strap holds catheter securely in place
- » Hook-and-loop locking system engages quickly and easily
- » Non-sterile

Item No.	Description	Sterility	HCPCS	Pkg.
DYND16800	Adjustable Foam*	29 x 1.5" (73.7 x 3.8 cm)	A4334	12/bx
DYND16900	Elastic*	22 x 2" (55.9 x 5.1 cm)	A4334	10/bx
DYND16905	Elastic**	51 x 2" (129.5 x 5.1 cm)	A4334	10/bx

*Fits Leg Circumferences up to 29-inch (73.7 cm) **Fits Leg Circumferences up to 65-inch (165.1 cm)



DYND16900



ALTERNATIVES TO FOLEY CATHETERS

- Intermittent Catheters 45
- Intermittent Catheterization T rays 47
- Male External Catheters 50
- Speci-Cath Kits 52

Intermittent Catheters

Medline intermittent catheters are made with care and precision. Catheters are offered in a wide selection of materials, configurations and French sizes. Available in female, pediatric and universal lengths.

- » Smooth tip for easy insertion and comfort
- » Staggered eyes allow for maximum flow
- » Funnel end for ease of handling
- » Color coded to indicate French size
- » All sizes are packaged sterile



Smooth eyes minimize urethral trauma. Clear vinyl intermittent catheter shown (item no. DYND10702).



6-inch (15.2 cm) Female Catheter available in 10, 12 and 14 Fr



10-inch (25.4 cm) Pediatric Catheter available in 6, 8 and 10 Fr



16-inch (40.6 cm) Universal Catheter available in 10, 12, 14, 16 and 18Fr

Clear Vinyl Intermittent Catheters

- » Not made with natural rubber latex

Item No.	Size	HCPCS	Pkg.
Female, 6" (15.2 cm)			
DYND10701	10 Fr	A4351	30/cs
DYND10702	12 Fr	A4351	30/cs
DYND10703	14 Fr	A4351	30/cs
Pediatric, 10" (25.4 cm)			
DYND10712	6 Fr	A4351	30/cs
DYND10713	8 Fr	A4351	30/cs
DYND10714	10 Fr	A4351	30/cs
Universal, 16" (40.6 cm)			
DYND10721	10 Fr	A4351	30/cs
DYND10722	12 Fr	A4351	30/cs
DYND10723	14 Fr	A4351	30/cs
DYND10700	14 Fr	A4351	100/cs
DYND10724	16 Fr	A4351	30/cs
DYND10725	18 Fr	A4351	30/cs

Vinyl Hydrophilic Intermittent Catheters

Hydrophilic coating minimizes friction during catheterization to help prevent discomfort and trauma.


- » Hydrophilic coating is activated with water to create a slick surface along the length of the catheter
- » Adhesive tab on packaging allows the user to attach packaging to a surface, keeping hands free during catheter prep
- » 16-inch (40.6 cm) length
- » Not made with natural rubber latex

Item No.	Size	HCPCS	Pkg.
Universal			
URO3PVC1016L	10 Fr	A4351	30/cs
URO3PVC1216L	12 Fr	A4351	30/cs
URO3PVC1416L	14 Fr	A4351	30/cs
URO3PVC1616L	16 Fr	A4351	30/cs
URO3PVC1816L	18 Fr	A4351	30/cs

Intermittent Catheters



URO3PVC1016C



Coude Tip

Vinyl Coude Intermittent Catheters

- » Featuring a curved, coude, tip for help in navigating past strictures or an enlarged prostate
- » 16-inch (40.6 cm) length
- » Not made with natural rubber latex

Item No.	Size	HCPCS	Pkg.
Universal			
URO3PVC1016C	10 Fr	A4352	30/cs
URO3PVC1216C	12 Fr	A4352	30/cs
URO3PVC1416C	14 Fr	A4352	30/cs
URO3PVC1616C	16 Fr	A4352	30/cs
URO3PVC1816C	18 Fr	A4352	30/cs



DYND13514



Universal



Coude Tip

Red Rubber Latex All-Purpose Intermittent Catheters

Urethral catheter used across the spectrum of health care for all patient types. High quality construction ensures consistency of use and patient comfort.

- » Smooth tip and funnel end made from a radiopaque material
- » Two large opposing eyes for optimum drainage
- » Suitable for male or female catheterization
- » Coude tip available
- » 16-inch (40.6 cm) length
- » Contains natural rubber latex

Item No.	Size	HCPCS	Pkg.
Universal			
DYND13508	8 Fr	A4351	12/cs
DYND13510	10 Fr	A4351	12/cs
DYND13512	12 Fr	A4351	12/cs
DYND13514	14 Fr	A4351	12/cs
DYND13515	15 Fr	A4351	50/cs
DYND13516	16 Fr	A4351	12/cs
DYND13518	18 Fr	A4351	12/cs
DYND13520	20 Fr	A4351	12/cs
DYND13522	22 Fr	A4351	12/cs
Coude Tip			
DYND13614	14 Fr	A4352	12/cs

Intermittent Catheterization Trays

Medline intermittent catheterization trays provide you with all of the necessary supplies required for catheterization.

All Sterile Kits Include:

- » Graduated basin
- » Gloves not made with natural rubber latex
- » Underpad
- » Fenestrated drape
- » Swabsticks, 3 pack
- » Lubricating jelly
- » Specimen container with label



Item No.	Peel-Lid Tray	CSR Wrapped Tray	Red Rubber Catheter	Vinyl Catheter	PVP Swabsticks 3/pk	BZK Swabsticks 3/pk	Aloetouch 3G Gloves	Vinyl Gloves	HCPCS	Pkg.
DYNC1820	•			14 Fr	•			•	A4353	20/cs
DYNC1816	•		15 Fr		•			•	A4353	20/cs
DYND10300		•		14 Fr	•		•		A4353	20/cs
DYND10305	•		No Cath	No Cath	•		•		A4310	20/cs
DYND10350		•	15 Fr		•		•		A4353	20/cs
DYND10500		•		14 Fr		•	•		A4353	20/cs

Total One Layer Intermittent Catheter Tray

This tray contains everything necessary to complete an intermittent catheterization, with the added benefit of a plush perineal cleansing wipe and pre-connected collection bag. The tray comes CSR wrapped with components layered in order of use, promoting proper insertion and minimizing risk of contamination.

Sterile Kit Includes:

- » Perineal care wipe
- » Hand sanitizer
- » Vinyl gloves
- » Underbuttocks drape
- » Fenestrated drape
- » PVP swabsticks, 3 pack
- » Lubricating jelly
- » Specimen container with label
- » 14 Fr 16" (40.6 cm) vinyl or poly catheter
- » Pre-attached urine collection bag with push-pull drain valve
- » Not made with natural rubber latex



URO4PVCPC14

Item No.	Description	Pkg.
URO4PVCPC14	Total One Layer 14 Fr PVC Catheter Tray	20/cs
URO4POLYPC14	Total One Layer 14 Fr Poly Catheter Tray	20/cs

Intermittent Catheterization Trays



Pre-Connected Intermittent Catheter Trays

This tray offers all the necessary components in order to complete an intermittent catheterization, with the additional benefit of a pre attached 1,000 mL drain bag.

All Sterile Kits Include:

- » Wrapped gloves
- » Underpad
- » Fenestrated drape
- » PVP swabsticks, 3 pack
- » Lubricating jelly
- » Specimen container with label

Item No.	Description	Size	HCPCS	Pkg.
DYND10400	Graduated Basin, Clear Vinyl Catheter*	14 Fr	A4353	20/cs
DYND10402	Compact Packaging, Clear Vinyl Catheter*	14 Fr	A4353	20/cs
DYND10407	Compact Packaging, Red Rubber Catheter	15 Fr	A4353	20/cs

*Not made with natural rubber latex.



DYND10440

My-Cath™ Touch-Free Self-Catheter System

The self-catheterization kit is ideal for active and independent patients. The closed system design reduces mess, eases disposal and reduces contamination risk by allowing a touch-free catheterization. The catheter and introducer are coated with a water soluble lubricant for ease of insertion.

All Sterile Kits Include:

- » Catheter guide – allows for simple, protected control of catheter during insertion
- » 1,500 mL collection bag – eliminates the need for a separate container, reducing the chance of spillage
- » Moisture-proof underpad
- » Wrapped gloves
- » PVP swabsticks, 3 pack
- » Not made with natural rubber latex

Item No.	Size	HCPCS	Pkg.
URO4PVCPC12MC	12 Fr	A4353	60/cs
DYND10440	14 Fr	A4353	60/cs
URO4PVCPC16MC	16 Fr	A4353	60/cs

Male External Catheters

Men's Liberty Acute by BioDerm

A revolutionary solution for acute male urinary management. With a completely external design, Men's Liberty Acute is a one size fits most. Men's Liberty Acute allows patients to avoid moisture related skin breakdown, keeping patients dry and comfortable.

- » Catheters have a completely external, noninvasive design that helps reduce the risk of infection associated with indwelling catheters
- » Direct urine away from the skin to avoid moisture-related skin breakdown while keeping patients dry and comfortable
- » Seven-petal faceplate and faceplate strip are easy to apply and form an occlusive seal around the urinary meatus for leak protection
- » Universal connection to any leg or bedside bag for convenience



BI02102610
BI0210268

Item No.	Description	Pkg.
Starter Kit		
BI02102610	Men's Liberty Acute with Cath Grip	35/bx
BI02102610H	Men's Liberty Acute with Cath Grip	1/ea
Replacement Kit		
BI0210268	Men's Liberty Acute	35/bx
BI0210268H	Men's Liberty Acute	1/ea

FreeDerm™ Adhesive Remover

FreeDerm is the proprietary adhesive remover formulated to instantly remove hydrocolloid, tapes and bandages without irritation and evaporates rapidly with no residue. It allows immediate re-application of an adhesive device without washing.

- » No-sting
- » Helps reduce skin trauma from tape and dressing removal including infection risk, inflammatory damage and edema
- » No harsh chemicals
- » Oil and lotion free



BOM522023

BOM52203BX

Item No.	Description	Pkg.
BOM52200	FreeDerm Adhesive Remover, 1-oz. Spray	36/bx
BOM522023	FreeDerm Adhesive Remover, 3-oz. Spray	63/bx
BOM52203BX	FreeDerm Adhesive Remover, Ampule	30/bx

BioPlus+™ Barrier Prep Wipes

BioPlus+ is the proprietary skin prep wipe that provides a barrier film layer directly on the skin to improve the adhesion of hydrocolloids, appliances, tapes and bandages while protecting skin from incontinence, bodily fluids and adhesive stripping.

- » No-sting
- » Forms transparent, protective barrier film
- » Improves adhesion by up to 50%
- » Dries in 15-20 seconds




Item No.	Description	Pkg.
BOM52300BX	BioPlus Wipe	25/bx

CathGrip® Securement Device

This universal securement device is made from proprietary hydrocolloid that is hypoallergenic. The skin-friendly adhesive helps prevent skin tearing and irritation while the patent pending soft straps are designed to reduce pressure injuries.

- » Available in single or dual straps in the large size
- » Hydrocolloid adhesive oval
- » Hypoallergenic
- » Soft straps with no-slip grip
- » Not made with natural rubber latex



BOM51300NSBX

Item No.	Description	Pkg.
BOM51250BX	Small, Single Strap, Sterile	30/bx
BOM51250NSBX	Small, Single Strap, Non-Sterile	30/bx
BOM51301BX	Medium, Double Strap, Sterile	30/bx
BOM51301NSBX	Medium, Double Strap, Non-Sterile	30/bx
BOM51350BX	Large, Single Strap, Sterile	30/bx
BOM51350NSBX	Large, Single Strap, Non-Sterile	30/bx
BOM51300BX	Large, Double Strap, Sterile	30/bx
BOM51300NSBX	Large, Double Strap, Non-Sterile	30/bx

Exo-Cath® Latex External Catheters

Flexible latex male external catheters include double-sided adhesive tape for proper securement.

- » Non-sterile



DYND12302

Item No.	Description	Size	HCPCS	Pkg.
DYND12303	Small	25 mm	A4349	25/cs
DYND12302	Medium	30 mm	A4349	25/cs
DYND12301	Large	35 mm	A4349	25/cs

Speci-Cath Kits

Medline's Speci-Cath kits provide an easy way to collect a sterile urine sample. The catheter is pre-inserted into the specimen tube to ensure sterile technique. Urine drains directly into the specimen tube.



Speci-Cath Kits

Kits available to fit your needs with a variety of options including gloves, PVP swabsticks and lubricating jelly.

Item No.	Size	Catheter Length	Vinyl 3G Gloves	PVP Swabsticks 3/pk	Lubricating Jelly	Drape	HCPCS	Pkg.
DYND10820	5 Fr	9" (23 cm)	•	•	•	•	A9270	25/cs
DYND10800	8 Fr	6.5" (16.5 cm)					A9270	25/cs
DYND10805	8 Fr	6.5" (16.5 cm)	•	•	•		A9270	25/cs
DYND10806	8 Fr	6.5" (16.5 cm)	•	•			A9270	25/cs
DYND10815	8 Fr	9" (23 cm)	•	•	•	•	A9270	25/cs
DYND10810	8 Fr	14" (35.6 cm)	•	•	•	•	A9270	25/cs



CKF105L

Centurion Urine Collection Kits

Centurion Urine Collection Kits create a completely closed system, from sampling to lab. The collection tube's flip-top cap features a vent that allows trapped air to escape without loosening the lid, enabling urine to flow into the tube quickly and easily, without spilling. Available in various sizes for female, male, and pediatric sampling.

Item No.	Size	Catheter Length	Vinyl 3G Gloves	PVP Swabsticks 3/pk	Lubricating Jelly	Drape	HCPCS	Pkg.
CKF155L	5 Fr	6" (15.2 cm)	•	•	•		A9270	25/cs
CKF5040L	5 Fr	9.5" (24.1 cm)	•	•	•		A9270	25/cs
CKF105	8 Fr	6" (15.2 cm)	•	•			A9270	25/cs
CKF105L	8 Fr	9" (23 cm)	•	•	•		A9270	25/cs



PEDIATRICS

Buddy the Brave® One-Layer Foley Catheter Trays **55**
 Foley Catheter Trays **56**
 Intermittent Catheters **60**
 Accessories **61**

Buddy the Brave ERASE CAUTI Trays

Introducing Medline's new patented Pediatric Catheter Tray* the latest addition to the innovative ERASE CAUTI Comprehensive Care Solution.

Sometimes, you just need a buddy. Buddy the Brave lion cub is here to help your youngest catheter patients. You'll find fun, upbeat activities and a bravery award sticker in every tray, along with some serious patient (and parent/guardian) educational resources.

But it's more than just fun. There's published evidence that distraction helps children tolerate unpleasant procedures more so than adult reassurance does; the enclosed bravery sticker and activity card are just the right distraction.

*For a current list of patent numbers, visit www.medline.com/patents.

buddy the brave®



Pediatric Catheter Tray



Bravery Sticker

Children's Activities



Foley Catheter Trays

Buddy the Brave is on the label, so you know it's kid-friendly. Sharp photography of the kit contents lets you know what's inside, and the illustrated instructions can be reviewed without opening the pack. The one-layer tray design holds all the equipment you need for a procedure in a logical, "best practices" order.

All Sterile Trays Include:

- » Aloetouch® 3G gloves
- » Foley InserTag®
- » Hand antiseptic packet
- » Underbuttocks drape
- » Fenestrated drape
- » Lubricating jelly syringe
- » Sterile water syringe
- » Specimen container
- » 3-pack PVP swabsticks
- » Pediatric securement device
- » Sureprep® pad
- » Alcohol prep pad



DYND160808

Pediatric One Layer SelectSilicone Foley Catheter Trays

- » Correct amount of sterile fluid required to fill the Foley included in the kit to help reduce the risk of over-inflation.
- » Features Medline's One Layer design to help ensure best practices are followed
- » Buddy the Brave is on the label, so you know it's kid-friendly
- » Pediatric 200 mL urine meter connected to a 2,000 mL drain bag with metal-free drainage tube
- » Not made with natural rubber latex

Item No.	Size	HCPCS	Pkg.
Buddy the Brave SelectSilicone Pediatric Trays			
DYND160806	6 Fr 3 mL, with 200 mL Urine Meter with 2,000 mL Drain Bag	A4315 + A4333	10/cs
DYND160808	8 Fr 5 mL, with 200 mL Urine Meter with 2,000 mL Drain Bag	A4315 + A4333	10/cs
DYND160810	10 Fr 5 mL, with 200 mL Urine Meter with 2,000 mL Drain Bag	A4315 + A4333	10/cs
Total One Layer Trays – Includes Peri-Care Wipe			
URO170806	6 Fr 3 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	A4315 + A4333	10/cs
URO170808	8 Fr 5 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	A4315 + A4333	10/cs
URO170810	10 Fr 5 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	A4315 + A4333	10/cs
Total One Layer Trays – Includes Peri-Care Wipe, Vented Drainage and SafeSecure			
URO180806S	6 Fr 3 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	No Code	10/cs
URO180808S	8 Fr 5 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	No Code	10/cs
URO180810S	10 Fr 5 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	No Code	10/cs



DYND160808T



DYND110510

100% Silicone catheter has an integrated thermistor which aids in monitoring the patient's core body temperature via the bladder.

Pediatric One Layer SelectSilicone Temperature Sensing Trays

Introducing the world's first temperature-sensing closed-system tray specifically designed for pediatric patients. Medline's revolutionary Buddy the Brave tray is now available with our 100% Silicone 400-series temperature-sensing Foley catheter. It can detect a patient's core body temperature with +/- 0.2° C accuracy while maintaining urine flow.

» Not made with natural rubber latex

Item No.	Size	HCPCS	Pkg.
SelectSilicone Temperature Sensing Trays			
DYND160806T*	6 Fr 1.5 mL, with 200 mL Urine Meter with 2,000 mL Drain Bag	No Code	10/cs
DYND160808T*	8 Fr 3 mL, with 200 mL Urine Meter with 2,000 mL Drain Bag	A4315 + A4333	10/cs
DYND160810T*	10 Fr 3 mL, with 200 mL Urine Meter with 2,000 mL Drain Bag	A4315 + A4333	10/cs
Total One Layer Pediatric Trays – Includes Peri-Care Wipe			
URO170806T*	6 Fr 3 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	No Code	10/cs
URO170808T*	8 Fr 5 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	A4315 + A4333	10/cs
URO170810T*	10 Fr 5 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	A4315 + A4333	10/cs
Total One Layer Pediatric Trays – Includes Peri-Care Wipe , Vented Drainage and SafeSecure			
URO180806T*	6 Fr 3 mL, with 400 mL Urine Meter with 2,500 mL Drain Bag	No Code	10/cs
SelectSilicone Temperature Sensing Foley Catheters			
DYND110508	8 Fr, 3 mL Balloon	A4344	10/cs
DYND110510	10 Fr, 3 mL Balloon	A4344	10/cs

*Conditionally MRI compatible

Foley Catheter Trays



Buddy the Brave Add-a-Cath Foley Catheterization System

Medline asked experts in pediatric care how we could help them reduce catheter-associated urinary tract infection (CAUTI), standardize best practices, and better serve their youngest catheter patients. Through this process, we developed the Buddy the Brave Add-a-Cath Foley Catheterization System.

- » One layer nurse-designed tray walks clinician through best practices
- » Add-a-cath style tray allows clinician to pull multiple catheters in case they are needed
- » Metal-free urine meter, Buddy the Brave securement device, and patient care card are all designed with your youngest patients in mind
- » Single sterile catheters are packaged with the correct amount of sterile water for optimal inflation, reducing the risk over over-inflating the balloon

Item No.	Size	Perineal Cleansing Wipe	Sterile Water Syringe	Pkg.
Buddy the Brave Add-a-Cath Trays				
UR0160300	200 mL Urine Meter with 2,000 mL Drain Bag			10/cs
UR0160350	400 mL Urine Meter with 2,500 mL Drain Bag		•	10/cs
UR0170350	400 mL Urine Meter with 2,500 mL Drain Bag	•		10/cs
UR0170350A	400 mL Urine Meter with 2,500 mL Drain Bag	•	•	10/cs



URO2S1003K

Item No.	Size	Pkg.
SelectSilicone Foley Catheter + Sterile Water Syringe Kits		
URO2S0615K	6 Fr 3 mL Foley, 3 mL Sterile Water Syringe	10/cs
URO2S0803K	8 Fr 5 mL Foley, 5 mL Sterile Water Syringe	10/cs
URO2S1003K	10 Fr 5 mL Foley, 5 mL Sterile Water Syringe	10/cs
URO2S0615TK	6 Fr 3 mL Temperature-Sensing Foley, 3 mL Sterile Water Syringe	10/cs
URO2S0803TK	8 Fr 5 mL Temperature-Sensing Foley, 5 mL Sterile Water Syringe	10/cs
URO2S1003TK	10 Fr 5 mL Temperature-Sensing Foley, 5 mL Sterile Water Syringe	10/cs

Pediatric 100% Silicone Foley Catheters

Medline's pediatric Foley catheters are available in 6, 8 and 10 Fr sizes with 3 mL and 5 mL balloons and short lengths for pediatric use. Completely inert catheters minimize tissue irritation during extended periods of indwelling use.

- » Pliable catheter with a firm tip to ease insertion
- » Large, smooth eyes for maximum comfort and drainage
- » Symmetrical balloon for added reliability
- » Not made with natural rubber latex



DYND11553

Item No.	Size	Balloon Size	HCPCS	Pkg.
DYND11552	6 Fr	3 mL	A4344	10/cs
DYND11553	8 Fr	3 mL	A4344	10/cs
DYND11554	10 Fr	3 mL	A4344	10/cs

Silicone-Elastomer Coated Latex Foley Catheters

Innovative manufacturing process provides the maximum amount of silicone coating available to increase patient comfort.

- » Shortened length, specifically for pediatric patients
- » Silicone-elastomer coating will resist encrustation
- » Symmetrical balloon for added reliability
- » Includes stylet



DYND11708

Item No.	Size	Balloon Size	HCPCS	Pkg.
DYND11706	6 Fr	3 mL	A4338	12/cs
DYND11708	8 Fr	5 mL	A4338	12/cs
DYND11710	10 Fr	5 mL	A4338	12/cs

Intermittent Catheters



Speci-Cath Kits

Medline's Speci-Cath kits provide an easy way to collect a sterile urine sample. The catheter is pre-inserted into the specimen tube to ensure sterile technique. All sizes are packaged sterile.

» Not made with natural rubber latex

Item No.	Size	Catheter Length	Vinyl Gloves	Pvp Swabsticks 3/pk	Lubricating Jelly	HCPCS	Pkg.
DYND10815	8 Fr	9" (23 cm)	•	•	•	A9270	25/cs
DYND10820	5 Fr	9" (23 cm)	•	•	•	A9270	25/cs

Centurion Pediatric Cath Kit

» Not made with natural rubber latex

Item No.	Size	Catheter Length	Aloetouch 3G Gloves	Pvp Swabsticks 3/pk	Lubricating Jelly	HCPCS	Pkg.
CKF155L	5 Fr	9" (23 cm)	•	•	•	A9270	25/cs
CKF5050L	5 Fr	9.5" (24.1 cm)	•	•	•	A9270	25/cs

Vinyl Intermittent Catheters

Medline's pediatric vinyl intermittent catheters are 10-inches (25.4 cm) long and suitable for both male and female patients. Smooth eyes provide increased patient comfort. All sizes are packaged sterile.

» Not made with natural rubber latex

Item No.	Size	Catheter Length	HCPCS	Pkg.
DYND10712	6 Fr	10" (25.4 cm)	A4351	30/cs
DYND10713	8 Fr	10" (25.4 cm)	A4351	30/cs
DYND10714	10 Fr	10" (25.4 cm)	A4351	30/cs



DYND10714

Pediatric Tube Securement Device

- » Made of a breathable, non-woven, child-friendly material
- » Adhesive under the hook-and-loop strap prevents pistoning movement of the catheter
- » Sterile

Item No.	Size	HCPCS	Pkg.
DYND7600P	Pediatric Foley Catheter Securement Device	A4333	100/cs



DYND7600P

Exact Measure Pediatric Urine Meter

Specifically designed for the pediatric patient population, the exact measure urine meter has 1 mL graduations for easy reading of output measurements.

- » 200 mL urine meter with 2,000 mL drain bag
- » Exclusive slide tap drainage port
- » 48- or 60-inches (122 or 152.4 cm) of tubing for direct flow into the bag
- » Luer lock needle-less access port for urine sampling
- » Packaged OR Sterile
- » Not made with natural rubber latex

Item No.	Description	Size	HCPCS	Pkg.
DYND17407	200 mL Urine Meter with 2,000 mL Drain Bag	48" (122 cm) Long Tubing	A4357	10/cs
DYND17407L	200 mL Urine Meter with 2,000 mL Drain Bag	60" (152.4 cm) Long Tubing	A4357	10/cs



DYND17407

Infant Urine Collector

The versatile infant urine collector can be used with both male and female infants.

- » Adhesive is gentle on skin and easy to remove, yet holds securely
- » Folds to form a spill-proof specimen transport
- » Available in both sterile and non-sterile packaging
- » Not made with natural rubber latex

Item No.	Description	Size	HCPCS	Pkg.
MDS190505	Non-Sterile Pediatric Urine Collector	200 mL	A4328	50/bx
MDS190510	Sterile Pediatric Urine Collector	200 mL	A4328	50/bx



MDS190510



IRRIGATION AND OTHER PRODUCTS

Irrigation Syringes **63**

Accessories **64**

Nursing Products **67**

Irrigation Syringes

Medline irrigation syringes have “ears” to minimize hand slippage and subsequent contamination of barrel or solution container. Flat bottom solution container is sturdy and does not tip when filled with liquid. The 60 mL syringe barrel features large, raised, easy-to-read graduations in both ounces and milliliters. Rigid packaging includes graduated tray basin.



Bulb Syringe

- » Soft, pliable, textured polymer bulb for maximum control of drawn or expelled fluid
- » Sterile
- » Not made with natural rubber latex

Item No.	Description	Package	Pad	Sq. Jar	Drape	Luer Tip	Tip Protector	HCPCS	Pkg.
DYND20100	Bulb	Rigid Base Wrapped in CSR	PVP	●	●		●	A4320	20/cs
DYND20102	Bulb	Rigid Base with Tyvek Lid	Alcohol	●	●		●	A4320	20/cs
DYND20125	Bulb	Soft with Tyvek Lid					●	A4322	50/cs
DYND20127*	Bulb	Bulk					●	A4322	50/cs

*Non-sterile



Piston Syringe

- » Thumb ring helps minimize hand slippage and prevent contamination
- » Sterile
- » Not made with natural rubber latex

Item No.	Description	Package	Pad	Sq. Jar	Drape	Luer Tip	Tip Protector	HCPCS	Pkg.
DYNC2303**	Piston	Rigid Base with Tyvek Lid	Alcohol	●	●		●	A4320	20/cs
DYND20300	Piston	Rigid Base Wrapped in CSR	PVP	●	●		●	A4320	20/cs
DYND20302	Piston	Rigid Base with Tyvek Lid	Alcohol	●	●		●	A4320	20/cs
DYND20325	Piston	Soft with Tyvek Lid				●	●	A4322	50/cs
DYND20327*	Piston	Bulk				●	●	A4322	50/cs
DYND20323	Toomey	Soft with Tyvek Lid, 70 mL				●	●	A4322	50/cs

*Non-sterile

**100 mL 0.9% NaCl

Accessories



DYND70280

DYND70277



URO70283



DYND70279



DYND2219

Ear Ulcer Syringes

Ideal for cleansing nasal passages of newborns and irrigating smaller wounds.





- » Made of PVC
- » Individually packaged in a soft package with peel lid
- » Sterile

Item No.	Size	Pkg.
DYND70279	1-oz. (29.6 mL) Sterile Ear Ulcer Syringe	50/cs
DYND70280	2-oz. (59.1 mL) Sterile Ear Ulcer Syringe	50/cs
DYND70277	3-oz. (88.7 mL) Sterile Ear Ulcer Syringe	50/cs
URO70283	2-oz. (59.1 mL) Red Sterile Ear Ulcer Syringe	24/cs
DYNDX48894	2-oz. (59.1 mL) Non-Sterile Ear Ulcer Syringe	200/cs

Luer Tip/Catheter Adapter

- » Single Luer tip easily attaches to a syringe
- » Attaches to the end of a Luer tip syringe to make it a catheter tip
- » Individually packaged
- » Sterile

Item No.	Size	Pkg.
DYND2219	Luer Tip Adapter	100/cs



Cysto/TUR Irrigation Sets

Ideal for use in cystoscopy procedures or for continuous bladder irrigation, these sets, which are not made with natural rubber latex, are 80- to 96-inches (203 to 244 cm) in length and compatible for use with solution bags.

- » Roller clamp
- » Flexible tubing
- » Drip chamber
- » Spike for solution bags
- » Packaged OR sterile
- » Not made with natural rubber latex

Item No.	Description	Length	HCPCS	Pkg.
Cysto/Tur Irrigation Sets				
DYND19120	Cysto/Tur Irrigation Set with 0.188" (4.8 mm) I.D. Tubing and Adjustable Roller Clamp	90"	A4355	20/cs
DYND19122	"Y" Type Cysto/Tur Irrigation Set with 0.28" (7.1 mm) I.D. Tubing, 2 Roberts Clamps Above Drip Chamber and Adjustable Roller Clamp Below	90"	A4355	20/cs
DYND19124	4 Lead Irrigation Set with 0.28" (7.1 mm) I.D. Tubing, 4 Roberts Clamps Above Drip Chamber and Adjustable Roller Clamp Below	90"	A4355	20/cs
Cysto/Tur Irrigation Sets (Tubing more flexible)				
DYND19140	Cysto/Tur Irrigation Set with 0.28" (7.1 mm) I.D. Tubing and Adjustable Roller Clamp	80"	A4355	20/cs
DYND19142	"Y" Type Cysto/Tur Irrigation Set with 0.28" (7.1 mm) I.D. Tubing, 2 Roberts Clamps Above Drip Chamber and Adjustable Roller Clamp Below	81"	A4355	20/cs
DYND19144	4 Lead Irrigation Set with 0.28" (7.1 mm) I.D. Tubing, Roberts Clamps Above Drip Chamber and Adjustable Roller Clamp Below	93"	A4355	20/cs

Accessories



Light-Source Adaptable Disposable Vaginal Speculum

Eliminates the need for positioning a separate light source.

- » Single use, reducing the risk of cross-contamination
- » Self-holding in a variety of positions
- » Made of clear, high-strength, non-toxic plastic
- » Easily adjusts with one hand
- » Non-sterile
- » Also available individually or in boxes of 25
- » Compatible with Welch Allyn corded illumination system

Item No.	Description	Pkg.
DYND70401L	Large	100/cs
DYND70401M	Medium	100/cs
DYND70401S	Small	100/cs



Disposable Vaginal Speculum

Help reduce the risk of cross-contamination with single-use speculums.

- » High-strength, clear non-toxic plastic
- » Self-holding in a variety of positions
- » Easily adjusts with one hand
- » Non-sterile
- » Also available in boxes of ten

Item No.	Description	Pkg.
DYND70420L	Large	100/cs
DYND70420M	Medium	100/cs
DYND70420S	Small	100/cs
VAGAP100	Centurion Vaginal Application	50/cs



Enema Bag Sets

Includes a spill-proof valve on the 1,500 cc enema bag.

- » Easy-to-use funnel top
- » 54-inch (137 cm) soft vinyl tube has a pre-lubricated, smooth, non-traumatic tip with rounded eyes for safety and easy insertion
- » Castile soap packet
- » Plastic shut-off clamp, designed for one-handed operation
- » Moisture-proof underpad
- » Individually packaged in a polybag or box; non-sterile

Item No.	Description	HCPCS	Pkg.
DYND70102	Enema Bag Set with Clamp, Polybag	A4458	48/cs
DYND70100	Enema Bag Set with Clamp, Boxed	A4458	48/cs



DYND70104

Enema Bucket Sets

The complete kit contains an enema bucket, 54-inch (137 cm) tube with plastic shut-off clamp, tube cap, castile soap packet and moisture-proof underpad.

- » 1,500 cc bucket is graduated in both ounces and milliliters
- » 54-inch (137 cm) soft vinyl tube has a pre-lubricated, smooth, non-traumatic tip with rounded eyes for safety and easy insertion
- » Clamp on tube is designed for one-handed operation
- » Individually wrapped in a polybag
- » Non-sterile

Item No.	Description	Pkg.
DYND70104	Enema Bucket Set with Clamp	50/cs

Perineal Cleansing Bottle

Regulate water flow by lifting cap and turning.

- » 8-oz. (236.5 mL) perineal cleansing bottle with screw top
- » Graduated
- » Not made with natural rubber latex
- » Non-sterile

Item No.	Description	Pkg.
DYND70125	Perineal Cleansing Bottle	50/cs



DYND70125

Trademarks

Accuryrn, Aloetouch, Buddy the Brave, ERASE CAUTI, Foley Insetag, SafeSecure, Select Silicone, Silvertouch, SimplyLatex and Sureprep are registered trademarks of Medline Industries, Inc.

BladderPoint, My-Cath, are trademarks of Medline Industries, Inc.

BioCon is a trademark of Mcubetechnology Co., Ltd.

BioPlus+ and FreeDerm are trademarks and CathGrip and Exo-Cath are registered trademarks of BioDerm, Inc.

GE is a registered trademark of General Electric Company

Mylar is a registered trademark of E. I. Du Pont De Nemours and Company

Philips is a trademark of Koninklijke Philips N.V.

Siemens Climate Solutions is a registered trademark of Siemens Aktiengesellschaft Joint Stock Company

Tyvek is a registered trademark of E. I. Du Pont De Nemours and Company

Alphabetical Index

B

- BioCon 750 Bladder Scanner, 24
- BioCon 900 Bladder Scanner, 25
- Bladder Scanner, 24
- Buddy the Brave
 - Foley Catheter, 55
 - Foley Catheter Trays, 56–58
 - Pediatric Foley Catheter Trays, 16

C

- Catheters
 - 3-way Foley Catheters, 28
- Catheter Securement Devices
 - Foley Catheter Securement, 43
- Catheter Trays, 8–9
- Cysto/TUR Irrigation Sets, 65

D

- Drain Bags, 36–38
 - Covers & Holder, 41
 - Hook, 41

E

- Ear Ulcer Syringes, 64
- Enema
 - Bag Sets, 66
 - Bucket Sets, 67

F

- Foley Catheter
 - Insertion Trays, 23
 - Leg Straps, 43
 - Pediatrics, 16, 56
 - Securement Devices, 43
 - SelectSilicone
 - Foley Catheters, 27–28
 - Silvertouch Foley Catheters, 29
 - Simply Latex
 - Foley Catheters, 30–31
 - Temperature-Sensing Foley Catheters, 32–33
 - Trays
 - 2-Layer, 18–21
 - Foley Catheter Insertion, 23
 - One Layer, 8
 - One Layer, 8–15
 - Pediatrics, 16, 56, 58

I

- Intermittent Catheterization Trays
 - My-Cath Touch-Free Self-Catheter System, 49
 - Pre-Connected Intermittent, 48
 - Total One Layer Intermittent, 47
- Intermittent Catheters, 45
 - Clear Vinyl, 45
 - Red Rubber Latex All-Purpose, 46
 - Vinyl Coude, 46
 - Vinyl Hydrophilic, 45
- Irrigation Syringes, 63

L

- Leg Bags, 39
- Luer Tip/Catheter Adapter, 64

M

- Male External Catheters
 - Exo-Cath Latex, 51
 - Men's Liberty Acute, 50

O

- One Layer Foley Catheter Trays
 - Add-a-Cath, 8
 - Coude, 11
 - Pediatric, 16
 - SelectSilicone, 12
 - SelectSilicone Temperature-Sensing, 14
 - Silvertouch 100% Silicone, 9
 - SimplyLatex, 10

P

- Pediatrics
 - 100% Silicone Foley Catheters, 59
 - Add-a-Cath Foley Catheterization System, 58
 - ERASE CAUTI 100% Silicone, 56
 - One Layer Foley Catheters, 56
 - Securement Device, 61
 - Silicone-Elastomer Coated Latex Foley, 59
 - Speci-Cath Kits, 60
 - Urine Collector, 61
 - Urine Meter, 61
 - Vinyl Intermittent Catheters, 60
- Perineal Cleansing Bottle, 67

S

- Silicone-Elastomer Coated Latex Foley Catheters
 - 3-Way, 31
 - Specialty Coude-Tip, 31
- Silicone Foley Catheters, 27
 - Pediatric, 28
 - Specialty Tiemann/Coude-Tip, 28
- Speci-Cath Kits, 52

T

- Temperature-Sensing Foley Catheters, 32
 - Cables, 33
 - Pediatric, 33
- Temperature-Sensing Trays, 14, 57
- 2-Layer Foley Catheter Trays
 - 100% Silicone, 21
 - Add-a-Cath, 18
 - Silicone-Elastomer Coated Latex, 20
 - Silvertouch 100% Silicone, 19

U

- Urine Collection
 - Accessories, 34–35
 - Drain Bags, 36–38
 - Leg Bags, 39
 - Urine Meters, 35

V

- Vaginal Speculum, 66







Medline Industries, Inc. Three Lakes Drive, Northfield, IL 60093

Medline United States
1-800-MEDLINE (633-5463)
medline.com | info@medline.com

Medline Canada
1-800-396-6996
medline.ca | canada@medline.com

Medline México
01-800-831-0898
medlinemexico.com | mexico@medline.com

FOLLOW US    

HCPCS DISCLAIMER: HCPCS codes and home health consolidated billing codes provided by Medline are intended as general guidelines only. Medline does not guarantee coverage or reimbursement of any products. The customer must address all coverage and reimbursement issues (including correctness and accuracy of codes) with their individual payers. It is the customer's responsibility to ensure the accuracy and appropriateness of each claim the customer submits, in accordance with all applicable payer requirements.

Some products may not be available for sale in Mexico or Canada. We reserve the right to correct any errors that may occur within this brochure.
© 2019 Medline Industries, Inc. Medline is a registered trademark of Medline Industries, Inc. MKT19101408 / LIT184 / 5M / SG / 3