

# F&P 850™ System

Adult and Infant Respiratory Care Catalogue



# TABLE OF CONTENTS

---

## F&P 850 System - Adult Care

---

Adult Invasive Ventilation	2
Adult Noninvasive Ventilation	4
Adult Optiflow™ and Oxygen Therapy	6

---

## F&P 850™ System - Infant Care

---

Infant Resuscitation	10
Infant Invasive Ventilation	12
Infant Noninvasive and F&P Bubble CPAP™	14
Infant Optiflow and Oxygen Therapy	16

---



## **F&P 850™ System** | Adult Care

The F&P 850 System has been designed to be simple to operate while answering a strong call for a high performing solution capable of operating in either Invasive or Noninvasive control modes.

Significant benefits flow from one system which can be configured to deliver optimal outcomes for all therapy needs.



# Adult Invasive Ventilation

The MR850 humidifier, in combination with the MR290 Autofill Chamber and RT series breathing circuits cohesively ensures patient care and therapy efficacy for adult and paediatric patients requiring intubation.

The F&P 850 System is designed to deliver optimal humidity and includes a range of breathing circuits with ISO standard connectors. These connectors allow seamless connection to ventilators and patient interfaces. Fisher & Paykel Healthcare's exclusive Evaqua™ expiratory limb technology minimizes mobile condensate in the expiratory limb.



## 1 Humidifier and Accessories

PART NO	DESCRIPTION	QUANTITY
MR850	MR850 Respiratory Humidifier	1/each
900MR805	Heater Wire Adapter for Dual Heated Breathing Circuits	1/each
900MR806	Heater Wire Adapter for Inspiratory Heated Breathing Circuits	1/each
900MR869	Temperature/Flow Probe for MR850 Compatible Breathing Circuits (1.5m)	1/each

## 2 Humidification Chambers

PART NO	DESCRIPTION	QUANTITY
MR290V	Auto-fill Adult/Infant Humidification Chamber with Vented Water Bag Spike	40/box
MR290HFV	Auto-fill Adult/Infant Humidification Chamber with Vented Water Bag Spike for High Frequency Ventilation Applications	40/box

## 3 Breathing Circuits and Kits

PART NO	DESCRIPTION	QUANTITY
RT380	Dual Limb Adult Circuit Kit with Evaqua 2 Technology and Inspiratory/Expiratory Filter	10/box
RT385	Dual Limb Adult Circuit Kit with Evaqua 2 Technology, Pressure Line and Inspiratory/Expiratory Filter	10/box
RT340	Dual Limb Adult Breathing Circuit Kit with Evaqua Technology and Inspiratory/Expiratory Filter	10/box
RT345	Dual Limb Adult Breathing Circuit Kit with Evaqua Technology, Pressure Line and Inspiratory/Expiratory Filter	10/box
RT200	Dual Limb Adult Breathing Circuit Kit	10/box
RT204	Dual Limb Adult Breathing Circuit Kit with Pressure Line	10/box
RT100	Dual Limb Adult Breathing Circuit Kit (excl. Chamber)	10/box
RT104	Dual Limb Adult Breathing Circuit (excl. Chamber) with Pressure Line	10/box
RT205	Dual Limb Adult Breathing Circuit Kit with Expiratory Water-trap	10/box
RT206	Dual Limb Adult Breathing Circuit Kit with Expiratory Water-trap and Pressure Line	10/box
RT105	Dual Limb Adult Breathing Circuit Kit (excl. Chamber) with Expiratory Water-trap	10/box
RT106	Dual Limb Adult Breathing Circuit (excl. Chamber) with Expiratory Water-trap and Pressure Line	10/box

## 4 Interfaces and Accessories

PART NO	DESCRIPTION	QUANTITY
RT021	F&P DirectConnect™ Dual Swivel Catheter Mount	20/box
RT019	Inspiratory/Expiratory Filter	20/box
RT020	Expiratory Filter	20/box
RT023	T-Piece Weaning Kit	20/box

## Adult Noninvasive Ventilation

The F&P 850 System for noninvasive ventilation is comprised of the MR850 humidifier with the MR290 Autofill Chamber and RT series breathing circuits with circuit step-down options. This system ensures exceptional performance for patients requiring noninvasive therapy.

Essential humidity is delivered to one of three types of FreeMotion™ mask. These are available as vented and non-vented full-face mask and vented nasal mask. The unique comfort and seal technologies minimize leak during patient movement.



## 1 Humidifier and Accessories

PART NO	DESCRIPTION	QUANTITY
MR850	MR850 Respiratory Humidifier	1/each
900MR805	Heater Wire Adapter for Dual Heated Breathing Circuits	1/each
900MR806	Heater Wire Adapter for Inspiratory Heated Breathing Circuits	1/each
900MR869	Temperature/Flow Probe for MR850 Compatible Breathing Circuits (1.5m)	1/each

## 2 Humidification Chambers

PART NO	DESCRIPTION	QUANTITY
MR290V	Auto-feed Adult/Infant Humidification Chamber with Vented Water Bag Spike	40/box
MR290HFV	Auto-feed Adult/Infant Humidification Chamber with Vented Water Bag Spike for High Frequency Applications	40/box

## 3 Breathing Circuits and Kits

PART NO	DESCRIPTION	QUANTITY
RT319	Single Limb Adult Breathing Circuit Kit with Low Resistance to Flow Inspiratory Limb, Flexible Extension, Exhalation Valve, Inspiratory Filter, and Pressure Line	10/box
RT201	Single Limb Adult Breathing Circuit Kit with Wye Piece and PEEP Valve Extension	10/box
RT101	Single Limb Adult Breathing Circuit (excl. Chamber) with Wye Piece and CPAP Extension	10/box
RT380	Dual Limb Adult Circuit Kit with Evaqua 2 Technology and Inspiratory/Expiratory Filter	10/box
RT385	Dual Limb Adult Circuit Kit with Evaqua 2 Technology, Pressure Line and Inspiratory/Expiratory Filter	10/box
RT340	Dual Limb Adult Breathing Circuit Kit with Evaqua Technology and Inspiratory/Expiratory Filter	10/box
RT345	Dual Limb Adult Breathing Circuit Kit with Evaqua Technology, Pressure Line and Inspiratory/Expiratory Filter	10/box
RT200	Dual Limb Adult Breathing Circuit Kit	10/box
RT204	Dual Limb Adult Breathing Circuit Kit with Pressure Line	10/box
RT100	Dual Limb Adult Breathing Circuit Kit (excl. Chamber)	10/box
RT104	Dual Limb Adult Breathing Circuit (excl. Chamber) with Pressure Line	10/box
RT205	Dual Limb Adult Breathing Circuit Kit with Expiratory Water-trap	10/box
RT206	Dual Limb Adult Breathing Circuit Kit with Expiratory Water-trap and Pressure Line	10/box
RT105	Dual Limb Adult Breathing Circuit Kit (excl. Chamber) with Expiratory Water-trap	10/box
RT106	Dual Limb Adult Breathing Circuit (excl. Chamber) with Expiratory Water-trap and Pressure Line	10/box
RT401 (S,M,L)	Adult CPAP Pack including Single Limb Adult Breathing Circuit Kit, Non-vented Face Mask with Anti-asphyxiation Valve, Adjustable PEEP Valve and Pressure Relief Valve	10/box
RT411 (S,M,L)	Adult CPAP Pack including Single Limb Adult Unheated Breathing Circuit, Non-vented Face Mask with Anti-asphyxiation Valve, Adjustable PEEP Valve and Pressure Relief Valve	10/box

## 4 Interfaces and Accessories


PART NO	DESCRIPTION	QUANTITY
RT040 (S,M,L)	Vented Face Mask with Anti-asphyxiation Valve	10/box
RT041 (S,M,L)	Non-vented Face Mask	10/box
RT042 (M,L)	Vented Nasal Mask	10/box
RT019	Inspiratory/Expiratory Filter	20/box
RT022	Unheated Flexible Extension	20/box
RT023	T-Piece Weaning Kit	20/box

## Adult Optiflow™ and Oxygen Therapy

The MR850 Humidifier together with the MR290 Autofill Chamber and RT series breathing circuits complete a system which is unmatched in terms of humidity performance and patient tolerance to oxygen therapy.

The soft, flexible nasal prongs of the Optiflow Nasal Cannula are designed for comfortable prolonged use replacing the need for a face mask.

Lightweight, flexible delivery tubing allows free movement during sleep and activity to make this a highly desirable, comfortable and unobtrusive therapy solution which delivers a wide range of flows to patients requiring oxygen therapy.





## 1 Humidifier and Accessories

PART NO	DESCRIPTION	QUANTITY
MR850	MR850 Respiratory Humidifier	1/each
900MR806	Heater Wire Adapter for Inspiratory Heated Breathing Circuits	1/each
900MR869	Temperature/Flow Probe for MR850 Compatible Breathing Circuits (1.5m)	1/each

## 2 Humidification Chambers

PART NO	DESCRIPTION	QUANTITY
MR290V	Auto-feed Adult/Infant Humidification Chamber with Vented Water Bag Spike	40/box
MR290HFV	Auto-feed Adult/Infant Humidification Chamber with Vented Water Bag Spike for High Frequency Applications	40/box

## 3 Breathing Circuits and Kits

PART NO	DESCRIPTION	QUANTITY
RT202	Single Limb Adult Breathing Circuit Kit	10/box
RT102	Single Limb Adult Breathing Circuit (excl. Chamber)	10/box
RT408	Single Limb Adult Breathing Circuit Kit with Flexible Extension and Venturi Air Entrainer*	10/box
RT307	Single Limb Adult Unheated Breathing Circuit Kit with Venturi Air Entrainer*	10/box

## 4 Interfaces and Accessories

PART NO	DESCRIPTION	QUANTITY
OPT542	Adult Optiflow Cannula - Small	20/box
OPT544	Adult Optiflow Cannula – Medium	20/box
OPT546	Adult Optiflow Cannula – Large	20/box
OPT570	Adult Optiflow Tracheostomy Interface	20/box
RT008	Venturi Air Entrainer*	20/box
900MR961	Dual Flow Adapter – For Connection of Gas Supply Tubing	20/box

\*Venturi Air Entrainer is not compatible with F&P Optiflow Nasal Cannula



## **F&P 850™ System** | Infant Care

Precious new life deserves the best possible start. From first breath to respiratory independence, the F&P 850 System helps to emulate the natural physiological balance in healthy, mature lungs.


As a result, caregivers can nurture life, confident they are using the best therapy solutions, delivered in the safest, most efficient ways.

# Infant Resuscitation

T-piece Resuscitation with Optimal Humidity (heated and humidified gas) can be delivered with the F&P Neopuff™ Infant T-Piece Resuscitator and MR850 Humidifier System.

The F&P Neopuff™ is a manually operated, gas-powered resuscitator designed to provide controlled Peak Inspiratory Pressure (PIP) and consistent Positive End Expiratory Pressure (PEEP).

The infant Resuscitation Masks are specially designed with soft and pliable edges to conform comfortably to an infant's face facilitating a seal around the infant's mouth and nose.





## 1 F&P Neopuff™, Humidifier and Accessories

PART NO	DESCRIPTION	QUANTITY
RD900	Neopuff Infant T-Piece Resuscitator	1/each
MR850	MR850 Respiratory Humidifier	1/each
900MR806	Heater Wire Adapter for Inspiratory Heated Breathing Circuits	1/each
900MR868	Temperature/Flow Probe for MR850 Compatible Breathing Circuits (1.1m)	1/each
900MR208	Reflective Shield for Temperature/Flow Probe	20/box

## 2 Humidification Chambers

PART NO	DESCRIPTION	QUANTITY
MR290V	Auto-feed Adult/Infant Humidification Chamber with Vented Water Bag Spike	40/box
MR225	Manual Feed Low Volume Infant Humidification Chamber	50/box

## 3 Breathing Circuits and Kits

PART NO	DESCRIPTION	QUANTITY
900RD110	Single limb Infant Resuscitation Circuit with T-Piece Attachment	10/box

## 4 Interfaces and Accessories

PART NO	DESCRIPTION	QUANTITY
RD800-EN	Infant Resuscitation Masks – starter kit (English text only)	5/box
RD800-SK	Infant Resuscitation Masks – starter kit (Universal)	5/box
RD803-10	Micro-Premature (XS) Mask	10/box
RD804-10	Premature (S) Mask	10/box
RD805-10	Neonatal (M) Mask	10/box
RD806-10	Infant (L) Mask	10/box
RD807-10	Pediatric (XL) Mask	10/box

## Infant Invasive Ventilation

Invasive Ventilation with Optimal Humidity is the best solution for supporting airway defence functions and ventilation in infant patients.

The F&P 850 system is designed to deliver Optimal Humidity, and includes a range of breathing circuits that ensures compatibility with contemporary infant ventilators.



## 1 Humidifier and Accessories

PART NO	DESCRIPTION	QUANTITY
MR850	MR850 Respiratory Humidifier	1/each
900MR805	Heater Wire Adapter for Dual Heated Breathing Circuits	1/each
900MR806	Heater Wire Adapter for Inspiratory Heated Breathing Circuits	1/each
900MR868	Temperature/Flow Probe for MR850 Compatible Breathing Circuits (1.1m)	1/each
900MR208	Reflective Shield for Temperature/Flow Probe	20/box

## 2 Humidification Chambers

PART NO	DESCRIPTION	QUANTITY
MR290V	Auto-feed Adult/Infant Humidification Chamber with Vented Water Bag Spike	40/box
MR290HFV	Auto-feed Adult/Infant Humidification Chamber with Vented Water Bag Spike for High Frequency Applications	40/box
MR225	Manual Feed Low Volume Infant Humidification Chamber	50/box
900MR191X	Waterfeed Extension Kit	10/box

## 3 Breathing Circuits and Kits

PART NO	DESCRIPTION	QUANTITY
RT265	Dual Limb Infant Breathing Circuit Kit with Evaqua 2 Technology and Pressure Line (Flow → 4L/min)	10/box
RT266	Dual Limb Infant Breathing Circuit Kit with Evaqua 2 Technology and Pressure Line (Flow → 0.3-4L/min)	10/box
RT267	Dual Limb Infant Breathing Circuit Kit with Evaqua 2 Technology and Pressure Line for SLE2000 Series Ventilators	10/box
RT268	Dual Limb Infant Breathing Circuit Kit with Evaqua 2 Technology and Pressure Line for SLE5000 Series Ventilators	10/box
RT235	Dual Limb Infant Breathing Circuit Kit with Evaqua Technology and Pressure Line (Flow → 4L/min)	10/box
RT236	Dual Limb Infant Breathing Circuit Kit with Evaqua Technology and Pressure Line (Flow → 0.3-4L/min)	10/box
RT225	Dual Limb Infant Breathing Circuit Kit with Expiratory Water-trap and Pressure Line (Flow → 4L/min)	10/box
RT226	Dual Limb Infant Breathing Circuit Kit with Expiratory Water-trap and Pressure Line (Flow → 0.3-4L/min)	10/box
RT227	Dual Limb Infant Breathing Circuit Kit with Expiratory Water-trap and Pressure Line for SLE2000 Series Ventilators	10/box
RT228	Dual Limb Infant Breathing Circuit Kit with Expiratory Water-trap and Pressure Line for SLE5000 Series Ventilators	10/box
RT125	Dual Limb Infant Breathing Circuit (excl. Chamber) with Expiratory Water-trap and Pressure Line (Flow → 4L/min)	10/box
RT126	Dual Limb Infant Breathing Circuit (excl. Chamber) with Expiratory Water-trap and Pressure Line (Flow → 0.3-4L/min)	10/box
RT127	Dual Limb Infant Breathing Circuit (excl. Chamber) with Expiratory Water-trap and Pressure Line for SLE2000 Series Ventilators	10/box

## 4 Interfaces and Accessories


PART NO	DESCRIPTION	QUANTITY
RT010	Inline Nebulizer Adapter Kit	20/box
RT030	Side-stream Nebulizer Adapter Kit	20/box
RT061	Adaptor Kit Drager Ventilators (2)	50/box
RT062	Luerlock Adaptor (2)	50/box
RT063	Pressure Port Adaptor (2)	100/box
RT064	Adaptor Kit Infant Oxygen/CPAP (2)	50/box
RT065	Adaptor Kit Infant Hamilton/EME (2)	10/box

# Infant Noninvasive and F&P Bubble CPAP™

Fisher & Paykel Healthcare developed the Bubble CPAP System for use with spontaneously breathing infants requiring respiratory support.

The F&P Bubble CPAP System incorporates unique features and world-leading humidification technology to provide a respiratory support system that is effective and easy to use.

The F&P Infant Interface complements the Bubble CPAP System, providing comfort and optimal CPAP support to the infant at risk using nasal prongs or nasal masks.





## 1 Humidifier and Accessories

PART NO	DESCRIPTION	QUANTITY
MR850	MR850 Respiratory Humidifier	1/each
900MR806	Heater Wire Adapter for Inspiratory Heated Breathing Circuits	1/each
900MR868	Temperature/Flow Probe for MR850 Compatible Breathing Circuits (1.1m)	1/each
900MR208	Reflective Shield for Temperature/Flow Probe	20/box


## 2 Breathing Circuits, Bubble CPAP™ Systems and Starter Kits

PART NO	DESCRIPTION	QUANTITY
RT324	Single Limb Infant Breathing Circuit Kit with Pressure Line, Expiratory Extension for Flow Driver CPAP, and Water-bag Extension	10/box
BC151-10	Single-heated Infant Breathing Circuit with MR290 Chamber, Pressure Manifold at 17cmH <sub>2</sub> O pop-off and CPAP Generator (for connection to T-bar prongs)	10/box
BC161-10	Single-heated Infant Breathing Circuit with MR290 Chamber, Pressure Manifold at 17cmH <sub>2</sub> O pop-off and CPAP Generator (for connection to F&P Interface)	10/box
BC171-10	Single-heated Infant Breathing Circuit with MR290 Chamber, Pressure Manifold, Pressure Line and CPAP Generator (for connection to Flow Drivers)	10/box
BC461-SK	Bubble CPAP Kit with BC161, 5 BC181 FlexiTrunk™, 10 Bonnets and 30 Prongs of various sizes	assorted/box
BC471-SK	Bubble CPAP Kit with BC171, 5 BC181 FlexiTrunk™, 10 Bonnets and 30 Prongs of various sizes	assorted/box

## 3 F&P Bubble CPAP/Ventilator CPAP Interface and Accessories

PART NO	DESCRIPTION	QUANTITY	PART NO	DESCRIPTION	QUANTITY
BC180-05	FlexiTrunk™ Nasal Tubing (50mm)	5/box	BC6070-10	Nasal Prong (6.0mm/7.0mm)	10/box
BC181-05	FlexiTrunk™ Nasal Tubing (70mm)	5/box	BC6570-10	Nasal Prong (6.5mm/7.0mm)	10/box
BC182-05	FlexiTrunk™ Nasal Tubing (100mm)	5/box	BC800-10	Nasal Mask Small	10/box
BC190-05	FlexiTrunk™ Nasal Tubing (50mm) universal fitting	5/box	BC801-10	Nasal Mask Medium	10/box
			BC802-10	Nasal Mask Large	10/box
BC191-05	FlexiTrunk™ Nasal Tubing (70mm) universal fitting	5/box	BC803-10	Nasal Mask Extra Large	10/box
			BC300-05	Infant Bonnet (17-22cm)	5/box
BC192-05	FlexiTrunk™ Nasal Tubing (100mm) universal fitting	5/box	BC303-05	Infant Bonnet (22-25cm)	5/box
			BC306-05	Infant Bonnet (25-29cm)	5/box
BC3020-10	Nasal Prong (3.0mm/2.0mm)	10/box	BC309-05	Infant Bonnet (29-36cm)	5/box
BC3520-10	Nasal Prong (3.5mm/2.0mm)	10/box	BC325-05	Head Gear (29-36cm)	5/box
BC4030-10	Nasal Prong (4.0mm/3.0mm)	10/box	BC328-05	Head Gear (35-40cm)	5/box
BC4540-10	Nasal Prong (4.5mm/4.0mm)	10/box	BC331-05	Head Gear (40-45cm)	5/box
BC5040-10	Nasal Prong (5.0mm/4.0mm)	10/box	BC351-10	Chin Strap (20-26cm)	10/box
BC5050-10	Nasal Prong (5.0mm/5.0mm)	10/box	BC353-10	Chin Strap (26-32cm)	10/box
BC5550-10	Nasal Prong (5.5mm/5.0mm)	10/box	BC355-10	Chin Strap (32-38cm)	10/box
BC5560-10	Nasal Prong (5.5mm/6.0mm)	10/box	BC357-10	Chin Strap (38-44cm)	10/box
BC6060-10	Nasal Prong (6.0mm/6.0mm)	10/box			

## Infant Optiflow and Oxygen Therapy



The combination of Optimal Humidity and nasal cannula provides comfortable, effective oxygen delivery of greater oxygen flows than previously possible.

Infant Optiflow is delivered with the MR850 System and a range of nasal cannula that can deliver flows up to 8L/min.



## 1 Humidifier and Accessories

PART NO	DESCRIPTION	QUANTITY
MR850	MR850 Respiratory Humidifier	1/each
900MR806	Heater Wire Adapter for Inspiratory Heated Breathing Circuits	1/each
900MR868	Temperature/Flow Probe for MR850 Compatible Breathing Circuits (1.1m)	1/each
900MR208	Reflective Shield for Temperature/Flow Probe	20/box

## 2 Humidification Chambers

PART NO	DESCRIPTION	QUANTITY
MR290V	Auto-feed Adult/Infant Humidification Chamber with Vented Water Bag Spike	40/box

## 3 Breathing Circuits and Kits

PART NO	DESCRIPTION	QUANTITY
RT330	Optiflow Junior Tubing Kit for MR850	10/box
RT329	Single Limb Infant Breathing Circuit Kit with Pressure Relief Valve	10/box
RT202	Single Limb Adult Breathing Circuit Kit for Headbox/Oxyhood	10/box
RT102	Single Limb Adult Breathing Circuit (excl. Chamber) for Headbox/Oxyhood	10/box

## 4 Interfaces and Accessories

PART NO	DESCRIPTION	QUANTITY
OPT312	Optiflow Junior Nasal Cannula - Premature (Flow up to 8L/min)	20/box
OPT314	Optiflow Junior Nasal Cannula - Neonatal (Flow up to 8L/min)	20/box
OPT316	Optiflow Junior Nasal Cannula - Infant (Flow up to 20L/min)	20/box
OPT318	Optiflow Junior Nasal Cannula - Pediatric (Flow up to 25L/min)	20/box
BC2425-20	Nasal Cannula – Premature (Flow up to 6L/min)	20/box
BC2435-20	Nasal Cannula – Neonatal (Flow up to 6L/min)	20/box
BC2745-20	Nasal Cannula – Infant (Flow up to 7L/min)	20/box
BC2755-20	Nasal Cannula – Intermediate (Flow up to 7L/min)	20/box
BC3780-20	Nasal Cannula – Pediatric (Flow up to 8L/min)	20/box
OPT010	Optiflow Junior Wigglepad Replacement - Premature	20/box
OPT012	Optiflow Junior Wigglepad Replacement – Neonatal, Infant, Pediatric	20/box
BC402-05	Gas Supply Line Kit for Dual Flow Meters	5/box
BC115-05	Pressure Relief Valve	5/box

For more information please contact  
your local Fisher & Paykel Healthcare representative

**Manufacturer** 

Fisher & Paykel Healthcare Ltd  
15 Maurice Paykel Place  
East Tamaki, Auckland 2013

PO Box 14 348 Panmure  
Auckland 1741  
New Zealand

Tel: +64 9 574 0100  
Fax: +64 9 574 0158  
Email: [info@fphcare.com](mailto:info@fphcare.com)  
Web: [www.fphcare.com](http://www.fphcare.com)

**Australia**

Tel: +61 3 9879 5022  
Fax: +61 3 9879 5232

**Austria**

Tel: 0800 29 31 23  
Fax: 0800 29 31 22

**Benelux**

Tel: +31 40 216 3555  
Fax: +31 40 216 3554

**China**

Tel: +86 20 3205 3486  
Fax: +86 20 3205 2132

**France**

Tel: +33 1 6446 5201  
Fax: +33 1 6446 5221

**Germany**

Tel: +49 7181 98599 0  
Fax: +49 7181 98599 66

**India**

Tel: +91 80 4284 4000  
Fax: +91 80 4123 6044

**Irish Republic**

Tel: 1800 409 011

**Italy**

Tel: +39 06 7839 2939  
Fax: +39 06 7814 7709

**Japan**

Tel: +81 3 3661 7205  
Fax: +81 3 3661 7206

**Northern Ireland**

Tel: 0800 132 189

**Spain**

Tel: +34 902 013 346  
Fax: +34 902 013 379

**Sweden**

Tel: +46 8 564 76 680  
Fax: +46 8 36 63 10

**Switzerland**

Tel: 0800 83 47 63  
Fax: 0800 83 47 54

**Taiwan**

Tel: +886 2 8751 1739  
Fax: +886 2 8751 5625

**Turkey**

Tel: +90 312 354 34 12  
Fax: +90 312 354 31 01

**UK  (EU Authorised Representative)**

Fisher & Paykel Healthcare Ltd  
Unit 16, Cordwallis Park  
Clivemont Road, Maidenhead  
Berkshire SL6 7BU, UK

Tel: +44 1628 626 136  
Fax: +44 1628 626 146

**USA/Canada**

Tel: +1 800 446 3908  
or +1 949 453 4000  
Fax: +1 949 453 4001