



## Deluxe Steel Drop-Arm Commode Knocked-Down Frame

- Easy-to-assemble frame
- Easy-to-release drop-arm mechanism allows for safe lateral patient transfers to and from commode
- Durable plastic snap-on seat and lid
- Easy-to-clean gray powder-coated steel finish
- Plastic armrests and back provide extra comfort and support
- Comes with 7 qt commode bucket with carry handle and splash shield

### PRODUCTS

#### Description

11125KD-1	Commode Drop Arm Knock Down, 1/cs
11125KD-2	CommodeDropArm, 2/cs

### SPECIFICATIONS

#### Description

Width	19.75" between arms   23.5" Outside
Seat Dimensions	13.75" (W) x 16.5" (D)
Height	20" - 24"
Weight	18 lbs.
Carton Shipping Weight (1   2)	20   39 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime

### ACCESSORIES

#### Description

RTL12085	Biodegradable Sanitary Commode Liner, 12/bx, 6 bx/cs
----------	--