



Deluxe Steel Drop-Arm Commode with Padded Seat Knocked-Down Frame



FIGURE A



FIGURE B



FIGURE C



FIGURE D

- Easy-to-assemble without tools in 4 easy steps (Figure A - D)
- Padded open front vinyl toilet seat
- Easy-to-release drop-arm mechanism allows for safe lateral patient transfers to and from commode
- Easy-to-clean gray powder-coated steel finish
- Padded armrests and back provide extra comfort and support
- Comes with 12 qt commode bucket with carry handle and splash shield

PRODUCTS	Description
11125PSKD-1	Commode DropArm, 1/cs

SPECIFICATIONS	Description
Width	19" between arms 23.75" Outside
Seat Dimensions	15.25" (W) x 16.75" (D)
Height	19" - 23"
Weight	21 lbs.
Carton Shipping Weight	23 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime

ACCESSORIES	Description
RTL12085	Biodegradable Sanitary Commode Liner, 12/bx, 6 bx/cs

PARTS	Description
11125PSKD-AL	Left Arm With Release Lever, 1/ea
11125PSKD-AR	Right Arm With Release Lever, 1/ea
11125-PS	Padded Seat, 1/ea
11125L	Leg, 1/ea
12002PB	Push Button, 1/ea
10386GB	Tips, 1/ea
10386G-8	Replacement Tip, 11/8" Grey, 1pr. per retail box, 8bx/cs
11108-1	Pail, 1/ea
11107-1	Splash Guard, 1/ea
11125PSKD-B	Bacl, 1/ea
11125KD-HDW	Hardware, 1/ea