

Extension Sets Macrobre CLAVE®

List No. **12544-01**


Macrobre Extension Set – with non-DEHP tubing, micro-CLAVE® spin lock T-connector, slide clamp and female luer lock. Latex-Free.



Nominal Length	Drops/mL	Approx. Priming Capacity (mL)	Tubing I.D.	Case Size
5" (13 cm)	—	0.6	0.090"	50

List No. **12312-01**

LifeShield® Macrobre Extension Set – with CLAVE® bonded onto T-connector, slide clamp and locking spin collar. Latex-Free. Sterile Pack.



Nominal Length	Drops/mL	Approx. Priming Capacity (mL)	Tubing I.D.	Case Size
5" (13 cm)	—	0.99	0.100"	50

List No. **19328-48**

LifeShield® Macrobre Extension Set – with CLAVE®, slide clamp, locking spin collar T and prepierced injection site.



Nominal Length	Drops/mL	Approx. Priming Capacity (mL)	Tubing I.D.	Case Size
6" (15 cm)	—	0.67	0.100"	50

List No. **12645-28**

LifeShield® Macrobre Extension Set – with CLAVE® and locking spin collar T. Latex-Free. Non-DEHP. Sterile Pack.



Nominal Length	Drops/mL	Approx. Priming Capacity (mL)	Tubing I.D.	Case Size
5.5" (13.97 cm)	—	1.01	0.100"	50

List No. **14929-48**

LifeShield® Extension Set – with removable CLAVE®, female adapter, slide clamp and Secure Lock with male adapter. Latex-Free.



Nominal Length	Drops/mL	Approx. Priming Capacity (mL)	Tubing I.D.	Case Size
7" (17.75 cm)	—	0.7	0.090"	48

List No. **12551-01**

Macrobre Extension Set – with non-DEHP tubing, micro-CLAVE®, female luer lock, slide clamp and Option-Lok®. Latex-Free.



Nominal Length	Drops/mL	Approx. Priming Capacity (mL)	Tubing I.D.	Case Size
7" (18 cm)	—	0.7	0.090"	50