

Premium Guard - Static Air Mattress Overlay



14427



FIGURE B

- Unique construction allows for efficient air flow throughout overlay for use in the prevention and treatment of pressure ulcers
- Two adapters for easy inflation and deflation
- Emergency CPR air release valve deflates overlay quickly and safely
- Air flow holes reduce heat and moisture build up
- Adjustable bed straps prevent overlay from moving
- Easy-to-clean with mild detergent
- Electric or Manual pumps sold separately

PRODUCTS

Description

14428	Air Mattress Overlay, E0197, 1/ea
-------	-----------------------------------

SPECIFICATIONS

Description

Weight Capacity	250 lbs.
-----------------	----------

Warranty	90 Days Limited on Seams & Valve
----------	----------------------------------