

## Balanced Aire Adjustable Skin Protection Cushion



RETAIL BOX



INSET B



- The Balanced Aire Adjustable Cushion is our premier clinically effective yet affordable, Adjustable Skin Protection cushion
- Flexible interconnected, independent air cells allow for deep immersion into the cushion
- Ideal for clients with limited seated tolerance, or those at risk for pressure injuries
- Low-shear stretch cover with adjustable locking straps and non-skid bottom is easily removed and washed
- Hand pump included to adjust the cushion to individual patient comfort
- Available in 4" and 2" adjustable air cell heights
- Repair kit included for quick fixes
- Meets flammability standard TB117

### PRODUCTS

#### Description

8047-16	Adjustable Air Cushion, E2622, 1/ea
8047-16-2	Adjustable Air Cushion, E2622, 1/ea
8047-18	Adjustable Air Cushion, E2622, 1/ea
8047-18-2	Adjustable Air Cushion, E2622, 1/ea
8047-1818	Adjustable Air Cushion, E2622, 1/ea
8047-1818-2	Adjustable Air Cushion, E2622, 1/ea
8047-20	Adjustable Air Cushion, E2622, 1/ea
8047-20-2	Adjustable Air Cushion, E2622, 1/ea
8047-2018	Adjustable Air Cushion, E2622, 1/ea
8047-2018-2	Adjustable Air Cushion, E2622, 1/ea
8047-2218	Adjustable Air Cushion, E2623, 1/ea
8047-2218-2	Adjustable Air Cushion, E2623, 1/ea
8047-2220	Adjustable Air Cushion, E2623, 1/ea
8047-2220-2	Adjustable Air Cushion, E2623, 1/ea

### SPECIFICATIONS

#### Description

Cover Material	Fluid-Resistant Stretch Top   Non-Skid Bottom
Weight Capacity	500 lbs.
Warranty	2 Year Limited

### ACCESSORIES

#### Description

8047-PUMP	Pump for Balanced Air Cushion, 1/ea
8047Patchkit	Balance Aire Cushion Patch Kit, 1/ea
8047-16-C	Replacement Cover, 1/ea
8047-16-2-C	Replacement Cover, 1/ea
8047-18-C	Replacement Cover, 1/ea
8047-1818-C	Replacement Cover, 1/ea
8047-18-2-C	Replacement Cover, 1/ea

**ACCESSORIES**

**Description**

8047-1818-2-C	Replacement Cover, 1/ea
8047-20-C	Replacement Cover, 1/ea
8047-2018-C	Replacement Cover, 1/ea
8047-20-2-C	Replacement Cover, 1/ea
8047-2018-2-C	Replacement Cover, 1/ea