

## Pressure-Sensitive Chair and Bed Patient Alarm



13606

- Alerts care giver with audio alarm when patient gets out of a chair or bed
- Pressure-sensitive pad connects to the audio alarm
- The alarm comes with a safety alert, On/Off switch and low battery warning
- Comes with 9V Battery (Included)

PRODUCTS	Description
13605	1/ea
13606	1/ea
13608	1/ea
13609	1/ea

SPECIFICATIONS	Description
Pad Dimensions: (chair   bed)	10" x 15"   11" x 30"
Warranty	(alarm) 6 months

PARTS	Description
13610	Alarm w/o Reset, 1/ea
13611	Alarm with reset, 1/ea
13610PL	Bed Pad, 1/ea
13610PS	Seat Pad, 1/ea
13610	Alarm w/o Reset, 1/ea
13610PS	Seat Pad, 1/ea
13610	Alarm w/o Reset, 1/ea
13610PL	Bed Pad, 1/ea
13610PS	Seat Pad, 1/ea
13610PL	Bed Pad, 1/ea