

Heavy Duty 50 PSI Compressor

HCPCS: E0565



FIGURE A



FIGURE B



FIGURE C



FIGURE D

- Newly designed compressor vibration mounts, clamps and lock washers, and wire harnesses improve upon an already durable design
- Ultra-powerful compressor for both humidification and high flow, high-pressure nebulization of viscous medications
- Designed to provide efficient, continuous high-pressure performance
- Lightweight and easy-to-clean

PRODUCTS

Description

| | |
|-------|---------------------------------|
| 18450 | Heavy Duty 50PSI Nebulizer,1/ea |
|-------|---------------------------------|

SPECIFICATIONS

Description

| | |
|-----------------------------|-----------------|
| Electrical Requirement | 115 VAC ~ 60 Hz |
| Flow Rate | 11~41 LPM |
| Operating Temperature Range | 50½ - 140½ F |
| Power Consumption | 160 W |
| Sound | ~ 60 dBA |
| Warranty | Limited 1 Year |

ACCESSORIES

Description

| | |
|-----------|------------------------------------|
| MASK 007A | Adult Tracheotomy Mask 50/case, CS |
| MASK 007P | PediTracheotomy Mask 50/CS, CS |
| MASK 001P | Pediatric Aerosol mask CS/50, CS |
| MASK 001A | Adult Aerosol mask CS/50, CS |

PARTS

Description

| | |
|---------------|--|
| 18450-NVALVE | Regulator Valve, 18450, 1/ea, CS |
| 18450-RF | Rear Filter, 18450, 1/ea, CS |
| 18450-TUBE | PETubingSet18450,1pkg, CS |
| 18450-SWITCH | PowerSwitch18450,1ea, CS |
| 18450-T | TeeConnector18450,1ea, CS |
| 18450-VALVE | ReliefValve18450,1ea, CS |
| 18450-FAN | Cooling Fan, 18450, 1/ea, CS |
| 18450-CB | CircuitBreaker18450,1ea, CS |
| 18450-HEADKIT | Compressor Head Rebuild Kit,18450, 2/pkg, CS |
| 18450-GAUGE | PressureGauge18450,1ea, CS |
| 18450-ADAPTER | OutletAdapter,50PSI Compressor,1/ea, CS |
| 18450-CAP | FilterCap18450,1ea, CS |

[Click here to view detailed product information and pricing.](#)

Call 800-397-5899 to order today!

PARTS

Description

18450-FILTER

Intake Filter (front) 18450, 3/pkg, CS