




Circuits Heated



RT110




RT111




RT114

Fisher & Paykel Healthcare High Performance Circuits


Fisher & Paykel Healthcare offer a line of high-performance circuits that have spiral heated wire technology that is unique in the industry. This technology increases the heating surface area in the circuit for superior condensate control and allows for delivery of optimal humidity to patients. Fisher & Paykel High Performance Circuits are part of the MR850 humidification system and may also be used with the MR730 humidifiers. Fisher & Paykel circuits also feature parallel Y, built-in MDI port, proximal pressure line, adapter kits and water traps, depending on the circuit configuration.




RT130



RT131



RT132



RT137

Code	Htd/Non Htd	Circuit Length	For Use With Ventilator	Humidifiers	Heater Wire Adapter	Y Style	Temp Probe	MDI Port	Pressure Line	Dry Gas Line	H ₂ O Trap	Pressure Pop Off	Air Entrainer
RT110	D	5'	1, 2, 3, 4, 6, 7, 9, 12, 13, 30, 32, 33	MR850 MR730	900MR805 900MR900	P	900MR869 900MR569	•		•			
RT111	D	5'	5, 11	MR850 MR730	900MR805 900MR900	P	900MR869 900MR569	•	•	•			
RT114	S	5'	2, 29, 30, 31, 33, 34	MR850 MR730	900MR805 900MR900		900MR869 900MR569			•			
RT130	DS	4'	25	MR850 MR730	900MR805 900MR900	SW	900MR869 900MR569		•	•			
RT131	DS	4'	18*, 19*, 20, 21, 22, 23*, 24, 26*, 27, 28	MR850 MR730	900MR805 900MR900	SW	900MR869 900MR569		•	•	•		
RT132	S	4'	17	MR850 MR730	900MR805 900MR900		900MR869 900MR569		•	•			
RT137	DS	5'	14	MR850 MR730	900MR805 900MR900	P	900MR869 900MR569	•		•	•		

* Verify Bias Flow Protocol is > 4 lpm otherwise use RT236

Key

Adult Ventilator

1. Draeger Evita 1, 2, Dura, Evita 4, Evita XL, Savina
2. eVent Inspiration
3. GE Engstrom Carestation, I vent
4. Hamilton Amadeus, Galileo, Raphael, Veolar
5. Infrasonics Adult Star 2000
6. Maquet 300A, 900C, Servo i
7. Newport E 100M, Wave E200, e500 Wave
9. PB 840, 7200
11. CareFusion Bear® 1, 2, 5, 1000 ventilators
12. CareFusion Bird® 8400 ventilator
13. CareFusion AVEA®, VELA®, TBIRD® ventilators
14. Respiroics Esprit, PB740, PB760

CPAP/NPPV Systems

17. CareFusion SiPAP, nCPAP and P driver systems

Neonatal Ventilators

18. PB840
19. Draeger Babylog 8000, Evita
20. Infrasonics Infant Star 500, 950
21. Newport Breeze® E150
22. Sechrist IV100B, IV200, Millenium
23. CareFusion AVEA®, Bear Cub 750 PSV
24. CareFusion Bird® VIP ventilator
25. Maquet 300A, Servo i
26. Hamilton G5
27. Biomed Devices
28. Impact Eagle

Homecare Ventilators

29. PB Achieva, LP10, LP6
30. PB 540
31. Newport HT50, HT70
32. CareFusion Legacy ventilator
33. CareFusion LTV® ventilator
34. Respiroics Trilogy100, PLV102

- DS Dual-Single Heat
S Single Limb
Y Standard "Y"
P Parallel "Y"
SW Swivel "Y"

Cat. No.	Description
High Performance	
RT110	Vent circuit, dual heated
RT111	Vent circuit, dual heated with pressure line
RT114	Vent circuit, single heated
RT130	Neonatal ventilator circuit < 4 LPM
RT131	Neonatal Ventilator Circuit > 4 LPM
RT137	Vent circuit, dual limb, single heated
Special Circuit	
RT132	Infant continuous flow breathing circuit

Accessories	
900MR805	Heater wire adapter-dual limb, MR850
900MR806	Heater wire adapter-single limb, MR850
900MR869	Temperature probe, MR850
900MR900	Heater wire adapter for dual heated wire circuits for a MR730
900MR901	Heater wire adapter for single heated wire circuits for a MR730
900MR569	Temperature probe, MR730