

## Med-Aire Plus 8" Alternating Pressure and Low Air Loss Mattress System with 10" Defined Perimeter



FIGURE A

- Active combination therapy system provides both alternating pressure and low air loss to optimize pressure redistribution, shear/friction reduction, and microclimate control
- Designed to prevent, treat and heal pressure ulcers in the home or long term care setting
- Digital control unit offers 4 cycle times and adjustable comfort weight settings for optimal immersion and patient satisfaction
- 3 diagnostic audio and visual alarms alert for issues with air pressure, the system, or power
- 10" deep static perimeter tubing surrounds the length of the mattress preventing patient falls or entrapment
- Pillow feature provides head stability & comfort
- Advanced Clinical and Care Giver Features: Auto-Firm, Seat Inflate, microprocessing feedback, lock-out and more!
- Cell-on-cell designed air cells prevent "bottoming out" and provide up to 24 hours of power outage protection
- 8" deep air cells are easily removed with quick disconnects for convenient servicing
- Removable 4 way stretch PU cover is fluid resistant, low shear and vapor permeable helping protect the skin from friction and moisture
- CPR valve allows for rapid deflation
- Tubing quickly and easily disconnects with EZ Lock/Quick Release design

| PRODUCTS        | Description                                  |
|-----------------|--|
| 14029DP         | MedAire8"LowAirLoss.w/Defined,1/cs           |
| SPECIFICATIONS  | Description                                  |
| Weight Capacity | 450 lbs                                      |
| ACCESSORIES     | Description                                  |
| 14029DPMN       | Mattress Only,14029DP,3Tube System,1/ea, EA  |
| 14029DPC        | Cover for 14029DP, 1/ea, EA                  |
| 14029B          | Carry bag for 14029,1/ea, EA                 |
| 14030XP         | 3-Outlet Pump,14030/48/54/60,1 ea, EA        |
| 14029DP-CELL    | Bladder Only,14029DP,1/ea, EA                |
| 14427           | Pump, Electric, for 14425,1/bx, BX           |
| 14029XP         | Digital 3-Outlet Pump 14029/14029DP,1/ea, EA |

[Click here to view detailed product information and pricing.](#)

**Call 800-397-5899 to order today!**