

## Natura® Unit Dose Kits



### Natura® Stomahesive® Flexible Skin Barrier and Drainable Pouch Unit Dose Kit – Non-Sterile

**Kit includes:**

- 1 Stomahesive® Flexible Skin Barrier with white tape collar with cut-to-fit opening
- 1 drainable pouch (30.5 cm (12”) transparent with 1-sided comfort panel)
- 1 tail closure

Product Code	Flange Size	Color Match	Box of
401920	45 mm (1 3/4”)	●	5 Kits
401921	70 mm (2 3/4”)	●	5 Kits



### Natura® Stomahesive® Skin Barrier and Drainable Pouch Unit Dose Kit – Non-Sterile

**Kit includes:**

- 1 Stomahesive® Skin Barrier with cut-to-fit opening
- 1 drainable pouch (35.5 cm (14”) transparent with 1-sided comfort panel)
- 1 tail closure

Product Code	Flange Size	Color Match	Box of
401922	100 mm (4”)	●	5 Kits



### Natura® Durahesive® Flexible Skin Barrier Urostomy Pouch Unit Dose Kit– Non-Sterile

**Kit includes:**

- 1 Durahesive® Flexible Skin Barrier with white tape collar with cut-to-fit opening
- 1 urostomy pouch (25.4 cm (10”) transparent with 1-sided comfort panel) with Accuseal® Tap with valve
- 1 Accuseal® night drainage adaptor

Product Code	Flange Size	Color Match	Box of
401923	45 mm (1 3/4”)	●	5 Kits