


SAFELINE® Needleless System

Description	Product Code	REF Number	Units Per cs
-------------	--------------	------------	--------------

Extension Sets with SAFELINE Split Septum Injection Sites (continued)


Components on this page do not contain DEHP or natural rubber latex.



Female adapter, large bore tubing, slide clamp, SAFELINE split septum injection sites 28 in. and 6 in. above distal end, slide clamp, SPIN-LOCK® connector. Priming volume: 6.0 mL, Length: 35 in. (89 cm)


NF1370 NF1370 50

Filterflow® Filtered Add-on Extension Set with SAFELINE Split Septum Injection Site



Female luer lock, 0.2 micron filter, SAFELINE split septum injection site 6 in. above distal end, on/off clamp, SPIN-LOCK connector. Priming volume: 5.3 mL, Length 16 in. (41 cm)


NF1380F 473993 50



Female luer lock, 1.2 micron filter, SAFELINE split septum injection site 6 in. above distal end, on/off clamp, SPIN-LOCK connector. Priming volume: 5.1 mL, Length 16 in. (40.0 cm)


NF1385F 473998 50

Basic Administration Sets with SAFELINE Split Septum Injection Site



(15 drops/mL) Universal spike, roller clamp, SAFELINE split septum injection sites 40 in. and 6 in. above distal end, slide clamp, SPIN-LOCK connector. Priming volume: 19.0 mL, Length: 90 in. (228 cm)

NF1150 NF1150 50



(60 drops/mL) Universal spike, roller clamp, SAFELINE split septum injection sites 38 in. and 6 in. above distal end, slide clamp, SPIN-LOCK connector. Priming volume: 18 mL, Length: 94 in. (239 cm)

NF1190 NF1190 50

Secondary Set for use with SAFELINE Split Septum Injection Site



(15 drops/mL) Universal spike, roller clamp, SPIN-LOCK connector, SAFELINE Clip Lock® cannula, plastic bag hanger. Priming volume: 10 mL, Length: 49 in. (125 cm)

NF1430 50