



Universal Aluminum Crutch with Accessories



FIGURE A



FIGURE B



FIGURE C



RETAIL BOX

- Universal adjustment for use by a child, adult or tall adult (Figure A)
- Euro-Style clips and pushpin adjustment make it easy-to-adjust height of hand grip (Figure B)
- Lighter than wood, the aluminum crutch provides stability and durability
- Crutch disassembles for convenient storage (Figure C)
- Non-skid, jumbo-size vinyl contoured tips provide excellent traction
- Comfortable, durable underarm pad and hand grip

PRODUCTS

Description

RTL10433	Knocked Down, 1/pr, 8 pr/cs
----------	-----------------------------

SPECIFICATIONS

Description

Height (armrest to floor)	38.5" - 62.5"
---------------------------	---------------

Patient Height	4'6" - 6'6"
----------------	-------------

Weight Limit	300 lbs.
--------------	----------

Warranty	Limited Lifetime
----------	------------------