

SiloLiner™



SiloLiner™ Pin Liner

- 6mm gel thickness
- Easy to Donn and Doff: Reinforced cuff prevents liner from rolling down and irritating the skin
- Medical grade mineral oil gel helps moisturize and condition skin grafts and scars
- Outer material allows for wide range of motion
- Not made with natural rubber latex & hypoallergenic

- 1460** - Length 12" (30.5cm), Flat Width 3" (7.6cm)
- 1461** - Length 12" (30.5cm), Flat Width 4" (10.2cm)
- 1462** - Length 12" (30.5cm), Flat Width 5" (12.7cm)
- 1463** - Length 12" (30.5cm), Flat Width 6" (15.2cm)
- 1464** - Length 14" (35.5cm), Flat Width 3" (7.6cm)
- 1465** - Length 14" (35.5cm), Flat Width 4" (10.2cm)
- 1466** - Length 14" (35.5cm), Flat Width 5" (12.7cm)
- 1467** - Length 14" (35.5cm), Flat Width 6" (15.2cm)



SILOLINER™ PIN LINER SIZING GUIDELINES

ITEM#		FLAT WIDTH	DISTAL CIRCUMFERENCE <i>at urethane attachment</i>	PROXIMAL CIRCUMFERENCE MPT
1460	1464	3" (7.6cm)	8"-12" (20cm-30.5cm)	8"-14" (20cm-35.5cm)
1461	1465	4" (10.2cm)	9.5"-13.5" (24cm-34cm)	9.5"-16" (24cm-40.5cm)
1462	1466	5" (12.7cm)	12"-15" (30.5cm-38cm)	12"-20" (30.5cm-51cm)
1463	1467	6" (15.2cm)	14"-16" (35.5cm-40.5cm)	14"-26" (35.5cm-66cm)

All sizes are approximate.