



### Quad-Support Cane Tip



RETAIL BOX



RTL10351

- Provides balance, stability and allows your cane to be self-standing
- Replaces existing cane tip
- Four-tier design allows for better traction
- For use with 3/4" tubing
- Fits most manufacturers' canes
- For use with Drive cane model #'s: 10302, 10303, 10304, 10306, 10307, 10309, 10317, 10331, 10335, 10339, 10342, 10345, 10347, 10350, 10351, 10370, 10372, 10381 series

**PRODUCTS**

**Description**

RTL10351

Black, 1/bx