

Shower Chair with Folding Back



FIGURE A



INSET B



RETAIL BOX



FIGURE B



FIGURE C

- Adds safety and support in the bath tub and shower
- Easy, tool-free assembly and height adjustable legs
- Folding back (Figure A) and removable legs make this shower chair great for storage and transport
- Lightweight design with stable legs and non-slip feet
- Add comfort and convenience with optional arms (RTL12606ARMS) and soap trays (RTL12606TRAY) (Sold Separately)
- Retail packaged

PRODUCTS

	Description
RTL12606	BATH SEAT WITH BACK, 1/ea

SPECIFICATIONS

	Description
Outside Legs	17.375" – 18.875"
Seat Dimensions	21.5" (W) x 13.125" (L)
Seat Height	17"-23"
Weight	8 lbs.
Carton Shipping Weight	8.6 lbs.
Weight Limit	300 lbs.
Warranty	Limited lifetime

ACCESSORIES

	Description
RTL12606TRAY	Soap Tray for RTL12606 (1-Right, 1-Left), 1/bx (R/L), 12bx/cs