


Wedge Cushion with Stretch Cover



- Wedge design helps prevent patient from sliding forward
- Firm foundation reduces pressure caused by sling-seat hammocks
- The top cover is fluid resistant, low shear stretch fabric
- The bottom cover is waterproof
- High-density foam increases stability and comfort

PRODUCTS

	Description
8134	Wedge Cushion with Stretch Cover, 1/ea, 4/cs

SPECIFICATIONS

	Description
Cover Material	Fluid-Resistant Stretch Top Waterproof Vinyl Bottom
Weight Capacity	250 lbs.
Warranty	1 Year Limited