



SUCTION PRODUCTS



AG Industries Filters & Accessories

Bacteria Filters



BF405
Suction Bacteria Filter
1/4" to 1/2"
Barb to Barb
BFE: 99.999+%
VFE: 99.999+%



BF201
Suction Bacteria Filter
High Flow
3/8" to 1/2"
Barb to Barb
BFE: 99.999+%
VFE: 99.999+%



BF401
Suction Bacteria Filter
1/8" Screw Mount
1/4" to 1/2" Barb
BFE: 99.999+%
VFE: 99.999+%



BF200S
Suction Bacteria Filter
1/4" to 3/8"
Barb to Barb
BFE: 99.999+%
VFE: 99.999+%



BF400
Suction Bacteria Filter
1/8" Screw Mount
1/4" to 3/8" Barb
BFE: 99.999+%
VFE: 99.999+%

Suction Tubing



AG615725
Blue Tip Suction
Tubing 72"



AG615730
Silicone Tubing 4"



AG615473
Blue Tip Suction
Tubing 16"

Connectors



AG7305D-023
Elbow Connector
(Drive/DeVilbiss®)



610E
Elbow Connector
Universal



PM1616
Elbow Connector
(Precision Medical®)

AG Industries Kits

Replacement Filter Kits



AG5604
Drive/DeVilbiss®
Replacement Filter
Kit Includes:
BF201, AG615730,
AG7305D-023



AG201SET
Replacement
Filter Kit Includes:
BF201, 5" & 6"
Connecting Hoses



AG400SET
Schuco® Replacement
Filter Kit Includes:
BF400, AG615473



AG502690
Precision Medical®
Power Vac Replacement
Filter Kit Includes:
BF200S,
2" & 14" Connecting Hoses



AG401SET
Schuco® Replacement
Filter Kit Includes:
BF401, AG615473



AG502012
Precision Medical®
Easy Vac Replacement
Filter Kit Includes:
BF200S, PM1616,
2" & 14" Connecting Hoses



AG402SET
Invacare® IRC1135
Replacement
Filter Kit Includes:
BF400, 610E,
16" Connecting Hose



AG503293
Precision Medical®
Easy Go Replacement
Filter Kit Includes:
BF400, PM1616,
4" Connecting Hose



AG405SET
Schuco® Replacement
Filter Kit Includes:
BF405, AG610E,
2 - 8" Connecting Hoses

Full Suction Kits



AG5603
Full Suction
Kit Includes:
S1160B,
AG615725,
AG5604
Filter Set



AGT5606
Full Suction
Kit Includes:
S1160B,
AG615725,
AG402SET
Filter Set



Suction Canisters (A7000)



Suction Canisters



S1160B
800 cc.
Disposable Canister,
Float Valve Shutoff
White Top
12/Case



CAN800-48
800 cc.
Disposable Canister,
Float Valve Shutoff
White Top
48/Case



AG434410
1200 cc. Hydrophobic
Disposable Canister,
Filter Shutoff
Red Top,
12/Case



AG784410
800 cc. Hydrophobic
Disposable Canister,
Filter Shutoff
Red Top
10/Case



AG484410
1200 cc. Hi-Flow
Disposable Canister,
Float Valve Shutoff
Green Top
12/Case



AG3695
1100 cc. Hydrophobic
Disposable Canister,
Filter Shutoff
Orange Top
12/Case



AG494410
2000 cc. Hi-Flow
Disposable Canister,
Float Valve Shutoff
Green Top
12/Case

Suction Catheters (A4605, A4624)

Suction Catheters & Tubing Kits



Q32001
Bulb Tip, Non-Vent
Yankauer Catheter,
Blue



Q320016
Bulb Tip, Non-Vent
Yankauer Catheter
with 6' Blue Tip Tubing



Q32002
Bulb Tip, with Vent
Yankauer Catheter,
Blue



Q320026
Bulb Tip, with Vent
Yankauer Catheter
with 6' Blue Tip Tubing



Q32013
Bulb Tip, with Vent
Yankauer Catheter,
Clear



Q320136
Bulb Tip, with Vent
Yankauer Catheter
with 6' Blue Tip Tubing



Q32014
Bulb Tip, Non-Vent
Yankauer Catheter,
Clear



Q320146
Bulb Tip, Non-Vent
Yankauer Catheter
with 6' Blue Tip Tubing



Tracheostomy Accessories

AG Industries Trach. Accessories



AG61075
Adult Tracheostomy Mask



AG64550
Drainage Bag



AG002002
Jet Nebulizer with Air Entrainment



AG8600
Corrugated Tubing, Segmented Every 6"
100Ft/Roll

Portex® HME



AG580011
ThermoVent® 600 HME
OEM# 580011

Applications: Heat and passive humidification for use between patient and circuit

Moisture Output: 27 mg H₂O/L at 34° C for 1000 ml tidal volume

Flow Resistance: 2.9 cm H₂O at 60 i/ml

Media: Corrugated Paper

Dead Space: 11.3 ml

Connectors: 15 mm I.D. x 22 mm O.D. on patient end, 15 mm O.D. on circuit end



AG570016
ThermoVent® T HME
OEM# 570016

Applications: Heat and passive humidification for use between patient and circuit

Moisture Output: 25 mg H₂O/L at 34° C for 1000 ml tidal volume

Flow Resistance: 2.0 cm H₂O at 60 i/ml

Media: Corrugated Paper

Connectors: 15 mm I.D. raised fenestration on patient end