

# **Product Guide**

Continence, Ostomy, Wound and Skin Care Effective June 2020

Keagan

# Making life easier.

"My stoma doesn't define me. My friends know about it and support me. I feel comfortable when we are hanging out." Keagan

#### About Keagan

"What will my friends think?" This was Keagan's first thought as an ostomate. Having been diagnosed with Crohn's disease at the age of 15, Keagan was admitted to the hospital with a perforated colon and quickly underwent ostomy surgery.

Like any teenager, being able to be 'normal' and to fit in is really important. Being able to do all the great stuff that life offers – just like everybody else, is also a real concern. Although these thoughts were a real worry at first, and even though it hasn't been an easy journey for her, circumstances have changed for the better for Keagan.

Today, Keagan has established a good routine. Importantly, those around her are aware of her condition.

Keagan is a SenSura<sup>®</sup> Mio user who has received compensation from Coloplast to provide this information. Each person's situation is unique and experiences may not be the same. Users should talk to their healthcare professional about which product may be right for them.

# Welcome to Coloplast

Coloplast is a leader in the market within intimate healthcare. For more than sixty years, we've worked closely with healthcare professionals and users to develop products and services that make life easier for people with intimate healthcare needs.

Our Mission Making life easier for people with intimate healthcare needs

# Making life **easier**

# Table of Contents

Continence Care	7-30
Straight Tip Intermittent Catheters	10-13
Coudé Tip Intermittent Catheters	14-15
Intermittent Catheters Closed Systems / Sets	16-19
Transanal Irrigation	20-21
Male External Catheters	22-24
Leg and Drain Bags	25-27
Indwelling Catheters	28-29
Accessories	

Wound	Care	31	-44
-------	------	----	-----

Foam Dressings	35-37
Alginate Dressings	
Hydrocolloid Dressings	
Hydrophilic Dressing / Wound Cleansers / Hydrog	els 42
Contact Layer	
Super Absorbents	

#### 

Cleansers	50-52
Moisturizers	
Skin Protectants	54-55
Antifungals	
Skin Fold Management	
Hand Hygiene	
Accessories	

Ostomy Care	61-184
Accessories	67-74

Irrigation	75-76
2-Piece Click	
2-Piece Flex	105-122
1-Piece	123-160
Post-op / Wound Management	161-172
Pediatric	173-184

Reference Tools	185-200
Education	
Continence Care	
Ostomy Care	188-189
Continence and Ostomy Education	
Wound & Skin Care	
Medicare Utilization	
Continence	
Ostomy	193-194
Wound Care	195
HCPCS Codes	196-200
Index	201

# **Product Guide**

Continence Care | Effective June 2020



SpeediCath<sup>®</sup> catheters are prescribed for use for patients who require bladder drainage due to urine retention or post void residual volume (PVR). Before use, carefully read all of the instructions. Call your doctor if you think you have a UTI or can't pass the catheter into the bladder. For more information regarding risks, potential complications and product support, call Coloplast Corp. at 1-866-226-6362 and/ or consult the company website at www.coloplast.us. Caution: Federal law restricts this device to sale by or on the order of a physician.

Self-Cath<sup>®</sup> Plus Intermittent Urinary Catheters are intended for use in male or female patients requiring bladder drainage as determined by their physician. This device is indicated for those individuals unable to promote a natural urine flow or for those who have a significant volume of residual urine following a natural bladder-voiding episode. The device can be used by the patient once appropriate training has taken place or by a trained healthcare professional. Single-use only. See the device manual for detailed information regarding the procedure, contraindications, warnings, precautions, and potential complications/adverse events. For further information, call Coloplast Corp. at 1-866-226-6362 and/or consult the company website at www.coloplast.us. Caution: Federal law restricts this device to sale by or on the order of a physician.

Self-Cath<sup>®</sup> Intermittent Catheters are intended for use in male, female, and pediatric patients (neonates, infants, children, adolescents, and transitional adolescents) requiring bladder drainage as determined by their physician. This device is indicated for those individuals unable to promote a natural urine flow or for those individuals who have a significant volume of residual urine following a natural bladder-voiding episode. Single-use only. See the device manual for detailed information regarding the procedure, contraindications, warnings, precautions, and potential complications/ adverse events. For further information, call Coloplast Corp. at 1-866-226-6362 and/ or consult the company website at www.coloplast.us. Caution: Federal law restricts this device to sale by or on the order of a physician.

SureCath® Catheters are indicated for use by patients for intermittent catheterization for the purpose of bladder drainage. See the device manual for detailed information regarding the instructions for use, contraindications, warnings, precautions, and potential complications/adverse events. For further information, call Coloplast Corp. at 1.866-226-6362 and/or consult the company website at www.coloplast.us. Caution: Federal law restricts this device to sale by or on the order of a physician.

Folysil<sup>®</sup> Indwelling Catheters empty urine from the bladder through the urethra via a small cut through the abdomen above the pubic bone (suprapubic) during a surgical procedure. The choice of the size, balloon volume, and catheter material is the responsibility of your healthcare professional based on your condition and the length of time to be used. Coloplast's indwelling catheters are available by prescription only. Coloplast's indwelling catheters are for single-use only; discard it after use. If you experience symptoms of a urinary tract infection, are unable to pass urine, or the balloon will not deflate, contact your healthcare professional. The risk information provided here is not comprehensive. To learn more, talk to your healthcare professional to understand the risks and benefits to determine if Folysil is right for you. Although many patients benefit from the use of this device, results may vary. For further questions, call Coloplast Corp. at 1-866-226-6362 and/or consult the company website at www.coloplast.us.

Indication: The Peristeen<sup>®</sup> Transanal Irrigation System is intended to instill water into the colon through a rectal catheter, which incorporates an inflatable balloon, inserted into the rectum to promote evacuation of the contents of the lower colon. The Peristeen Transanal Irrigation System is indicated for use by children (2 years - <12 years old), adolescent (12 years - <18 years old), transitional adolescent (18 - <21 years old) and adult patients with neurogenic bowel dysfunction who suffer from fecal incontinence, chronic constipation, and/or time-consuming bowel management procedures.

**Contraindications:** Peristeen Transanal Irrigation must not be used in the following situations: known anal or colorectal stenosis, colorectal cancer, radiotherapy to the pelvis, and recent abdomino-perineal surgery, active inflammatory bowel disease, acute diverticulitis, severe diverticulitis, previous diverticulitis and diverticular abscess, chronic symptomatic diverticular disease, within 3 months of abdominal, anal or colorectal surgery, within 4 weeks of endoscopic polypectomy, recent colonic biopsy, recent endoscopic mucosal resection (EMR) and recent endoscopic sub-mucosal dissection (ESD), severe autonomic dysreflexia, ischemic colitis, during spinal cord shock phase, complex diverticular disease, in patients who are pregnant and have not used the system before\*. Since the list is not exhaustive, the healthcare professional should always consider individual patient factors as well.

\*If the patient is pregnant and has never used transanal irrigation before, they should not start the irrigation procedure during pregnancy.

Warning: Peristeen Transanal Irrigation procedure should always be carried out with care. Bowel perforation is an extremely rare, but serious and potentially lethal complication to transanal irrigation and will require immediate admission to a hospital, often requiring surgery. See the device manual for complete user instructions, contraindications, warnings, precautions, and potential complications/adverse events. For further information, call Coloplast Corp at 1-855-605-7594 or consult the company website at www.coloplast.us. Caution: Federal law restricts this device to sale by or on the order of a physician.

## Ensure your patients get the product you recommend. Ask for SpeediCath® and Self-Cath® by brand name.

Intermittent catheter users who have insurance can pay the same price when buying any catheter within the same category: straight tip,<sup>1</sup> coudé tip,<sup>2</sup> or catheter sets.<sup>3</sup>



\* Dispense as written; exact language to specify brand dispensed varies by state. Contact your Coloplast Sales Rep or National Association of Boards of Pharmacy for more information.

Advise your patients to find a supplier who accepts assignment for the intermittent catheter you prefer by contacting Coloplast Consumer Care at 1-866-226-6362.

#### Steps to ensure patients get the catheter they want:

1. Include the following information on each prescription:

Patient Nam Address	ie		Type of catheter (straight, coudé, kit / closed system)	-
Dx: Urir	Cath Compact Se ne retention ters / day x 30 do		0 catheters / month	
Brand of catheter			atheterizations / day and needed / month	
Refill	Times	NPI Number		
Date		Print Last No	ame	

#### 2. Enroll your patients in Coloplast® Care

Call - contact Coloplast Care at 1-866-226-6362 Email - send request to care-us@coloplast.com Fax - send patient enrollment form to 1-855-676-2594 Website - www.ColoplastCareEnroll.com

1. HCPCS A4351 (straight)\*\*

- 2. HCPCS A4352 (coudé)\*\*
- 3. HCPCS A4353 (catheter set)\*\*

#### SpeediCath® Standard

The world's first ready-to-use hydrophilic-coated intermittent catheter

Easy to carry, store and use. The catheter is pre-hydrated and ready-to-use, eliminating wait time. The hydrophilic coating is smooth and uniform, designed to minimize urethral friction. The eyelets are polished for comfort. The uncoated Grip Zone™ area is created to offer control in handling. **This catheter is not made with natural rubber latex and is DEHP and phthalate free.** 

#### SpeediCath® Standard

SpeediCath catheter only. No insertion supplies.

#### Straight Tip

Code	Size	Units	Description	HCPCS*
28408	• 8 FR	30	14" length	A4351
28410	• 10 FR	30	14" length	A4351
28412	O 12 FR	30	14" length	A4351
28414	• 14 FR	30	14" length	A4351
28416	🛑 16 FR	30	14" length	A4351
28418	• 18 FR	30	14" length	A4351
Female S	Straight Tip			
28506	• 6 FR	30	6" length	A4351
28508	• 8 FR	30	6" length	A4351
28510	• 10 FR	30	6" length	A4351
28512	O 12 FR	30	6" length	A4351
28514	• 14 FR	30	6" length	A4351
28516	🗕 16 FR	30	6" length	A4351
Pediatri	c Straight Ti	р		
28706	• 6 FR	30	6" length	A4351
28708	• 8 FR	30	6" length	A4351
28710	• 10 FR	30	6" length	A4351
Boy Stro	aight Tip			
28606	• 6 FR	30	10" length	A4351
28608	• 8 FR	30	10" length	A4351
28610	• 10 FR	30	10" length	A4351
28612	O 12 FR	30	10" length	A4351

#### SpeediCath<sup>®</sup> Soft

SpeediCath® Soft is designed to make catheterization simple and easy with fast preparation and a confident grip.

- Pre-lubricated intermittent catheter that is instantly ready to use with a unique hydrophilic coating, eliminating the need for additional lubricants or waiting for a catheter to hydrate.
- Insertion grip made from a soft, flexible material designed to ensure a confident and hygienic grip of the catheter, without touching the tube or tip.
- Soft catheter that is easy to use the first time and every time
- Not made with natural rubber latex and is DEHP and phthalate free.

#### SpeediCath<sup>®</sup> Soft

SpeediCath Soft straight tip intermittent catheter.

Code	Size	Units	Description	HCPCS*
27010	• 10 FR	30	13" length	A4351
27012	O 12 FR	30	13" length	A4351
27014	• 14 FR	30	13" length	A4351
27016	🗕 16 FR	30	13" length	A4351

#### SpeediCath<sup>®</sup> Compact Female and Compact Plus

The most compact, ready-to-use hydrophilic-coated intermittent catheters for women.

SpeediCath Compact Female and Compact Plus offer the benefits of our ready-to-use hydrophilic-coated intermittent catheters in a small, discreet package. Designed for women, the length of SpeediCath Compact Female matches the length of the female urethra, and when unopened it is no longer than a tube of lipstick. This catheter is not made with natural rubber latex and is DEHP and phthalate free.

#### SpeediCath<sup>®</sup> Compact Female

SpeediCath Compact straight tip intermittent catheter.

Code	Size	Units	Description	HCPCS*
28578	• 8 FR	30	2.75" length	A4351
28580	• 10 FR	30	2.75" length	A4351
28582	O 12 FR	30	2.75" length	A4351
28584	• 14 FR	30	2.75" length	A4351

#### SpeediCath<sup>®</sup> Compact Plus

## SpeediCath Compact Female with additional catheter length for women who want discretion but require a longer catheter.

Code	Size	Units	Description	HCPCS*
28810	• 10 FR	30	3.5" length	A4351
28812	O 12 FR	30	3.5" length	A4351
28814	• 14 FR	30	3.5" length	A4351

\*For important information regarding the use of HCPCS codes listed and code definitions see page 192 | For important safety information see page 8

#### Self-Cath<sup>®</sup> Plus

Self-Cath<sup>®</sup> Plus intermittent catheters have fire-polished eyelets designed for comfort and a siliconized surface created for smooth insertion. The catheter has a hydrophilic coating that activates immediately upon exposure to water. **This catheter is not made with natural rubber latex.** 

#### Self-Cath<sup>®</sup> Plus

Hydrophilic-coated intermittent catheter.

#### Straight Tip

Code	Size	Units	Description	HCPCS*
4408	• 8 FR	30	16" length, funnel end	A4351
4410	• 10 FR	30	16" length, funnel end	A4351
4412	O 12 FR	30	16" length, funnel end	A4351
4414	• 14 FR	30	16" length, funnel end	A4351
4416	🗕 16 FR	30	16" length, funnel end	A4351
4418	🖲 18 FR	30	16" length, funnel end	A4351
Female	Straight Tip			
4210	• 10 FR	30	6" length, funnel end	A4351
4212	O 12 FR	30	6" length, funnel end	A4351
4214	• 14 FR	30	6" length, funnel end	A4351
Pediatri	c Straight Ti	р		
4306	• 6 FR	30	10" length, funnel end	A4351
4308	• 8 FR	30	10" length, luer end	A4351
4310	• 10 FR	30	10" length, luer end	A4351

#### Self-Cath<sup>®</sup> Plus Soft

The Self-Cath® Plus Soft is the same as Self-Cath® Plus but designed for individuals that prefer a softer catheter.

#### Straight Tip

Code	Size	Units	Description	HCPCS*
4110	• 10 FR	30	16" length, funnel end	A4351
4112	O 12 FR	30	16" length, funnel end	A4351
4114	• 14 FR	30	16" length, funnel end	A4351
4116	🗕 16 FR	30	16" length, funnel end	A4351

WARNING: This product can expose you to Di(2-ethylhexyl)phthalate
(DEHP), which is known to the State of California to cause cancer and birth
defects or other reproductive harm. For more information go to www.
P65Warnings.ca.gov

#### Self-Cath®

Self-Cath intermittent catheters have fire-polished eyelets designed for comfort, a siliconized surface created for smooth insertion, and are available in a variety of options. This catheter is not made with natural rubber latex.

#### Self-Cath®

#### Straight Tip

Code	Size	Units	Description	HCPCS*
408	• 8 FR	30	16" length, funnel end	A4351
410	• 10 FR	30	16" length, funnel end	A4351
412	O 12 FR	30	16" length, funnel end	A4351
414	• 14 FR	30	16" length, funnel end	A4351
416	🗕 16 FR	30	16" length, funnel end	A4351
418	🛑 18 FR	30	16" length, funnel end	A4351
450	• 14 FR	50	16" length, funnel end, curved packaging	A4351
460	O 12 FR	50	16" length, funnel end, curved packaging	A4351

#### Female Straight Tip

Code	Size	Units	Description	HCPCS*
208	• 8 FR	30	6" length, funnel end	A4351
210	• 10 FR	30	6" length, funnel end	A4351
212	O 12 FR	30	6" length, funnel end	A4351
214	• 14 FR	30	6" length, funnel end	A4351
240	• 14 FR	30	6" length, luer end	A4351

#### Pediatric Straight Tip

Code	Size	Units	Description	HCPCS*
305	O 5 FR	30	10" length, funnel end	A4351
306	• 6 FR	30	10" length, funnel end	A4351
308	• 8 FR	30	10" length, luer end	A4351
310	• 10 FR	30	10" length, luer end	A4351

#### Self-Cath<sup>®</sup> Soft

The Self-Cath $^{\otimes}$  Soft is the same as Self-Cath $^{\otimes}$  but designed for individuals that prefer a softer catheter.

<b>.</b>			
Stra		ht	lin
<b>J</b> U 4			- <b>P</b>
	-		

Code	Size	Units	Description	HCPCS*
110	• 10 FR	30	16" length, funnel end	A4351
112	O 12 FR	30	16" length, funnel end	A4351
114	• 14 FR	30	16" length, funnel end	A4351
116	🛑 16 FR	30	16" length, funnel end	A4351

A WARNING: This product can expose you to Di(2-ethylhexyl)phthalate (DEHP), which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www. P65Warnings.ca.gov

\*For important information regarding the use of HCPCS codes listed and code definitions see page 192 | For important safety information see page 8

000

#### SpeediCath<sup>®</sup> Flex Coudé Pro

Every detail makes it easy. Hydrophilic coated intermittent catheter that is instantly ready-touse. Flexible, soft catheter that is easy to handle with a curved tip which can bend to follow the curves of the urethra. Hygienic with a protective dry sleeve and soft-squeeze grip. Discreet standard length catheter packaged in two compact sizes. This catheter is not made with natural rubber latex and is DEHP and phthalate free.



Code	Size	Units	Description	HCPCS*
20010	• 10 FR	30	13" length	A4352
20012	O 12 FR	30	13" length	A4352
20014	• 14 FR	30	13" length	A4352
20016	🗕 16 FR	30	13" length	A4352
20021	• 10 FR	30	13" length Pocket	A4352
20022	O 12 FR	30	13" length Pocket	A4352
20024	• 14 FR	30	13" length Pocket	A4352

#### SpeediCath® Standard Coudé Tip

The world's first ready-to-use hydrophilic-coated intermittent catheter

This tapered coudé tip catheter is easy to carry, store and use. It is pre-hydrated and ready-to-use, eliminating wait time. The hydrophilic coating is smooth and uniform, designed to minimize urethral friction. The eyelets are polished and designed for comfort. The uncoated Grip Zone<sup>™</sup> area is designed to offer control in handling. **This catheter is not made with natural rubber latex and is DEHP and phthalate free**.

Speedic	atn° Stan	aara Tapere	a Coude Tip	
Code	Size	Units	Description	

Code	Size	Units	Description	HCPCS*
28490	• 10 FR	30	14" length	A4352
28492	O 12 FR	30	14" length	A4352
28494	• 14 FR	30	14" length	A4352
28496	🗕 16 FR	30	14" length	A4352

\*For important information regarding the use of HCPCS codes listed and code definitions see page 192 | For important safety information see page 8

#### Self-Cath<sup>®</sup> Plus Coudé Tip

This coudé tip intermittent catheter has fire-polished eyelets designed for comfort and a siliconized surface created for smooth insertion. The catheter has a hydrophilic coating that activates immediately upon exposure to water. The blue Guide Stripe® along the top of the catheter allows monitoring of the coudé tip throughout the catheterization process. This catheter is not made with natural rubber latex.

#### Self-Cath<sup>®</sup> Plus Tapered Coudé Tip with Guide Stripe<sup>®</sup>

Code	Size	Units	Description	HCPCS*
4608	• 8 FR	30	16" length, funnel end	A4352
4610	• 10 FR	30	16" length, funnel end	A4352
4612	O 12 FR	30	16" length, funnel end	A4352
4614	• 14 FR	30	16" length, funnel end	A4352

0 - 10

10

#### Self-Cath<sup>®</sup> Plus Olive Coudé Tip with Guide Stripe<sup>®</sup>

Code	Size	Units	Description	HCPCS*
4808	• 8 FR	30	16" length, funnel end	A4352
4810	• 10 FR	30	16" length, funnel end	A4352
4812	O 12 FR	30	16" length, funnel end	A4352
4814	• 14 FR	30	16" length, funnel end	A4352
4816	🗕 16 FR	30	16" length, funnel end	A4352

#### Self-Cath<sup>®</sup> Coudé Tip

Self-Cath coudé tip intermittent catheters have fire-polished eyelets designed for comfort, a siliconized surface created for smooth insertion. The blue Guide Stripe® along the top of the catheter allows monitoring of the coudé tip throughout the catheterization process. They are available in a variety of options. **This catheter is not made with natural rubber latex.** 

#### Self-Cath<sup>®</sup> Olive Coudé Tip with Guide Stripe<sup>®</sup>

Code	Size	Units	Description	HCPCS*
806	• 6 FR	30	16" length, funnel end	A4352
808	• 8 FR	30	16" length, funnel end	A4352
810	• 10 FR	30	16" length, funnel end	A4352
812	O 12 FR	30	16" length, funnel end	A4352
814	• 14 FR	30	16" length, funnel end	A4352
816	🗕 16 FR	30	16" length, funnel end	A4352
818	• 18 FR	30	16" length, funnel end	A4352

#### Self-Cath<sup>®</sup> Tapered Coudé Tip with Guide Stripe<sup>®</sup>

Code	Size	Units	Description	HCPCS*
608	• 8 FR	30	16" length, funnel end	A4352
610	• 10 FR	30	16" length, funnel end	A4352
612	O 12 FR	30	16" length, funnel end	A4352
614	• 14 FR	30	16" length, funnel end	A4352

 WARNING: This product can expose you to Di(2-ethylhexyl)phthalate (DEHP), which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov Code

28702

#### SpeediCath<sup>®</sup> Compact Male with SpeediBag<sup>®</sup>

The most compact catheter for men.

Hydrophilic coated catheter that is pre-hydrated and readyto use. This catheter is hygienic with a no-touch insertion. Its telescoping mechanism allows expansion from its small modest tube into a full length intermittent catheter for males. A 750 mL SpeediBag is included. **This catheter is not made with natural rubber latex and is DEHP and phthalate free**.

#### SpeediCath<sup>®</sup> Compact Male with SpeediBag<sup>®</sup>

O 12/18 FR

Size

		9.00 ° ° ° ° ° ° ° ° ° ° ° ° ° ° ° ° ° °
Units	Description	HCPCS*
20	13.2" length	A4353

includes funnel

#### SpeediCath® Compact Set

The most compact set solution. An instantly ready-to-use hydrophiliccoated intermittent catheter in a compact design with an integrated



#### reddot design award winner 2013

bag available for both males and females. Provides ease of use and discretion and thanks to its award winning simple and compact design it is small enough to keep in a handbag or pocket. This catheter is not made with natural rubber latex and is DEHP and phthalate free.

#### SpeediCath<sup>®</sup> Compact Set



#### Female

Code	Size	Units	Description	HCPCS*
28520	• 10 FR	20	3.5" length	A4353
28522	O 12 FR	20	3.5" length	A4353
28524	• 14 FR	20	3.5" length	A4353
Male				
28422	O12/18 FR	20	13.2" length includes funnel	A4353

#### SpeediCath® Standard with Insertion Supplies

The world's first ready-to-use hydrophilic-coated intermittent catheter.

Easy to carry, store and use. The catheter is pre-hydrated and ready-to-use, eliminating wait time. The hydrophilic coating is smooth and uniform, designed to minimize urethral friction. The eyelets are polished for comfort. The uncoated Grip Zone<sup>TM</sup> area is designed to offer control in handling. This catheter is not made with natural rubber latex and is DEHP and phthalate free.

#### SpeediCath® Standard with Insertion Supplies

SpeediCath Standard intermittent catheter with one pair of non-latex exam gloves, one waterproof drape, two iodine prep pads and urine collection bag.

#### Straight Tip

Code	Size	Units	Description	HCPCS*
28481	• 8 FR	100	14" length	A4353
28482	• 10 FR	100	14" length	A4353
28483	O 12 FR	100	14" length	A4353
28484	• 14 FR	100	14" length	A4353
28485	🗕 16 FR	100	14" length	A4353
28486	• 18 FR	100	14" length	A4353
Female S	Straight Tip			
28499	• 6 FR	100	6" length	A4353
28500	• 8 FR	100	6" length	A4353
28501	• 10 FR	100	6" length	A4353
28502	O 12 FR	100	6" length	A4353
28503	• 14 FR	100	6" length	A4353
28504	🗕 16 FR	100	6" length	A4353

#### SureCath<sup>®</sup> Set

For intermittent self-catheterization.

Hydrophilic-coated intermittent catheter set with sterile solution enclosed in a single-use, collection bag. The catheter is available with and without insertion supplies. **This catheter is not made with natural rubber latex**.

#### SureCath<sup>®</sup> Set with Insertion Supplies

SureCath Set intermittent catheter and collection bag, one pair of exam gloves, one CSR wrap and two iodine prep pads.

#### Straight Tip 700mL

Code	Size	Units	Description	HCPCS*
28022S	• 8 FR	100	14" length	A4353
28023S	• 10 FR	100	14" length	A4353
28024S	O 12 FR	100	14" length	A4353
28032	• 14 FR	20	14" length	A4353
28032S	• 14 FR	100	14" length	A4353
28033S	🗕 16 FR	100	14" length	A4353

#### Straight Tip 1200 mL

Code	Size	Units	Description	HCPCS*
28036	• 14 FR	20	14" length	A4353
28037	🗕 16 FR	20	14" length	A4353
28036S	• 14 FR	100	14" length	A4353
28037S	🗕 16 FR	100	14" length	A4353

#### SureCath<sup>®</sup> Set without Insertion Supplies

SureCath hydrophilic-coated intermittent catheter and collection bag only. No additional insertion supplies.

#### Straight Tip 700 mL

Code	Size	Units	Description	HCPCS*
28005S	• 8 FR	100	14" length	A4353
28006S	• 10 FR	100	14" length	A4353
28007S	O 12 FR	100	14" length	A4353
28008S	• 14 FR	100	14" length	A4353
28009	🛑 16 FR	20	14" length	A4353
28009S	🖲 16 FR	100	14" length	A4353
Straight T	<b>ip</b> 1200 mL			
28026S	• 10 FR	100	14" length	A4353
28027	O 12 FR	20	14" length	A4353
28027S	O 12 FR	100	14" length	A4353
28028	• 14 FR	20	14" length	A4353
28028S	• 14 FR	100	14" length	A4353
28029S	🗕 16 FR	100	14" length	A4353

WARNING: This product can expose you to chemicals including Di(2ethylhexyl) phthalate (DEHP) and Methyl isobutyl ketone (MIBK), which are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov

\*For important information regarding the use of HCPCS codes listed and code definitions see page 192 | For important safety information see page 8

#### Self-Cath<sup>®</sup> Closed System

The Self-Cath Closed System has a Self-Cath pre-lubricated catheter enclosed in a single-use, collection bag. The catheter is available with and without insertion supplies. This catheter is not made with natural rubber latex.

#### Self-Cath<sup>®</sup> Closed Systems with Insertion Supplies

Self-Cath pre-lubricated intermittent catheter, gloves, 1100 mL seamless collection bag, povidone-iodine swabsticks, BZK wipe (towelette), underpad and instructions for use.

#### Straight Tip

Code	Size	Units	Description	HCPCS*
1008	• 8 FR	50	16" length	A4353
1010	• 10 FR	50	16" length	A4353
1012	O 12 FR	50	16" length	A4353
1014	• 14 FR	50	16" length	A4353
1016	🛑 16 FR	50	16" length	A4353
Female	Straight Tip			
Code	Size	Units	Description	HCPCS*
2214		ГО	C" longth	A 42F2

3214	14 FR	50	6" length	A4353
Soft Stro	light Tip			
3114	• 14 FR	50	16" length	A4353
Olive Co	udé Tip with Gu	ide Stripe®		
C3814	• 14 FR	50	16" length	A4353
C3816	• 16 FR	50	16" length	A4353
Tapered	Coudé Tip with	Guide Strip	0e <sup>®</sup>	
3614	• 14 FR	50	16" length	A4353

#### Self-Cath<sup>®</sup> Closed System Single Units

Self-Cath pre-lubricated intermittent catheter and 1100 mL collection bag only. No additional insertion supplies.

#### Straight Tip

Code	Size	Units	Description	HCPCS*
1108	• 8 FR	50	16" length	A4353
C1110	• 10 FR	50	16" length	A4353
1112	O 12 FR	50	16" length	A4353
1114	• 14 FR	50	16" length	A4353
1116	🛑 16 FR	50	16" length	A4353
Female S	Straight Tip			
2214	• 14 FR	50	6" length	A4353
Soft Stra	ight Tip			
2114	• 14 FR	50	16" length	A4353
Olive Co	udé Tip with	Guide Stri	pe®	
2814	• 14 FR	50	16" length	A4353
2816	🗕 16 FR	50	16" length	A4353
Tapered	Coudé Tip wi	ith Guide S	Stripe®	
2614	• 14 FR	50	16" length	A4353

 WARNING: This product can expose you to Di(2-ethylhexyl)phthalate (DEHP), which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov



## Peristeen<sup>®</sup> Transanal Irrigation System A better path to bowel control

Transanal irrigation with Peristeen is a well-documented bowel emptying technique that is effective at preventing fecal incontinence and chronic constipation for patients with neurogenic bowel dysfunction.

To learn more visit www.coloplastcare.us/bowel.



# Transanal Irrigation

#### Peristeen<sup>®</sup> Transanal Irrigation

Peristeen is a FDA Class II cleared transanal irrigation system intended to instill water into the colon through a rectal catheter to promote evacuation of the contents of the lower colon. It is indicated for use by children (2 and up), adolescents, transitional adolescents and adults with neurogenic bowel dysfunction.

Peristeen is specially designed for portability and ease of use. It offers independence and a sense of security to the user. In clinical trials with patients with neurogenic bowel dysfunction, Peristeen system reduced fecal incontinence and constipation and improved quality of life.<sup>1,2</sup>

#### Peristeen® Transanal Irrigation

Code			
29129	<b>Ն</b> Րջ	System	1 control unit 1 bag (with lid) 2 leg straps 1 toiletries bag
29122	ŢĴ	Accessory unit	15 standard size rectal catheters 1 bag (no lid)
29127	ĴĴ	Accessory unit	15 small size rectal catheters 1 bag (no lid)

#### HCPCS code A4459

"Manual pump enema system, includes balloon, catheter and all accessories, reusable, any type" A4459 would include 29129 (1 each) and either 29122 or 29127. Available via prior authorization process for Medicaid and private insurance<sup>\*</sup>.

 Christensen P et al. A randomized, controlled trial of transanal irrigation versus conservative bowel management in spinal cord injured patients. Gastroenterology 2006; 131:738-747

 Mosiello G., Marshall D., Rolle U., Crétolle C., Santacruz B.G., Frischer J. and Benninga M. Consensus Review of Best Practice of Transanal Irrigation in Children. JPGN; Volume 64, Number 3, March 2017

WARNING: This product can expose you to Di-n-butyl phthalate (DBP), which is known to the State of California to cause birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov

This warning does not apply to catalog number 29129 but does apply to catalog numbers 29122 and 29127  $\,$ 

#### Male External Catheters (Condom Catheters)

Coloplast offers a complete line of male external catheters. A variety of styles and sheath lengths are designed to suit end-user needs. A complete incontinence system can be achieved by combining a Coloplast male external catheter and a Coloplast leg bag.

#### Conveen® Optima

A self-adhering, 100% silicone clear male external catheter packaged in a discreet, recyclable smart pack.

- Easy on the double strip grip ensures an easy and accurate application
- Triple action system with anti-kink bellow, sure grip ribs and push ring provide maximum security
- The PSX silicone is soft, flexible and breathable
- The GX adhesive is reliable and skin friendly
- Not made with natural rubber latex



2007 WINNER
-------------

Cha		a se al la	Len	-
STO	nac	Ira.	Len	arn
				g

Code	Size	Units	HCPCS*
22025	25 mm	30	A4349
22028	28 mm	30	A4349
22030	30 mm	30	A4349
22035	35 mm	30	A4349
22040	40 mm	30	A4349

#### Sport Length

	5		
22121	21 mm	30	A4349
22125	25 mm	30	A4349
22130	30 mm	30	A4349
22135	35 mm	30	A4349

#### Conveen<sup>®</sup> Security+ Self-Sealing Male External Catheter

A self-adhering male external catheter features a placement ring to assist with individual preference and an anti-kink design to ensure security. Easy to apply with the applicator and roll-out strip. **Not made with natural rubber latex.** 



Code	Size	Units	HCPCS*
5221	21 mm	35	A4349
5225	25 mm	35	A4349
5230	30 mm	35	A4349
5235	35 mm	35	A4349
5240	40 mm	35	A4349
5225H	25 mm	100	A4349
5230H	30 mm	100	A4349
5235H	35 mm	100	A4349
5240H	40 mm	100	A4349

#### Freedom Clear Advantage® with Aloe

The Clear Advantage male external catheter is a self-adhering, silicone catheter. It features aloe in the adhesive and an anti-kink reservoir. Not made with natural rubber latex.

Size	Units
23 mm	100
28 mm	100

6200	28 mm	100	A4349
6300	31 mm	100	A4349
6400	35 mm	100	A4349
6500	40 mm	100	A4349
6130	23 mm	30	A4349
6230	28 mm	30	A4349
6330	31 mm	30	A4349
6430	35 mm	30	A4349
6530	40 mm	30	A4349

#### Freedom Clear®

Code

6100

The Freedom Clear male external catheter offers a silicone sheath with an anti-kink reservoir. Not made with natural rubber latex.

Code	Size	Units	HCPCS*
5100	23 mm	100	A4349
5200	28 mm	100	A4349
5300	31 mm	100	A4349
5400	35 mm	100	A4349
5500	40 mm	100	A4349

#### Freedom Clear<sup>®</sup> LS

The Freedom Clear LS (Long Seal) male external catheter is designed with a longer-length adhesive seal on a silicone sheath featuring an anti-kink reservoir. Not made with natural rubber latex.

Code	Size	Units	HCPCS*
5190	23 mm	100	A4349
5290	28 mm	100	A4349
5390	31 mm	100	A4349
5490	35 mm	100	A4349
5590	40 mm	100	A4349

#### Freedom Clear<sup>®</sup> SS


The Freedom Clear SS (Sport Sheath) male external catheter offers a shortened silicone sheath for less roll-up of the shaft of the catheter. This product features an anti-kink reservoir. Not made with natural rubber latex.

Code	Size	Units	HCPCS*
5110	23 mm	100	A4349
5210	28 mm	100	A4349
5310	31 mm	100	A4349
5410	35 mm	100	A4349

HCPCS\*

A4349





#### Freedom Cath®

Freedom Cath male external catheter is a selfadhering, one-piece latex catheter. Combines a latex sheath with a wide, watertight adhesive seal featuring an anti-kink reservoir.



#### Standard Length

Code	Size	Units	HCPCS*
8000	23 mm	100	A4349
8200	28 mm	100	A4349
8205	31 mm	100	A4349
8400	35 mm	100	A4349
8030	23 mm	30	A4349
8230	28 mm	30	A4349
8235	31 mm	30	A4349
8430	35 mm	30	A4349

#### Freedom Active Cath®

Freedom Active Cath male external catheter is a self-adhering, extended wear, one-piece latex catheter. Combines a shortened latex sheath with a wide, watertight adhesive seal. Features an anti-kink reservoir.



#### Sport Length

Code	Size	Units	HCPCS*
8100	23 mm	100	A4349
8300	28 mm	100	A4349
8305	31 mm	100	A4349
8500	35 mm	100	A4349

#### Manhood<sup>®</sup> Absorbent Pouch

Manhood Absorbent Pouch is designed for light to moderate incontinence. Using a super-absorbent polymer, the pouch can hold 250 mL of liquid.



Code	Size	Units	Description
27-4200-В	250 mL	30	One size fits most

Features a discreet design and color scheme to be less visible under clothing. Two types are available, clamp or lever tap. Anti-kink tubing, soft cloth backing, straps and extended outlet tubing for urine emptying process.

#### Conveen<sup>®</sup> Security+ Leg Bag

Secure and reliable urine collection bag with anti-kink tubing designed to prevent accidental leakage. Not made with natural rubber latex and is DEHP and phthalate free.

Code	Size	Units	Description	HCPCS*
Non-Ster	ile			
21026	17 oz. (500 mL)	10	20" adjustable tubing, clamp closure, flocked backing, 1" leg bag straps	A4358
21027	34 oz. (1000 mL)	10	20" adjustable tubing, clamp closure, flocked backing, 1" leg bag straps	A4358
Sterile				
21029	17 oz. (500 mL)	10	3" non-adjustable tubing, lever tap, 1" leg bag straps	A4358
21034	17 oz. (500 mL)	10	20" non-adjustable tubing, lever tap, 1" leg bag straps	A4358
21053	34 oz. (1000 mL)	10	3" non-adjustable tubing, lever tap, 1" leg bag straps	A4358
21054	34 oz. (1000 mL)	10	20" non-adjustable tubing, lever tap, 1" leg bag straps	A4358

#### Conveen® Contour Leg Bag\*

Deluxe leg bag with pre-attached, adjustable, anti-kink tubing with fabric straps. Secure, user-friendly closure, soft cloth backing, and chambered panels for discretion and low profile. **Not made with natural rubber latex.** 

#### Non-Sterile

Code	Size	Units	Description	HCPCS*
5170	21 oz. (600 mL)	10	Medium with 18" adjustable tubing, clamp closure, 1" leg bag straps	A4358
5174 Sterile	28 oz. (800 mL)	10	Large with 18" adjustable tubing, clamp closure, 1" leg bag straps	A4358
5171	21 oz. (600 mL)	10	Medium with 18" non- adjustable tubing, clamp closure, 1" leg bag straps	A4358

\* (A WARNING: This product can expose you to Di(2-ethylhexyl)phthalate (DEHP), which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www. P65Warnings.ca.gov

#### Conveen<sup>®</sup> Standard Drainage Bag

Large drainage bag with soft cloth backing. Secure, user-friendly closure. **Not made with natural rubber latex and is DEHP and phthalate free.** 



Code	Size	Units	Description	HCPCS*
5062	51 oz. (1500 mL)	10	Non-sterile, extra large with 36" tubing, clamp closure, straps not included	A4357

#### Conveen<sup>®</sup> Security+ Drainage Bag (formerly Moveen<sup>®</sup>)

Large volume (2 liters) bedside drainage bag. Nonsterile package includes 30 bags and tubes and 3 hangers. Sterile package includes 10 bags and 2 hangers. Not made with natural rubber latex and is DEHP and phthalate free.

#### Non-Sterile

Code	Size	Units	Description	HCPCS*
21346	2 liters	30	52" adjustable tubing, slide closure	A4357
Sterile				
21356	2 liters	10	Non-adjustable 52" tubing, slide closure	A4357

#### Conveen<sup>®</sup> Active Leg Bag

The Conveen Active Leg bag provides a reliable and discreet solution for patients. It wraps around the thigh with the bag facing inward for minimal visibility, even when wearing shorts. Developed to fit a normal active lifestyle and easy to use. Not made with natural rubber latex.

Code	Size	Units	Description	HCPCS*
25502	8.5 oz. (250 mL)	30	Non-sterile, 3" adjustable tubing. Wrap 13" long and 4.5" wide	A4358

Neugraphic Editor Ledg Dugg Only the Coloplast Alpine® Leg Bag has the patented KRAYLEX® odor barrier, offset port options and integrated loops providing total confidence and comfort.						
				-		
Alpine® Reusable Leg Bag						
Code	Capacity	Port	Size	Units	HCPCS*	
Sport/Sm	all					
68006 R or L	15 oz.	Offset R or L	11 x 4" (27.9 x 10.2 cm)	1	A5112	
Short Wic	le					
68004	26 oz.	Straight	13 x 5" (33 x 12.7 cm)	1	A5112	
Large						
68001	32 oz.	Straight	17 x 5" (43.2 x 12.7 cm)	1	A5112	
Extra Lar	ge					
68007	44 oz.	Straight	20 x 5" (50.8 x 12.7 cm)	1	A5112	

**Reusable Latex Leg Bags** 

Cysto-Care® Folysil® 100% Silicone Indwelling Catheters.

#### Folysil<sup>®</sup> 2-way Pediatric

100,950 2 0					
Code	Balloon Size	Size	Description	Units	KCPCS*
Straight Tip					
AA 6106	1.5 cc	6 FR	12" length	5	A4344
AA 6108	3 cc	8 FR	12" length	5	A4344
AA 6110	3 cc	10 FR	12" length	5	A4344
Coudé Tip					
AA 6308	3 cc	8 FR	12" length	5	A4340
AA 6310	3 cc	10 FR	12" length	5	A4340
Open-tip					
AA 6408	3 cc	8 FR	10" length	5	A4344
AA 6410	3 cc	10 FR	12" length	5	A4344

#### Folysil<sup>®</sup> 2-way Straight Tip

Code	Balloon Size	Size	Description	Units	HCPCS*
AA 6112	10 cc	12 FR	16" length	5	A4344
AA 6114	10 cc	14 FR	16" length	5	A4344
AA 6116	15 cc	16 FR	16" length	5	A4344
AA 6118	15 cc	18 FR	16" length	5	A4344
AA 6120	15 cc	20 FR	16" length	5	A4344
AA 6122	15 cc	22 FR	16" length	5	A4344
AA 6124	15 cc	24 FR	16" length	5	A4344
AA 6C18	30 cc	18 FR	16" length	5	A4344
AA 6C20	30 cc	20 FR	16" length	5	A4344
AA 6C22	30 cc	22 FR	16" length	5	A4344
AA 6C24	30 cc	24 FR	16" length	5	A4344



#### Folysil<sup>®</sup> 2-way Coudé Tip

	Balloon				
Code	Size	Size	Description	Units	HCPCS*
AA 6312	10 cc	12 FR	16" length	5	A4340
AA 6314	10 cc	14 FR	16" length	5	A4340
AA 6316	15 cc	16 FR	16" length	5	A4340
AA 6318	15 cc	18 FR	16" length	5	A4340
AA 6320	15 cc	20 FR	16" length	5	A4340
AA 6322	15 cc	22 FR	16" length	5	A4340
AA 6324	15 cc	24 FR	16" length	5	A4340

# Folysil<sup>®</sup> 2-way Open Tip

	Balloon				
Code	Size	Size	Description	Units	HCPCS*
AA 6412	10 cc	12 FR	16" length	5	A4344
AA 6414	10 cc	14 FR	16" length	5	A4344
AA 6416	15 cc	16 FR	16" length	5	A4344
AA 6418	15 cc	18 FR	16" length	5	A4344
AA 6422	15 cc	22 FR	16" length	5	A4344
AA 6424	15 cc	24 FR	16" length	5	A4344

#### Intermittent Catheter Accessories

This tubing attaches to the end of an intermittent catheter. It is designed to allow the individual to catheterize when unable to reach the toilet from a slight distance.

#### Self-Cath® Intermittent Catheter Extension Tube

Code	Size	Units	HCPCS*
475	24"	12	A4331

WARNING: This product can expose you to Di(2-ethylhexyl)phthalate (DEHP), which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov

#### Male External Catheter/Bag Accessories

#### Bag Hanger

Reusable, sturdy, plastic drainage bag hanger secures tubing and ensures bag stability when attached to a bed or wheelchair. Not made with natural rubber latex



Code	Size	Units	HCPCS*
5070	One size	1	Medicare A9900
			Other A9999

#### Conveen<sup>®</sup> Fabric Leg Bag Straps

Soft, reusabl	e, elastic fabric	, button straps with	
hook and loc	p closure. Each	n pair includes a	
1.5" width ar	id 1" width stra	p.	-
Not made w	th natural rub	ber latex	e
Code	Sizo	Units	

Code	Size	Units	HCPCS*
50501	23" length	10 pair	A5114

			Elener .
_			
Alpine®	Fresh Cleaner		
For use w	ith Alpine reusable leg	bags.	
Code	Size/Description	Units	HCPCS*
90000	16 fl. oz. bottle	12	A5131

<sup>\*</sup>For important information regarding the use of HCPCS codes listed and code definitions see page 192

# Product Guide

Wound Care | Effective June 2020



#### **Brief statements**

Biatain® Ag Foam Adhesive and Non-Adhesive Dressings are indicated for use in patients with moderately to highly exuding leg ulcers and pressure sores, second degree burns, donor sites, post-operative wounds and skin abrasions. Biatain® Ag Foam Non-Adhesive Dressings are also indicated for diabetic foot ulcers. Patients being treated for wounds should follow the advice of, and direct questions about use of the product to, their medical professional. Before using the device, carefully read the product labels and information accompanying the device including the instructions for use which contain additional safety information. Coloplast's wound dressings are for single-use only; discard it after use. The use of Biatain Ag dressings during pregnancy, lactation and on children has not been demonstrated. Biatain Ag foam dressings should not be used on patients with a known sensitivity to silver. If you experience symptoms of an allergic reaction, contact your healthcare professional. The risk information provided here is not comprehensive. To learn more, talk to your healthcare provider. For further information, call Coloplast Corp. at 1-800-533-0464 and/ or consult the company website at www.coloplast.us.

Biatain<sup>®</sup> Alginate Ag Dressings are indicated for wounds with moderately to heavily exuding partial to full thickness wounds such as post-operative wounds, trauma wounds, leg ulcers, pressure ulcers, and diabetic ulcers. Patients being treated for wounds should follow the advice of, and direct questions about use of the product to, their medical professional. Before using the device, carefully read the product labels and information accompanying the device including the instructions for use which contain additional safety information. Coloplast's wound dressings are for single-use only; discard it after use. The use of Biatain Ag dressings during pregnancy, lactation and on children has not been demonstrated. Biatain Ag foam dressings should not be used on patients with a known sensitivity to silver. If you experience symptoms of an allergic reaction, contact your healthcare professional. The risk information provided here is not comprehensive. To learn more, talk to your healthcare provider. For further information, call Coloplast.cus.

Triad<sup>®</sup> hydrophilic wound dressing is indicated for the local management of pressure and venous stasis ulcers, dermal lesions/injuries, superficial wounds, scrapes, first-and second-degree burns, partial-and full-thickness wounds. Suitable for use on neonates. Triad is contraindicated for third degree burns or infected wounds. See device label for complete information including warnings and precautions. For further information call Coloplast Corp. at 1-800-533-0464 and/or consult the company website at www.trytriad.com.

# **Dressing Selection Chart**

	Minimal Drainage	rainage	Moderate Drainage	Drainage	Heavy Drainage	rainage
	Comfeel <sup>®</sup> Plus Transparent	Triad®	Comfeel <sup>®</sup> Plus	Biatain <sup>®</sup> Silicone Lite	Biatain <sup>®</sup> Silicone	Biatain <sup>®</sup> Super
		#1 2				0
Deep Wound (Less than 2 cm deep)	Triad+ Comfeel® Plus Transparent	Triad+ Comfeel® Plus	Biatain <sup>®</sup> Silicone	Biatain° Silicone Lite	Biatain <sup>®</sup> Silicone	Biatain <sup>®</sup> Non-Adhesive
Deep Wound (Greater than 2 cm	Triad+ Gauze + Comfeel® Plus	Triad+ Gauze + Comfeel® Plus	Biatain® Alginate + Biatain®	Biatain° Alginate + Biatain° Silicone	Biatain® Alginate + Biatain® Silicone	Biatain° Alginate + Biatain° Super
aeep or sinus tracts/ undermining)	Iransparent		Silcone Lite	0	0	0
						0



# Discover the difference of Biatain<sup>®</sup> dressings with 3DFit Technology™

🚔 Coloplast

Biatain foam dressings feature 3DFit Technology for optimal moist wound healing



#### Conforms to the wound bed

Upon contact with exudate, the unique foam conforms to the wound bed and reduces exudate pooling



Absorbs vertically Microcapillaries go into action to absorb the exudate vertically and create optimal healing conditions



Retains exudate Exudate is lockedaway and retained even under compression, reducing risk of leakage and maceration

#### Visit www.biatainsiliconeus.com to learn more

#### Biatain<sup>®</sup> Silicone Foam Dressing

- Features 3DFit Technology<sup>™</sup> for optimal moist wound healing
- Gentle and secure fit
- Soft and flexible
- Aseptic 3-piece non-touch opening



Code	Size	Units	HCPCS*
33434	3 x 3" (7.5 x 7.5 cm)	10	A6212
33435	4 × 4" (10 × 10 cm)	10	A6212
33436	5 x 5" (12.5 x 12.5 cm)	10	A6212
33437	6 x 6" (15 x 15 cm)	5	A6213
33438	7 x 7" (17.5 x 17.5 cm)	5	A6213
33400	4 x 8" (10 x 20 cm)	5	A6245
33401	4 x 12" (10 x 30 cm)	5	A6246
33406	Heel 7 x 7" (18 x 18 cm)	5	A6245
39000**	Sacral 9.8 x 9.8" (25 x 25 cm)	5	A6213
39001**	Sacral 6 x 7.5" (15 x 19 cm)	5	A6210
39010**	Multishape 5.5 x 7.6" (14 x 19.5 cm)	5	A6212

#### 00-000 Exudate Level

\*\*indicated for pressure injury prevention

#### Biatain<sup>®</sup> Silicone Lite Foam Dressing

- Features 3DFit Technology<sup>™</sup> for optimal moist wound healing
- Gentle and secure fit
- Thin, soft, and ultra flexibleAseptic 3-piece non-touch opening



Code	Size	Units	HCPCS*
33452	2 x 2" (5 x 5 cm)	5	A6413
33453	2 x 5" (5 x 12.5 cm)	5	A6212
33444	3 x 3" (7.5 x 7.5 cm)	10	A6212
33445	4 x 4" (10 x 10 cm)	10	A6212
33446	5 x 5" (12.5 x 12.5 cm)	10	A6212

0-00 Exudate Level

#### Biatain® Non-Adhesive Foam Dressing

- Features 3DFit Technology™ for optimal moist wound healing
- Soft and flexible dressing with beveled edges
- · Suitable for compression wraps and fragile skin

Code	Size	Units	HCPCS*
6105	2 x 2.75" (5 x 7 cm)	10	A6209
3410	4 x 4" (10 x 10 cm)	10	A6209
3413	6 x 6" (15 x 15 cm)	5	A6210
3416	8 x 8" (20 x 20 cm)	5	A6211

00-000 Exudate Level



#### Biatain<sup>®</sup> Soft-Hold Foam Dressing

- Features 3DFit Technology™ for optimal moist wound healing
- Skin-friendly contact adherent provides a "third hand" during dressing changes
- Helps minimizes adhesive trauma, easy to apply and remove



Code	Size	Units	HCPCS*
3473	2 x 2.75" (5 x 7 cm)	5	A6209
3470	4 x 4" (10 x 10 cm)	5	A6209
3475	6 x 6" (15 x 15 cm)	5	A6210

00-000 Exudate Level

#### Biatain® Adhesive Foam Dressing

- Features 3DFit Technology<sup>™</sup> for optimal moist wound healing
- Absorptive and retentive foam with a hydrocolloid-based adhesive border
- Designed for wounds with healthy periwound skin that need extra dressing adhesion

Code	Size	Units	HCPCS*
3430	4 x 4" (10 x 10 cm) circle pad	10	A6212
3420	5 x 5" (12.5 x 12.5 cm)	10	A6212
3423	7 x 7" (18 x 18 cm)	5	A6213
3485	Sacral 9 x 9" (23 x 23 cm)	5	A6213
3488	Heel 7.5 x 8" (19 x 20 cm)	5	A6212

00-000 Exudate Level



Biatain<sup>®</sup> Ag Foam Dressings combine 3DFit Technology<sup>™</sup> with controlled ionic silver-release. The dressings demonstrate in-vitro antibacterial activity in certain strains known to be detrimental to wound healing.<sup>1</sup> Clinical evidence shows a 45-56% ulcer size reduction within 4 weeks.<sup>2-5</sup> Not made with natural rubber latex; sterile.

## Biatain® Ag Non-Adhesive Foam Dressing

- Features 3DFit Technology<sup>™</sup> for optimal moist wound healing
- Combines effective antimicrobial technology with exudate management
- · Soft and flexible dressing with beveled edges
- · Suitable for compression wraps and fragile skin

Code	Size	Units	HCPCS*
9622	4 x 4" (10 x 10 cm)	5	See PDAC
9625	6 x 6" (15 x 15 cm)	5	See PDAC

00-000 Exudate Level

## Biatain<sup>®</sup> Ag Adhesive Foam Dressing

- Features 3DFit Technology<sup>™</sup> for optimal moist wound healing
- Combines effective antimicrobial technology
   with exudate management
- Absorptive and retentive foam with a hydrocolloid-based
   adhesive border
- Designed for wounds with healthy periwound skin that need extra dressing adhesion

Code	Size	Units	HCPCS*
9632	5 x 5" (12.5 x 12.5 cm)	5	See PDAC
9635	7 x 7" (18 x 18 cm)	5	See PDAC
9641	Sacral 9 x 9" (23 x 23 cm)	5	See PDAC
9643	Heel 7.5 x 8" (19 x 20 cm)	5	See PDAC


## 00-000 Exudate Level

- 1. Data on file at Coloplast Corp.
- Münter K-C et al. Effect of a sustained silver-releasing dressing on ulcers with delayed healing: the CONTOP study. Journal of Wound Care 2006, 15(5), 199-206.
- Jørgensen B et al. The silver-releasing foam dressing, Contreet Foam, promotes faster healing of critically colonised venous leg ulcers: a randomised, controlled trial. International Wound Journal 2005, 2(1), 64-73.
- 4. Rayman et al. Sustained silver-releasing dressing in the treatment of diabetic foot ulcers. British Journal of Nursing 2005, 14(2), 109-114.
- 5. Karlsmark T. et al. Clinical performance of a new silver dressing, Contreet Foam, for chronic exuding venous leg ulcers. Journal of Wound Care 2003, 12(9), 351-354.



## **Biatain®** Alginate

Biatain<sup>®</sup> Alginate is a highly absorbent wound dressing composed of alginate and hydrocolloid (sodium carboxymethylcellulose). On contact with wound exudate it converts to a soft cohesive gel which provides an optimal moist wound healing environment. Not made with natural rubber latex; sterile.



## **Biatain® Alginate Dressing**

- Fast-gelling alginate
- One-piece removal
- · Recommended for cavity filling

Code	Size	Units	HCPCS*
3705	2 x 2" (5 x 5 cm)	30	A6196
3710	4 x 4" (10 x 10 cm)	10	A6196
3715	6 x 6" (15 x 15 cm)	10	A6197

00-000 Exudate Level

## **Biatain®** Alginate Rope

- Fast-gelling alginate
- One-piece removal
- Recommended for cavity filling



Code	Size	Units	HCPCS*
3740	1 x 17.5" rope (3 x 44 cm)	6	A6199
00-00	Exudate Level		



## Biatain<sup>®</sup> Alginate Ag Dressing

- Fast-gelling alginate containing an ionic silver complex with broad antimicrobial profile
- One-piece removal
- Recommended for cavity filling

Code	Size	Units	HCPCS*
3755	2 x 2" (5 x 5 cm)	30	A6196
3760	4 x 4" (10 x 10 cm)	10	A6196
3765	6 x 6" (15 x 15 cm)	10	A6197

00-000 Exudate Level

## Biatain® Alginate Ag Rope

- Fast-gelling alginate rope containing an ionic silver complex with broad antimicrobial profile
- One-piece removal
- Recommended for cavity filling

Code	Size	Units	HCPCS*
3780	1 x 17.5" rope (3 x 44 cm)	10	A6199

♦♦−♦♦♦ Exudate Level

The Comfeel® Plus portfolio offers a range of dressings to manage and protect non- to moderately-exuding wounds.

- Exudate management with fewer dressing changes and less disruption to wound bed
- Specialty shapes and sizes for difficult-to-dress areas
- Color guide indicates best time to change dressing
- Non-touch application system for ease of use
- Not made with natural rubber latex; sterile

#### Comfeel® Plus Transparent Dressing (Thin)

- Comfortable, flexible hydrocolloid adaptable to wound and body contours
- wound and body contoursClear dressing and mapping grid allow for
- visual inspection without disturbing wound
- Does not contain alginate

Code	Size	Units	HCPCS*
33530	2 x 2.75" (5 x 7 cm)	10	A6234
33533	4 x 4" (10 x 10 cm)	10	A6234
33536	3.5 x 5.5" (9 x 14 cm)	10	A6235
33539	6 x 6" (15 x 15 cm)	5	A6235

Exudate Level

## Comfeel<sup>®</sup> Plus Dressing

· Visual wound monitoring with mapping grid

Contains alginate for enhanced absorption
 capacity



Code	Size	Units	HCPCS*
33146	1.5 x 2.5" (4 x 6 cm)	30	A6234
33110	4 x 4" (10 x 10 cm)	10	A6234
33115	6 x 6" (15 x 15 cm)	5	A6235
33120	8 x 8" (20 x 20 cm)	5	A6236

Comfeel® Plus Contour Dressing • Butterfly-shaped hydrocolloid dressing for better conformity and adhesion • Specifically designed for difficult-to-dress areas • Contains alginate for enhanced absorption					
Code	Size	Units	HCPCS*		
33280	24 sq. in. (6 x 8 cm)	5	A6235		
33283	42 sq. in. (9 x 11 cm)	5	A6235		
	idate Level				
Comfeel <sup>®</sup> Plus Sacral Dressing • Triangular-shaped hydrocolloid dressing designed for the sacral area • Contains alginate for enhanced absorption					
Code	Size	Units	HCPCS*		
33285	7 x 7" (17 x 17 cm)	5	A6235		

#### Triad<sup>®</sup> Hydrophilic Wound Dressing

- Zinc oxide-based hydrophilic dressing absorbs moderate levels of wound exudate
- Facilitates autolytic debridement
- Ideal alternative for difficult-to-dress areas



Not made with natural rubber latex; sterile

Code	Size	Units	HCPCS*
1964	2.5 oz. tube (71 g)	12	A6240
1967	6 oz. tube (170 g)	12	A6240

## $\wedge - \wedge \wedge$ Exudate Level

#### Sea-Clens<sup>®</sup> Wound Cleanser

- · Saline-based solution for cleansing and irrigating acute and chronic wounds
- No-rinse, isotonic, pH-balanced formula
- · Easy-to-use spray bottle with SPRAY and STREAM setting with irrigation pressure

Code	Size	Units	HCPCS*
1063	6 fl. oz. (178 mL)	12	A6260
1061	12 fl. oz. (355 mL)	12	A6260

0-000 Exudate Level

## Woun'Dres<sup>®</sup> Collagen Hydrogel


- · Rehydrates dry wounds and eschar
- · Promotes moist wound environment and autolytic debridement
- Not made with natural rubber latex

Code	Size	Units	HCPCS*
1166	1 oz. tube (28 g)	36	A6248
7690	3 oz. tube (84 g)	12	A6248

Exudate Level

## Purilon<sup>®</sup> Gel

· Sterile hydrogel with hydrating and absorbing properties to promote autolytic debridement



- Contains alginate
- · Not made with natural rubber latex; sterile

Code	Size	Units	HCPCS*
3906	0.28 oz. (8 g)	10	A6248
3900	0.5 oz. (15 g)	10	A6248
3903	0.88 oz. (25 g)	10	A6248

∆-∆∆ Exudate Level



## Biatain<sup>®</sup> Contact

Biatain Contact is a one-sided silicone wound contact layer designed to protect wounds from secondary dressing changes.

- · Easy application and removal
- Protects against in-growth
- Bigger pores for enhanced exudate transfer and reduced risk of pooling
- Improved transparency to facilitate wound assessment without removing the dressing

Code	Size	Units	HCPCS*
33560	2 x 3" (5 x 7.5 cm)	10	A6206
33561	3 x 4" (7.5 x 10 cm)	10	A6206
33562	4 x 7" (10 x 18 cm)	10	A6207
33563	6 x 10" (15 x 25 cm)	5	A6208
33564	7 x 12" (18 x 30 cm)	5	A6208

0-000 Exudate Level



## Biatain<sup>®</sup> Super, Adhesive

- Fast absorption, distribution and retention of exudate
- Suitable for highly exuding wounds surrounded by healthy skin



Code	Size	Units	HCPCS*
46100	4 x 4" (10 x 10 cm)	10	A6203
46120	5 x 5" (12.5 x 12.5 cm)	10	A6203
46250	5 x 8" (12.5 x 20 cm)	10	A6204
46150	6 x 6" (15 x 15 cm)	10	A6203
46200	8 x 8" (20 x 20 cm)	10	A6204

#### Biatain® Super, Non-Adhesive

- Fast absorption, distribution and retention of exudate
- Suitable for highly exuding wounds surrounded by fragile skin



463004 x 4" (10 x 10 cm)10A6251463205 x 5" (12.5 x 12.5 cm)10A6251464505 x 8" (12.5 x 20 cm)10A6252463506 x 6" (15 x 15 cm)10A6251463908 x 8" (20 x 20 cm)10A6252	Code	Size	Units	HCPCS*
46450         5 x 8" (12.5 x 20 cm)         10         A6252           46350         6 x 6" (15 x 15 cm)         10         A6251	46300	4 x 4" (10 x 10 cm)	10	A6251
<b>46350</b> 6 x 6" (15 x 15 cm) 10 A6251	46320	5 x 5" (12.5 x 12.5 cm)	10	A6251
	46450	5 x 8" (12.5 x 20 cm)	10	A6252
<b>46390</b> 8 x 8" (20 x 20 cm) 10 A6252	46350	6 x 6" (15 x 15 cm)	10	A6251
	46390	8 x 8" (20 x 20 cm)	10	A6252

000-0000 Exudate Level

# Product Guide

Skin Care | Effective June 2020



#### **Brief statement**

InterDry<sup>®</sup> is a skin protectant for the treatment of skin folds. The silver complex protects against microbial colonization in the fabric. This single-use fabric provides moisture translocation to keep skin dry while the antimicrobial in the fabric reduces odor. The safe use of InterDry fabric during pregnancy, lactation and on children has not been demonstrated. Although many patients benefit from the use of this device, results may vary. See the device labeling for detailed information including directions for use, indications, description, precautions, warnings, and other information. For further questions, call Coloplast Corp. at 1-800-533-0464 and/or consult the company website at www.coloplast.us.



# We believe.

At Coloplast, healthy skin is our purpose, our mission, our cause, and our promise. And we deliver on our promise every day by focusing on three key areas: breadth, simplicity, and partnership.

## Color-coded labels simplify product selection



## Brand name

A brand name you've come to trust

2 Category Defines the need

3 Description Defines the product



## 4 Where to use Defines the

product use



## 5 Symbol

Graphic representation of skin care needs

## 6 Cap color

Further differentiates product function





## Easy product identification:

Symbols and colors represent product function

CLEANSERS

MOISTURIZERS







SKIN PROTECTANTS



SKIN PROTECTANTS WITH ANTIFUNGAL



Barrier with Antifungal

ANTIFUNGALS



SKIN FOLD MANAGEMENT







Rinse



No-Rinse



## Bedside-Care® Foam

No-rinse foam body wash, shampoo and incontinence cleanser

- Cleanses, moisturizes and conditions the skin
- pH-balanced to help maintain the natural acid mantle of the skin
- Alcohol-free
- CHG compatible

Product	Code	Size	Units
Scented	67143	4.1 fl. oz. (120mL)	36
	67145	8.1 fl. oz. (240mL)	12
Unscented	67147	4.1 fl. oz. (120mL)	36
	67146	8.1 fl. oz. (240mL)	12

## **Bedside-Care® Spray**

No-rinse body wash, shampoo and incontinence cleanser

- Tear-free, hypoallergenic
- · Cleanses and conditions the skin
- Alcohol-free
- pH-balanced to help maintain the natural acid mantle of the skin
- CHG compatible

Product	Code	Size	Units
Scented	61768	4.1 fl. oz. (120mL)	36
	61762	8.1 fl. oz. (240mL)	12
Unscented	61760	4.1 fl. oz. (120mL)	36
	61761	8.1 fl. oz. (240mL)	12

## Bedside-Care® Sensitive Skin Foam

No-rinse foam body wash, shampoo and incontinence cleanser for all ages including neonates<sup>1</sup>

- Tear-free, hypoallergenic
- · Cleanses and conditions the skin
- Alcohol-free
- pH-balanced to help maintain the natural acid mantle of the skin
- CHG compatible

Code	Size	Units
7300	4 fl. oz. (118 mL)	36
7301	8 fl. oz. (237 mL)	12

## Bedside-Care® Sensitive Skin Spray

No-rinse body wash, shampoo and incontinence cleanser for all ages including neonates<sup>1</sup>

- Tear-free, hypoallergenic
- · Cleanses and conditions the skin
- Alcohol-free
- pH-balanced to help maintain the natural acid mantle of the skin
- · CHG compatible

Code	Size	Units
1300	4 fl. oz. (118 mL)	36





## EasiCleanse®Bath



Soft, sudsing disposable washcloths for a pleasant patient bathing experience

#### EasiCleanse® Bath

No-rinse body wash, shampoo, meatal, perineal and incontinence cleanser in a self-sudsing disposable washcloth

- · 3 ingredients with no preservatives
- No filmy residue left on skin
- No basin needed, to help avoid the risk of basin-associated contamination
- Suitable for all ages including neonates<sup>1</sup>

Code	Size	Units
7055	30 washcloths/pk	30 pks/cs
7056	5 washcloths/pk	100 pks/cs
7058	8 washcloths/pk	80 pks/cs





## Now available in 8pk for total body bathing Visit www.easicleansebath.com to learn more

1. Association of Women's Health, Obstetric and Neonatal Nurses. Neonatal Skin Care, Third Edition. Evidence-Based Clinical Practice Guideline, 2013.



## Gentle Rain® Extra Mild

Moisturizing body wash, shampoo and hand wash

• Contains mild ingredients and added emollients for use on fragile skin



- Dermatologist tested, hypoallergenic and dye-free
- pH-balanced to help maintain the natural acid mantle of the skin

Code	Size	Units
7231	5 mL single application packet	300
7230	2 fl. oz. (59 mL)	36
7229	4 fl. oz. (118 mL)	36
7235	8 fl. oz. (237 mL)	36
7233	21 fl. oz. (621 mL)†	12
7234	34 fl. oz. (1000 mL)†	12

- <sup>†</sup> See dispensers on page 59.
- 1. Association of Women's Health, Obstetric and Neonatal Nurses. Neonatal Skin Care, Third Edition. Evidence-Based Clinical Practice Guideline, 2013.

52

## Sween<sup>®</sup> 24

Once-a-day moisturizing body cream

- Once-a-day application
- · Contains 6% dimethicone
- Fragrance, petrolatum, lanolin, alcohol and preservative-free
- CHG compatible

. . .. • 5

<ul> <li>Suitable for all ag</li> </ul>	es including neonates <sup>1</sup>	
Code	Size	Units
7090	4 g single application packet	300
7091	2 oz. (57 g)	12
7092	5 oz. (142 g)	12
7095	9 oz. (256 g)	12


## Sween Cream®

## Moisturizing cream

- Non-occlusive formula
- · Will not interfere with tape adhesion
- · CHG compatible

Code	Size	Units
0283	2 g single application packet	300
7067	3 oz. (85 g)	12
7068	6.5 oz. (184 g)	12
7069	12 oz. (340 g) jar	12





- Non-occlusive formula
- CHG compatible



## Atrac-Tain® Cream

Moisturizing cream with urea & alpha hydroxy acid (AHA)

- 10% urea and 4% AHA
- · Moisturizing cream for extremely dry, scaly cracked skin and fissured heels
- · Fragrance, dye and preservative-free
- CHG compatible

0	0	100 100

Code	Size	Units
1843	2 g single application packet	300
1802	2 oz. (57 g)	12
1814	5 oz. (142 g)	12



#### Critic-Aid<sup>®</sup> Clear

Clear moisture barrier ointment

- Clear moisture barrier
- · Minimal ingredients and preservative-free
- Adheres to intact and eroded skin
- Easy to apply and remove
- CHG compatible
- Suitable for all ages including neonates1



	5 5	
Code	Size	Units
7565	4 g single application packet	300
7566	2.5 oz. (71 g)	12
7567	6 oz. (170 g)	12

## Critic-Aid<sup>®</sup> Skin Paste

Thick moisture barrier paste

Thick, zinc oxide-based skin protectant
Adheres to intact and moist, denuded skin
Contains only 4 ingredients to help minimize potential skin irritation
CHG compatible
Suitable for all ages including neonates<sup>1</sup>

Code	Size	Units
1944	2.5 oz. (71 g)	12
1947	6 oz. (170 g)	12

## Baza® Clear Ointment

Moisture barrier ointment

- Clear, petrolatum-based moisture barrier
- Enriched with skin conditioners
- Contains 96% petrolatum
- Clear formulation makes it easy to visually
   assess the skin
- CHG compatible

Code	Size	Units
7015	4 g single application packet	300
1005	1.75 oz. (50 g)	12
1006	5 oz. (142 g)	12

<sup>1.</sup> Association of Women's Health, Obstetric and Neonatal Nurses. Neonatal Skin Care, Third Edition. Evidence-Based Clinical Practice Guideline, 2013.

## Baza® Protect Cream

Moisture barrier cream

- Thick, zinc oxide-based skin protectant
- Enriched with skin conditioners
- CHG compatible



Code	Size	Units
1873	4 g single application packet	300
1877	2 oz. (57 g)	12
1880	5 oz. (142 g)	12

## Baza Cleanse & Protect<sup>®</sup> Lotion

All-in-one perineal lotion

- One step provides mild cleansing,
- a breathable barrier and moisturization
- Available with odor control
- pH-balanced to help maintain the natural acid mantle of the skin
- CHG compatible



Product	Code	Size	Units
Unscented	7710	4 g single application packet	300
	7712	8 fl. oz. (237 mL)	12
Scented	7725	8 fl. oz. (237 mL)	12

## Sween® Body Powder

All-body powder

- Absorbs moisture and perspiration
- Non-talcum based formula
- Helps reduce friction and prevent shearing



Code	Size	Units
0506	3 oz. (85 g)	36
0505	8 oz. (227 g)	36

## Critic-Aid® Clear Antifungal Ointment

Clear moisture barrier ointment with antifungal

- 2% miconazole nitrate
- Contains only 5 ingredients to help minimize potential skin irritation
- Adheres to intact and moist, denuded skin
  Clear formulation makes it easy to visually inspect the skin
- CHG compatible



Code	Size	Units
7570	4 g single application packet	300
7571	2 oz. (57 g)	12
7572	5 oz. (142 g)	12

## Baza® Antifungal Cream

Moisture barrier antifungal cream

- 2% miconazole nitrate
- Enriched with skin conditioners
- CHG compatible



Code	Size	Units
1622	4 g single application packet	300
1611	2 oz. (57 g)	12
1607	5 oz. (142 g)	12

## Micro-Guard<sup>®</sup> Powder

#### Antifungal powder

- Non-talcum based formula with 2% miconazole nitrate
- Effective treatment of topical fungal infections
- Helps reduce the effects of friction and shear

Code	Size	Units
1337	3 oz. (85 g)	12

## InterDry<sup>®</sup>

Moisture-wicking fabric with antimicrobial silver

- Manages all four causes of intertrigo; moisture, friction and bacterial or fungal organisms in skin folds<sup>1</sup> and other skin-to-skin and skin-to-device contact areas
- Manages the symptoms associated with intertrigo such as erythema, painful rash, denudement, satellite lesions and odor<sup>2</sup>
- Silver within the fabric provides effective antimicrobial action for up to 5 days

By targeting all four causes of Intertrigo, InterDry is clinically proven to provide complete, cost-effective relief within 5 days for as little as \$5<sup>3</sup>.

#### InterDry<sup>®</sup>

Code	Size	Units
7910	10 x 144"	10 rolls/cs
7912	10 x 36"	10 pks/box
7915	10 x 18"	40 pks/box



## Visit www.interdry.com to learn more

- 1. Vorbeck, E. Evaluation of a Skin Fold Management Textile with Antimicrobial Silver Complex in a Variety of Case Studies
- Freyberg J, et al. Moisture Management Challenges for the WOC Nurse presented at the WOCN Society 39th Annual Conference, June 2007
- 3. Kennedy-Evans KL, Viggiano B, Henn T, Smith D. Multi-site feasibility study using a new textile with silver for management of skin conditions located in skin folds. Poster presented at: 20th Annual Symposium on Advanced Wound Care; April 28 May 1, 2007; Tampa, FL and 39th WOCN<sup>®</sup> Society Annual Conference; June 9-13, 2007; Salt Lake City, UT.

## Isagel®

- No-rinse, antiseptic hand cleansing gel
- Formulated with emollients for extra mildness
- Contains 60% ethanol alcohol
- CHG and latex glove compatible



Code	Size	Units
1644	4 fl. oz. (118 mL)	36
1645	21 fl. oz. (621 mL)†	12
7041	27 fl. oz. (800 mL)†	12
7041	27 fl. oz. (800 mL) <sup>1</sup>	

 $^{\dagger}$  See dispensers on page 59.

## Hex-On®

Product

Fresh Linen

Fragrance Free

Odor antagonist

- Helps control airborne odors
- Non-aerosol spray



CCessori

Size	Units
2 fl. oz. (59 mL)	12
2 fl. oz. (59 mL)	12

## Bedside-Care<sup>®</sup> Gentle Touch Cloths

Disposable washcloths for use with all-body bathing and perineal cleansing products.

Code

7583

67589

- Rayon, for softness and absorbency
- · Polyester, for strength and tear resistance



Code	Size	Units
7057	25 washcloths/pack	24

## Coloplast<sup>®</sup> Wall Dispenser

- Holds 27 fl. oz. and 34 fl. oz. (800 mL and 1000 mL) refill bags
- · Completely closed to prevent possible contamination
- · Metered 1 mL pump dispenses appropriate amount
- No-drip nozzle leaves facility clean and safe
- · Easy-to-operate push pad meets Americans with Disabilities Act guidelines
- Keyless opening for simple cartridge replacement
- · Installs quickly with choice of screws or tape

Code	Size	Units
7251	Holds 27 fl. oz. and 34 fl. oz. (800 mL and 1000 mL)	12
7262	Drip tray	60

## Embrace<sup>™</sup> Bracket


- · Holds 21 fl. oz. (621 mL) Coloplast containers securely in place
- · Improved tamper and theft-resistant design
- · Optional locking collar conceals cap
- Installs quickly with choice of screws or tape

Code	Size	Units
0089	Holds 21 fl. oz. (621 mL)	12



59

## **Product Guide** Ostomy Care | Effective June 2020



## SenSura®Mio



## Through BodyFit Technology<sup>®</sup>, we offer a secure fit to real bodies to help prevent leakage and skin issues

Leakage and skin issues are closely related challenges. The right fit helps prevent both.



Skin issues

## 76%

of people with a stoma have experienced leakage within the last six months.<sup>1</sup>

## Leakage

is often caused by a poor fit to the body.<sup>2</sup>

## Skin issues

There is a significant correlation between leakage frequency and the severity of the peristomal skin condition.<sup>2</sup>

## 61%

of users attributed their most recent skin issue to leakage under the barrier.<sup>3</sup>

Ostomy

Results from the 2014 Ostomy Life Study - an international survey aimed to better understand the challenges people living with an ostomy face (n=4000 ostomates). Published: Claessens I, Probert R, et. al. The Ostomy Life Study: the everyday challenges faced by people living with a stoma in a snapshot. Gastrointestinal Nursing. 2015;13(5), 18-25.

Results from a multinational, randomized controlled cross-over trial comparing SenSura Mio Convex to the study participants' own appliance (n=129). Published: Walker H, Hopkins G, et. al. Raising the bar: New flexible convex ostomy appliance - a randomised controlled trial. World Council of Enterostomal Therapists Journal. 2016;36(1) S6-S11.

Results from the 2018-2019 Ostomy Life Study - an international survey (n=3112 ostomates); Coloplast market research, Data on file.

## The fit is the difference between sitting around and

hitting the town

## SenSura®Mio

SenSura<sup>®</sup> Mio with BodyFit Technology<sup>®</sup> is designed to provide a *secure fit for all body shapes*. The BodyFit Technology in all Mio products



allows the barrier to attain and maintain a secure fit. BodyFit Technology is founded on an elastic double-layer adhesive and combined with innovative barrier design to secure a better fit.

## The fit is the difference with SenSura® Mio



## **Products of Distinction**



SenSura® Mio Convex Flip

Designed to fit curves, bulges and hernias with less creasing and folding



SenSura® Mio Convex Three different types of convexity with proven leakage reduction



## SenSura® Mio Baby Long wear time without compromising skin friendliness



## Brava® Ostomy Support Belt

Relieves the heavy sensation of hernias and may help prevent hernias post-operatively



SenSura® Mio High Output

The only 1-piece convex option with a spout outlet for fecal output



Coloplast Fistula and Wound Management System (FWMS)

Available in three sizes to accommodate a large surgical incision or management of type II fistulas

## **Brava**<sup>®</sup>

# Security to live everyday life

"I was determined to play sports again. The Brava Protective Seal gave me confidence."

- Kevin\*

\*Kevin is a Brava user who has received compensation from Coloplast to provide this information.

Accessories



# Brava<sup>®</sup> supporting products

For customizing pouching system fit and for issues such as leakage and skin irritation.



## Brava<sup>®</sup> Protective Seal

- Easily shaped
- Leaves minimal residue behind
- Follows body movements
- Significantly reduces leakage compared to competitor rings

The **Brava® Ostomy Support Belt** brings relief by relieving the heavy sensation of hernias.

 Comfortable and breathable fabric



- Stretchable to fit body contours
- Easy-to-use pocket closure
- Silicone grip to prevent sliding and rolling

## What's right for your patients?



Accessories

## Think again. Think Brava®

Your complete solution to address leakage and care for your skin.

## Brava<sup>®</sup> Elastic Barrier Strips

Skin-friendly strips are used to "picture frame" around an ostomy barrier, providing extra security.

- Elastic (follow the contours of the body)
- Ensure the position of the barrier and minimize roll-up
- Skin-friendly and secure alternative to tape
- Not made with natural rubber latex

Brava® Elastic Barrier Strips- Curved For round or oval barrier shapes

Code	Size	Units	HCPCS*
120700**	5½" strips	20 strips	A4362
**2 strips =	1 billable unit of A4362		

Brava® Elastic Barrier Strips- Y shape Designed for users with challenging body shapes, e.g. hernias or skin folds

Code	Size	Units	HCPCS
120721**	5¾" strips	30 strips	A5121
** 2	1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1		

\*\*3 strips = 1 billable unit of A5121

Brava<sup>®</sup> Elastic Barrier Strips- Straight For square or rectangular barrier shapes



Code	Size	Units	HCPCS*
120740**	5¼" strips	20 strips	A4362

\*\*2 strips = 1 billable unit of A4362



Brava® Elastic Barrier Strip XL Designed for increased security & support around the barrier. Wider shape accommodates the petals of Convex Flip.

Code Size		Units HCPCS*
12076** 6 <sup>3</sup> /4" str	ips 30	strips A5121
**3 strips = 1 billable	unit of A5121	

## Brava<sup>®</sup> Protective Seal

Moldable barrier ring designed to adhere securely and protect the skin around the stoma.

- Soft, flexible material adapts to body shape as the wearer bends and stretches
- Can be easily and cleanly removed
- Alcohol-free and sting-free
- Multiple sizes for a custom fit
- Not made with natural rubber latex



Code	Stoma size	Adhesive diameter	Thickness	Seal width	Units	HCPCS*
12035	3⁄4″	1 7/8"	2.5 mm	Standard (5/8")	10	A4385
12036	3⁄4″	2 ¼"	2.5 mm	Wide (¾")	10	A4385
12038	3⁄4″	2 1⁄2″	2.5 mm	Wide (7/8")	10	A4385
12037	1 1/8″	2 ¼"	2.5 mm	Standard (5/8")	10	A4385
12039	1 ³/8″	2 1⁄2″	2.5 mm	Standard (5/8")	10	A4385
12045	3⁄4″	1 7/8"	4.2 mm	Standard (5/8")	10	A4385
12046	3⁄4″	2 ¼"	4.2 mm	Wide (¾")	10	A4385
12048	3⁄4″	2 1⁄2″	4.2 mm	Wide (7/8")	10	A4385
12047	1 1/8″	2 ¼"	4.2 mm	Standard (5/8")	10	A4385
12049	1 ³/8″	2 1⁄2″	4.2 mm	Standard (5/8")	10	A4385
Coming	soon					
12032	3⁄4″	3"	2.5 mm	Extra wide (1¹/՛ଃ")	10	A4385
12033	1 1/8″	3"	2.5 mm	Extra wide ( <sup>15</sup> / <sub>16</sub> ")	10	A4385

## Brava® Moldable Ring

Designed for dual-action leakage reduction - used to fill in uneven skin surfaces and protect the skin.

- Molds and adheres securely to dips, creases, and folds to help create a tight seal around the stoma
- Helps keep skin healthy by protecting it from output and absorbing moisture from the skin
- Easy to handle, mold and apply
- Alcohol-free and sting-free

#### Not made with natural rubber latex

Code	Size	Units	HCPCS*
120307	2 mm (thin)	10	A4385
120427	4.2 mm (thick)	10	A4385

#### **Brava® Strip Paste**

Pectin-based product that protects peristomal skin.

- Protects peristomal skin from stoma output
- Used to fill uneven skin surfaces to create a flat pouching surface
- Alcohol-free and sting-free
- Not made with natural rubber latex

Code	Size	Units	HCPCS*
26555	2 oz. total package (.2 oz. per strip)	10 strips	A4406



## Brava<sup>®</sup> Paste

A sting-free alternative to conventional alcohol-based paste.

- Creates a tight, effective seal between peristomal
- skin and the ostomy/fistula appliance • No alcohol = no need to "set"; contains pectin
- Not made with natural rubber latex

Code	Size	Units	HCPCS*
12050	2 oz.	1	A4406

## Coloplast<sup>®</sup> Paste

Used to fill in uneven skin surfaces to create a flat pouching surface.

- Provides tight sealing around the stoma
- Protects against stoma output, minimizing risk of maceration
- · Low alcohol content; no pectin
- Not made with natural rubber latex

Code	Size	Units	HCPCS*
2650	2 oz.	12	A4405

## Brava<sup>®</sup> Skin Cleanser Wipes

The wipe that cleans and cares. A convenient way to cleanse and maintain healthy skin when changing the ostomy pouching system.

- Contains ProVitamin B5 to provide a soothing and moisturizing effect
- Does not contain alcohol or perfumes
- No residue left on the skin that could affect barrier adhesion
- Not made with natural rubber latex

Code	Size	Units	HCPCS*
12080	1 wipe	15 wipes/pack	A9286

## Brava<sup>®</sup> Lubricating Deodorant

A pouch lubricant + deodorizer in one. Powder-fresh liquid neutralizes odor and allows output to slide to the bottom of the pouch for easier emptying. Not made with natural rubber latex

Code	Size	Units	HCPCS*
12060	.25 fl. oz. (7.5 mL) sachet	20 sachets	A4394
12061	8 oz. bottle	16	A4394

## Brava<sup>®</sup> Adhesive Remover Spray

- Designed to lift the edges of the barrier for ease in removal
- Dry within seconds without leaving an oily residue behind
- Alcohol-free, silicone-based, and sting-free

## Not made with natural rubber latex

Code	Size	Units	HCPCS*
120105	1.7 oz. (50 mL)	1	A4455





CCessori



#### Brava<sup>®</sup> Adhesive Remover Wipes

- Wipe away adhesive residue
- Dry within seconds without leaving an oily residue behind
- Alcohol-free, silicone-based, and sting-free
- Not made with natural rubber latex

Code	Size	Units	HCPCS*
120115	1 wipe	30 wipes/box	A4456

## Brava<sup>®</sup> Skin Barrier Spray

- · Protects against stoma output and adhesives
- Dry within seconds without building a thick or tacky layer
- Alcohol-free, silicone-based, and sting-free
- Not made with natural rubber latex

Code	Size	Units	HCPCS*
120205	1.7 oz. (50 mL)	1	A4369

## Brava® Skin Barrier Wipes

- Protects against stoma output and adhesives
- Dry within seconds without building a thick or tacky layer
- Tip: wring the wipe like a wet rag to expel the liquid, then wipe to disperse on skin
- Alcohol-free, silicone-based, and sting-free
- Not made with natural rubber latex

Code	Size	Units	HCPCS*
120215	1 wipe	30 wipes/box	A5120

## **Brava® Protective Sheet**

Hydrocolloid skin barrier that protects the skin from stoma output.

- Offers effective protection for peristomal skin and "3-way stretch"
- Provides a ready surface for good adhesion
- Adheres to moist skin
- Not made with natural rubber latex

Code	Size	Units	HCPCS*
32105	4 x 4" (10 x 10 cm)	10	A4362
32155	6 x 6" (15 x 15 cm)	5	A5121
32205	8 x 8" (20 x 20 cm)	5	A5122

## Brava<sup>®</sup> Powder

Non-sterile protective powder for moisture absorption.

- May use in conjunction with Brava<sup>®</sup> Skin Barrier Spray or PREP<sup>™</sup> to "crust", thereby protecting moist skin
- Not made with natural rubber latex

Code	Size	Units	HCPCS*
19075	1 oz. (25 g)	16	A4371







## Hex-On®

#### Odor Antagonist

- Helps control airborne odors
- Non-aerosol spray
- Not made with natural rubber latex



Accessorie

Product	Code	Size	Units
Fresh Linen	7583	2 fl. oz./59 mL	12
Fragrance Free	67589	2 fl. oz./59 mL	12

## **Ostomy Scissors**

Curved ostomy scissors.



Code	Size	Units	HCPCS*
95050	One size	30/bx	N/A

## Coloplast<sup>®</sup> Skin Barrier Rings

Not made with natural rubber latex

Soft and flexible barrier rings provide skin protection for stomas, fistulas, and draining wounds.

- Non-pectin based
- Not made with natural rubber latex



Code	Size	Units	HCPCS*
2310	³⁄ଃ" (10 mm)	30	A4404
2315	³∕₅" (15 mm)	30	A4404
2320	¾″ (20 mm)	30	A4404
2325	1" (25 mm)	30	A4404
2330	11/8" (30 mm)	30	A4404
2340	1³⁄₅" (40 mm)	30	A4404
2350	2" (50 mm)	30	A4404

## Prep<sup>™</sup> Protective Skin Barrier

Clear film forms a breathable barrier protecting against skin irritants such as adhesives, urine and feces. No need to remove between applications.



Not made with natural rubber latex

Wipes			
Code	Size	Units	HCPCS*
2041	54 wipes/pack	12 packs/box	A5120
Dabber			
Code	Size	Units	HCPCS*
0925	2 fl. oz. (59 mL)	12	A4369

## Pouch Filters (Filtrodor®)

Charcoal-activated, deodorizing filter with self-sealing foam center; adheres to pouch.



• Not made with natural rubber latex

Code	Units	HCPCS*	
0509	50	A4368	
	each box contains 5 packs of 10		

## Drainable Pouch Clamp

Simple closure clamp is easy to manipulate, even for those with low dexterity. • Not made with natural rubber latex



Code	Units	HCPCS*
9500	20	A4363

## Brava<sup>®</sup> Ostomy Support Belt

This belt relieves the heavy sensation of hernias and may help prevent hernias post-operatively. Comfortable and breathable fabric, stretchable to fit body G

contours, intuitive pocket closure, and includes a silicone grip to prevent sliding and rolling.

## • Not made with natural rubber latex

#### • All sizes are an 8" width

Code	Size	Units	HCPCS*
12003	S (29-32")	1	A4396
12004	M (31-35")	1	A4396
12005	L (34-39")	1	A4396
12006	XL (37-42")	1	A4396
12007	2XL (40-46")	1	A4396
12008	3XL (44-51")	1	A4396
Coming :	soon		
12009	4XL (50-58")	1	A4396

#### Brava<sup>®</sup> Belt for SenSura<sup>®</sup> Mio

Neutral gray color was designed to blend in under clothing. Four-tab connection provides extra stability across the entire peristomal area. Designed specifically for use with SenSura® Mio products with belt tabs.

## Not made with natural rubber latex

Code	Size	Units	HCPCS*
4237	40" length	15	A4367
4247	61" length	15	A4367

# Brava<sup>®</sup> Belt

Comfortable, adjustable belt offers extra security for people with a flush or retracted stoma by supporting a convex barrier and allowing the stoma to protrude.

- Neutral color won't show under clothing
- Not made with natural rubber latex

Code	Size	Units	HCPCS*
4215	43¹∕₃" length	15	A4367
4220	49" length	15	A4367

# **Urostomy Night Bag**


Includes a fixed connector that fits directly to SenSura® Mio, SenSura® and Assura® urostomy pouches. The bag holds up to 2 liters and has a 120 cm long flexible tube with anti-reflux valve. Bag hanger (#5070) not included. Order separately.




Not made with natural rubber latex

Code	Units	HCPCS*
21365	10	A4357

# **Bed Drainage Bag**





# Not made with natural rubber latex

Code	Units	HCPCS*
14010	6	A4357

# **Bag Hanger**

Reusable, sturdy, plastic drainage bag hanger secures tubing and ensures bag stability. Not made with natural rubber latex



Code	Size	Units	HCPCS*
5070	One Size	1	Medicare A9900
			Other A4421

# Extra Drain Port

Code

14015

For use with Fistula Wound Management Systems (page 167). Drain port can be easily applied to the lid to allow for catheter insertion. For catheters ranging from 18FR - 22FR, the drain port is ready for use. For catheters ranging from 24FR-30FR, pull off the top of the drain port - DO NOT CUT. Not made with natural rubber latex



Units

12

# Coloplast® Transparent Fistula Lids

Replacement lids match the transparent lid included with each fistula pouch. • Not made with natural rubber latex

Size	Units
Mini (matches 14050)	10
Midi (matches 14060)	10
Maxi (matches 14070)	10
	Mini (matches 14050) Midi (matches 14060) Maxi (matches 14070)

# Irrigation Sets

# Coloplast<sup>®</sup> Economy Version

Durable water bag with integrated thermometer, bold print, large opening. Can be used with 1- or 2-piece ostomy systems. • Not made with natural rubber latex

Code	Units	Includes	HCPCS*
1500	1	1-Water Bag/Regulator	A4398
		1-Cone (#1110)	A4399

# Assura<sup>®</sup> Hospital Version

Durable water bag with integrated thermometer, bold print, large opening. Can be used with 1- or 2-piece ostomy systems. • Not made with natural rubber latex

Code	Units	Includes	HCPCS*
1500H	1	1-Water Bag/Regulator 1-Cone (#1110) 4-Adhesive Irrigation Sleeves (#1003) 2-Closure Clips (#9500)	A4398 A4399 A4397

# SenSura® or Assura® Deluxe Version

Components supplied in a travel bag.

Not made with natural rubber latex



Code	Units	Includes	HCPCS*
12830	1	1-Water Bag/Regulator	A4398
		1-Cone (#1110)	A4399
		1-Irrigation Faceplate (#12820)	A4361
		1-Belt (#4215)	A4367
		2-Irrigation Sleeves (#12836)	A4397

# Irrigation Sleeves

# Irrigation Sleeves (with adhesive)

Disposable irrigation sleeve (adhesive-backed).

# Not made with natural rubber latex

Transparent		
Code	Units	HCPCS*
1003	10	A4397

# Irrigation Sleeves (without adhesive)

Irrigation sleeves for use with SenSura® Click or Assura® barriers or faceplate. Consult chart below to order the correct transparent sleeves.

# Not made with natural rubber latex

Code	Size	Units	HCPCS*
12834		5	A4397
12835		5	A4397
12836		5	A4397

Irrigation Sleeve #	Use with SenSura® Click Product Code	Units
12834	- 10011, 10012, 10013, 10015- 10017, 11011, 11012, 11015, 11016	5
12835	- 10021, 10022, 10023, 10025- 10027, 11021-11025	5
12836	- 10031, 10032, 10035, 10036, 11031, 11032, 11035, 11036	5



# Irrigation Accessories

# Assura® Irrigation Faceplate

For use with ostomy belt 4215 or 4220. (Use with #12836 Irrigation Sleeve) Note: not compatible with Brava® Belt for SenSura® Mio



# • Not made with natural rubber latex

Code	Units	HCPCS*
12820	1	A4361

# **Coloplast Stoma Cone**

Soft and flexible stoma cone with tubing.	
<ul> <li>Not made with natural rubber latex</li> </ul>	

Code	Units	HCPCS*
1110	1	A4399



# 2-piece Click (Mechanical coupling)



Pouch can be repositioned while it's still on the body



Tactile experience of pressing two pieces together can be reassuring



Lock ring gives an audible cue that it's securely locked



Coupling is strong and secure!

# SenSura®Mio

# The fit is the difference between sitting on the edge and diving right in

Mike is a SenSura® Mio Click user.

# SenSura<sup>®</sup> Mio Click

The BodyFit Technology<sup>®</sup> in all SenSura Mio barriers combines an elastic adhesive, that bends and stretches with the body, with an innovative barrier design to provide a secure and comfortable fit.



NOTE: SenSura^ Mio Click is NOT compatible with SenSura^ Click or Assura systems.

# SenSura® Mio Click Standard Wear Barrier

- Flat
- Belt tabs for use with 4237, 4247
- Not made with natural rubber latex

Cut-to-fit			Units: 5
Code	Stoma Size	Color	HCPCS*
10502	³⁄8-1³∕8″ (10-35 mm)		A4414
10512	<sup>3</sup> / <sub>8</sub> -1 <sup>3</sup> / <sub>4</sub> " (10-45 mm)		A4414
10522	<sup>3</sup> / <sub>8</sub> -2 <sup>1</sup> / <sub>8</sub> " (10-55 mm)		A4414
10531	<sup>3</sup> / <sub>8</sub> -2 <sup>9</sup> /16" (10-65 mm)	•	A4415
Pre-cut			Units: 5
Code	Stoma Size	Color	HCPCS*
10503	<sup>3</sup> ⁄4" (20 mm)		A4414
10506	∛8″ (22 mm)		A4414
10504	1" (25 mm)		A4414
10505	1¼″ (28 mm)		A4414
10513	1³∕16″ (30 mm)		A4414
10514	1¼" (32 mm)		A4414
10515	1¾" (35 mm)		A4414
10516	1½" (38 mm)		A4414
10523	1%16″ (40 mm)		A4414
10524	1¾" (45 mm)		A4414
10525	2" (50 mm)		A4414

# SenSura® Mio Click Convex Light Barrier 🦟

# Standard wear

- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility

# • Not made with natural rubber latex

Cut-to-fit			Units: 5
Code	Stoma Size	Color	HCPCS*
16901	3⁄8 - 3⁄4" (10-20 mm)		A4373
16911	5⁄8 - 1³∕16″ (15-30 mm)		A4373
16921	5⁄8 - 19∕16″ (15-40 mm)		A4373
16931	<sup>5</sup> / <sub>8</sub> - 2 <sup>1</sup> / <sub>16</sub> " (15-53 mm)		A4373

Indicates convexity



# SenSura® Mio Click Convex Light Barrier

- Standard wear
- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility





Pre-cut			Units: 5
Code	Stoma Size	Color	HCPCS*
16903	<sup>11</sup> /16" (18 mm)		A4373
16913	<sup>13</sup> / <sub>16</sub> " (21 mm)		A4373
16914	1" (25 mm)		A4373
16915	11⁄8" (28 mm)		A4373
16923	11⁄4" (31 mm)		A4373
16924	1¾" (35 mm)		A4373
16933	11⁄2" (38 mm)		A4373
16934	15⁄8" (41 mm)		A4373

# SenSura® Mio Click Deep Convex Barrier

- Standard wear
- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- Not made with natural rubber latex

Cut-to-fit			Units: 5
Code	Stoma Size	Color	HCPCS*
16941	¾ - ¾" (10-20 mm)		A4373
16951	5⁄8 - 1³∕16″ (15-30 mm)		A4373
16961	5⁄8 - 1⅔16″ (15-40 mm)		A4373
16971	<sup>5</sup> / <sub>8</sub> - 2 <sup>1</sup> / <sub>16</sub> " (15-53 mm)	-	A4373
Pre-cut			Units: 5
Code	Stoma Size	Color	HCPCS*
16942	<sup>11</sup> / <sub>16</sub> " (18 mm)		A4373
16953	<sup>13</sup> / <sub>16</sub> " (21 mm)		A4373
16954	1" (25 mm)		A4373
16955	1½″ (28 mm)		A4373
16963	1¼″ (31 mm)		A4373
16964	1³⁄₃" (35 mm)		A4373
16973	11⁄2" (38 mm)		A4373
16974	15⁄8″ (41 mm)		A4373



NOTE: SenSura® Mio Click is NOT compatible with SenSura® Click or Assura® systems.

# SenSura<sup>®</sup> Mio Click Convex Flip 2-Piece Barrier **(**

- Designed for curves, bulges and hernias
- Curved, star-shaped barrier conforms to outward areas with less creasing and folding
- "Flip to fit" design allows for easy application
- Belt tabs for use with 4237, 4247
- Not made with natural rubber latex

### Cut-to-fit Units: 5 Code Stoma Size Barrier Size HCPCS\* Color s 18500 <sup>3</sup>/8 - 1<sup>3</sup>/8" (10-35 mm) A4373 3/8 - 13/4" (10-45 mm) 18510 A4373 18520 3/8 - 21/8" (10-55 mm) A4373 18530 <sup>3</sup>/8 - 2<sup>9</sup>/16" (10-66 mm) A4373

Pre-cut					Units: 5
Code	Stoma Size	Max Cut	Barrier Size	Color	HCPCS*
18501	¾" (20 mm)	1³/8" (35 mm)	S		A4373
18511	1" (25 mm)	1¾" (45 mm)	M		A4373
18512	1 ³/16" (30 mm)	1¾" (45 mm)	M		A4373
18513	1 ³/8" (35 mm)	1¾" (45 mm)	M		A4373
18521	1 ³/16" (30 mm)	2¹/8″ (55mm)			A4373
18522	1 ³/8" (35 mm)	2¹/8″ (55 mm)			A4373
18523	1 <sup>9</sup> /16" (40 mm)	2¹/8" (55 mm)			A4373
18524	1 ¾" (45 mm)	2¹/8" (55 mm)			A4373
18531	1 ¾" (45 mm)	2º/16" (66 mm)	<u>×</u>		A4373
18532	2" (50 mm)	2º/16" (66 mm)	×-		A4373





# SenSura<sup>®</sup> Mio Click MAXI Drainable Pouch

- EasiClose<sup>™</sup> WIDE Outlet
- Full-circle filter
- Neutral gray, water-resistant textile
- Opaque pouches feature peek-a-boo inspection
   window for easy application
- Not made with natural rubber latex

# Transparent 11" (28 cm) / 655 mL

Code	Color	HCPCS*
11463		A4427
11473		A4427
11483		A4427
11487	-	A4427
Opaque 11" (28	3 cm) / 655 mL	Units: 10
Code	Color	HCPCS*
11462		A4427
11472		A4427
11482		A4427
11486	-	A4427

	uch VIDE Outlet	
Transparent 1	1" (28 cm) / 655 mL	Units: 10
Code	Color	HCPCS*
11474		A4426
11484		A4426
11488		A4426

### SenSura<sup>®</sup> Mio Click MIDI **Drainable Pouch** EasiClose<sup>™</sup> WIDE Outlet Full-circle filter Neutral gray, water-resistant textile · Peek-a-boo inspection window for easy application Not made with natural rubber latex Opaque 10<sup>1</sup>/4" (26 cm) / 485 mL Units: 10 Code Color HCPCS\* 11461 A4427 11471 A4427

Units: 10

# 2-Piece SenSura® Mio Click

# SenSura<sup>®</sup> Mio Click High Output Pouch with Tap Outlet

- Special full-circle filter designed specifically
  for ileostomy output
- Neutral gray, water-resistant textile
- Transparent pouches feature a printed scale for visual indication of output volume
- Opaque pouches feature a peek-a-boo inspection window and tuckaway option to keep the outlet from rubbing on the wearer's leg
- Wider-lumen outlet provides effective drainage for patients with thicker output; connects to Coloplast fistula drain bag #14010
- Not made with natural rubber latex

Transparent 12 ½" (32 cm) / 1080 mL		Units: 10
Code	Color	HCPCS*
18610		A4413
18611		A4413
18612		A4413
18613	•	A4413

Opaque 12 ½" (32 cm) / 1080 mL		Units: 10
Code	Color	HCPCS*
18600		A4413
18601		A4413
18602		A4413
18603		A4413

# SenSura<sup>®</sup> Mio Click High Output Pouch with Soft Outlet

- Special full-circle filter designed specifically for ileostomy output
- Neutral gray, water-resistant textile
- Transparent pouches feature a printed scale for visual indication of output volume
- Opaque pouches feature a peek-a-boo
   inspection window and tuck-away option to keep the outlet from
   rubbing on the wearer's leg
- Soft outlet provides effective drainage; connects to a standard drain bag or Coloplast fistula drain bag #14010
- Not made with natural rubber latex

Transparent 12 ½" (32 cm) / 1080 mL		Units: 10
Code	Color	HCPCS*
18640		A4413
18641		A4413
18642		A4413
Opaque 12 ½" (32 cm) / 1080 mL		Units: 10
Opaque 12 ½" (32 cm) / 1080 mL Code	Color	Units: 10 HCPCS*
	Color	
Code	Color	HCPCS*

\*For important information regarding the use of HCPCS codes listed and code definitions see page 192



# SenSura® Mio Click Drainable Pouch with Tap Outlet

 Special full-circle filter designed specifically for ileostomy output



Neutral gray, water-resistant textile

# Opaque pouches feature a peek-a-boo inspection window and tuck-away option to keep the outlet from rubbing on the wearer's leg

- Wider-lumen outlet provides effective drainage for patients with thicker output; connects to Coloplast fistula drain bag #14010
- Not made with natural rubber latex

# Opaque 9" (23 cm) / 650 ml

Opaque 9" (23 cm) / 650 mL		Units: 10
Code	Color	HCPCS*
18620		A4427
18621		A4427
18622		A4427

# SenSura<sup>®</sup> Mio Click MAXI Closed Pouch

- Full-circle filter
- Neutral gray, water-resistant textile
- · Opaque pouches feature peek-a-boo inspection window for easy application
- Not made with natural rubber latex



Transparent 8¼" (20½ cm) / 570 mL		Units: 30
Code	Color	HCPCS*
11426		A4423
<mark>Opaque</mark> 8¼" (20½ cm) / 570 mL		Units: 30
Code	Color	HCPCS*
11402		A4423
11413		A4423
11422		A4423
11425		A4423

# SenSura<sup>®</sup> Mio Click MIDI **Closed Pouch**

- Full-circle filter
- Neutral gray, water-resistant textile
- Peek-a-boo inspection window for easy application
- Not made with natural rubber latex



Units: 30

### Opaque 71/4" (181/2 cm) / 275 mL

Code	Color	HCPCS*
11401		A4423
11412		A4423
11421		A4423

# Urostomy

# SenSura® Mio Click MAXI Urostomy Pouch

- Multi-chambers minimize sloshing
- Wider anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Neutral gray, water-resistant textile
- Not made with natural rubber latex



Transparent 10½" (26½ cm) / 480 ml	L	Units: 10
Code	Color	HCPCS*
11492		A4433
11498		A4433
11500		A4433
Opaque 10½" (26½ cm) / 480 mL		Units: 10
Code	Color	HCPCS*
11491		A4433
11497		A4433
11499		A4433

NOTE: SenSura® Mio Click is NOT compatible with SenSura® Click or Assura<sup>®</sup> systems. SenSura<sup>®</sup> Click and Assura<sup>®</sup> are compatible (pages 88-93).

# SenSura<sup>®</sup> Click Flat Barrier

- Standard wear
- Belt tabs for use with 4215, 4220
- Not made with natural rubber latex



Cut-to-fit			Units: 5
Code	Stoma Size	Color	HCPCS*
10011	³⁄ଃ - 1³∕ଃ" (10-35 mm)		A4414
10021	<sup>3</sup> / <sub>8</sub> - 1 <sup>3</sup> / <sub>4</sub> " (10-45 mm)		A4414
10031	³⁄8 - 2¹⁄4″ (10-55 mm)		A4414
10041	<sup>3</sup> / <sub>8</sub> - 2 <sup>1</sup> / <sub>2</sub> " (10-65 mm)		A4415
Pre-cut			Units: 5
Code	Stoma Size	Color	HCPCS*
10012	<sup>13</sup> ⁄ <sub>16</sub> " (20 mm)		A4414
10013	1" (25 mm)		A4414
10022	1 <sup>13</sup> ⁄16" (30 mm)		A4414
10023	1¾" (35 mm)		A4414
10032	1%16" (40 mm)		A4414

# SenSura® Click Xpro Flat Barrier

- Extended wear
- Belt tabs
- Not made with natural rubber latex

Cut-to-fit			Units: 5
Code	Stoma Size	Color	HCPCS*
10015	¾ - 1¾" (10-35 mm)		A4409
10025	¾ - 1¾" (10-45 mm)		A4409
10035	¾ - 2¼″ (10-55 mm)		A4409
10045	<sup>3</sup> / <sub>8</sub> - 2 <sup>1</sup> / <sub>2</sub> " (10-65 mm)		A4410
Pre-cut			Units: 5
Code	Stoma Size	Color	HCPCS*
10016	<sup>13</sup> /16″ (20 mm)		A4409
10017	1" (25 mm)		A4409
10026	1¾16″ (30 mm)		A4409
10027	1¾″ (35 mm)		A4409
10036	1%16″ (40 mm)		A4409



# SenSura<sup>®</sup> Click Convex Light Barrier 🖊

- Standard wear
- Belt tabs
- Not made with natural rubber latex



Units: 5

HCPCS\*

Cut-to-fit			Units: 5
Code	Stoma Size	Color	HCPCS*
11011	5⁄8 - 7⁄8″ (15-22 mm)		A4373
11021	⁵⁄8 - 1⁵⁄16″ <b>(15-33 mm)</b>		A4373
11031	5/8 - 1 <sup>11</sup> /16" (15-43 mm)		A4373
11041	<sup>5</sup> /8 - 2 <sup>1</sup> /16" (15-53 mm)		A4373

Pre-cut	
Code	Stoma Size
11012	7∕8″ (22 mm)
11022	1" (25 mm)

11012	7∕8″ (22 mm)	A4373
11022	1" (25 mm)	A4373
11023	11/8" (28 mm)	A4373
11024	11⁄4" (31 mm)	A4373
11032	1¾" (35 mm)	A4373

Color



Cut-to-fit			Units: 5
Code	Stoma Size	Color	HCPCS*
11015	⁵⁄8 - ⅔" (15-22 mm)		A4407
11025	5⁄8 - 15∕16″ (15-33 mm)		A4407
11035	<sup>5</sup> /8 - 1 <sup>11</sup> /16" (15-43 mm)		A4407
11045	5⁄8 - 21⁄16" (15-53 mm)		A4408

# Pre-cut

Pre-cut		Units: 5	
Code	Stoma Size	Color	HCPCS*
11016	⅔″ (22 mm)		A4407
11026	1" (25 mm)		A4407
11027	11/8" (28 mm)		A4407
11028	11⁄4" (31 mm)		A4407
11036	1³∕≋" (35 mm)		A4407

For Assura® CONVEX barriers, see pages 93-95.

Indicates convexity

# Xpro vs. Standard Wear: how do I choose? Extended wear doesn't necessarily indicate longer wear time!



SenSura® Xpro (Extended Wear) is engineered to offer extra protection for liquid ileostomies and urostomies.

SenSura® Standard Wear focuses on moisture absorption and may offer better adhesion for those in more humid climates.

SenSura <sup>®</sup> Click MAXI Drainable Pouch • EasiClose <sup>®</sup> WIDE Outlet • Filter • Soft cloth on back (and front of opaque) • Not made with natural rubber latex		
Transparent 11½" (29½ cm) / 655 mL		Units: 20
Code	Color	HCPCS*
11194		A4427
11195		A4427
11196		A4427
11197		A4427
Opaque 11½" (29½ cm) / 655 mL		Units: 20
Code	Color	HCPCS*
11124		A4427
11125		A4427
11126		A4427
11127		A4427
SenSura <sup>®</sup> Click MIDI Drainable Pouch • EasiClose <sup>®</sup> WIDE Outlet • Filter • Soft cloth on front and back • Not made with natural rubber latex		
<mark>Opaque</mark> 10½" (26½ cm) / 485 mL		Units: 20
Code	Color	HCPCS*
11114		A4427
11115		A4427

SenSura® Click MAXI Drainable Pouch • EasiClose" WIDE Outlet • Without filter • Soft cloth on back (and front of opaque) • Not made with natural rubber latex		
Transparent 111/2" (291/2 cm) / 655 mL		Units: 20
Code	Color	HCPCS*
11185		A4426
11186		A4426
Opaque 11½" (29½ cm) / 655 mL		Units: 20
Code	Color	HCPCS*
11135		A4426
11136		A4426

# SenSura® Click Magnum High Output Drainable Pouch

- Wider-lumen outlet provides effective drainage for patients with thicker output; connects to Coloplast fistula drain bag #14010
- Filter
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex

Transparent 12½" (32 cm) / 1160 mL		Units: 10
Code	Color	HCPCS*
19050		A4413
19051		A4413
19052		A4413
<mark>Opaque</mark> 12½" (32 cm) / 1160 mL		Units: 10
Code	Color	HCPCS*
19040		A4413
19041		A4413
19041 19042		A4413 A4413

For Assura<sup>®</sup> Ileo Night Bag, see page 99.

SenSura® Click MAXI Closed Pouch • Filter • Soft cloth on back (and front of opaque) • Not made with natural rubber latex		Ô
Transparent 81/2" (21 cm) / 570 mL		Units: 30
Code	Color	HCPCS*
10187		A4423
Opaque 81/2" (21 cm) / 570 mL		Units: 30
Code	Color	HCPCS*
10164		A4423
10165		A4423
10166		A4423
10167	-	A4423

# SenSura<sup>®</sup> Click MIDI **Closed Pouch**

- Filter
- Soft cloth on front and back
- Not made with natural rubber latex

6	J	~		
	6		1	١
		2	6/	
1				y

# Opaque 7" (18 cm) / 275 mL

Opaque 7" (18 cm) / 275 mL		Units: 30
Code	Color	HCPCS*
10154		A4423
10155		A4423
10156		A4423

For Assura<sup>®</sup> stoma caps, see page 100.

# Urostomy


# SenSura® Uro

# Wide anti-reflux valves support varying ways of managing stents:

- Stents threaded into a pouch connected to a night drainage bag
- Rolling the stents in the top of the pouch, making a loop
- Stents that have been cut

# Welded "multi-chambers" minimize ballooning

Wide anti-reflux valves



Units: 10

# SenSura® Click Uro MAXI Urostomy Pouch

- Multi-chambers minimize sloshing
- Wider anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex

# Transparent 103/8" (26 cm) / 480 mL

Code	Color	HCPCS*
11854		A4433
11855		A4433
11856		A4433
<mark>Opaque</mark> 10⅔" (26 cm) / 480 mL		Units: 10
Code	Color	HCPCS*
11844		A4433
11845		A4433
11846		A4433



# 2-Piece SenSura® Click

# SenSura<sup>®</sup> Click Uro MIDI **Urostomy Pouch**

- Multi-chambers minimize sloshing
- · Wider anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex




Transparent 91/2" (24 cm) / 430 mL		Units: 10
Code	Color	HCPCS*
11851		A4433

Opaque 91/2" (24 cm) / 430 mL		Units: 10
Code	Color	HCPCS*
11841		A4433

For Assura® urostomy micro-pouches and minicaps, see page 100.

# Urostomy Night Bag and Bag Hanger

See Accessories page 73 for more information.



# Assura<sup>®</sup> Flat Barrier

- Standard wear
- Belt tabs for use with 4215, 4220
- Not made with natural rubber latex

### Cut-to-fit

Cut-to-fit			Units: 5
Code	Stoma Size	Color	HCPCS*
2881	³⁄ଃ - 1³⁄ଃ" (10-35 mm)		A4414
2882	¾ - 1¾" (10-45 mm)		A4414
2883	<sup>3</sup> ⁄ <sub>8</sub> - 2 <sup>1</sup> ⁄ <sub>4</sub> " (10-55 mm)		A4414
Pre-cut			Units: 5
Code	Stoma Size	Color	HCPCS*
12841	<sup>13</sup> / <sub>16</sub> " (20 mm)		A4414
12842	1" (25 mm)		A4414
12844	1³⁄16" (30 mm)		A4414
12845	1³⁄8" (35 mm)		A4414
12847	1%16″ (40 mm)		A4414
12848	1¾″ (45 mm)		A4414



NOTE: Assura® is compatible with SenSura® Click , but NOT with SenSura® Mio Click.

# Assura<sup>®</sup> Flat Barrier

- Extended wear
- Belt tabs for use with 4215, 4220
- Not made with natural rubber latex

Units: 5

Cut-to-fit			Units: 5
Code	Stoma Size	Color	HCPCS*
2831	∛8 - 1¾8″ (10-35 mm)		A4409
2832	¾ - 1¾" (10-45 mm)		A4409
2833	<sup>3</sup> / <sub>8</sub> - 2 <sup>1</sup> / <sub>4</sub> " (10-55 mm)		A4409
Pre-cut			Units: 5
Code	Stoma Size	Color	HCPCS*
14250	5⁄⁄8″ (15 mm)		A4409
14251	<sup>13</sup> ⁄16″ (20 mm)		A4409
14252	1" (25 mm)		A4409
14254	1³⁄16" (30 mm)		A4409
14255	1³⁄ଃ" (35 mm)		A4409
14257	1%16″ (40 mm)		A4409
14258	1¾" (45 mm)		A4409

# Assura<sup>®</sup> Convex Light Barrier 🦟

- Standard wear
- Belt tabs for use with 4215, 4220
- Not made with natural rubber latex

# Cut-to-fit

Code	Stoma Size	Color	HCPCS*
14261	<sup>5</sup> /8 - <sup>15</sup> /16" (15-23 mm)		A4373
14262	<sup>5</sup> /8 - 1 <sup>5</sup> /16" (15-33 mm)		A4373
14263	<sup>5</sup> /8 - 1 <sup>11</sup> /16" (15-43 mm)		A4373

Pre-cut			Units: 5
Code	Stoma Size	Color	HCPCS*
14271	<sup>3</sup> ⁄4" (18 mm)		A4373
14272	<sup>7</sup> /8" (22 mm)		A4373
14273	1" (25 mm)		A4373
14274	1 <sup>1</sup> /8" (28 mm)		A4373
14275	11⁄4" (31 mm)		A4373
14276	1¾″ (35 mm)		A4373
14277	1½″ (38 mm)		A4373
14278	1 <sup>5</sup> /8" (41 mm)		A4373



Indicates convexity

NOTE: Assura® is compatible with SenSura® Click, but NOT with SenSura® Mio Click.

# Assura<sup>®</sup> Deep Convex Light Barrier

- Extended wearBelt tabs for use with 4215, 4220
- Not made with natural rubber latex



Cut-to-fit			Units: 5
Code	Stoma Size	Color	HCPCS*
14281	<sup>5</sup> /8 - <sup>15</sup> /16" (15-23 mm)		A4407
14282	<sup>5</sup> /8 - 1 <sup>5</sup> /16" (15-33 mm)		A4407
14283	<sup>5</sup> /8 - 1 <sup>11</sup> / <sub>16</sub> " (15-43 mm)		A4407
Pre-cut			Units: 5
Code	Stoma Size	Color	HCPCS*
14291	<sup>3</sup> ⁄4" (18 mm)		A4407
14292	7∕₃″ (22 mm)		A4407
14293	1" (25 mm)		A4407
14294	11/8" (28 mm)		A4407
14295	1¼4" (31 mm)		A4407
14296	1¾" (35 mm)		A4407
14297	11/2" (38 mm)		A4407
14298	15⁄8" (41 mm)		A4407
<ul><li>Standard</li><li>Belt tabs</li></ul>	Deep Convex Barrier		Ò

( ) i		

Cut-to-fit			Units: 5
Code	Stoma Size	Color	HCPCS*
12716	5⁄8 - 15∕16″ (15-33 mm)		A4373
12719	5⁄8 - 111⁄16" (15-43 mm)		A4373



# Assura<sup>®</sup> Convex Barrier *—*

Standard wear

Pre-cut

- Belt tabs for use with 4215, 4220
- Not made with natural rubber latex



Units: 5

Code	Stoma Size	Color	HCPCS*
12701	⁵⁄₃" (15 mm)		A4373
12702	¾" (18 mm)		A4373
12703	⅔″ (22 mm)		A4373
12704	1" (25 mm)		A4373
12705	1½″ (28 mm)		A4373
12706	1¼″ (31 mm)		A4373
12707	1¾" (35 mm)		A4373
12708	1½″ (38 mm)		A4373
12709	15⁄%" (41 mm)		A4373

Assura <sup>®</sup>	Deep	Convex	Barrier	

- Extended wear
- Belt tabs for use with 4215, 4220
- Not made with natural rubber latex



Cut-to-fit			Units: 5
Code	Stoma Size	Color	HCPCS*
14243	<sup>5</sup> /8 - <sup>15</sup> /16" (15-23 mm)		A4407
14246	<sup>5</sup> /8 - 1 <sup>5</sup> /16" (15-33 mm)		A4407
14249	<sup>5</sup> /8 - 1 <sup>11</sup> /16" (15-43 mm)		A4407
Pre-cut			Units: 5
Code	Stoma Size	Color	HCPCS*
14232	¾" (18 mm)		A4407
14233	<sup>7</sup> /8" (22 mm)		A4407
14234	1" (25 mm)		A4407
14235	1 <sup>1</sup> /8" (28 mm)		A4407
14236	11⁄4" (31 mm)		A4407
14237	1¾″ (35 mm)		A4407
14238	11⁄2" (38 mm)		A4407
14239	1 <sup>5</sup> /8" (41 mm)		A4407

Indicates convexity

NOTE: Assura® is compatible with SenSura® Click, but NOT with SenSura® Mio Click.

# Assura<sup>®</sup> MAXI Drainable Pouch

- EasiClose<sup>™</sup> WIDE Outlet
- Filter
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex



# Transparent 111/2" (291/2 cm) / 655 mL

PCS* 1427
427
427
1427
: 10
PCS*
427
427

Assura <sup>®</sup> MAXI	
Drainable Dauch	

•

Drainable Pouch	1
EasiClose™ WIDE Outlet	1000
Without filter	
Soft cloth on front and back	
Not made with natural rubber latex	A DECK

Opaque 11 <sup>1</sup> /2" (29 <sup>1</sup> /2 cm) / 655 mL		Units: 10
Code	Color	HCPCS*
15982		A4426
15983		A4426

A	ssura <sup>®</sup> MIDI
D	rainable Pouch
•	EasiClose" WIDE Outlet
•	Filter
•	Soft cloth on front and back
•	Not made with natural rubber latex

Opaque 101/2" (261/2 cm) / 485 mL		Units: 10
Code	Color	HCPCS*
15954		A4427
15955		A4427
15956		A4427

Assura® MAXI Drainable Pouch • EasiClose" Outlet • Filter • Soft cloth on back (and front of opaque) • Not made with natural rubber latex		
Transparent 111½" (29½ cm) / 655 mL		Units: 10
Code	Color	HCPCS*
13974		A4427
13975		A4427
13976		A4427
Opaque 111/2" (291/2 cm) / 655 mL		Units: 10
Code	Color	HCPCS*
13984		A4427
13985		A4427
13986		A4427
Assura® MIDI		

# Drainable Pouch

_		
	EasiClose <sup>™</sup>	Outlet

- Filter
- Soft cloth on front and back
- Graphic pattern
- Not made with natural rubber latex

# Opaque 10<sup>1</sup>/<sub>2</sub>" (26<sup>1</sup>/<sub>2</sub> cm) / 485 mL

Code	Color	HCPCS*
13964		A4427
13965		A4427
13966		A4427

Assura® MINI Drainable Pouch • EasiClose <sup>®</sup> Outlet • Filter • Soft cloth on front and back • Not made with natural rubber latex		
Opaque 9¼" (23½ cm) / 280 mL		Units: 10
Code	Color	HCPCS*
13924		A4427
13925		A4427

Units: 10

NOTE: Assura® is compatible with SenSura® Click, but NOT with SenSura® Mio Click.

# Assura<sup>®</sup> MAXI Drainable Pouch

- Clamp closure
- Without filter
- Soft cloth on front and back
- Not made with natural rubber latex



Opaque 11¼" (28½ cm) / 760 mL		Units: 10
Code	Color	HCPCS*
12577		A4426
12578		A4426
12579		A4426

Assura® MIDI Drainable Pouch	C
Clamp closure	
Without filter	
<ul> <li>Soft cloth on front and back</li> </ul>	C. T. V
Not made with natural rubber latex	
$O_{1000} = 03/2/2(241/com)/(400 m)$	Linita 10

Opaque 93/4" (241/2 cm) / 400 mL		Units: 10
Code	Color	HCPCS*
12527		A4426
12528		A4426
12529		A4426

# Assura<sup>®</sup> High Output Drainable Pouch

- Soft outlet can connect to Coloplast's fistula drain bag #14010, or a drain bag with a standard connection
- Filter
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex

Transparent 121/2" (32 cm) / 750 mL		Units: 10
Code	Color	HCPCS*
8114		A4413
8115		A4413
8116		A4413

Opaque 121⁄2" (32 cm) / 750 mL		Units: 10
Code	Color	HCPCS*
2846		A4413
2847		A4413

For bed drainage bag and bag hanger, see page 73.

2837	
Assura® MAXI Closed Pouch • Filter • Soft cloth on back (and front of opaque) • Not made with natural rubber latex	(
Transparent 81/2" (21 cm) / 570 mL	
Code	Color
12374	
12375	
12376	
Opaque 81/2" (21 cm) / 570 mL	
Code	Color
12384	

# Assura<sup>®</sup> Ileo Night Bag

- 2 liter capacity for very high output stomas
- Clamp closure
- Comes with bag hanger
- Filter
- Not made with natural rubber latex

Opaque 2 L		Units: 5
Code	Color	HCPCS*
2836		A4413
2837		A4413

# Assu Close

# Filt

Code	Color	HCPCS*
12374		A4423
12375		A4423
12376		A4423
Opaque 81⁄2" (21 cm) / 570 mL		Units: 30
Code	Color	HCPCS*
12384		A4423
12385		A4423
12386		A4423

# Assura<sup>®</sup> MIDI

### **Closed Pouch** Filter

- Soft cloth on back (and front of opaque) •
- Not made with natural rubber latex

Transparent 7" (18 cm) / 275 mL		Units: 30
Code	Color	HCPCS*
12344		A4423
12345		A4423
12346		A4423
Opaque 7" (18 cm) / 275 mL		Units: 30
Code	Color	HCPCS*
12354		A4423
10055		4 4 4 2 2
12355		A4423
12356		A4423 A4423



Units: 30

NOTE: Assura® is compatible with SenSura® Click, but NOT with SenSura® Mio Click.

# Ileostomy / Colostomy

# Assura<sup>®</sup> Stoma Cap

### Filter

- Small size is optimal for intimate moments or colostomates who irrigate
- Not made with natural rubber latex

0	-	1	Y.
	Ju	)	

Opaque		Units: 30
Code	Color	HCPCS*
2801		A5055
2802		A5055
2803		A5055

# Urostomy

# Assura<sup>®</sup> Urostomy Micro-Pouch

Not made with natural rubber latex

- · Miniature size prevents twisting
- · Soft outlet can be connected to a drain bag for nighttime use, or a leg bag for patients in a wheelchair



# Transparent 6" (15 cm)

Soft cloth on front

Transparent 6" (15 cm)		Units: 10
Code	Color	HCPCS*
14204		A4433
14205		A4433
14206		A4433

Try a leg bag for discreet day wear — available on page 25-27. For nighttime, use a urostomy night bag found on page 73.

# Assura<sup>®</sup> Urostomy Minicap

- Ultra-absorbent material inside
- Small size is optimal for intimate use
- Not made with natural rubber latex



Opaque 100 mL		Units: 30
Code	Color	HCPCS*
2807		A5055
2808		A5055

<ul> <li>Assura<sup>®</sup> MAXI</li> <li>Urostomy Pouch</li> <li>Multi-chambers minimize sloshing</li> <li>Soft outlet can be connected to urostomy night bag #21365 or any standard night bag</li> <li>Soft cloth on back (and front of opaque)</li> <li>Not made with natural rubber latex</li> </ul>		Ô
Transparent 10 <sup>1</sup> / <sub>2</sub> " (26 <sup>1</sup> / <sub>2</sub> cm) / 488 mL		Units: 10
Code	Color	HCPCS*
14227		A4433
14228		A4433
14229		A4433
Opaque 10 <sup>1</sup> /2" (26 <sup>1</sup> /2 cm) / 488 mL		Units: 10
Code	Color	HCPCS*
14224		A4433
14225		A4433
14226		A4433
Assura <sup>®</sup> MAXI Urostomy Pouch • Soft outlet can be connected to urostomy night bag #21365 or any standard night bag • Soft cloth on back (and front of opaque) • Not made with natural rubber latex		
Transparent 10¾" (27¼ cm) / 375 mL	<b>C</b>	Units: 10
Code	Color	HCPCS*
1757		A4433
1758		A4433
1759		A4433
Opaque 10¾" (27¼ cm) / 375 mL		Units: 10
Code	Color	HCPCS*
1754		A4433
1755		A4433
1756		A4433

NOTE: Assura® is compatible with SenSura® Click, but NOT with SenSura® Mio Click.

# Assura<sup>®</sup> MIDI Urostomy Pouch

- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Soft cloth on back
- Not made with natural rubber latex



Transparent 91/2" (24 cm) / 150 mL		Units: 10
Code	Color	HCPCS*
1751		A4433
1752		A4433

# Urostomy Night Bag and Bag Hanger


See Accessories page 73 for more information.



# SenSura® Mio Click Standard Wear Barrier + Drainable Pouch • Barrier #10522 Flat with belt tabs Pouch #11483 MAXI / EasiClose<sup>™</sup> WIDE Outlet Instructions for use Stoma measuring guide

Not made with natural rubber latex

# Cut-to-fit / Transparent



<ul> <li>+ Drainab</li> <li>Barrier #</li> <li>Pouch # WIDE Out</li> <li>Instruction</li> <li>Stoma mage</li> </ul>	E10031 Flat with belt tabs 11196 MAXI / EasiClose™	Ô	HO CO
Cut-to-fit /	Transparent		Units: 5
Code	Stoma Size	Color	HCPCS*
11009	³⁄8 - 2¹⁄4″ (10-55 mm)		A4414 A4427

<ul> <li>Drainab</li> <li>Barrier #</li> <li>Pouch # WIDE Out</li> <li>Instruction</li> <li>Stoma model</li> </ul>			Ó
Cut-to-fit /	Transparent		Units: 5
Code	Stoma Size	Color	HCPCS*
11010	5/ 111/ 1/15 <b>(</b> 2) mm)		A4373



5/8 - 111/16" (15-43 mm)



A4427

# SenSura® Click Xpro Extended Wear Barrier + Urostomy Pouch • Barrier #10035 Flat with belt tabs Pouch #11856 MAXI / soft outlet Urostomy adapter

- Instructions for use
- Stoma measuring guide
- Not made with natural rubber latex



### Cut-to-fit / Transparent Units: 5 Code HCPCS\* Stoma Size Color A4409 11004 3/8 - 21/4" (10-55 mm) A4433

See page 73 for night drainage bag and hanger, sold separately.

# Assura<sup>®</sup> Standard Wear Barrier + Drainable Pouch

- Barrier #2883 Flat with belt tabs
- Pouch #15976 MAXI / EasiClose™ WIDE Outlet
- Instructions for use
- · Stoma measuring guide
- Not made with natural rubber latex

# Cut-to-fit / Transparent



0000000	,		0
Code	Stoma Size	Color	HCPCS*
11002	¾ - 2¼″ (10-55 mm)		A4414

# Assura<sup>®</sup> Standard Wear Barrier + Drainable Pouch 🖊

- · Barrier #14263 Convex Light with belt tabs
- Pouch #15976 MAXI / EasiClose<sup>™</sup> WIDE Outlet
- Instructions for use
- · Stoma measuring guide
- Not made with natural rubber latex

Cut-to-fit / Transparent			Units: 5
Code	Stoma Size	Color	HCPCS*
11001	5⁄8 - 111⁄16" (15-43 mm)		A4373 A4427





# **2-piece Flex** (Adhesive coupling)



Ultra-flexible system bends with body movements



The low profile of a 1-piece – the convenience of a 2-piece



Floating flange = minimal pressure to apply



Built-in security channel + strong adhesive = security against leakage SenSura®Mio

# The fit is the difference between staying in and working out

Deirdre is a SenSura® Mio user.

# SenSura® Mio Flex

The BodyFit Technology<sup>®</sup> in all SenSura Mio barriers combines an elastic adhesive, that bends and stretches with the body, with an innovative barrier design to provide a secure and comfortable fit.



**2-Piece Flex** 

NOTE: SenSura<sup>®</sup> Mio Flex, SenSura<sup>®</sup> Flex, and Assura<sup>®</sup> AC **are compatible** with one another

# SenSura® Mio Flex Flat Barrier

- Standard wear
- Not made with natural rubber latex

Cut-to-fit			Units: 5
Code	Stoma Size	Color	HCPCS*
10551	³⁄8- 1⁵∕16" (10-33 mm)		A4414
10561	³⁄ଃ- 1⅔" (10-48 mm)		A4414
10571	<sup>3</sup> / <sub>8</sub> - 2 <sup>11</sup> / <sub>16</sub> " (10-68 mm)		A4414
Pre-cut			Units: 5
Code	Stoma Size	Color	HCPCS*
10562	1" (25 mm)		A4414
10566	1½" (28 mm)		A4414
10563	1³∕16″ (30 mm)		A4414
10567	1¼" (32 mm)		A4414
10564	1¾″ (35 mm)		A4414
10568	1½" (38 mm)		A4414
10565	1%16" (40 mm)		A4414
10572	1¾" (45 mm)		A4415
10573	2" (50 mm)		A4415

- SenSura<sup>®</sup> Mio Flex Convex Light Barrier
- Standard wear
- Belt tabs for use with 4237, 4247FlexShell is designed for stability and flexibility
- -0-
- Not made with natural rubber latex

Cut-to-fit			Units: 5
Code	Stoma Size	Color	HCPCS*
16471	<sup>3</sup> / <sub>8</sub> - <sup>3</sup> / <sub>4</sub> " (10-20 mm)		A4373
16481	5⁄8-1³∕16" (15-30 mm)		A4373
16483	5⁄8 - 19∕16″ (15-40 mm)		A4373
16491	<sup>5</sup> /8 - 2 <sup>1</sup> /16" (15-53 mm)	-	A4373
Pre-cut			Units: 5
Code	Stoma Size	Color	HCPCS*
16472	<sup>11</sup> /16" (18 mm)		A4373
16484	<sup>13</sup> / <sub>16</sub> " (21 mm)		A4373
16485	1" (25 mm)		A4373
16486	11/8" (28 mm)		A4373
16487	1¼″ (31 mm)		A4373
16488	1³⁄ଃ" (35 mm)		A4373
16493	11⁄2" (38 mm)	-	A4373
16494	15⁄8″ (41 mm)		A4373



Indicates convexity

NOTE: SenSura® Mio pouches are compatible with SenSura Xpo extended wear barriers

# SenSura<sup>®</sup> Mio Flex MAXI Drainable Pouch

- EasiClose™ WIDE Outlet
- Full-circle filter
- Neutral gray, water-resistant textile



Units: 10

- Opaque pouches feature peek-a-boo inspection
   window for easy application
- Not made with natural rubber latex

### Transparent 11" (28 cm) / 655 mL

	- (	0111001 20
Code	Color	HCPCS*
12273		A4425
12283	-	A4425

Opaque 11" (28 cm) / 655 mL		Units: 10
Code	Color	HCPCS*
12262		A4425
12272		A4425
12282		A4425

<ul> <li>Drainable P</li> <li>EasiClose™</li> <li>Without filt</li> <li>Neutral grade</li> </ul>	WIDE Outlet	
Transparent	11" (28 cm) / 655 mL	Units: 10
Code	Color	HCPCS*
12274		A5063
12284		A5063

## SenSura<sup>®</sup> Mio Flex MIDI Drainable Pouch

- EasiClose<sup>™</sup> WIDE Outlet
- Full-circle filter
- Neutral gray, water-resistant textile
- · Peek-a-boo inspection window for easy application
- Not made with natural rubber latex

Opaque 10	0¼″ (26 cm) / 485 mL	Units: 10
Code	Color	HCPCS*
12281	-	A4425

### SenSura® Mio Flex High Output Pouch with Tap Outlet

- Special full-circle filter designed specifically for ileostomy output
- Transparent pouches feature a printed scale for visual indication of output volume
- Opaque pouches feature a peek-a-boo inspection window and tuck-away option to keep the outlet from rubbing on the wearer's leg
- · Wider-lumen outlet provides effective drainage for patients with thicker output; connects to Coloplast fistula drain bag #14010
- Not made with natural rubber latex

### Transparent 12 1/2" (32 cm) / 1080 mL Units: 10 Code HCPCS\* Color 18655 A4413 18656 A4413 Opaque 12 1/2" (32 cm) / 1080 mL Units: 10 HCPCS\* Code Color Color 18650 A4413 A4413

18651

### SenSura<sup>®</sup> Mio Flex High Output Pouch with Soft Outlet

- Special full-circle filter designed specifically for ileostomy output
- Transparent pouches feature a printed scale for visual indication of output volume
- Opaque pouches feature a peek-a-boo inspection window and tuck-away option to keep the outlet from rubbing on the wearer's leg
- Wider-lumen outlet provides effective drainage for patients with thicker output; connects to Coloplast fistula drain bag #14010
- Not made with natural rubber latex

Transparent 12 ½" (32 cm) / 1080 mL		Units: 10
Code	Color	HCPCS*
18657		A4413
18658		A4413
Opaque 12 ½" (32 cm) / 1080 mL		Units: 10
Code	Color	HCPCS*
18652		A4413
18653		A4413

\*For important information regarding the use of HCPCS codes listed and code

definitions see page 192





NOTE: SenSura® Mio Flex, SenSura® Flex, and Assura® AC are compatible with one another

### SenSura<sup>®</sup> Mio Flex MAXI **Closed Pouch** Full-circle filter

- Neutral gray, water-resistant textile
- Opaque pouches feature peek-a-boo inspection window for easy application
- Not made with natural rubber latex



### Transparent 81/4" (201/2 cm) / 570 mL Units: 30 HCPCS\* Code Color 12203 A4419 12216 A4419 12223 A4419 Opaque 81/4" (201/2 cm) / 570 mL Units: 30 Code HCPCS\* Color 12202 A4419 12213 A4419 12222 A4419

### SenSura<sup>®</sup> Mio Flex MAXI **Urostomy Pouch**

- · Multi-chambers minimize sloshing
- · Anti-reflux valves prevent urine from washing over the stoma



- · Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Neutral gray, water-resistant textile
- Not made with natural rubber latex

Transparent 10 <sup>1</sup> /2" (27 cm) / 400 mL		Units: 10
Code	Color	HCPCS*
12292		A5073
12298		A5073

<mark>Opaque</mark> 10½" (27 cm) / 400 mL		Units: 10
Code	Color	HCPCS*
12291		A5073
12297		A5073

### SenSura® Mio Flex MIDI Urostomy Pouch

- Multi-chambers minimize sloshing
- · Anti-reflux valves prevent urine from washing over the stoma
- · Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- · Neutral gray, water-resistant textile
- Not made with natural rubber latex

### Opaque 91/2" (24 cm) / 330 mL Units: 10 HCPCS\* Code Color 12296 A5073

S

NOTE: SenSura® Mio Flex, SenSura® Flex, and Assura® AC are compatible with one another

### SenSura<sup>®</sup> Flex Flat Barrier

- Standard wear
- Not made with natural rubber latex

Cut-to-fit		Units: 10 (**Units fo	or 10108: 5)
Code	Stoma Size	Color	HCPCS*
10101	<sup>3</sup> ⁄8 - 1 <sup>5</sup> ⁄16" (10-33 mm)		A4414
10102	³⁄ଃ - 1⅔″ (10-48 mm)		A4414
10103	<sup>3</sup> / <sub>8</sub> - 2 <sup>3</sup> / <sub>4</sub> " (10-68 mm)		A4415
10108**	<sup>3</sup> ⁄ <sub>8</sub> - 3 <sup>1</sup> ⁄ <sub>2</sub> " (10-88 mm)		A4415

### SenSura<sup>®</sup> Flex Xpro Flat Barrier

• Extended wear

~

- Compatible with SenSura Mio pouches
- Not made with natural rubber latex

Cut-to-fit			Units: 10
Code	Stoma Size	Color	HCPCS*
10104	³⁄8 - 1⁵⁄16″ (10-33 mm)		A4409
10105	³⁄₃ - 1⅔" (10-48 mm)		A4409
10106	<sup>3</sup> ⁄ <sub>8</sub> - 2 <sup>3</sup> ⁄ <sub>4</sub> " (10-68 mm)		A4410

### SenSura<sup>®</sup> Flex Convex Light Barrier

- Standard wear
- Belt tabs for use with 4215, 4220
- Not made with natural rubber latex

Cut-to-fit			Units: 5
Code	Stoma Size	Color	HCPCS*
11301	5⁄8 - 15⁄16" (15-23 mm)		A4373
11302	5⁄8 - 15⁄16" (15-33 mm)		A4373
11303	5⁄8 - 1¹¹∕₁6″ (15-43 mm)		A4373
11307	<sup>5</sup> / <sub>8</sub> - 2 <sup>1</sup> / <sub>16</sub> " (15-53 mm)		A4373

### SenSura<sup>®</sup> Flex Convex Light Barrier

- Standard wear
- Belt tabs for use with 4215, 4220
- Not made with natural rubber latex

Pre-cut			Units: 5
Code	Stoma Size	Color	HCPCS*
11312	‰" (22 mm)		A4373
11313	1" (25 mm)		A4373
11314	11/8" (28 mm)		A4373
11315	11⁄4" (31 mm)		A4373
ndi	icates convexity		

\*For important information regarding the use of HCPCS codes listed and code definitions see page 192

### SenSura<sup>®</sup> Flex Xpro Convex Light Barrier

- Extended wear
- · Compatible with SenSura Mio pouches
- Belt tabs for use with 4215, 4220
- Not made with natural rubber latex



Cut-to-fit			Units: 5
Code	Stoma Size	Color	HCPCS*
11304	5⁄8 - 15⁄16" (15-23 mm)		A4407
11305	5⁄8 - 15∕16″ (15-33 mm)		A4407
11306	5⁄8 - 111⁄16" (15-43 mm)		A4407
11308	5⁄8 - 2¹∕16″ (15-53 mm)		A4408

SenSura®	Flex MAXI	Drainable Pouch
	TM	

- EasiClose<sup>™</sup> WIDE Outlet
- Filter
- · Soft cloth on back (and front of opaque)
- Not made with natural rubber latex

Units: 20
HCPCS*
A4425
A4425
A4425
Units: 20
HCPCS*
A4425
A4425
A4425
A4425

<ul> <li>EasiClose<sup>™</sup> I</li> <li>Without filte</li> <li>Soft cloth o</li> </ul>	er		9
Transparent 1	.1½" (29½ cm) / 655 mL		Units: 20
Code		Color	HCPCS*
11612			A5063
11690	750 mL XXL pouch		A5063

NOTE: SenSura® Mio Flex, SenSura® Flex, and Assura® AC are compatible with one another

- Filter
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex



### Assura<sup>®</sup> AC Flat Barrier

- Standard wear
- Not made with natural rubber latex

Cut-to-fit		Units: 10 (**Units fo	or 14309: 5)
Code	Stoma Size	Color	HCPCS*
14301	<sup>3</sup> /8 - 1 <sup>5</sup> /16" (10-33 mm)		A4414
14302	<sup>3</sup> /8 - 1 <sup>7</sup> /8" (10-48 mm)		A4414
14303	<sup>3</sup> / <sub>8</sub> - 2 <sup>3</sup> / <sub>4</sub> " (10-68 mm)		A4415
14309**	³⁄8 - 3¹∕2″ (10-88 mm)		A4415

### Assura<sup>®</sup> AC Flat Barrier

- Extended wear
- Compatible with SenSura Mio pouches
- Not made with natural rubber latex

Cut-to-fit			Units: 5
Code	Stoma Size	Color	HCPCS*
14304	³⁄8 - 1⁵⁄16″ (10-33 mm)		A4409
14305	³⁄ଃ - 1⅔″ (10-48 mm)		A4409
14306	¾ - 2¾" (10-68 mm)		A4410

For Assura® AC urostomy pouches see page 118. For XXL closed pouch, see page 116.

Indicates convexity



### Assura<sup>®</sup> AC Convex Light Barrier 🖊

- Standard wear
- Belt tabs for use with 4215, 4220
- Not made with natural rubber latex



Cut-to-fit			Units: 5
Code	Stoma Size	Color	HCPCS*
14401	<sup>5</sup> / <sub>8</sub> - <sup>15</sup> / <sub>16</sub> " (15-23 mm)		A4373
14403	5⁄8 - 111⁄16" (15-43 mm)		A4373
14406	5⁄8 - 2³⁄8" (15-61 mm)		A4373
Pre-cut			Units: 5
Code	Stoma Size	Color	HCPCS*
Code 14642	Stoma Size 7⁄8" (22 mm)	Color	HCPCS* A4373
		Color	
14642	‰″ (22 mm)	Color	A4373
14642 14643	%" (22 mm) 1″ (25 mm)	Color	A4373 A4373

### Assura® AC Convex Light Barrier 🦟

- Extended wear
- Compatible with SenSura<sup>®</sup> Mio pouches
- Belt tabs for use with 4215, 4220
- Not made with natural rubber latex



Cut-to-fit			Units: 5
Code	Stoma Size	Color	HCPCS*
14601	<sup>5</sup> / <sub>8</sub> - <sup>15</sup> / <sub>16</sub> " (15-23 mm)		A4407
14602	5⁄8 - 15∕16″ (15-33 mm)		A4407
14603	<sup>5</sup> / <sub>8</sub> - 1 <sup>11</sup> / <sub>16</sub> " (15-43 mm)		A4407
14408	5⁄8 - 23∕8″ (15-61 mm)		A4408

### Assura<sup>®</sup> AC XXL Drainable Pouch

- EasiClose<sup>™</sup> Outlet
- Filter
- Soft cloth on front and back
- Not made with natural rubber latex

F	-
C	J
	. 1
-	
Units	5: 20

Opaque	121⁄2″	(32	cm) /	750	mL
--------	--------	-----	-------	-----	----

Code	Color	HCPCS*
14379		A4425

Assura® AC XXL Drainable Pouch • EasiClose" Outlet • Without filter • Soft cloth on back • Not made with natural rubber latex		
Transparent 121/2" (32 cm) / 750 mL		Units: 20
Code	Color	HCPCS*
14527		A5063
Assura <sup>®</sup> AC MAXI Drainable Pouch • EasiClose <sup>™</sup> Outlet • Filter • Soft cloth on back • Not made with natural rubber latex		9
Transparent 111/4" (281/2 cm) / 600 mL		Units: 20
Code	Color	HCPCS*
14356		A4425
14357		A4425
14359		A4425
Assura® AC MAXI Drainable Pouch • EasiClose" Outlet • Filter • Soft cloth on front and back • Not made with natural rubber latex		
Opaque 11¼" (28½ cm) / 600 mL		Units: 20
Opaque 11¼" (28½ cm) / 600 mL Code	Color	Units: 20 HCPCS*
	Color	
Code	Color	HCPCS*

### Assura<sup>®</sup> AC MIDI Drainable Pouch

EasiClose<sup>™</sup> Outlet

- Filter
- · Soft cloth on back (and front of opaque)
- Not made with natural rubber latex



Transparent 10¼" (26 cm) / 470 mL		Units: 20
Code	Color	HCPCS*
14342**		A4425
14344		A4425

\*\*Custom Peel. Includes perforated rings providing the option of an exact fit around the stoma. Rings are simply removed for the best fit.

Opaque 10¼" (26 cm) / 470 mL		Units: 20
Code	Color	HCPCS*
14346		A4425
14347		A4425
14349		A4425

Assura® AC MAXI Drainable Pouch • EasiClose" Outlet • Without filter • Soft cloth on back (and front of opaque) • Not made with natural rubber latex		
Transparent 11¼" (28½ cm) / 600 mL		Units: 20
Code	Color	HCPCS*
14521		A5063
Opaque 11¼" (28½ cm) / 600 mL		Units: 20
Code	Color	HCPCS*
14525		A5063

Assura® AC XXL Closed Pouch • Filter • Soft cloth on front and back • Not made with natural rubber latex		<b>O</b>
Opaque 101/2" (26 1/2 cm) / 850 mL		Units: 30
Code	Color	HCPCS*
14355		A4419

### Assura<sup>®</sup> AC MAXI Closed Pouch

Filter

14317 14319

- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex



A4419

A4419

Transparent 81/2" (21 cm) / 500 mL		Units: 30
Code	Color	HCPCS*
14326		A4419
14327		A4419
14328**		A4419
14329		A4419

\*\*Custom Peel. Includes perforated rings providing the option of an exact fit around the stoma. Rings are simply removed for the best fit.

Opaque 81⁄2" (21 cm) / 500 mL		Units: 30
Code	Color	HCPCS*
14331		A4419
14332		A4419
14333**		A4419
14334		A4419

\*\*Custom Peel. Includes perforated rings providing the option of an exact fit around the stoma. Rings are simply removed for the best fit.

Assura® AC MIDI Closed Pouch • Filter • Soft cloth on front and back • Not made with natural rubber latex		Q
Opaque 7" (18 cm) / 250 mL		Units: 30
Code	Color	HCPCS*
14316		A4419

*For important information regarding the use of HCPCS codes listed and code
definitions see page 192

### Urostomy

Assura® AC MAXI Urostomy Pouch • Multi-chambers minimize sloshing • Soft outlet can be connected to urostom night bag #21365 or any standard night • Soft cloth on back (and front of opaque) • Not made with natural rubber latex Transparent 10 <sup>1</sup> /2" (26 <sup>1</sup> /2 cm) / 475 mL	,	Units: 20
Code	Color	HCPCS*
	COIOI	
14552		A5073
14554		A5073
Opaque 101⁄2" (261⁄2 cm) / 475 mL		Units: 20
Code	Color	HCPCS*
14553		A5073
14555		A5073

# **Urostomy Night Bag and Bag Hanger** See Accessories page 73 for more information.



# Single Use Kits

### SenSura<sup>®</sup> Mio Flex Standard Wear Barrier + MAXI Drainable Pouch

- Instructions for use
- Stoma measuring guide
- Not made with natural rubber latex
- Barrier #10561 Flat
- Pouch #12273 EasiClose<sup>™</sup> WIDE Outlet / filter

Cut-to-fit / Transparent		Units: 5	
Code	Stoma Size	Color	HCPCS*
21009	³⁄ଃ - 1⅔" (10-48 mm)		A4414 A4425

- Barrier #10571 Flat
- Pouch #12283 EasiClose<sup>™</sup> WIDE Outlet / filter

Cut-to-fit / Transparent			Units: 5
Code	Stoma Size	Color	HCPCS*
21010	∛8 - 2¾″ (10-68 mm)	-	A4414 A4425

- Barrier #10561 Flat
- Pouch #12274 EasiClose<sup>™</sup> WIDE Outlet / without filter

Cut-to-fit / Transparent			Units: 5
Code	Stoma Size	Color	HCPCS*
21012	¾ - 1%″ (10-48 mm)		A4414 A5063

- Barrier #10571 Flat
- Pouch #12284 EasiClose™ WIDE Outlet / without filter

Cut-to-fit / Transparent			Units: 5
Code	Stoma Size	Color	HCPCS*
21013	³⁄ଃ - 2³⁄₄″ (10-68 mm)	-	A4414 A5063

<ul> <li>HAXI Dra</li> <li>Instruction</li> <li>Stoma me</li> </ul>	<b>Elex Standard Wear Barrier</b> <b>inable Pouch</b> is for use asuring guide <b>with natural rubber latex</b>	C	
<ul><li>Barrier #1</li><li>Pouch #11</li></ul>	0102 Flat 1.512 EasiClose <sup>™</sup> WIDE Outlet / filter		- 🕑
Cut-to-fit / Tr	ansparent		Units: 5
Code	Stoma Size	Color	HCPCS*
11011K	∛8 - 1%" (10-48 mm)		A4414 A4425

\*For important information regarding the use of HCPCS codes listed and code definitions see page 192

<ul> <li>MAXI I</li> <li>Instruct</li> <li>Stoma</li> <li>Not mo</li> <li>Barrier</li> </ul>	Flex Standard Wear Barrier Drainable Pouch ions for use measuring guide ide with natural rubber latex #10103 Flat #11517 EasiClose <sup>®</sup> WIDE Outlet / filter	C	20
	/ Transparent		Units: 5
Code	Stoma Size	Color	HCPCS*
21001	∛8 - 2¾″ (10-68 mm)		A4415

Barrier #10108 Flat

Pouch #11690 EasiClose<sup>™</sup> WIDE Outlet / without filter

Cut-to-fit	/ Transparent		Units: 5
Code	Stoma Size	Color	HCPCS*
21011	³⁄₃ - 3½″ (10-88 mm)		A4415 A5063

### SenSura<sup>®</sup> Flex Xpro (Extended Wear) Barrier

### + MAXI Drainable Pouch

- Instructions for use
- Stoma measuring guide
- Not made with natural rubber latex
- Barrier #10105 Flat
- Pouch #11512 EasiClose<sup>™</sup> WIDE Outlet / filter

### Cut-to-fit / Transparent

Units: 5 Code Stoma Size Color HCPCS\* A4409 21002 3/8 - 17/8" (10-48 mm) A4425

### • Barrier #10106 Flat

Pouch #11517 EasiClose<sup>™</sup> WIDE Outlet / filter

Cut-to-fit / Transparent			Units: 5
Code	Stoma Size	Color	HCPCS*
21003	∛8 - 2¾″ (10-68 mm)	-	A4410 A4425

A4425

<ul> <li>Instruct</li> <li>Stoma</li> <li>Not mo</li> <li>Barrier</li> </ul>	Drainable Pouch tions for use measuring guide ade with natural rubber latex #14302 Flat #14521 EasiClose <sup>®</sup> Outlet / with	Que filter	0
	/ Transparent		Units: 5
Code	Stoma Size	Color	HCPCS*
11006	¾ - 1%″ (10-48 mm)		A4414 A5063

### Assura<sup>®</sup> AC Standard Wear Barrier

Assura<sup>®</sup> AC Standard Wear Barrier

### + MAXI Drainable Pouch

- Instructions for use
- Stoma measuring guide
- Not made with natural rubber latex
- Barrier #14303 Flat
- Pouch #14344 EasiClose<sup>™</sup> Outlet / filter

Cut-to-fit / Transparent			Units: 5
Code	Stoma Size	Color	HCPCS*
11003	³⁄8 - 2³⁄4″ (10-68 mm)	-	A4415 A4425

<ul> <li>MAXI E</li> <li>Instruct</li> <li>Stoman</li> <li>Not ma</li> <li>Barrier</li> </ul>	Flex Xpro (Extended Wear) Ba Drainable Pouch ions for use measuring guide de with natural rubber latex #11308 Convex Light #11517 EasiClose <sup>®</sup> WIDE Outlet / filte		2
Cut-to-fit /	' Transparent		Units: 5
Code	Stoma Size	Color	HCPCS*
21004	5⁄8 - 2¹∕16″ (15-53 mm)		A4410

A4425

### SenSura<sup>®</sup> Mio Flex Standard Wear + MAXI Urostomy Pouch

- Urostomy adapter
- Instructions for use
- Stoma measuring guide
- Not made with natural rubber latex
- Barrier #10561 Flat
- Pouch #12298 Soft outlet / multi-chamber

Cut-to-fit / Transparent			Units: 5
Code	Stoma Size	Color	HCPCS*
110072	³⁄ଃ - 1⅛" (10-48 mm)		A4414 A5073

S

### Urostomy Night Bag and Bag Hanger

See Accessories page 73 for more information.







The most flexible choice: bends with body movements



No coupling = pouch and barrier are welded together and apply as one piece L-piece



Low-profile

# The fit is the difference between asking for help and doing it yourself

Walton is a SenSura® Mio Convex Flip user.

SenSura®Mic

# SenSura<sup>®</sup> Mio 1-piece Pouching System

The BodyFit Technology<sup>®</sup> in all SenSura Mio barriers combines an elastic adhesive, that bends and stretches with the body, with an innovative barrier design to provide a secure and comfortable fit.



### SenSura<sup>®</sup> Mio Flat MAXI **Drainable Pouch**

- Standard wear
- EasiClose<sup>™</sup> WIDE Outlet
- Full-circle filter
- Neutral gray, water-resistant textile
- Opaque pouches feature peek-a-boo inspection window for easy application • Not made with natural rubber latex
- Transparent 11" (28 cm) / 655 mL



	( ))		
Code	Cut	Stoma Size	HCPCS*
10481	Cut-to-fit	<sup>3</sup> / <sub>8</sub> -2 <sup>1</sup> / <sub>8</sub> " (10-55 mm)	A4424
10482	Pre-cut	1" (25 mm)	A4424
10488	Pre-cut	1 <sup>1</sup> /8" (28 mm)	A4424
10483	Pre-cut	1 <sup>3</sup> / <sub>16</sub> " (30 mm)	A4424
10490	Pre-cut	1 <sup>1</sup> ⁄4" (32 mm)	A4424
10484	Pre-cut	1³⁄₃" (35 mm)	A4424

10484	Pre-cut	1³⁄8″ (35 mm)	A4424
10491	Pre-cut	1 <sup>1</sup> ⁄2" (38 mm)	A4424
10487	Pre-cut	1 <sup>9</sup> ⁄16″ (40 mm)	A4424
10485	Pre-cut	1¾″ (45 mm)	A4424
10486	Pre-cut	2" (50 mm)	A4424

### Opaque 11" (28 cm) / 655 mL

Opaque 11" (28 cm) / 655 mL			Units: 10
Code	Cut	Stoma Size	HCPCS*
10471	Cut-to-fit	<sup>3</sup> / <sub>8</sub> -2 <sup>1</sup> / <sub>8</sub> " (10-55 mm)	A4424
10472	Pre-cut	1" (25 mm)	A4424
10473	Pre-cut	1 <sup>1</sup> /8" (28 mm)	A4424
10474	Pre-cut	1 <sup>3</sup> ⁄16″ (30 mm)	A4424
10475	Pre-cut	1 <sup>1</sup> ⁄4" (32 mm)	A4424
10476	Pre-cut	1¾" (35 mm)	A4424
10480	Pre-cut	1½" (38 mm)	A4424
10479	Pre-cut	1 <sup>9</sup> /16" (40 mm)	A4424
10477	Pre-cut	1¾" (45 mm)	A4424
10478	Pre-cut	2" (50 mm)	A4424

### SenSura<sup>®</sup> Mio Flat MIDI Drainable Pouch

- Standard wear
- EasiClose<sup>™</sup> WIDE Outlet
- Full-circle filter
- Peek-a-boo inspection window allows easy application
- Neutral gray, water-resistant textile
- Not made with natural rubber latex



Units: 10

### <mark>Opaque</mark> 10<sup>1</sup>⁄4" (26 cm) / 485 mL

Code	Cut	Stoma Size	HCPCS*
10461	Cut-to-fit	<sup>3</sup> / <sub>8</sub> - 1 <sup>3</sup> / <sub>4</sub> " (10-45 mm)	A4424

### SenSura® Mio Flat MAXI Drainable Pouch

- Standard wear
- EasiClose<sup>™</sup> WIDE Outlet
- Without filter
- Neutral gray, water-resistant textile
- Not made with natural rubber latex

### Transparent 11" (28 cm) / 655 mL



U	nı	ts:	10	)
-				-

Code	Cut	Stoma Size	HCPCS*
10489	Cut-to-fit	<sup>3</sup> / <sub>8</sub> - 2 <sup>1</sup> / <sub>8</sub> " (10-55 mm)	A5061
10492	Pre-cut	1" (25 mm)	A5061
10493	Pre-cut	1 <sup>1</sup> /8" (28 mm)	A5061
10494	Pre-cut	1¼" (32 mm)	A5061
10495	Pre-cut	1³⁄/₃" (35 mm)	A5061
10496	Pre-cut	1 <sup>1</sup> ⁄2" (38 mm)	A5061

### SenSura<sup>®</sup> Mio Soft Convex MAXI Drainable Pouch

- Standard wear
- Standard wear
- Belt tabs for use with 4237, 4247
- EasiClose<sup>™</sup> WIDE Outlet
- FlexShell is designed for stability and flexibility
- Full-circle filter
- Neutral gray, water-resistant textile
- Opaque pouches feature peek-a-boo-inspection
   window for easy application
- Not made with natural rubber latex

### Transparent 11" (28 cm) / 655 mL

Code	Cut	Stoma Size	HCPCS*
16715	Cut-to-fit	<sup>5</sup> ⁄8-1⁵∕16″ (15-33 mm)	A4389
16716	Cut-to-fit	<sup>3</sup> / <sub>8</sub> - 2" (10-50 mm)	A4389

### Opaque 11" (28 cm) / 655 mL

Code	Cut	Stoma Size	HCPCS*
16705	Cut-to-fit	<sup>5</sup> ∕8 - 1 <sup>5</sup> ∕16″ (15-33 mm)	A4389
16706	Cut-to-fit	<sup>3</sup> / <sub>8</sub> - 2" (10-50 mm)	A4389



Units: 10

### SenSura® Mio Convex Light MAXI Drainable Pouch

- Standard wear
- Belt tabs for use with 4237, 4247
- EasiClose™ WIDE Outlet
- FlexShell is designed for stability and flexibility
- Full-circle filter
- Neutral gray, water-resistant textile
- Opaque pouches feature peek-a-boo-inspection window for easy application
- Not made with natural rubber latex

### Transparent 11" (28 cm) / 655 mL

Code	Cut	Stoma Size	HCPCS*
16737	Cut-to-fit	<sup>3</sup> / <sub>8</sub> - <sup>7</sup> / <sub>8</sub> " (10-23 mm)	A4389
16738	Cut-to-fit	<sup>5</sup> / <sub>8</sub> - 1 <sup>5</sup> / <sub>16</sub> " (15-33 mm)	A4389
16739	Cut-to-fit	<sup>3</sup> /8 - 1 <sup>11</sup> /16" (10-43 mm)	A4389
16740	Pre-cut	<sup>13</sup> / <sub>16</sub> " (21 mm)	A4389
16741	Pre-cut	1" (25 mm)	A4389
16742	Pre-cut	1 <sup>1</sup> /8" (28 mm)	A4389

### Opaque 11" (28 cm) / 655 mL

	. ,,		
Code	Cut	Stoma Size	HCPCS*
16725	Cut-to-fit	<sup>3</sup> / <sub>8</sub> - <sup>7</sup> / <sub>8</sub> " (10-23 mm)	A4389
16726	Cut-to-fit	⁵⁄₃ - 1⁵⁄ı₅" (15-33 mm)	A4389
16733	Cut-to-fit	<sup>3</sup> / <sub>8</sub> - 1 <sup>11</sup> / <sub>16</sub> " (10-43mm)	A4389
16730	Pre-cut	<sup>13</sup> /16" (21 mm)	A4389
16731	Pre-cut	1" (25 mm)	A4389
16732	Pre-cut	1 <sup>1</sup> /8" (28 mm)	A4389
16727	Pre-cut	1 <sup>1</sup> ⁄4" (31 mm)	A4389
16734	Pre-cut	1³⁄8″ (35 mm)	A4389
16735	Pre-cut	1 <sup>1</sup> ⁄2" (38 mm)	A4389

### SenSura® Mio Convex Light MAXI Drainable Pouch

- Standard wear
- Belt tabs for use with 4237, 4247
- EasiClose<sup>™</sup> WIDE Outlet
- FlexShell is designed for stability and flexibility
- Without filter
- Neutral gray, water-resistant textile
- Opaque pouches feature peek-a-boo-inspection
   window for easy application

### Not made with natural rubber latex

Transparent 11" (28 cm) / 655 mL		Units: 10	
Code	Cut	ut Stoma Size	
16748	Cut-to-fit	<sup>3</sup> /8 - 1 <sup>11</sup> /16" (10-43 mm)	A4389



Units: 10

Units: 10

Indicates convexity

### SenSura<sup>®</sup> Mio Deep Convex MAXI Drainable Pouch

- Standard wear
- Belt tabs for use with 4237, 4247
- EasiClose<sup>™</sup> WIDE Outlet
- FlexShell is designed for stability and flexibility
- Full-circle filter
- Neutral gray, water-resistant textile
- Opaque pouches feature peek-a-boo-inspection window for easy application
- Not made with natural rubber latex

### Transparent 11" (28cm) / 655 ml

Transparent 11" (28cm) / 655 mL			Units: 10
Code	Cut	Stoma Size	HCPCS*
16765	Cut-to-fit	<sup>3</sup> / <sub>8</sub> - <sup>7</sup> / <sub>8</sub> " (10-23 mm)	A4389
16766	Cut-to-fit	5⁄8-15∕16″ (15-33 mm)	A4389
16767	Cut-to-fit	<sup>3</sup> /8 - 1 <sup>11</sup> /16" (10-43 mm)	A4389
16770	Pre-cut	<sup>13</sup> / <sub>16</sub> " (21 mm)	A4389
16771	Pre-cut	1" (25 mm)	A4389
16772	Pre-cut	1¹⁄ଃ" (28 mm)	A4389
16773	Pre-cut	1 <sup>1</sup> ⁄4" (31 mm)	A4389

### Opaque 11" (28 cm ) / 655 mL

Opaque 11" (28 cm ) / 655 mL			Units: 10
Code	Cut	Stoma Size	HCPCS*
16755	Cut-to-fit	<sup>3</sup> / <sub>8</sub> - <sup>7</sup> / <sub>8</sub> " (10-23 mm)	A4389
16756	Cut-to-fit	<sup>5</sup> /8 - 1 <sup>5</sup> /16" (15-33 mm)	A4389
16768	Cut-to-fit	<sup>3</sup> / <sub>8</sub> - 1 <sup>11</sup> / <sub>16</sub> " (10-43 mm)	A4389
16760	Pre-cut	<sup>13</sup> ⁄16" (21 mm)	A4389
16761	Pre-cut	1" (25 mm)	A4389
16762	Pre-cut	1 <sup>1</sup> /8" (28 mm)	A4389
16757	Pre-cut	1¼" (31 mm)	A4389
16763	Pre-cut	1 <sup>3</sup> ⁄8" (35 mm)	A4389
16764	Pre-cut	1 <sup>1</sup> /2" (38 mm)	A4389

### SenSura® Mio Deep Convex MAXI Drainable Pouch

- Standard wear
- Belt tabs for use with 4237, 4247
- EasiClose<sup>™</sup> WIDE Outlet
- · FlexShell is designed for stability and flexibility
- Without filter
- · Neutral gray, water-resistant textile
- Opaque pouches feature peek-a-boo-inspection window for easy application
- Not made with natural rubber latex

### Transparent 11" (28 cm) / 655 ml

			011103. ±0
Code	Cut	Stoma Size	HCPCS*
16778	Cut-to-fit	<sup>3</sup> / <sub>8</sub> - 1 <sup>11</sup> / <sub>16</sub> " (10-43 mm)	A4389



Linits: 10



# L-Piece SenSura® Mio

# SenSura® Mio MAXI Drainable Pouch with Soft Outlet

- Standard wear
- Special full-circle filter designed specifically for ileostomy output
- Transparent pouches feature a printed scale for visual indication of output volume
- Opaque pouches feature a peek-a-boo inspection window and tuck-away option to keep the outlet from rubbing on the wearer's leg
- Soft outlet can connect to Coloplast's fistula drain bag #14010, or most standard night bags
- Not made with natural rubber latex

Transparent 11" (28 cm) / 630 mL			Units: 10
Code	Cut	Stoma Size	HCPCS*
18665	Cut-to-fit	<sup>3</sup> /8 - 2 <sup>3</sup> ⁄4" (10-70 mm)	A4424
18666	Cut-to-fit	<sup>3</sup> / <sub>8</sub> -4" (10-100 mm)	A4424
Ondaue		and a second s	United 10
opuque -	11" (28 cm) / 630	/ mL	Units: 10
Code	Cut	Stoma Size	HCPCS*

### SenSura® Mio Soft Convex MAXI Drainable Pouch with Soft Outlet 🗡

- Standard wear
- Belt tabs for use with 4237, 4247
- Special full-circle filter designed specifically for ileostomy output
- Transparent pouches feature a printed scale for visual indication of output volume
- Opaque pouches feature a peek-a-boo inspection window and tuckaway option to keep the outlet from rubbing on the wearer's leg
- Soft outlet can connect to Coloplast's fistula drain bag #14010, or any standard night bag
- Not made with natural rubber latex

Transparent 11" (28 cm) / 630 mL			Units: 10
Code	Cut	Stoma Size	HCPCS*
18671	Cut-to-fit	<sup>3</sup> /8 - 2" (10-50 mm)	A4389
Opggue 1			
Opuque -	L1" (28 cm) / 630	mL	Units: 10
Code	Cut (28 cm) / 630	ML Stoma Size	HCPCS*



Indicates convexity

### SenSura® Mio Convex Light MAXI Drainable Pouch with Soft Outlet 🖊

- Standard wear
- Belt tabs for use with 4237, 4247
- Special full-circle filter designed specifically for ileostomy output
- Transparent pouches feature a printed scale for visual indication of output volume
- Opaque pouches feature a peek-a-boo inspection window and tuck-away option to keep the outlet from rubbing on the wearer's leg
- Soft outlet can connect to Coloplast's fistula drain bag #14010, or any standard night bag
- Not made with natural rubber latex

Transpare	Units: 10		
Code	Cut	Stoma Size	HCPCS*
18676	Cut-to-fit	<sup>3</sup> / <sub>8</sub> - 1 <sup>11</sup> / <sub>16</sub> " (10-43 mm)	A4389
Opaque 1	L1" (28 cm) / 630	mL	Units: 10
Opaque 1 Code	L1" (28 cm) / 630 Cut	mL Stoma Size	Units: 10 HCPCS*





hdicates Convex Flip

### SenSura® Mio Convex Flip MAXI Drainable Pouch 🎓

- Designed for curves, bulges and hernias
- Curved, star-shaped barrier conforms to outward areas with less creasing and folding
- "Flip to fit" design allows for easy application
- Opaque pouches feature a peek-a-boo inspection window for easy application
- Full-circle filter
- Not made with natural rubber latex

Transparent 11" (28 cm) / 655 mL

Code	Cut	Stoma Size	Max Cut	Barrier Size	HCPCS*
18305	Cut-to-fit	<sup>3</sup> / <sub>8</sub> - 1 <sup>9</sup> / <sub>16</sub> " (10-40 mm)	1º/16" (40 mm)	S	A4389
18321	Cut-to-fit	<sup>3</sup> ⁄8 - 2" (10-50 mm)	2" (50 mm)	M	A4389
18322	Pre-Cut	1" (25 mm)	2" (50 mm)		A4389
18323	Pre-Cut	1¾16" (30 mm)	2" (50 mm)	M	A4389
18324	Pre-Cut	1¾" (35 mm)	2" (50 mm)		A4389
18341	Cut-to-fit	<sup>3</sup> ⁄8 - 2 <sup>3</sup> ⁄8" (10-60 mm)	2¾" (60 mm)		A4389
18342	Pre-Cut	1 <sup>9</sup> ⁄16" (40 mm)	2¾" (60 mm)		A4389
18343	Pre-Cut	1¾" (45 mm)	2¾" (60 mm)		A4389
18344	Pre-Cut	2" (50 mm)	2 <sup>3</sup> ⁄8″ (60 mm)		A4389

Орадие 11" (28 ст) / 655 mL Ur					
Code Cut	Stoma Size	Max Cut	Barrier Size	HCPCS*	
18301 Pre-Cut	1" (25 mm)	1 <sup>9</sup> / <sub>16</sub> " (40 mm)	S	A4389	
18302 Pre-Cut	1³∕₁₀" (30 mm)	1 <sup>9</sup> /16" (40 mm)	S	A4389	
<b>18312</b> Cut-to-fit	<sup>3</sup> ⁄8- 2" (10-50 mm)	2" (50 mm)		A4389	
18313 Pre-cut	1" (25 mm)	2" (50 mm)	M	A4389	
18314 Pre-cut	1³⁄16" (30 mm)	2" (50 mm)	M	A4389	
18315 Pre-cut	1¾" (35 mm)	2" (50 mm)	M	A4389	
<b>18330</b> Cut-to-fit	<sup>3</sup> ⁄8- 2 <sup>3</sup> ⁄8" (10-60 mm)	2 <sup>3</sup> ⁄8" (60 mm)		A4389	
18331 Pre-Cut	1 <sup>9</sup> /16" (40 mm)	2³⁄8″ (60 mm)		A4389	
18332 Pre-Cut	1¾" (45 mm)	2 <sup>3</sup> /8" (60 mm)		A4389	
18333 Pre-Cut	2" (50 mm)	2 <sup>3</sup> ⁄8″ (60 mm)		A4389	





### SenSura<sup>®</sup> Mio Convex Flip MAXI Drainable Pouch 🖠

### Without filter

- Designed for curves, bulges and hernias
- · Curved, star-shaped barrier conforms to outward areas with less creasing and folding
- · "Flip to fit" design allows for easy application
- Not made with natural rubber latex

### sparont 11" (28 cm) / 655 ml

Transparent 11"		Units: 10		
Code Cut	Stoma Size	Max Cut	Barrier	HCPCS*
			Size	
<b>18320</b> Cut-to-fit	<sup>3</sup> ⁄8- 2" (10-50 mm)	2" (50 mm)	M	A4389
18340 Cut-to-fit	<sup>3</sup> / <sub>8</sub> - 2 <sup>3</sup> / <sub>8</sub> " (10–60 mm)	2 <sup>3</sup> ⁄8″ (60 mm)		A4389

### SenSura<sup>®</sup> Mio Flat MAXI **Closed Pouch**

- Standard wear
- Full-circle filter
- Neutral gray, water-resistant textile
- · Peek-a-boo inspection window allows easy application
- Not made with natural rubber latex

### Opaque 8<sup>1</sup>/<sub>4</sub>" (21 cm) / 550 mL

<mark>Οραque</mark> 8 <sup>1</sup> ⁄4" (21 cm) / 550 mL			Units: 30
Code	Cut	Stoma Size	HCPCS*
10881	Cut-to-fit	<sup>5</sup> / <sub>8</sub> - 1 <sup>3</sup> / <sub>4</sub> " (15-45 mm)	A4416
10885***	Cut-to-fit	1 <sup>9</sup> / <sub>16</sub> - 2 <sup>1</sup> /8" (40-55 mm)	A4416
10883	Pre-cut	1³⁄16″ (30 mm)	A4416
10884	Pre-cut	1³⁄₃" (35 mm)	A4416
10885***	Pre-cut	1 <sup>9</sup> ⁄16" (40 mm)	A4416

\*\*\*Because of its larger starter hole, 10885 functions as both a pre-cut and cut-to-fit solution.

### SenSura<sup>®</sup> Mio Flat MIDI **Closed Pouch**

- Standard wear
- Full-circle filter
- Neutral gray, water-resistant textile
- Peek-a-boo inspection window allows easy application
- Not made with natural rubber latex

Opaque 7¼" (	Units: 30		
Code	Cut	Stoma Size	HCPCS*
10871	Cut-to-fit	<sup>5</sup> / <sub>8</sub> - 1 <sup>3</sup> / <sub>4</sub> " (15-45 mm)	A4416
10872	Pre-cut	1" (25 mm)	A4416
10873	Pre-cut	1³⁄16″ (30 mm)	A4416
10874	Pre-cut	1³⁄8" (35 mm)	A4416







### SenSura® Mio Standard Wear MINI Closed Pouch

- Full-circle filter
- · Neutral gray, water-resistant textile
- Peek-a-boo inspection window allows easy application
- · Comfortable and discreet stoma cap alternative
- Not made with natural rubber latex

### Opaque 5<sup>1</sup>/<sub>4</sub>" (13 cm) / 110 mL

Code	Cut	Stoma Size	HCPCS*
10870	Cut-to-fit	<sup>5</sup> /8 - 1 <sup>3</sup> /4" (15-45 mm)	A4416

### SenSura® Mio Soft Convex MAXI Closed Pouch

- Standard wear
- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- Full-circle filter
- · Neutral gray, water-resistant textile
- Opaque pouches feature peek-a-boo-inspection
   window for easy application

### Not made with natural rubber latex

Transpare	ent 8 <sup>1</sup> /4" (21 cm)	/ 550 mL	Units: 10
Code	Cut	Stoma Size	HCPCS*
16315	Cut-to-fit	<sup>3</sup> /8-2" (10-50 mm)	A4417

Opaque 8	Units: 10		
Code	Cut	Stoma Size	HCPCS*
16305	Cut-to-fit	<sup>5</sup> ∕8 - 1 <sup>5</sup> ∕16″ (15-33 mm)	A4417
16306	Cut-to-fit	<sup>3</sup> / <sub>8</sub> - 2" (10-50 mm)	A4417
16310	Pre-cut	1" (25 mm)	A4417
16311	Pre-cut	1 <sup>1</sup> /8" (28 mm)	A4417
16312	Pre-cut	1 <sup>1</sup> ⁄4" (31 mm)	A4417
16313	Pre-cut	1 <sup>3</sup> ⁄8" (35 mm)	A4417

### SenSura® Mio Soft Convex MIDI Closed Pouch

- Standard wear
- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- Full-circle filter
- · Neutral gray, water-resistant textile
- Opaque pouches feature peek-a-boo-inspection window for easy application
- Not made with natural rubber latex

Indicates convexity

Opaque	7 <sup>1</sup> ⁄4" (19 cm) / 34	0 mL	Units: 10
Code	Cut	Stoma Size	HCPCS*
16301	Cut-to-fit	<sup>5</sup> /8-1 <sup>5</sup> /16" (15-33 mm)	A4417





### SenSura<sup>®</sup> Mio Convex Light MAXI Closed Pouch

- Standard wear
- Belt tabs for use with 4237, 4247
- · FlexShell is designed for stability and flexibility
- Full-circle filter
- Neutral gray, water-resistant textile
- Opaque pouches feature peek-a-boo-inspection window for easy application
- Not made with natural rubber latex

	Transparent	<b>8</b> <sup>1</sup> /4″	(21	cm) /	550	mL
--	-------------	---------------------------	-----	-------	-----	----

	()		
Code	Cut	Stoma Size	HCPCS*
16335	Cut-to-fit	<sup>3</sup> / <sub>8</sub> - 1 <sup>11</sup> / <sub>16</sub> " (10-43 mm)	A4417

### Opaque 8<sup>1</sup>/4" (21 cm) / 550 mL

Opaque 8	Units: 10		
Code	Cut	Stoma Size	HCPCS*
16325	Cut-to-fit	<sup>3</sup> / <sub>8</sub> - <sup>7</sup> / <sub>8</sub> " (10-23 mm)	A4417
16326	Cut-to-fit	⁵⁄ଃ-1⁵∕ı₀" (15-33 mm)	A4417
16330	Pre-cut	1" (25 mm)	A4417
16331	Pre-cut	1 <sup>1</sup> ⁄/8" (28 mm)	A4417
16327	Pre-cut	1 <sup>1</sup> ⁄4" (31 mm)	A4417
16333	Pre-cut	1 <sup>3</sup> /8" (35 mm)	A4417
16334	Pre-cut	1½" (38 mm)	A4417

### SenSura® Mio Convex Light MIDI Closed Pouch

- Standard wear
- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- Full-circle filter
- Neutral gray, water-resistant textile
- Opaque pouches feature peek-a-boo-inspection window for easy application
- Not made with natural rubber latex

### Opaque 7<sup>1</sup>/4" (19 cm) / 340 mL

a provide a	( ) /		
Code	Cut	Stoma Size	HCPCS*
16321	Cut-to-fit	<sup>3</sup> /8 - <sup>7</sup> /8" (10-23 mm)	A4417
16322	Cut-to-fit	<sup>5</sup> /8 - 1 <sup>5</sup> /16" (15-33 mm)	A4417
16323**	Cut-to-fit / Pre-cut	1 <sup>1</sup> / <sub>4</sub> -1 <sup>11</sup> / <sub>16</sub> " (31-43 mm)	A4417

\*\*Because of its larger starter hole, 16323 functions as both a pre-cut and cut-to-fit solution.



Units: 10

### Transparent 8<sup>1</sup>/4" (21 cm) / 550 mL

window for easy applicationNot made with natural rubber latex

Full-circle filter

Code	Cut	Stoma Size	HCPCS*
16355	Cut-to-fit	<sup>3</sup> / <sub>8</sub> - 1 <sup>11</sup> / <sub>16</sub> " (10-43 mm)	A4417

### Opaque 8<sup>1</sup>/4" (21 cm) / 550 mL

Code	Cut	Stoma Size	HCPCS*
16345	Cut-to-fit	<sup>3</sup> / <sub>8</sub> - <sup>7</sup> / <sub>8</sub> " (10-23 mm)	A4417
16346	Cut-to-fit	⁵⁄8-1⁵∕16" (15-33 mm)	A4417
16350	Pre-cut	1" (25 mm)	A4417
16351	Pre-cut	1 <sup>1</sup> /s" (28 mm)	A4417
16347	Pre-cut	1 <sup>1</sup> ⁄4" (31 mm)	A4417
16353	Pre-cut	1³⁄8″ (35 mm)	A4417
16354	Pre-cut	1 <sup>1</sup> ⁄2" (38 mm)	A4417

### SenSura® Mio Deep Convex MIDI

- Closed Pouch
- Standard wear
- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- Full-circle filter
- Neutral gray, water-resistant textile
- Opaque pouches feature peek-a-boo-inspection
   window for easy application
- Not made with natural rubber latex

### Opaque 7<sup>1</sup>/<sub>4</sub>" (19 cm) / 340 mL

Code	Cut	Stoma Size	HCPCS*
16341	Cut-to-fit	<sup>3</sup> ⁄8 - <sup>7</sup> ⁄8" (10-23 mm)	A4417
16342	Cut-to-fit	<sup>5</sup> /8 - 1 <sup>5</sup> /16" (15-33 mm)	A4417
16343**	Cut-to-fit / Pre-cut	1 <sup>1</sup> / <sub>4</sub> - 1 <sup>11</sup> / <sub>16</sub> " (31-43 mm)	A4417

\*\*Because of its larger starter hole, 16343 functions as both a pre-cut and cut-to-fit solution.

Indicates convexity

Belt tabs for use with 4237, 4247FlexShell is designed for stability and flexibility

Neutral gray, water-resistant textile

Opaque pouches feature peek-a-boo-inspection



Units: 10

Units: 10

.-Piece SenSura® Mi



### SenSura® Mio Convex Flip MAXI Closed Pouch

- Designed for curves, bulges and hernias
- Curved, star-shaped barrier conforms to outward areas with less creasing and folding
- "Flip to fit" design allows for easy application
- Peek-a-boo inspection window for easy application
- Not made with natural rubber latex



Opaque 8 ¼" (21	Opaque 8 ¼" (21 cm) / 550 mL Units: 10				
Code Cut	Stoma Size	Max Cut	Barrier Size	HCPCS*	
18152 Cut-to-Fit	<sup>3</sup> / <sub>8</sub> - 1 <sup>9</sup> / <sub>16</sub> " (10-40 mm)	1 <sup>9</sup> / <sub>16</sub> " (40 mm)	S	A4417	
18170 Cut-to-Fit	<sup>3</sup> ⁄8 - 2" (10-50 mm)	2" (50 mm)	M	A4417	
<b>18171</b> Pre-Cut	1" (25 mm)	2" (50 mm)		A4417	
18172 Pre-Cut	1¾16" (30 mm)	2" (50 mm)	M	A4417	
18173 Pre-Cut	1¾" (35 mm)	2" (50 mm)		A4417	
18180 Cut-to-Fit	<sup>3</sup> / <sub>8</sub> - 2 <sup>3</sup> / <sub>8</sub> " (10-60 mm)	2¾" (60 mm)		A4417	
<b>18181</b> Pre-Cut	1 <sup>9</sup> /16" (40 mm)	2 <sup>3</sup> /8" (60 mm)		A4417	
18182 Pre-Cut	1¾" (45 mm)	2¾" (60 mm)		A4417	

### SenSura<sup>®</sup> Mio Original Flat MAXI Closed Pouch

- Standard wear
- Filter
- Soft cloth on front and back
- Peek-a-boo inspection window allows easy application
- Not made with natural rubber latex


### Opaque 8<sup>1</sup>/2" (21 cm) / 550 mL



Units: 30

	( ))		
Code	Cut	Stoma Size	HCPCS*
10831	Cut-to-fit	⁵⁄8-1³⁄4″ (15-45 mm)	A4416
10835**	Cut-to-fit	1 <sup>9</sup> / <sub>16</sub> -2 <sup>1</sup> / <sub>4</sub> " (40-55 mm)	A4416
10833	Pre-cut	1 <sup>3</sup> / <sub>16</sub> " (30 mm)	A4416
10834	Pre-cut	1³⁄8" (35 mm)	A4416
10835**	Pre-cut	1 <sup>9</sup> /16" (40 mm)	A4416

\*\*Because of its larger starter hole, 10835 functions as both a pre-cut and cut-to-fit solution.



### SenSura® Mio Original Flat MIDI Closed Pouch

- Standard wear
- Filter
- Soft cloth on front and back
- Peek-a-boo inspection window allows easy application
- Not made with natural rubber latex



# 1-Piece SenSura® Mio

Opaque	/" (18 cm) / 340 mL	-	Units: 30
Code	Cut	Stoma Size	HCPCS*
10821	Cut-to-fit	<sup>5</sup> /8-1 <sup>3</sup> /4" (15-45 mm)	A4416
10822	Pre-cut	1" (25 mm)	A4416
10823	Pre-cut	1³⁄16" (30 mm)	A4416
10824	Pre-cut	1 <sup>3</sup> /8" (35 mm)	A4416

### SenSura<sup>®</sup> Mio Original Flat MINI Closed Pouch

- Standard wear
- Filter
- Soft cloth on front and back
- Peek-a-boo inspection window allows easy application
- Not made with natural rubber latex

### Opaque 5<sup>1</sup>/2" (14 cm) / 110 mL



Code	Cut	Stoma Size	HCPCS*
10811	Cut-to-fit	<sup>5</sup> / <sub>8</sub> - 1 <sup>3</sup> / <sub>4</sub> " (15-45 mm)	A4416

## Urostomy

### SenSura® Mio Flat MAXI Urostomy Pouch

- Standard wear
- Multi-chambers minimize sloshing
- Anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to urostomy night bag #21365 or most standard night bags
- Neutral gray, water-resistant textile on back
- Not made with natural rubber latex

### Transparent 10<sup>1</sup>/2" (27 cm) / 400 mL

Code	Cut	Stoma Size	HCPCS*
10586	Cut-to-fit	<sup>3</sup> / <sub>8</sub> - 1 <sup>3</sup> / <sub>4</sub> " (10-45 mm)	A5071

### SenSura<sup>®</sup> Mio Soft Convex MAXI Urostomy Pouch

- Standard wear
- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- Multi-chambers minimize sloshing
- Anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Neutral gray, water-resistant textile
- Not made with natural rubber latex

### Transparent 10<sup>1</sup>/2" (27 cm) / 400 mL

nunspure		/ 400 IIIL	Units. 10
Code	Cut	Stoma Size	HCPCS*
16810	Cut-to-fit	<sup>3</sup> /8 - 2" (10-50 mm)	A4392

<mark>Opaque 10<sup>1</sup>⁄2" (27 cm) / 400 mL</mark>			Units: 10
Code	Cut	Stoma Size	HCPCS*
16805	Cut-to-fit	<sup>3</sup> /8 - 1 <sup>5</sup> /16" (10-33 mm)	A4392

### SenSura® Mio Soft Convex MIDI Urostomy Pouch

- Standard wear
- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- Multi-chambers minimize sloshing
- Anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Neutral gray, water-resistant textile
- Not made with natural rubber latex

### Opaque 9<sup>1</sup>/2" (24 cm) / 330 mL

Code	Cut	Stoma Size	HCPCS*
16801	Cut-to-fit	<sup>3</sup> /8 - 1 <sup>5</sup> /16" (10-33 mm)	A4392







Units: 10



Linite: 10

### SenSura<sup>®</sup> Mio Convex Light MAXI Urostomy Pouch >

- Standard wear
- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- Multi-chambers minimize sloshing
- · Anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Neutral gray, water-resistant textile
- Not made with natural rubber latex

### Transparent $10^{1}/_{7}$ " (27 cm) / 400 ml



Units: 10

Transparent 10 <sup>1</sup> /2" (27 cm) / 400 mL			Units: 10
Code	Cut	Stoma Size	HCPCS*
16835	Cut-to-fit	<sup>3</sup> /8- <sup>7</sup> /8" (10-23 mm)	A4392
16836	Cut-to-fit	⁵⁄8-1⁵∕16″ (15-33 mm)	A4392
16837	Cut-to-fit	<sup>3</sup> / <sub>8</sub> - 1 <sup>11</sup> / <sub>16</sub> " (10-43 mm)	A4392
16840	Pre-cut	<sup>13</sup> /16" (21 mm)	A4392
16841	Pre-cut	1" (25 mm)	A4392
16842	Pre-cut	1 <sup>1</sup> /8" (28 mm)	A4392

### Opaque 10<sup>1</sup>/2" (27 cm) / 400 mL

Code	Cut	Stoma Size	HCPCS*
16825	Cut-to-fit	³⁄₃- ⅔″ (10-23 mm)	A4392
16826	Cut-to-fit	<sup>5</sup> /8-1 <sup>5</sup> /16" (15-33 mm)	A4392
16830	Pre-cut	<sup>11</sup> / <sub>16</sub> " (18 mm)	A4392
16831	Pre-cut	<sup>13</sup> / <sub>16</sub> " (21 mm)	A4392
16832	Pre-cut	1" (25 mm)	A4392
16833	Pre-cut	1 <sup>1</sup> /8" (28 mm)	A4392
16827	Pre-cut	1 <sup>1</sup> ⁄4" (31 mm)	A4392

### SenSura<sup>®</sup> Mio Convex Light MIDI Urostomy Pouch

- Standard wear
- Belt tabs for use with 4237, 4247
- · FlexShell is designed for stability and flexibility
- Multi-chambers minimize sloshing
- Anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Neutral gray, water-resistant textile
- Not made with natural rubber latex

### Opeque $9^{1}/2^{"}$ (24 cm) / 330 ml

Opaque 9 <sup>1</sup> /2" (24 cm) / 330 mL			Units: 10
Code	Cut	Stoma Size	HCPCS*
16821	Cut-to-fit	<sup>3</sup> / <sub>8</sub> - <sup>7</sup> / <sub>8</sub> " (10-23 mm)	A4392
16822	Cut-to-fit	<sup>5</sup> ∕8 - 1 <sup>5</sup> ∕16″ (15-33 mm)	A4392
16823	Pre-cut	1 <sup>1</sup> ⁄4" (31 mm)	A4392



### SenSura® Mio Deep Convex MAXI Urostomy Pouch >

- Standard wear
- Belt tabs for use with 4237, 4247
- · FlexShell is designed for stability and flexibility
- Multi-chambers minimize sloshing
- · Anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- · Neutral gray, water-resistant textile
- · Not made with natural rubber latex



Transparent 10 <sup>1</sup> /2" (27 cm) / 400 mL			Units: 10
Code	Cut	Stoma Size	HCPCS*
16865	Cut-to-fit	<sup>3</sup> /8 - <sup>7</sup> /8" (10-23 mm)	A4392
16866	Cut-to-fit	5⁄8-15∕16″ (15-33 mm)	A4392
16867	Cut-to-fit	<sup>3</sup> / <sub>8</sub> - 1 <sup>11</sup> / <sub>16</sub> " (10-43 mm)	A4392
16870	Pre-cut	<sup>13</sup> / <sub>16</sub> " (21 mm)	A4392
16871	Pre-cut	1" (25 mm)	A4392
16872	Pre-cut	1 <sup>1</sup> /8" (28 mm)	A4392

### Opaque 10<sup>1</sup>/2" (27 cm) / 400 mL

	- ( - ))		
Code	Cut	Stoma Size	HCPCS*
16855	Cut-to-fit	<sup>3</sup> /8 - <sup>7</sup> /8" (10-23 mm)	A4392
16856	Cut-to-fit	<sup>5</sup> /8 - 1 <sup>5</sup> /16" (15-33 mm)	A4392
16860	Pre-cut	<sup>11</sup> / <sub>16</sub> " (18 mm)	A4392
16861	Pre-cut	<sup>13</sup> / <sub>16</sub> " (21 mm)	A4392
16862	Pre-cut	1" (25 mm)	A4392
16863	Pre-cut	1 <sup>1</sup> /8" (28 mm)	A4392
16857	Pre-cut	1 <sup>1</sup> ⁄4" (31 mm)	A4392

### SenSura<sup>®</sup> Mio Convex Flip MAXI Urostomy Pouch 🌒

- Designed for curves, bulges and hernias
- · Curved, star-shaped barrier conforms to outward areas with less creasing and folding
- "Flip to fit" design allows for easy application
- · Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Not made with natural rubber latex

### Opaque 10 ½" (27 cm) / 400 ml



Units: 10

Opaque 10 22 (27 cm) / 400 mL					Units: 10
Code	Cut	Stoma Size	Max Cut	Barrier Size	HCPCS*
18401	Cut-to-Fit	<sup>3</sup> /8 - 1 <sup>9</sup> /16'' (10-40 mm)	1º/16" (40 mm)	S	A4392
18402	Pre-Cut	¾" (20 mm)	1º/16" (40 mm)	S	A4392
18403	Pre-Cut	1" (25 mm)	1º/16" (40 mm)	S	A4392
18404	Pre-Cut	1³⁄16" (30 mm)	1º/16" (40 mm)	S	A4392
Transparent 10 ½" (27 cm) / 400 mL					Units: 10

nunsp		2 (27 CH)/ 400	/ III E		011113. 10
Code	Cut	Stoma Size	Max Cut	Barrier Size	HCPCS*
18405	Cut-to-Fit	<sup>3</sup> / <sub>8</sub> - 1 <sup>9</sup> / <sub>16</sub> " (10-40 mm)	1º/16″ (40 mm)	S	A4392

Indicates convexity

Indicates Convex Flip

www.coloplast.us



SenSura® Flat MAXI Drainable Pouch

- Standard wear
- EasiClose<sup>™</sup> WIDE Outlet
- Filter
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex Transparent  $11\frac{1}{2}$ " (29<sup>1</sup>/<sub>2</sub> cm) / 655 mL



1	
1	Pie
F	ece
20	Se
CS*	S
24	5
24	D ®
2	

Code	Cut	Stoma Size	HCPCS*
15521	Cut-to-fit	<sup>3</sup> / <sub>8</sub> -3" (10-76 mm)	A4424
15522	Pre-cut	1" (25 mm)	A4424
15523	Pre-cut	1³⁄16″ (30 mm)	A4424
15524	Pre-cut	1³⁄8″ (35 mm)	A4424
Opaque 11	<sup>1</sup> ⁄2" <b>(29</b> <sup>1</sup> ⁄2 <b>cm) / 655 mL</b>		Units: 20
Code	Cut	Stoma Size	HCPCS*
15531	Cut-to-fit	<sup>3</sup> ⁄8 - 3" (10-76 mm)	A4424
15532	Pre-cut	1" (25 mm)	A4424
15533	Pre-cut	1³⁄16″ (30 mm)	A4424

### SenSura® Flat MAXI Drainable Pouch

Pre-cut

Standard wear

15534

- EasiClose™ WIDE Outlet
- Without filter
- Soft cloth on back
- Not made with natural rubber latex

### Transparent 11<sup>1</sup>/<sub>2</sub>" (29<sup>1</sup>/<sub>2</sub> cm) / 655 mL

Code	Cut	Stoma Size	HCPCS*
15821	Cut-to-fit	<sup>3</sup> / <sub>8</sub> -3" (10-76 mm)	A5061





### SenSura<sup>®</sup> Xpro Flat MAXI Drainable Pouch

Extended wear

- EasiClose<sup>™</sup> WIDE Outlet
- Filter
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex



Transparer	nt 11 <sup>1</sup> ⁄2" (29 <sup>1</sup> ⁄2 cm) /	655 mL	Units: 10
Code	Cut	Stoma Size	HCPCS*
15981	Cut-to-fit	<sup>3</sup> / <sub>8</sub> - 3" (10-76 mm)	A5056

Opaque 1	Units: 10		
Code	Cut	Stoma Size	HCPCS*
15991	Pre-cut	1" (25 mm)	A5056
15992	Pre-cut	1³⁄16" (30 mm)	A5056
15993	Pre-cut	11⁄4" (32 mm)	A5056
15994	Pre-cut	1¹⁄₂" (38 mm)	A5056

### SenSura<sup>®</sup> Convex Light MAXI Drainable Pouch

- Standard wear
- EasiClose<sup>™</sup> WIDE Outlet
- Filter
- Belt tabs for use with 4215, 4220
- · Soft cloth on back (and front of opaque)
- Not made with natural rubber latex

C	10	
U		
-	-	

Units: 10

# Transparent 11<sup>1</sup>⁄<sub>2</sub>" (29<sup>1</sup>⁄<sub>2</sub> cm) / 655 mL

Code	Cut	Stoma Size	HCPCS*
15606	Cut-to-fit	<sup>5</sup> /8-1 <sup>11</sup> /16" (15-43 mm)	A4389
15623	Pre-cut	1" (25 mm)	A4389
15624	Pre-cut	1 <sup>1</sup> /8" (28 mm)	A4389

Opaque 11 <sup>1</sup> ⁄2" (29 <sup>1</sup> ⁄2 cm) / 655 mL		Units: 10	
Code	Cut	Stoma Size	HCPCS*
15608	Cut-to-fit	<sup>5</sup> /8 - 1 <sup>11</sup> /16" (15-43 mm)	A4389
15633	Pre-cut	1" (25 mm)	A4389
15634	Pre-cut	1 <sup>1</sup> / <sub>8</sub> " (28 mm)	A4389

Indicates convexity

### SenSura<sup>®</sup> Xpro Convex Light MAXI Drainable Pouch 🖊

- Extended wear
- EasiClose<sup>™</sup> WIDE Outlet
- Filter
- Belt tabs for use with 4215, 4220
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex



Transparent 11 <sup>1</sup> / <sub>2</sub> " (29 <sup>1</sup> / <sub>2</sub> cm) / 655 mL		Units: 10	
Code	Cut	Stoma Size	HCPCS*
15695	Cut-to-fit	⁵⁄8 - 1⁵⁄16″ (15-33 mm)	A5057
15696	Cut-to-fit	<sup>5</sup> /8 - 1 <sup>11</sup> /16" (15-43 mm)	A5057

<mark>Opaque</mark> 11 <sup>1</sup> ⁄2" (29 <sup>1</sup> ⁄2 cm) / 655 mL			Units: 10
Code	Cut	Stoma Size	HCPCS*
15697	Pre-cut	⅔" (22 mm)	A5057
15698	Pre-cut	1" (25 mm)	A5057

### SenSura<sup>®</sup> Flat MAXI **Drainable Pouch**

- Standard wear
- Soft outlet provides effective drainage; connects to a standard drain bag or Coloplast fistula drain bag #14010
- Filter
- Soft cloth on front and back
- Peek-a-boo inspection window allows easy application
- Not made with natural rubber latex

### Opeque 12<sup>1</sup>/<sub>2</sub>" (32 cm) / 630 mL

Opaque 1	2 <sup>1</sup> ⁄2" (32 cm) / 630 m	L	Units: 30
Code	Cut	Stoma Size	HCPCS*
19004	Cut-to-fit	<sup>3</sup> / <sub>8</sub> -3" (10-76 mm)	A4424

For bed drainage bag, see page 73.

### SenSura<sup>®</sup> Flat MAXI **Closed Pouch**

- Standard wear
- Filter
- Soft cloth on front and back
- Not made with natural rubber latex



### Opque 8<sup>1</sup>/2" (21 cm) / 550 mL

opuque ovr (Er eni)/ 550 me			011103. 50	
Code	Cut	Stoma Size	HCPCS*	
15480	Cut-to-fit	<sup>3</sup> ⁄8-3" (10-76 mm)	A4416	
15484	Pre-cut	1¾" (35 mm)	A4416	





## Urostomy

### SenSura<sup>®</sup> Flat MAXI Urostomy Pouch

- Standard wear
- Multi-chambers minimize sloshing
- Anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
  Soft cloth on back
- Not made with natural rubber latex





Transparer	nt 10 3/8" (26 cm) / 400 mL	-	Units: 10
Code	Cut	Stoma Size	HCPCS*
11804	Cut-to-fit	<sup>3</sup> ⁄ <sub>8</sub> - 3" (10-76 mm)	A5071

### SenSura® Flat MIDI Urostomy Pouch

- Standard wear
- Multi-chambers minimize sloshing
- Anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Soft cloth on back
- Not made with natural rubber latex

Transparent 9 <sup>1</sup> / <sub>2</sub> " (24 cm) / 330 mL			Units: 10
Code	Cut	Stoma Size	HCPCS*
11803	Cut-to-fit	<sup>3</sup> / <sub>8</sub> - 2 <sup>1</sup> / <sub>2</sub> " (10-66 mm)	A5071

#### SenSura® Xpro Flat MAXI Urostomy Pouch

- Extended wear
- Multi-chambers minimize sloshing
- Anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex



Transparer	nt 10 <sup>3</sup> /8" (26 cm)	/ 400 mL	Units: 10
Code	Cut	Stoma Size	HCPCS*
11809	Cut-to-fit	<sup>3</sup> / <sub>8</sub> - 3" (10-76 mm)	A4391
Opaque 1	0³⁄ଃ" (26 cm) / 40	0 mL	Units: 10
Code	Cut	Stoma Size	HCPCS*
11805	Pre-cut	⅔" (22 mm)	A4391
11806	Pre-cut	1" (25 mm)	A4391
11807	Pre-cut	1³⁄16″ (30 mm)	A4391
11808	Pre-cut	1 <sup>1</sup> ⁄4" (32 mm)	A4391

#### SenSura<sup>®</sup> Convex Light MAXI Urostomy Pouch

- Standard wear
- Multi-chambers minimize sloshing
- Anti-reflux valves prevent urine from washing
   over the stoma
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Belt tabs for use with 4215, 4220
- Soft cloth on back
- Not made with natural rubber latex

#### Transparent 10<sup>3</sup>/<sub>8</sub>" (26 cm) / 400 mL

	===== (=====,,,,		
Code	Cut	Stoma Size	HCPCS*
11814	Cut-to-fit	<sup>5</sup> ∕8-1 <sup>5</sup> ∕16″ (15-33 mm)	A4392



Units: 10

Indicates convexity

#### SenSura® Xpro Convex Light MAXI Urostomy Pouch

- Extended wear
- Belt tabs for use with 4215, 4220
- Multi-chambers minimize sloshing
- · Anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex



Transpare	nt 10 <sup>3</sup> /8" (26 cm)	/ 400 mL	Units: 10
Code	Cut	Stoma Size	HCPCS*
11819	Cut-to-fit	⁵⁄8-1⁵∕16″ (15-33 mm)	A4393

Opaque 103/8" (26 cm) / 400		00 mL	Units: 10
Code	Cut	Stoma Size	HCPCS*
11816	Pre-cut	<sup>7</sup> ⁄₃" (22 mm)	A4393
11817	Pre-cut	1" (25 mm)	A4393
11818	Pre-cut	1¼" (32 mm)	A4393

#### Urostomy Night Bag and Bag Hanger

See Accessories page 73 for more information.



#### Assura<sup>®</sup> Flat MAXI **Drainable Pouch**

- Standard wear
- EasiClose<sup>™</sup> WIDE Outlet
- Filter
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex

Transparent 11 <sup>1</sup> / <sub>2</sub> " (29 cm) / 530 mL
---

	. ,,		
Code	Cut	Stoma Size	HCPCS*
15860	Cut-to-fit	<sup>3</sup> ⁄ <sub>8</sub> -3" (10-75 mm)	A4424

Opaque 11½" (29 cm) / 530 mL			Units: 10	
Code	Cut	Stoma Size	HCPCS*	
15870	Cut-to-fit	<sup>3</sup> ⁄8-3" (10-75 mm)	A4424	
15874	Pre-cut	1" (25 mm)	A4424	
15875	Pre-cut	1 <sup>3</sup> /16" (30 mm)	A4424	
15876	Pre-cut	1³⁄8" (35 mm)	A4424	

#### Assura<sup>®</sup> Flat MAXI **Drainable Pouch**

- Standard wear
- EasiClose<sup>™</sup> WIDE Outlet
- Without filter
- Soft cloth on front and back
- Not made with natural rubber latex

Opaque 11 <sup>1</sup> / <sub>2</sub> " (29 cm) / 530 mL	Opaque	11½"	(29	cm)	/	530	mL
--	--------	------	-----	-----	---	-----	----

Opaque 11 <sup>1</sup> /	2" (29 cm) / 530	0 mL	Units: 10
Code	Cut	Stoma Size	HCPCS*
15871	Cut-to-fit	<sup>3</sup> /8-3" (10-75 mm)	A5061

#### Assura<sup>®</sup> Flat MAXI **Drainable Pouch**

- Extended wear
- EasiClose<sup>™</sup> WIDE Outlet
- Without filter
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex

Transparent	11¼"	(28½	cm) ,	/ 530	mL
-------------	------	------	-------	-------	----

Code	Cut	Stoma Size	HCPCS*
15869	Cut-to-fit	<sup>3</sup> / <sub>8</sub> -3" (10-75 mm)	A4388
Opaque 1	1¼" (28½ cm) / 5	30 mL	Units: 10
Code	Cut	Stoma Size	HCPCS*
15879	Cut-to-fit	<sup>3</sup> /8-3" (10-75 mm)	A4388



Units: 10



Units: 10

# 1-Piece Assura®

#### Assura<sup>®</sup> Convex Light MAXI Drainable Pouch

- Standard wear
- EasiClose<sup>™</sup> WIDE Outlet
- Filter
- Belt tabs for use with 4215, 4220
- Soft cloth on back
- Not made with natural rubber latex

#### Transparent 11½" (29 cm) / 530 mL



mansparer	IL 11/2 (29 Cm)/	530 ML	Units: 10
Code	Cut	Stoma Size	HCPCS*
14511	Cut-to-fit	5⁄8 - 15∕16″ (15-33 mm)	A4389
14512	Cut-to-fit	<sup>5</sup> /8 - 1 <sup>11</sup> /16" (15-43 mm)	A4389

#### Assura® Flat MAXI Drainable Pouch

- Standard wear
- EasiClose<sup>™</sup> Outlet
- Filter
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex



Units: 10

#### Transparent 11¼" (28<sup>1</sup>/<sub>2</sub> cm) / 530 mL

Code	Cut	Stoma Size	HCPCS*
13860	Cut-to-fit	<sup>3</sup> ⁄ <sub>8</sub> - 3" (10-75 mm)	A4424
Opaque 1	1¼" (28½ cm) / 5	30 mL	Units: 10
Code	Cut	Stoma Size	HCPCS*
13870	Cut-to-fit	<sup>3</sup> ⁄ <sub>8</sub> - 3" (10-75 mm)	A4424
13874	Pre-cut	1" (25 mm)	A4424
13875	Pre-cut	1 <sup>3</sup> /16" (30 mm)	A4424
13876	Pre-cut	1¾" (35 mm)	A4424

#### Assura<sup>®</sup> Flat MIDI Drainable Pouch

- Standard wear
- EasiClose<sup>™</sup> Outlet
- Filter
- Soft cloth on front and back
- Not made with natural rubber latex

#### Opaque 101/4" (26 cm) / 390 mL



Units: 10

Code	Cut	Stoma Size	HCPCS*
13840	Cut-to-fit	<sup>3</sup> / <sub>8</sub> - 2 <sup>1</sup> / <sub>2</sub> " (10-65 mm)	A4424
13844	Pre-cut	1" (25 mm)	A4424
13845	Pre-cut	1³⁄16″ (30 mm)	A4424
13846	Pre-cut	1 <sup>3</sup> /8" (35 mm)	A4424

Indicates convexity

#### Assura<sup>®</sup> Convex Light MAXI Drainable Pouch

- Standard wear
- EasiClose<sup>™</sup> Outlet
- Filter
- Belt tabs for use with 4215, 4220
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex

		100	
		0	
			1
F	F	-12	q
ŕ	ŕ	1	9

Units: 10

Transparent 11¼" (28½ cm) / 530 mL			Units: 10
Code	Cut	Stoma Size	HCPCS*
14411	Cut-to-fit	<sup>5</sup> / <sub>8</sub> -1 <sup>5</sup> / <sub>16</sub> " (15-33 mm)	A4389
14412	Cut-to-fit	<sup>5</sup> /8 - 1 <sup>11</sup> /16" (15-43 mm)	A4389

#### Opaque 11¼" (28½ cm) / 530 mL

opuque 11/4	(20)2 cm)/ 550	1112	011103. 10
Code	Cut	Stoma Size	HCPCS*
14414	Cut-to-fit	<sup>5</sup> ∕ <sub>8</sub> -1 <sup>5</sup> ∕ <sub>16</sub> " (15-33 mm)	A4389
14415	Cut-to-fit	<sup>5</sup> / <sub>8</sub> - 1 <sup>11</sup> / <sub>16</sub> " (15-43 mm)	A4389
14421	Pre-cut	<sup>7</sup> ∕₃" (22 mm)	A4389
14422	Pre-cut	1" (25 mm)	A4389
14423	Pre-cut	1 <sup>1</sup> /8" (28 mm)	A4389
14424	Pre-cut	1 <sup>1</sup> ⁄4" (31 mm)	A4389

#### Assura<sup>®</sup> Convex Light MIDI Drainable Pouch

- Standard wear
- EasiClose<sup>™</sup> Outlet
- Filter

Opaque

- Belt tabs for use with 4215, 4220
- Soft cloth on front and back
- Not made with natural rubber latex



Units: 10

10¼″	(26 cm) / 39	0 ml
TO 14	(200011) / 33	

Code	Cut	Stoma Size	HCPCS*
14404	Cut-to-fit	<sup>5</sup> / <sub>8</sub> -1 <sup>5</sup> / <sub>16</sub> " (15-33 mm)	A4389
14405	Cut-to-fit	<sup>5</sup> / <sub>8</sub> - 1 <sup>11</sup> / <sub>16</sub> " (15-43 mm)	A4389

Indicates convexity

#### Assura® Deep Convex MAXI Drainable Pouch 🖊

- Standard wear
- EasiClose<sup>™</sup> Outlet
- Filter
- Belt tabs for use with 4215, 4220
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex



Units: 10

Transparent 11¼" (28½ cm) / 530 mL

Code	Cut	Stoma Size	HCPCS*
14103	Cut-to-fit	<sup>5</sup> / <sub>8</sub> - 1 <sup>11</sup> / <sub>16</sub> " (15-43 mm)	A4389
14163	Pre-cut	<sup>7</sup> ⁄₀" (22 mm)	A4389
14164	Pre-cut	1" (25 mm)	A4389
14165	Pre-cut	1 <sup>1</sup> ⁄8" (28 mm)	A4389
14166	Pre-cut	1¼" (31 mm)	A4389

#### Opaque 11¼" (28½ cm) / 530 mL

Opaque 11¼" (28½ cm) / 530 mL			Units: 10
Code	Cut	Stoma Size	HCPCS*
14106	Cut-to-fit	<sup>5</sup> / <sub>8</sub> - 1 <sup>11</sup> / <sub>16</sub> " (15-43 mm)	A4389
14173	Pre-cut	<sup>7</sup> ⁄⁄₀" (22 mm)	A4389
14174	Pre-cut	1" (25 mm)	A4389
14175	Pre-cut	1 <sup>1</sup> ⁄⁄8" (28 mm)	A4389
14176	Pre-cut	1 <sup>1</sup> ⁄4" (31 mm)	A4389



#### Assura<sup>®</sup> Deep Convex MIDI Drainable Pouch

- Standard wear
- EasiClose<sup>™</sup> Outlet
- Filter
- Belt tabs for use with 4215, 4220
- Soft cloth on front and back
- Not made with natural rubber latex



1 I I
D
$\mathbf{O}$
D
$\rightarrow$
5
Š.
Q
Ē
6

Opaque	10¼"	(26 cm) /	390 mL
--------	------	-----------	--------

Code	Cut	Stoma Size	HCPCS*
14196	Cut-to-fit	<sup>5</sup> /8 - 1 <sup>11</sup> /16" (15-43 mm)	A4389

#### Assura<sup>®</sup> Convex Light MAXI Drainable Pouch

- Extended wear
- EasiClose<sup>™</sup> Outlet
- Without filter
- Belt tabs for use with 4215, 4220
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex

#### t of opaque) **per latex**

#### Transparent 11¼" (28½ cm) / 530 mL



Units: 10

Code	Cut	Stoma Size	HCPCS*
14413	Cut-to-fit	5⁄% -15∕16″ (15-33 mm)	A4390
14416	Cut-to-fit	<sup>5</sup> / <sub>8</sub> - 1 <sup>11</sup> / <sub>16</sub> " (15-43 mm)	A4390
<mark>Орадие</mark> 11¼" (28½ ст) / 530 mL			Units: 10
Code	Cut	Stoma Size	HCPCS*
14419	Cut-to-fit	<sup>5</sup> /8 - 1 <sup>11</sup> /16" (15-43 mm)	A4390

#### Assura<sup>®</sup> Flat MAXI Drainable Pouch

- Standard wear
- Clamp closure
- Without filter
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex



Transparent 11¼" (28½ cm) / 530 mL			Units: 10
Code	Cut	Stoma Size	HCPCS*
12660	Cut-to-fit	<sup>3</sup> / <sub>8</sub> - 2 <sup>3</sup> / <sub>4</sub> " (10-70 mm)	A5061
Ораque 11¼" (28½ ст) / 530 mL			Units: 10
Code	Cut	Stoma Size	HCPCS*
12670	Cut-to-fit	<sup>3</sup> / <sub>8</sub> - 2 <sup>3</sup> / <sub>4</sub> " (10-70 mm)	A5061

Indicates convexity

#### Assura<sup>®</sup> Flat MINI Drainable Pouch

- Standard wear
- Clamp closure
- Without filter
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex



Units: 10

Transparent 934" (24½ cm) / 230 mL			Units: 10
Code	Cut	Stoma Size	HCPCS*
12630	Cut-to-fit	<sup>3</sup> / <sub>8</sub> - 2¼″ (10-55 mm)	A5061

#### Opaque 9¾" (24½ cm) / 230 mL

Code	Cut	Stoma Size	HCPCS*
12620	Cut-to-fit	<sup>3</sup> / <sub>8</sub> - 2 <sup>1</sup> / <sub>4</sub> " (10-55 mm)	A5061

#### Assura® Deep Convex MAXI Drainable Pouch

- Standard wear
- Clamp closure
- Without filter
- Belt tabs for use with 4215, 4220
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex

Units: 10

Linits: 10

#### Transparent 11¼" (28½ cm) / 530 mL

Code	Cut	Stoma Size	HCPCS*
13706	Cut-to-fit	<sup>5</sup> / <sub>8</sub> - 1 <sup>11</sup> / <sub>16</sub> " (15-43 mm)	A4389

#### Opaque 11¼" (28½ cm) / 530 mL

			011103. 10
Code	Cut	Stoma Size	HCPCS*
12533	Cut-to-fit	<sup>5</sup> / <sub>8</sub> - 1 <sup>11</sup> / <sub>16</sub> " (15-43 mm)	A4389


#### Assura® Deep Convex MINI Drainable Pouch

- Standard wear
- Clamp closure
- Without filter
- Belt tabs for use with 4215, 4220
- Soft cloth on back
- Not made with natural rubber latex

#### Transparent 9¾" (24½ cm) / 230 mL



Code	Cut	Stoma Size	HCPCS*
13703	Cut-to-fit	<sup>5</sup> / <sub>8</sub> - 1 <sup>11</sup> / <sub>16</sub> " (15-43 mm)	A4389



#### Assura<sup>®</sup> Deep Convex MAXI Drainable Pouch 🖊

- Extended wear
- Clamp closure
- Without filter
- Belt tabs for use with 4215, 4220
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex

Transparent 11¼" (28½ cm) / 530 mL			Units: 10
Code Cut Stoma Size			HCPCS*
12416	Cut-to-fit	<sup>5</sup> / <sub>8</sub> - 1 <sup>11</sup> / <sub>16</sub> " (15-43 mm)	A4390
Opaque 1	Units: 10		

Code	Cut	Stoma Size	HCPCS*
12419	Cut-to-fit	<sup>5</sup> ∕8 - 1 <sup>11</sup> ∕16″ (15-43 mm)	A4390

#### Assura<sup>®</sup> Flat MAXI **Closed Pouch**

- Standard wear
- Filter

12177

- Soft cloth on front and back
- Not made with natural rubber latex

Pre-cut



A4416

Opaque 8½" (22 cm) / 530 mL			Units: 30
Code	Cut	Stoma Size	HCPCS*
12170	Cut-to-fit	<sup>13</sup> ⁄ <sub>16</sub> - 3" (20-75 mm)	A4416
12174	Pre-cut	1" (25 mm)	A4416
12175	Pre-cut	1³⁄16" (30 mm)	A4416
12176	Pre-cut	1³⁄8" (35 mm)	A4416

1<sup>9</sup>/16" (40 mm)

#### Assura<sup>®</sup> Flat MIDI Closed Pouch

- Standard wear
- Filter
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex

#### Transparent 7" (18 cm) / 350 mL

	( ))		
Code	Cut	Stoma Size	HCPCS*
12130	Cut-to-fit	<sup>13</sup> / <sub>16</sub> - 2½" (20-65 mm)	A4416
12134	Pre-cut	1" (25 mm)	A4416
12135	Pre-cut	1³⁄16″ (30 mm)	A4416
12136	Pre-cut	1 <sup>3</sup> / <sub>8</sub> " (35 mm)	A4416

Opaque 7	" (18 cm) / 350 mL		Units: 30
Code	Cut	Stoma Size	HCPCS*
12140	Cut-to-fit	<sup>13</sup> / <sub>16</sub> - 2 <sup>1</sup> / <sub>2</sub> " (20-65 mm)	A4416
12144	Pre-cut	1" (25 mm)	A4416
12145	Pre-cut	1³⁄16" (30 mm)	A4416
12146	Pre-cut	1¾" (35 mm)	A4416
12147	Pre-cut	1 <sup>9</sup> ⁄16" (40 mm)	A4416
12149	Pre-cut	2" (50 mm)	A4416

#### Assura<sup>®</sup> Flat MINI Closed Pouch

- Standard wear
- Filter
- Soft cloth on front and back
- Not made with natural rubber latex

1			-	6
1	6		1	A
	W-			
	1	~		
	-	-		y
	~		1	

Units: 30

#### Opaque 5" (12 cm) / 170 mL

Code	Cut	Stoma Size	HCPCS*
12110	Cut-to-fit	<sup>13</sup> / <sub>16</sub> - 2 <sup>1</sup> / <sub>2</sub> " (20-65 mm)	A4416



#### Assura<sup>®</sup> Convex Light MAXI Closed Pouch

- Standard wear
- Filter
- Belt tabs for use with 4215, 4220
- Soft cloth on back (and front of opaque)

Not made with natural rubber latex



Transparer	Units: 10		
Code	Cut	Stoma Size	HCPCS*
14441	Cut-to-fit	<sup>5</sup> / <sub>8</sub> - 1 <sup>5</sup> / <sub>16</sub> " (15-33 mm)	A4417
14442	Cut-to-fit	<sup>5</sup> /8 - 1 <sup>11</sup> /16" (15-43 mm)	A4417

#### Opaque 8½" (22 cm) / 560 mL

Code	Cut	Stoma Size	HCPCS*
14444	Cut-to-fit	⁵⁄≋ - 1⁵⁄₁₀″ (15-33 mm)	A4417
14445	Cut-to-fit	<sup>5</sup> / <sub>8</sub> - 1 <sup>11</sup> / <sub>16</sub> " (15-43 mm)	A4417

#### Assura® Convex Light MIDI Closed Pouch

- Standard wear
- Filter
- Belt tabs for use with 4215, 4220
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex

#### Transparent 7" (18 cm) / 350 mL

Code	Cut	Stoma Size	HCPCS*
14431	Cut-to-fit	⁵⁄8 - 1⁵⁄16" (15-33 mm)	A4417
14432	Cut-to-fit	⁵⁄8 - 1 <sup>11</sup> ∕16″ (15-43 mm)	A4417

#### Opaque 7" (18 cm) / 350 mL

Code	Cut	Stoma Size	HCPCS*
14434	Cut-to-fit	<sup>5</sup> / <sub>8</sub> - 1 <sup>5</sup> / <sub>16</sub> " (15-33 mm)	A4417
14435	Cut-to-fit	<sup>5</sup> / <sub>8</sub> - 1 <sup>11</sup> / <sub>16</sub> " (15-43 mm)	A4417
14451	Pre-cut	<sup>7</sup> ⁄₃″ (22 mm)	A4417
14452	Pre-cut	1" (25 mm)	A4417
14453	Pre-cut	1 <sup>1</sup> /8" (28 mm)	A4417
14454	Pre-cut	1 <sup>1</sup> ⁄4" (31 mm)	A4417

Indicates convexity



Units: 10

#### Assura<sup>®</sup> Flat Stoma Cap

- Standard Wear
- Filter
- Belt tabs
- · Soft cloth on back
- Small pouch size for intimate moments or colostomates who irrigate



Units: 30

• Not made with natural rubber latex

#### Opaque

opuque			011103. 50	
Code	Cut	Stoma Size	HCPCS*	
2501	Cut-to-fit	<sup>13</sup> / <sub>16</sub> - 2 <sup>1</sup> / <sub>4</sub> " (20-55 mm)	A5055	

#### Urostomy

#### Assura<sup>®</sup> Flat MAXI Urostomy Pouch

- Standard wear
- Multi-chambers minimize sloshing
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Soft cloth on back
- Not made with natural rubber latex



Transparent 10¾" (27½ cm) / 470 mL			Units: 10
Code	Cut	Stoma Size	HCPCS*
14222	Cut-to-fit	<sup>3</sup> / <sub>8</sub> - 2 <sup>1</sup> / <sub>4</sub> " (10-55 mm)	A5071

#### Assura<sup>®</sup> Convex Light MAXI Urostomy Pouch

- Standard wear
- Multi-chambers minimize sloshing
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Belt tabs for use with 4215, 4220
- Soft cloth on back
- Not made with natural rubber latex

Transparent 10<sup>3</sup>/<sub>4</sub>" (27<sup>1</sup>/<sub>2</sub> cm) / 470 mL



Code	Cut	Stoma Size	HCPCS*
14717	Cut-to-fit	<sup>5</sup> / <sub>8</sub> - 1 <sup>11</sup> / <sub>16</sub> " (15-43 mm)	A4392



#### Assura<sup>®</sup> Flat MAXI **Urostomy Pouch**

- Standard wear
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex Transparent 10<sup>3</sup>/<sub>4</sub>" (27 cm) / 465 mL



Units: 10

A5071

	( //		
Code	Cut	Stoma Size	HCPCS*
5585	Cut-to-fit	<sup>3</sup> / <sub>8</sub> - 2 <sup>1</sup> / <sub>4</sub> " (10-55 mm)	A5071
Opaque 1	Units: 10		
Code	Cut	Stoma Size	HCPCS*

<sup>3</sup>/<sub>8</sub> - 2<sup>1</sup>/<sub>4</sub>" (10-55 mm)

#### Assura<sup>®</sup> Flat MIDI **Urostomy Pouch**

#### Standard wear

5580

- · Soft outlet can be connected to urostomy
- night bag #21365 or any standard night bag
- · Soft cloth on back (and front of opaque) Not made with natural rubber latex

Cut-to-fit

	i.			2	
1					
1					
16		-	÷.		7)
		3		D	1
			27		

Units: 10

Transparent	7" (18 cm) /	265 mL
Cada	Cut	Stores Size

Code	Cut	Stoma Size	HCPCS*
5575	Cut-to-fit	<sup>3</sup> / <sub>8</sub> - 2 <sup>1</sup> / <sub>4</sub> " (10-55 mm)	A5071
Opaque 7	" (18 cm) / 265 mL		Units: 10
Code	Cut	Stoma Size	HCPCS*
5570	Cut-to-fit	<sup>3</sup> / <sub>8</sub> - 2 <sup>1</sup> / <sub>4</sub> " (10-55 mm)	A5071

#### Assura<sup>®</sup> Flat MAXI Urostomy Pouch

#### Extended wear

- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Soft cloth on back
- Not made with natural rubber latex

#### Transparent 91/3" (24 cm) / 375 ml

Transparent 91/2" (24 cm) / 375 mL			Units: 10
Code	Cut	Stoma Size	HCPCS*
12474	Cut-to-fit	<sup>3</sup> / <sub>8</sub> - 2 <sup>1</sup> / <sub>4</sub> " (10-55 mm)	A4391

#### Assura<sup>®</sup> Flat MIDI **Urostomy Pouch**

- Extended wear
- Soft outlet can be connected to urostomy
- night bag #21365 or any standard night bag
- Soft cloth on back (and front of opaque)
- Not made with natural rubber latex



Transparent 9½" (24 cm) / 265 mL			Units: 10
Code Cut Stoma Size			HCPCS*
12472	Cut-to-fit	<sup>3</sup> ∕8 - 2 <sup>1</sup> ⁄4″ (10-55 mm)	A4391
Opaque 91/2" (24 cm) / 265 mL			Units: 10

Code	Cut	Stoma Size	HCPCS*
12471	Cut-to-fit	<sup>3</sup> / <sub>8</sub> - 2 <sup>1</sup> / <sub>4</sub> " (10-55 mm)	A4391

#### Assura<sup>®</sup> Deep Convex MAXI Urostomy Pouch 🖊

- Standard wear
- · Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Belt tabs for use with 4215, 4220
- Soft cloth on back
- Not made with natural rubber latex

#### Transparent 10 ¾" (27 cm) / 465 mL



Units: 10

			011103. 10
Code	Cut	Stoma Size	HCPCS*
12595	Cut-to-fit	⁵⁄₃ - 1⁵∕ı₅" (15-33 mm)	A4392
12596	Cut-to-fit	<sup>5</sup> / <sub>8</sub> - 1 <sup>11</sup> / <sub>16</sub> " (15-43 mm)	A4392
12991	Pre-cut	⁵⁄‰" (15 mm)	A4392
12992	Pre-cut	<sup>3</sup> ⁄4″ (18 mm)	A4392
12993	Pre-cut	<sup>7</sup> ⁄₀" (22 mm)	A4392
12994	Pre-cut	1" (25 mm)	A4392
12995	Pre-cut	1 <sup>1</sup> /8" (28 mm)	A4392
12996	Pre-cut	1 <sup>1</sup> ⁄4" (31 mm)	A4392
12997	Pre-cut	1³⁄ଃ" (35 mm)	A4392

#### Assura<sup>®</sup> Deep Convex MAXI Urostomy Pouch /

- Extended wear
- · Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Belt tabs for use with 4215, 4220
- Soft cloth on back
- Not made with natural rubber latex

#### Transparent 10<sup>3</sup>/<sub>4</sub>" (27 cm) / 465 mL

Transparent 10 <sup>3</sup> /4" (27 cm) / 465 mL			Units: 10
Code	Cut	Stoma Size	HCPCS*
12432	Cut-to-fit	<sup>5</sup> ∕ <sub>8</sub> - 1 <sup>5</sup> ∕ <sub>16</sub> " (15-33 mm)	A4393
12433	Cut-to-fit	<sup>5</sup> / <sub>8</sub> - 1 <sup>11</sup> / <sub>16</sub> " (15-43 mm)	A4393
ndice	ates convexity		



#### SenSura® Mio 1-piece Single Use Kit

- Pouch #10481 MAXI / EasiClose<sup>™</sup> WIDE Outlet
- Instructions for use
- Stoma measuring guide
- Not made with natural rubber latex

#### Transparent

Transparent			Units: 5
Code Cut Stoma Size			HCPCS*
21007	Cut-to-fit	<sup>3</sup> / <sub>8</sub> - 2 <sup>1</sup> / <sub>8</sub> " (10-55 mm)	A4424

#### SenSura® 1-piece Single Use Kit

- Pouch #15521 MAXI / EasiClose<sup>™</sup> WIDE Outlet
- Instructions for use
- Stoma measuring guide
- Not made with natural rubber latex

#### Transparent

Code	Cut	Stoma Size	HCPCS*
11008	Cut-to-fit	<sup>3</sup> / <sub>8</sub> - 3" (10-76 mm)	A4424

#### Assura<sup>®</sup> 1-piece Single Use Kit

- Pouch #15860 MAXI / EasiClose<sup>™</sup> WIDE Outlet
- Instructions for use
- Stoma measuring guide
- Not made with natural rubber latex

Transparent	

Code	Cut	Stoma Size	HCPCS*
11005	Cut-to-fit	<sup>3</sup> ⁄ <sub>8</sub> - 3" (10-76 mm)	A4424

#### SenSura® Xpro 1-piece Urostomy Single Use Kit

- Pouch #11809 MAXI / soft outlet
- Instructions for use
- Stoma measuring guide
- Urostomy adapter
- Not made with natural rubber latex

#### Transparent 10<sup>3</sup>/<sub>8</sub>" (26 cm) / 400 mL

Code	Cut	Stoma Size	HCPCS*
21005	Cut-to-fit	<sup>3</sup> ⁄ <sub>8</sub> - 3" (10-76 mm)	A4391

#### SenSura® Xpro 1-piece Convex Light Urostomy Single Use Kit

- · Pouch #11819 MAXI / soft outlet
- Instructions for use
- Stoma measuring quide
- Urostomy adapter
- Not made with natural rubber latex

Transparent 10 <sup>3</sup> /s" (26 cm) / 400 mL		Units:	
Code	Cut	Stoma Size	HCPCS

Code	Cut	Stoma Size	HCPCS*
21006	Cut-to-fit	<sup>5</sup> ∕8-1 <sup>5</sup> ∕16″ (15-33 mm)	A4393

Indicates convexity

\*For important information regarding the use of HCPCS codes listed and code definitions see page 192



Units: 5

5


# Download the Coloplast Stoma Care App



Designed with nurses in mind

#### Features:

- 1. Easy access to Coloplast® Care
- 2. BodyCheck<sup>™</sup> Assessment
- 3. Step-by-step guides with how-to-use videos
- 4. Common problems and recommended actions
- 5. General information on urostomies, ileostomies and colostomies



Download the Coloplast Stoma Care app now! Search for "Stoma Care" in the App Store and Google Play to download for free today.





# Post op / fistulas



Elastic adhesive fits over folds and contours, yet remains secure against effluent



Two outlet types for different options:

- Soft (for liquid stool; fits standard drain bags)
- Tap (for thicker output; fits fecal drain bag with wide-lumen tubing)

Post-op



Printed measurement scale on post op pouches helps visually track output volume



# Security and comfort from hospital to home.

Your patients deserve all the benefits of SenSura<sup>®</sup> Mio from day one following their operation.

Getting used to a stoma is challenging enough without the added complexity of having to learn a different pouching system. Why not give your patients the security and comfort of Mio - right from the start? Options are available in different sizes, outlet types, and coupling systems to meet patient needs.



2-piece Click Soft | Tap \*Also available in smaller size! pg. 83



2-piece Flex Soft | Tap pg. 109



1-piece Soft pg. 129



Post Op with window Soft Sterile | Non-Sterile pg. 163

#### SenSura<sup>®</sup> Mio 1-piece Post Op Pouch with Window

- Neutral gray, water-resistant textile
- Printed scale provides visual indication of output volume
- Easy-to-open and close post-op window attached to pouch
- Soft outlet can connect to Coloplast's fistula drain bag #14010, or a drain bag with a standard connection
- Not made with natural rubber latex

#### Transnarent

in an op al offe		01110010	
Code	Cut	Volume	Stoma Size
18690	Cut-to-fit	680mL	<sup>3</sup> ⁄ <sub>8</sub> - 2 <sup>3</sup> ⁄ <sub>4</sub> " (10-70 mm)
18691	Cut-to-fit	725mL	³⁄8 - 4″ (10-100 mm)

#### SenSura<sup>®</sup> Mio 1-piece STERILE Post-Op Pouch with Window

- Neutral gray, water-resistant textile Printed scale provides visual indication of output volume
- Easy-to-open and close access window attached to pouch
- · Soft outlet can connect to Coloplast's fistula drain bag #14010, or a drain bag with a standard connection
- Not made with natural rubber latex

#### Transparent

Code	Cut	Volume	Stoma Size
18680	Cut-to-fit	680mL	<sup>3</sup> / <sub>8</sub> - 2 <sup>3</sup> / <sub>4</sub> " (10-70 mm)
18681	Cut-to-fit	725mL	<sup>3</sup> / <sub>8</sub> - 4" (10-100 mm)

#### SenSura<sup>®</sup> Standard Wear 1-piece Post-Op Pouch w/ Window

- Flat
- Soft outlet can connect to Coloplast's fistula drain bag #14010, or a drain bag with a standard connection
- · Easy-to-open and close access window attached to pouch
- Soft cloth on back
- Not made with natural rubber latex

Code	Cut	Stoma Size	HCPCS*
19020	Cut-to-fit	³⁄ଃ - 3″ (10-76 mm)	A6154
19021	Cut-to-fit	<sup>3</sup> / <sub>8</sub> - 4 <sup>1</sup> / <sub>2</sub> " (10-115 mm)	A6154

#### Transparent 13" (33 cm) / 830 mL





Units: 5

Units: 5



Units: 5

-Piece Post-Op / Wound Pouche

#### SenSura<sup>®</sup> Standard Wear 1-piece STERILE Post-Op Pouch w/ Window

- Flat
- Soft outlet can connect to Coloplast's fistula drain bag #14010, or a drain bag with a standard connection
- · Easy-to-open and -close access window attached to pouch
- · Soft cloth on back
- Not made with natural rubber latex

#### ansparent 13" (33 cm) / 830 ml



Transparent 13" (33 cm) / 830 mL			Units: 5
Code	Cut	Stoma Size	HCPCS*
19010	Cut-to-fit	³⁄ଃ - 3″ (10-76 mm)	A6154
19011	Cut-to-fit	<sup>3</sup> / <sub>8</sub> - 4 <sup>1</sup> / <sub>2</sub> " (10-115 mm)	A6154

#### SenSura® Standard Wear 1-piece Post-Op Pouch

#### • Flat

- Soft outlet can connect to Coloplast's fistula drain bag #14010, or a drain bag with a standard connection
- · Effective pouching solution for small fistulas
- Soft cloth on back
- Not made with natural rubber latex .

#### 1 21/// (21 --> / coo



Units: 5

Transparent 1244" (31 cm) / 690 mL			Units: 5
Code	Cut	Stoma Size	HCPCS*
19030	Cut-to-fit	³⁄ଃ - 3″ (10-76 mm)	A6154

#### SenSura<sup>®</sup> Standard Wear 1-piece STERILE **Post-Op Pouch**

- Flat
- · Soft outlet can connect to Coloplast's fistula drain bag #14010, or a drain bag with a standard connection
- Soft cloth on back
- Not made with natural rubber latex

#### Transparent 121/4" (31 cm) / 690 mL

Code	Cut	Stoma Size	HCPCS*
19001	Cut-to-fit	¾ - 3″ (10-76 mm)	A6154

For SenSura<sup>®</sup> Magnum high output pouches, see page 89.

#### Assura<sup>®</sup> Standard Wear 1-piece Post-Op Pouch w/ Window

- Flat
- Soft outlet can connect to Coloplast's fistula drain bag #14010, or a drain bag with a standard connection
- Easy-to-open and -close access window
- Soft cloth on back
- Not made with natural rubber latex

#### Transparent 13" (33 cm) / 830 mL

Code	Cut	Stoma Size	HCPCS*
12800	Cut-to-fit	<sup>1</sup> / <sub>2</sub> - 2 <sup>3</sup> / <sub>4</sub> " (13-70 mm)	A6154
12802	Cut-to-fit	<sup>1</sup> / <sub>2</sub> - 4" (13-100 mm)	A6154

#### Assura<sup>®</sup> Standard Wear 1-piece STERILE Post-Op Pouch w/ Window

- Flat
- Soft outlet can connect to Coloplast's fistula drain bag #14010, or a drain bag with a standard connection
- Easy-to-open and -close access window
- Soft cloth on back
- Not made with natural rubber latex

#### Transparent 13" (33 cm) / 830 mL



#### Assura<sup>®</sup> Standard Wear 1-piece Post-Op Pouch

- Flat
- Soft outlet can connect to Coloplast's fistula drain bag #14010, or a drain bag with a standard connection
- Soft cloth on back
- Not made with natural rubber latex

#### Transparent 121/4" (31 cm) / 690 mL

Code	Cut	Stoma Size	HCPCS*
12805	Cut-to-fit	¹⁄₂ - 2³⁄₄″ (13-70 mm)	A6154

#### Assura® Standard Wear 1-piece STERILE Post-Op Pouch

- Flat
- Soft outlet can connect to Coloplast's fistula drain bag #14010, or a drain bag with a standard connection
- Soft cloth on back
- Not made with natural rubber latex

#### Transparent 121/4" (31 cm) / 690 mL

Units: 5

Units: 5

Code	Cut	Stoma Size	HCPCS*
12818	Cut-to-fit	<sup>1</sup> / <sub>2</sub> -2 <sup>3</sup> / <sub>4</sub> " (13-70 mm)	A6154

For Assura® high output pouches and ileo night pouches, see page 98-99.



#### Coloplast® Fistula and Wound Management System

Each system includes:

- 1 pouch and tracing guide
- 1 drain port (extras sold separately on next page)
- 1 lid (extras sold separately on next page)
- 1 clear drape
- 1 clamp
- Not made with natural rubber latex

Each box includes one pump.



90



14050 Mini

14060 Midi

14070 Maxi

Units: 3 Units for 14050: 6

Code	Cut	Stoma Size	HCPCS*
14050	Mini	41/8 x 61/4" (104-159 mm)	A6154
14060	Midi	6 <sup>1</sup> / <sub>8</sub> x 9" (156-228 mm)	A6154
14070	Maxi	8¼8 x 11¾" (208-297 mm)	A6154

#### Fistula / Wound Pouch Accessories

#### Coloplast<sup>®</sup> Fistula Drain Port

Can be easily applied to the lid to allow for catheter insertion. For catheters ranging from 18FR - 22FR, the drain port is ready for use. For catheters ranging from 24FR - 30FR, pull off the top of the drain port - DO NOT CUT. • Not made with natural rubber latex



Code	Size	Units
14015	One size	12

#### Coloplast<sup>®</sup> Bed Drainage Bag

Bag tubing specifically designed for fecal output. Holds up to 2000 mL. Long tubing can be cut easily, allowing nurses to find the right length for bed drainage.

Not made with natural rubber latex



Code	Size	Units	HCPCS*
14010	2000 mL	6	A4357

For matching bag hanger, see page 73.

#### Coloplast® Transparent Fistula Lids

Replacement lids match the transparent lid included with each fistula pouch. • Not made with natural rubber latex

Code	Size	Units
14021	Mini (matches 14050)	10
14031	Midi (matches 14060)	10
14041	Maxi (matches 14070)	10

#### Coloplast<sup>®</sup> Pouch Clamp

Simple closure clamp is easy to manipulate, even for those with low dexterity. • Not made with natural rubber latex



Code	Size	HCPCS*
9500	One size	A4363

#### Extra Protection

#### **Brava® Elastic Barrier Strips**

- Elastic (follows the contours of the body)
- Helps ensure the position of the appliance and minimize roll-up
- · Skin-friendly and secure alternative to tape
- Not made with natural rubber latex

Brava® Elastic Barrier Strips- Curved For round or oval barrier shapes

Code	Size	Units	HCPCS*
120700"	5 <sup>1</sup> / <sub>2</sub> " strips	20 strips	A4362
"2 strips =	1 billable unit of A4362		

Brava® Elastic Barrier Strips- Y-shape Extra-wide shape is designed to extend the barrier edges at '3' and '9' for users with challenging body shapes - like an outward area / hernia or a deep skin fold.

Code	Size	Units	HCPCS*
120721"	5 <sup>3</sup> /4" strips	30 strips	A5121
** 2	1		

"3 strips = 1 billable unit of A5121

Brava<sup>®</sup> Elastic Barrier Strips- Straight Rectangular strip designed to fit the straight edges of a square ostomy barrier.

Code	Size	Units	HCPCS*
120740"	5¼" strips	20 strips	A4362

"2 strips = 1 billable unit of A4362

#### Brava® Elastic Barrier Strip XL Designed for increased security & support

around the barrier.

Code	Size	Units	HCPCS*
12076**	6 <sup>3</sup> /4" strips	30 strips	A5121
**3 strips - 1	billable unit of A5121		

3 strips = 1 billable unit of A5121

#### Brava<sup>®</sup> Strip Paste

- · Pectin-based product that protects peristomal skin
- · Protects peristomal skin from stoma output
- · Used to fill uneven skin surfaces to create a flat pouching surface
- Alcohol-free and sting-free
- Not made with natural rubber latex

Code	Size	Units	HCPCS*
26555	2 oz. package (.2 oz / strip)	10 strips	A4406

\*For important information regarding the use of HCPCS codes listed and code definitions see page 192









#### Brava® Protective Seal

Moldable barrier ring designed to adhere securely and protect the skin around the stoma.

- Soft, flexible material adapts to body shape as the wearer bends and stretches
- · Can be easily and cleanly removed
- Alcohol-free and sting-free
- Multiple sizes for a custom fit
- Not made with natural rubber latex



Code	Stoma size	Adhesive diameter	Thickness	Seal width	Units	HCPCS*
12035	3⁄4″	1 7/8"	2.5 mm	Standard (5/8")	10	A4385
12036	3⁄4″	2 ¼"	2.5 mm	Wide (¾")	10	A4385
12038	3/4"	2 1⁄2″	2.5 mm	Wide (7/8")	10	A4385
12037	1 1/8″	2 ¼"	2.5 mm	Standard (5/8")	10	A4385
12039	1 ³/8″	2 1⁄2″	2.5 mm	Standard (5/8")	10	A4385
12045	3⁄4″	1 7/8"	4.2 mm	Standard (5/8")	10	A4385
12046	3⁄4″	2 ¼"	4.2 mm	Wide (¾")	10	A4385
12048	3⁄4″	2 1⁄2″	4.2 mm	Wide (7/8")	10	A4385
12047	1 1/8"	2 ¼"	4.2 mm	Standard (5/8")	10	A4385
12049	1 ³/8″	2 1⁄2″	4.2 mm	Standard (5/8")	10	A4385
Coming	soon					
12032	3⁄4″	3″	2.5 mm	Extra wide (1¹/՛ଃ")	10	A4385
12033	1 1/8"	3″	2.5 mm	Extra wide ( <sup>15</sup> / <sub>16</sub> ")	10	A4385

#### Brava<sup>®</sup> Moldable Ring

Designed for dual-action leakage reduction - used to fill in uneven skin surfaces and protect the skin.

- fill in uneven skin surfaces and protect the skin.
  Molds and adheres securely to dips, creases, and folds to help create a tight seal around the stoma
- Helps keep skin healthy by protecting it from output
   and absorbing moisture from the skin
- Easy to handle, mold and apply
- Alcohol-free and sting-free
- Not made with natural rubber latex

Code	Size	Units	HCPCS*
120307	2 mm (thin)	10	A4385
120427	4.2 mm (thick)	10	A4385

#### Brava<sup>®</sup> Paste

Pectin-based product that protects peristomal skin

- · Protects peristomal skin from stoma output
- Used to fill uneven skin surfaces to create a flat pouching surface
- · Alcohol-free and sting-free
- Not made with natural rubber latex

Code	Size	Units	HCPCS*
12050	2 oz.	1	A4406



<sup>\*</sup>For important information regarding the use of HCPCS codes listed and code definitions see page 192

#### Coloplast<sup>®</sup> Paste (non-pectin)

Used to fill in uneven skin surfaces to create a flat pouching surface.

- Provides tight sealing around the stoma
- Protects against stoma output, minimizing risk of maceration
- Low alcohol content; no pectin
- Not made with natural rubber latex

Code	Size	Units	HCPCS*
2650	2 oz.	12	A4405

#### Brava<sup>®</sup> Powder

- Provides effective moisture absorption and protection
- Not made with natural rubber latex

Code	Size	Units	HCPCS*
19075	1 oz.	16	A4371

#### Brava® Adhesive Remover Spray

- Designed to lift the edges of the barrier for ease in removal
- Dry within seconds without leaving an oily residue behind
- Alcohol-free, silicone-based, and sting-free
- Not made with natural rubber latex

Code	Size	Units	HCPCS*
120105	1.7 oz. (50 mL)	1	A4455

#### Brava® Adhesive Remover Wipes

- Wipe away adhesive residue
- Dry within seconds without leaving an oily residue behind
- Alcohol-free, silicone-based, and sting-free

Tip: wring the wipe like a wet rag to expel the liquid, then wipe to disperse on skin and wipe away residues

#### Not made with natural rubber latex

Code	Size	Units	HCPCS*
120115	1 wipe	30 wipes/box	A4456

#### Brava<sup>®</sup> Skin Barrier Spray

- Protect against effluent and adhesives
- Dry within seconds without building a thick or tacky layer
- Alcohol-free, silicone-based, and sting-free
- Not made with natural rubber latex

Code	Size	Units	HCPCS*
120205	1.7 oz. (50 mL)	1	A4369

#### Brava<sup>®</sup> Skin Barrier Wipes

- Protect against effluent and adhesives
- Dry within seconds without building a thick or tacky layer
- Alcohol-free, silicone-based, and sting-free

Tip: wring the wipe like a wet rag to expel the liquid, then wipe to disperse on skin for extra protection

#### Not made with natural rubber latex

Code	Size	Units	HCPCS*
120215	1 wipe	30 wipes/box	A5120











# **BodyCheck**<sup>™</sup>

Identifying ostomy solutions based on individual needs

Together with 300 Wound, Ostomy, Continence (WOC) nurses, we've developed BodyCheck<sup>™</sup> which uses Body Profile Terminology, a common language to describe the peristomal body profile.

Assess your patient's overall parastomal body profile as a Regular, Inward, or Outward shape



Visit **www.bodycheck.coloplast.us** to complete the assessment

Room American
and the second se
Comments of the local division of the local
Reduct to

Receive personalized ostomy solutions (may include accessory options)

### Visit www.bodycheck.coloplast.us



# Coloplast Pediatric Pouching Systems

## For little fighters

Babies and kids present different pouching challenges than adults - so we've applied **BodyFit Technology**<sup>®</sup> to develop specific solutions for the different stages in a child's life.

In 2019, we launched a new range of **SenSura® Mio** pouches specifically designed for babies and kids.



SenSura<sup>®</sup> Mio Baby Long wear time. Minimal disturbance.



SenSura<sup>®</sup> Mio Kids So kids can be kids.

### Ask us about our pediatric education!











#### Pediatric FAQ

A quick tipsheet with answers to frequently asked questions, like "what do I do if friends and family - or other kids - have questions?" and "what do I do if my child pulls the pouch off?"

Caregivers can request this via our Coloplast<sup>®</sup> Care program.

#### • My Child Has An Ostomy: A Resource Sheet for Caregivers

A fill-in-the-blank sheet to help parents communicate about their child's needs to day care providers, school nurses, and other caregivers.

Caregivers can request this via our Coloplast<sup>®</sup> Care program.

#### Life After Your Child's Ostomy

A high-level overview of anatomy, how to change a pouch, pouching system options, where to get supplies, dietary notes, and other information. Caregivers can request this via our Coloplast® Care program.

## • When I Met Tipster: A coloring book for children with ostomies

A simple coloring book to help kids engage with their ostomy surgery and "new normal" - or learn more about a sibling's ostomy.

To receive a book, please contact your Coloplast sales representative. Patients can also receive via the Coloplast<sup>®</sup> Care program.

#### Jake Gets a Stoma

A children's book written by Jake's mom, who wanted a better way to explain the surgery and stoma management to her son and other kids going through the same thing.

To receive a book, please contact your Coloplast sales representative. Patients can also receive via the Coloplast® Care program.

#### SenSura<sup>®</sup> Mio Bunny

A furry friend for babies and kids with a stoma (or their siblings!). Hang tag includes the number for Coloplast<sup>®</sup> Care in case parents or caregivers have any product or lifestyle questions - and the belly is sized to fit a SenSura<sup>®</sup> Mio Baby barrier (thus providing a good opportunity for practicing management skills). To receive a bunny, please contact your Coloplast sales representative.

Pediatric

#### SenSura<sup>®</sup> Mio Baby

2-piece pouching system designed with BodyFit Technology<sup>®</sup> to tackle the challenges of preemies and babies from 30 weeks of gestation to 6 months of age.

#### SenSura® Mio Baby Flex 2-Piece Barrier

- Oval barrier accommodates one or multiple stomas
- No starter hole custom cut
- · Gentle adhesive sticks securely and removes gently with Brava<sup>®</sup> Baby Adhesive Remover
- · Soft flange enables the pouch to be removed and replaced
- Sterile, single-packed
- Not made with natural rubber latex

#### Cut-to-fit

Cut-to-fit			Units: 5
Code	Stoma Size	Color	HCPCS*
18705	0 - 1 <sup>9</sup> /16'' (0-40 mm)		A4414

#### SenSura® Mio Baby Flex 2-Piece Pouch

- · Adhesive coupling is easy to apply and remove
- Access window
- Soft outlet fits to enteral feeding systems
- Sterile, single-packed
- Not made with natural rubber latex

#### Transparent, No Filter Units: 10 HCPCS\* Code Color 18700 A5063

#### Brava® Baby Adhesive Remover

- · Oil-based adhesive remover designed specifically for babies
- · Liquid sachet (no wipe) can be used with cotton swabs or gauze
- Gently removes adhesive / adhesive residues from SenSura<sup>®</sup> Mio **Baby Barrier**
- Not recommended for use with SenSura Mio Kids or adult products
- Not made with natural rubber latex

		Units: 30
Code	Size	HCPCS*
12014	3 mL sachet	A4455



#### SenSura<sup>®</sup> Mio Kids

1-piece and 2-piece pouching systems designed with BodyFit Technology<sup>®</sup> to fit kids' active movements and rounded bellies. Designed with options to step up as kids grow, the range covers 6 months to 4 years of age (larger options for older kids are found in the SenSura Mio adult line).

<ul> <li>Star-sh flexibilit</li> <li>No star</li> <li>Not ma</li> </ul>	<sup>®</sup> Mio Kids Flex 2-Piece Bar aped, elastic barrier can be rota y in placement ter hole – custom cut ide with natural rubber latex		
Cut-to-fit			Units: 5
Code	Stoma Size	Color	HCPCS*
18735	0 - 1¹/16″ (0-27 mm)		A4414
18747	0 - 1 <sup>3</sup> /8" (0-35 mm)		A4414
Drainabl • EasiClo • Full-circ • Neutral • Small si over tin	se <sup>™</sup> Outlet <mark>cle filter</mark> gray, water-resistant textile ize enables kids to step up in pou	-	
Transpare	nt 61/8" folded		Units: 10
Code		Color	HCPCS*
18721			A4425
Neutral gr	ay with inspection window 6 <sup>1</sup> /8"	folded	Units: 10
Code		Color	HCPCS*
18720			A4425

Transparent 6¾" folded Code Colo	
Code Colo	Units: 10
	or HCPCS*
18723	A4425
18725	A4425
Neutral gray with inspection window 634" folded	Units: 10
Code Colo	or HCPCS*
18722	A4425
18724	A4425
<ul> <li>SenSura<sup>®</sup> Mio Kids Flex 2-piece MIDI Urostomy Pouch</li> <li>Neutral gray, water-resistant textile</li> <li>Soft outlet can be connected to a drain bag for nighttime use, or a leg bag for kids in a wheelchai</li> <li>Not made with natural rubber latex</li> </ul>	r Linits: 10

Opaque 6¾" folded		Units: 10
Code	Color	HCPCS*
18730		A5073
18731		A5073

2-Piece SenSura® Mio Kids

#### SenSura<sup>®</sup> Mio Kids 1-Piece MIDI Drainable Pouch

- Star-shaped, elastic barrier with no starter hole custom cut
- EasiClose<sup>™</sup> Outlet
- Full-circle filter
- · Neutral gray, water-resistant textile
- Not made with natural rubber latex

#### Transparent 6¾" folded | Cut-to-fit

Code	Stoma size	HCPCS*
18711	0 – 1³/ɛ" (0-35 mm)	A4424
Neutral gray	y with inspection window 6¾" folded   Cut-to-fit	Units: 10
Code	Stoma size	HCPCS*
18710	0 – 1³/8" (0-35 mm)	A4474

Units: 10

#### SenSura<sup>®</sup> Mio Kids 1-Piece MIDI Urostomy Pouch

- Star-shaped, elastic barrier with starter hole
- · Neutral gray, water-resistant textile
- Soft outlet can be connected to a drain bag for nighttime use, or a leg bag for kids in a wheelchair
- Not made with natural rubber latex

Transparen	Units: 10	
Code	Stoma size	HCPCS*
18712	<sup>3</sup> /8-1 <sup>3</sup> /8" (10-35 mm)	A5071
Neutral gray with inspection window		Units: 10
Code	Stoma size	HCPCS*
18713	<sup>3</sup> /8-1 <sup>3</sup> /8" (10-35 mm)	A5071

#### Assura® AC (Adhesive Coupling) Pediatric

A simple, soft, secure system for neonates! This ostomy pouching system is based on adhesives that make it easy to apply the pouch to the barrier. The pouching system is easy to clean, and simple to handle with the EasiClose<sup>™</sup> standard integrated outlet.

#### Assura<sup>®</sup> AC Pediatric Flat Barrier

#### Standard wear

- No starter hole
- Not made with natural rubber latex

Cut-to-fit			Units: 10
Code	Stoma Size	Color	HCPCS*
14307	0 - <b>5⁄8</b> ″ (0-17 mm)		A4414
14308	0 - 1" (0-27 mm)		A4414

# Assura<sup>®</sup> AC Pediatric Drainable Pouch EasiClose<sup>~</sup> Outlet Filter

Soft cloth on back

14692

Not made with natural rubber latex

6<sup>1</sup>/2" (16 cm)

Transparent			Units: 10
Code	Length	Color	HCPCS*
14681	5 <sup>1</sup> ⁄2" (14 cm)		A4425
14682	6 <sup>1</sup> ⁄2" (16 cm)		A4425

# Assura® AC Pediatric Drainable Pouch • EasiClose" Outlet • Filter • Soft cloth on front and back • Not made with natural rubber latex Opaque (Teddy Bear) Color HCPCS\* 14691 5½" (14 cm) A4425

A4425

#### Assura® 2-piece Pediatric Pouching Systems

Specifically developed adhesive and a smaller barrier size are ideal for pediatrics, while a mix-and-match "green" flange size offers versatility. The SenSura® and Assura® pouches and barrier may be used interchangeably with "green"-sized adult SenSura® Click or Assura® products for a broad array of pouching options.

Assura <sup>®</sup> Pediatric Standard Wear Barrier		and and a second
<ul> <li>Flat</li> </ul>		0
Not made with natural rubber latex		2
Cut-to-fit		Units: 5
Code	Color	HCPCS*
<b>2182</b> 3/8 - 13/8" (10-35 mm)		A4414
Assura <sup>®</sup> Pediatric Drainable Pouch <ul> <li>Clamp closure</li> <li>Soft cloth on back</li> <li>Not made with natural rubber latex</li> </ul>		
Transparent 81/2" (21 cm) / 150 mL		Units: 10
Code	Color	HCPCS*
2156		A4426
Assura <sup>®</sup> Pediatric Closed Pouch <ul> <li>Soft cloth on front and back</li> <li>Not made with natural rubber latex</li> </ul>		<b>*</b>
Opaque 5¾" (14½ cm)		Units: 30
Code	Color	HCPCS*
2161		A4423
Assura <sup>®</sup> Stoma Cap • With filter • Small size is optimal for swimming at the poo • Not made with natural rubber latex Opaque	l	Units: 30
Code	Color	HCPCS*
2801		A5055
<ul> <li>Assura<sup>®</sup> Pediatric Urostomy Pouch</li> <li>Soft outlet can be connected to a drain bag f nighttime use, or a leg bag for kids in a wheelchair - see page 73 for options</li> <li>Soft cloth on back</li> <li>Not made with natural rubber latex</li> </ul>	or	O I
Transparent 81⁄2" (21 cm) / 150 mL		Units: 10
Code	Color	HCPCS*
1789		A4433

<ul> <li>Assura<sup>®</sup> Urostomy Micro-Pouch</li> <li>Miniature size prevents twisting</li> <li>Soft outlet can be connected to a drain b nighttime use, or a leg bag for people in a wheelchair - see page 73 for options</li> <li>Soft cloth on front</li> <li>Not made with natural rubber latex</li> </ul>	0	Ø
Transparent 6" (15 cm)		Units: 10
Code	Color	HCPCS*
14204		A4433
Assura <sup>®</sup> Urostomy Minicap <ul> <li>Ultra-absorbent material inside</li> <li>Small size is optimal for swimming at the</li> <li>Not made with natural rubber latex</li> </ul> Opague 100 mL	pool	Units: 30
Code	Color	HCPCS*
2807		A5055

\*For important information regarding the use of HCPCS codes listed and code definitions see page 192

## ColoKids<sup>™</sup> 1-piece Pediatric Pouching Systems

Specially developed Assura® "swiss roll" technology provides moisture absorption and balance while offering consistent tack and adhesion to the skin. Small pouch sizes provide a pediatric-friendly pouching option.

## ColoKids<sup>™</sup> Flat 1-piece **Drainable Pouch**

- Standard wear
- Clamp closure
- Soft cloth on back
- Not made with natural rubber latex

### Transparent



Transparent			Units: 10
Code	Cut	Stoma Size	HCPCS*
2115	Cut-to-fit	<b>¾</b> - 1 <b>¾</b> ″ (10-35 mm)	A5061

## ColoKids<sup>™</sup> Flat 1-piece **Closed Pouch**

- Standard wear
- Soft cloth on back
- Not made with natural rubber latex



Units: 10


Transpare	nt		Units: 30
Code	Cut	Stoma Size	HCPCS*
2125	Cut-to-fit	³⁄8 - 1³∕8″ (10-35 mm)	A4416

## ColoKids<sup>™</sup> Flat 1-piece **Urostomy Pouch**

- Standard wear
- Soft outlet can be connected to a drain bag for night-. time use, or a leg bag for people in a wheelchair
- Soft cloth on back
- Not made with natural rubber latex

### Transparent

Code	Cut	Stoma Size	HCPCS*
8009	Cut-to-fit	5∕8 - 1³∕8″ (15-35 mm)	A5071



Pediatric Accessories: Caring for skir

Bedside-Care<sup>®</sup> EasiCleanse<sup>™</sup> (page 51)

**Moisturizers** 

Bathing

Sween Cream<sup>®</sup> (page 53) Sween® 24 Cream (page 53) Sween<sup>®</sup> Lotion (page 53)

**Skin Protectants** 

Critic-Aid<sup>®</sup> Clear (page 54) Critic-Aid<sup>®</sup> Skin Paste (page 54)

Wound Dressings

## Triad<sup>™</sup> Hydrophilic Wound Dressing (page 42)

All product recommendations on this page are based on: Association of Women's Health, Obstetric and Neonatal Nurses. Neonatal Skin Care, Fourth Edition. Evidence-Based Clinical Practice Guideline, 2018.

## Product Guide

Reference tools | Effective June 2020





Coloplast<sup>®</sup> Care is the only manufacturer-provided patient support program with better patient outcomes supported by peer-reviewed data.



## About the team

Coloplast<sup>®</sup> Care Advisors assist patients with their intimate healthcare needs. They are available Monday through Friday from 7 a.m. to 7 p.m. (CT). The team offers one-on-one assistance, as well as:

## Find the right product solution

Coloplast Care can help locate a medical supplies dealer in-network with your patient's insurance to ensure that they have access to the correct products.

## **Answer questions**

Our team can answer questions regarding product lines and educational materials, as well as insurance and reimbursement issues. Coloplast Care will provide the information you and your patients need.

## Provide product samples

Care Advisors are trained support staff, offering free samples from the complete Coloplast ostomy and continence product lines. Call toll free at 1-888-726-7872 or email care-us@coloplast.com to have trial products shipped today.

## 1-888-726-7872

# **Continence Care Education**

## Coloplast<sup>®</sup> Care Program

The program provides patients with product samples\*, educational resources, and consumer care support.

Please visit **bladder.coloplastcare.us** or contact your local Coloplast Care Advisor for further information: **1-866-226-6362**.

For Peristeen<sup>®</sup> inquiries, contact your Coloplast Care Advisor at **1-855-605-7594** or **peristeen@coloplast.com**. Visit **www. coloplastcare.us/bowel** for more information.

## DVDs/USBs

Intermittent Self-Catheterization and Bowel Management

- Intermittent Self-Catheterization, Adult Female, English/Spanish
- Intermittent Self-Catheterization, Adult Male, English/Spanish
- SpeediCath<sup>®</sup> Instructional DVD, Animated, English/Spanish
- Visit www.coloplastcare.us/bowel to view interactive how-to-use videos for Peristeen.
- Visit bladderguides.coloplastcare.us for interactive how-to-use videos for intermittent catheters.

## **Measuring Guides**

- Conveen<sup>®</sup> Optima Male External Catheter Selection and Fitting Guide
- Freedom Male External Catheter Selection and Fitting Guide

## Educational Guides

- $\bullet$  Making Peristeen  $^{\otimes}$  Work for Your Child, Pediatric Patient Information Kit, English
- $\bullet$  Making Peristeen  $^{\otimes}$  Work for Your Child, Pediatric Patient Information Kit, Spanish
- Welcome to Peristeen®, Patient Information Kit Adult, English
- Welcome to Peristeen®, Patient Information Kit Adult, Spanish
- Peristeen® Training Guide for Healthcare Professionals
- Male IC Care Kit, Patient Information Kit- Adult, English & Spanish
- Female IC Care Kit, Patient Information Kit- Adult, English & Spanish

### Wellness Education

- · Socializing and Peristeen®, English
- Socializing and Peristeen®, Spanish
- Traveling and Peristeen®, English
- Traveling and Peristeen<sup>®</sup>, Spanish
- Hand Hygiene, English & Spanish
- Basics of Intermittent Self Catheterization
- Bladder Hygiene, Better Bladder Health
- Medicare Reimbursement Guidelines
- Understanding UTIs
- Cathing Reference Guide
- Basics of MECs
- Why Intermittent Self-Catheterization
- Bladder Diary



Clinician-validated education A dedicated website and emails, with articles and reliable advice, relevant for your situation



My Continence Check Keep track of your bladder issues through regular self-assessment





Product access and guidance Access to product samples and support to help find a supplier that meets your insurance needs

\*SpeediCath®, Self-Cath®, and SureCath® catheters are available by prescription only. Samples are not available for Peristeen

## Coloplast<sup>®</sup> Care Program

The program provides patients with individualized support and product access coordination from an ostomy Care Advisor, nurse-validated education, selfassessment tools, and product samples. Please visit www.ostomy.coloplastcare.us or contact your local Care Advisor for further information: 1-877-858-2656.

## With Coloplast<sup>®</sup> Care, your patients receive:



Individualized education and connections to community support or product access coordination led by patient-focused Ostomy Advisors



Routine selfassessments to identify and address concerns via the My Ostomy Check



Personalized, **nurse**validated education • Coloplast Care Kit • Website and Emails

## **On-Site Education**

- Product in-services on appropriate use of Coloplast products
- Customized clinical education programs conducted by certified wound, ostomy, continence (WOC) nurses

## **Ostomy Tools for Nurses**

Evidence-based nursing tools provide a standardized language to assess peristomal skin conditions and quality of life

- · Stoma QoL (Quality of Life) Questionnaire
- Ostomy Skin Tool
- DET score
- Discharge Plan
- Ostomy Educational Table/Pocket Guide a diagram for staff nurses covering types of ostomies, pouching selection, and peristomal and parastomal challenges.

## **Ostomy Self-Assessment Tools**

- BodyCheck<sup>™</sup>
- Ostomy Check
- Troubleshooter
- Apply-Remove-Check (ARC)

## **Educational Guides**

- Life After: Ileostomy
- Life After: Colostomy
- Life After: Urostomy
- Life After: Your Child's Ostomy
- Teen Life with an Ostomy

## Wellness Education

- · Colostomy: Self-Care and Dietary Guidelines
- Urostomy: Self-Care and Dietary Guidelines
- Ileostomy: Self-Care and Dietary Guidelines

## Wellness Education (cont.)

- Ostomy Pouching Systems
- What to Wear
- Socializing
- Keeping Your Skin Healthy
- Intimacy
- Ostomy Accessory Products
- Getting into Exercise
- Tips for Dining Out
- Tips for Traveling
- Reach Out and Get Connected
- General Medicare Guidelines Related to Ostomy
- Pregnancy with an Ostomy
- My Child has an Ostomy
- Pediatric FAQ

## Product Pictorial Guide

## NOTE: These are abbreviated instructions designed to visually show how to use the product. Always refer to the product labeling (including the package insert) before using.

- SenSura<sup>®</sup> Mio 1-piece Pouching System
- SenSura<sup>®</sup> Mio Click 2-piece Pouching System
- SenSura<sup>®</sup> Mio Flex 2-piece Pouching System
- SenSura<sup>®</sup> 1-piece Pouching System
- SenSura<sup>®</sup> Click 2-piece Pouching System
- SenSura® Flex 2-piece Pouching System
- Fistula Wound Management System
- SenSura<sup>®</sup> Mio 1-piece Urostomy Pouching System
- SenSura® Mio Click 2-piece Urostomy Pouching System
- SenSura<sup>®</sup> Mio Flex 2-piece Urostomy Pouching System

## Spanish Educational Tools Educational Guides

- La Vida Después de la Colostomía Life After: Colostomy
- La Vida Después de la Ileostomía Life After: Ileostomy
- La Vida Después de la Urostomía Life After: Urostomy

## **Educational Sheets**

- Cómo Cuidar Su Nueva Ostomía Caring for your New Ostomy
- Directriz Alimenticio de Colostomía Colostomy Dietary Guidelines
- Directriz Alimenticio de lleostomía Ileostomy Dietary Guidelines
- Directriz Alimenticio de Urostomía Urostomy Dietary Guidelines

For additional educational resources, visit **www.coloplast.us** 

## Professional



## **Online Learning** with CE Contact Hours

Courses for Healthcare Professionals working in continence, ostomy and wound clinical spaces.

## **Expert** Education Programs

Live or virtual education on a variety of clinical topics presented by content experts, providing interactive learning experiences.

## Facility-led Education

Train the Trainer courses include slide presentations, training resources, toolkits and more to help you facilitate a meaningful training experience.

## Education Resources

Videos, podcasts and printables designed to support your patients in everyday situations and help you promote positive outcomes for your patients.



## Register by using one of the following methods:

Website - www.coloplast.us/professional Email - CProfessional\_US@coloplast.us





## Professional

Coloplast Professional is a

comprehensive program that suits a variety of learning needs and provides



ongoing support to deliver best-in-class education and training for your staff.

## Product Training

This hands-on training program is tailored to your staff and your schedule to establish proper product usage and compliance to protocols. Coloplast provides training at implementation and as needed to meet your needs.

## **Expert** Education Programs

A Coloplast Clinical Consultant or Key Opinion Leader in the field of wound and skin care will provide

a live learning experience at your facility, offsite, or via webinar. Learn best practices while sharing experiences with colleagues.

## Online Learning

Coloplast HEAL is an online education program that offers CE credits for healthcare professionals. The program aims to expand knowledge of modern wound healing principles and improve the standard of care for patients. To learn more and register, visit www.coloplast.us/professional

## Facility-led Education

This train-the-trainer program will provide expertise and materials

to key individuals in your organization so they can drive education across the organization. Materials provided include presentations, training tools, and more to help you facilitate a meaningful training experience.

To learn more about the education programs offered by Coloplast, contact Customer Service at 1-800-533-0464.





The information in this document is informational only, general in nature, and does not cover all payers' rules or policies. This information was obtained from third party sources and is subject to change without notice as a result of changes in reimbursement regulations and payer policies. This document represents no promise or guarantee by Coloplast Corp. regarding coverage or payment for products or procedures by CMS or other payers. Providers are responsible for reporting the codes that most accurately describe the patient's medical condition, procedures performed and products used. Providers should check Medicare bulletins, manuals, program memoranda and Medicare guidelines to ensure compliance with Medicare requirements. Inquires should be directed to the appropriate other payer for non-Medicare coverage situations.

## Medicare Utilization - Continence Care

Usual maximum limits for some products from the Medicare Local Coverage Determinations (LCD) are shown below. Suppliers should examine their applicable Medicare contractor policies for full examination of policy requirements, including medical necessity and documentation.

HCPCS*	Long Description	Usual Maximum
A4310	Insert tray without drainage bag and without catheter (accessories only)	1 insertion tray per episode of indwelling catheter insertion
A4332	Lubricant, individual sterile packet, each	200/month
A4349	Male external catheter, with or without adhesive, disposable, each	35/month
A4351	Intermittent urinary catheter; straight tip, with or without coating (Teflon, silicone, silicone elastomer, or hydrophilic, etc.), each	200/month
A4352	Intermittent urinary catheter; Coudé (curved) tip, with or without coating (Teflon, silicone, silicone elastomer, or hydrophilic, etc.) each	200/month
A4353	Intermittent urinary catheter, with insertion supplies	200/month
A4357	Bedside drainage bag, day or night, with or without anti-reflux device, with or without tube, each	2/month
A4358	Urinary drainage bag, leg or abdomen, vinyl, with or without tube, with straps, each	2/month

## Medicare Utilization - Ostomy

Usual maximum limits for some products from the Medicare Local Coverage Determinations (LCD) are shown below. The amount of ostomy supplies varies by patient based on the type of ostomy, its location, its construction, and the condition of the peristomal skin. The table below shows the maximum number of items that are usually medically necessary for patients with ostomies. It may be more or less depending on the factors that affect the frequency of barrier and pouch change. The medical necessity for a greater quantity of supplies must be clearly documented in the patient's medical record (physician's signed and dated order).

HCPCS*	Short Description	Usual Maximum
A4357	Bedside drainage bag	2 ea/month
A4361	Ostomy faceplate, each	3 ea/6 months
A4362	Skin barrier, solid, 4x4 or equivalent; each*	20 ea/month
A4363	Ostomy clamp replacement	none*
A4367	Ostomy belt, each	1 ea/month
A4368	Ostomy filter	none*
A4369	Ostomy skin barrier; liquid, per ounce	2 oz./month
A4371	Ostomy skin barrier powder, per ounce	10 oz/6 months
A4372	Ostomy skin barrier, 4 x 4, standard wear, CONVEX	none*
A4373	Ostomy skin barrier, with flange, CONVEX	none*
A4385	Ostomy skin barrier, 4 x 4, extended wear, FLAT	none*
A4388	One-piece ostomy drainable pouch and barrier, extended wear, FLAT	none*
A4389	One-piece ostomy drainable pouch and barrier, CONVEX	none*
A4390	One-piece ostomy drainable pouch, extended wear barrier, CONVEX	none*
A4391	One-piece urinary pouch, extended wear barrier, FLAT	none*
A4392	One-piece urinary pouch, standard wear barrier, FLAT	none*
A4394	Ostomy deodorant per ounce	none*
A4397	Ostomy irrigation supply; sleeve, each	4 ea/month
A4398	Ostomy irrigation supply; bag, each	2 ea/6 months
A4399	Ostomy irrigation supply; cone/catheter, with or without brush	2 ea/6 months
A4404	Ostomy ring, each	10 ea/month
A4405	Ostomy skin barrier, non-pectin based, paste, per ounce	4 oz/month
A4406	Ostomy skin barrier, pectin-based, per ounce	4 oz/month
A4407	Two-piece ostomy skin barrier, with flange, extended wear, CONVEX, (4 x 4 inches or smaller)	none*
A4408	Two-piece ostomy skin barrier, with flange extended wear, CONVEX, (larger than 4 x 4 inches)	none*

\* No published Medicare usual maximum. Monthly quantity subject to medical necessary documentation.

\*\*Elastic Barrier Strips - Curved or Straight: 40 strips/month

HCPCS*	Short Description	Usual Maximum
A4409	Two-piece ostomy skin barrier, with flange, extended wear, FLAT, (4 x 4 inches or smaller)	none*
A4410	Two-piece ostomy skin barrier, with flange, extended wear, FLAT, (larger than 4 x 4 inches)	none*
A4411	Two-piece ostomy solid skin barrier, extended wear, CONVEX	none*
A4412	Two-piece high output ostomy drainable pouch without filter	none*
A4413	Two-piece high output ostomy drainable pouch with filter	none*
A4414	Two-piece ostomy skin barrier with flange, FLAT (4x4 inches or smaller)	20 ea/month
A4415	Two-piece ostomy skin barrier, with flange, FLAT, greater than 4 x 4 inches	20 ea/month
A4416	One-piece ostomy closed pouch and barrier with filter, FLAT	60 ea/month
A4417	One-piece ostomy closed pouch and barrier with filter, CONVEX	60 ea/month
A4418	One-piece ostomy closed pouch with filter and barrier	60 ea/month
A4419	Two-piece ostomy closed pouch with filter for non-locking system	60 ea/month
A4423	Two-piece ostomy closed pouch with filter for locking system	60 ea/month
A4424	One-piece ostomy drainable pouch and barrier with filter, FLAT	20 ea/month
A4425	Two-piece ostomy drainable pouch with filter for non-locking system	20 ea/month
A4426	Two-piece ostomy drainable pouch without filter for locking system	20 ea/month
A4427	Two-piece ostomy drainable pouch with filter for locking system	20 ea/month
A4433	Two-piece ostomy urinary pouch for locking system	20 ea/month
A4455	Adhesive remover or solvent, per ounce	16 oz/6 months
A4456	Adhesive remover, wipes	none*
A5055	Stoma cap	31 ea/month
A5056	One-piece drainable ostomy pouch, extended wear with filter	40 ea/month
A5057	One-piece drainable ostomy pouch extended wear with filter, CONVEX	40 ea/month
A5061	One-piece ostomy drainable pouch and barrier, FLAT	20 ea/month
A5063	Two-piece ostomy drainable pouch for barrier with flange	20 ea/month
A5071	One-piece ostomy urinary pouch	20 ea/month
A5073	Two-piece ostomy urinary pouch	20 ea/month
A5081	Stoma plug or seal	31 ea/month
A5120	Skin barrier, wipes or swabs, each	150 ea/6 months
A5121	Solid skin barrier, 6x6 or equivalent, each**	20 ea/month
A5122	Solid skin barrier, 8x8 or equivalent, each	20 ea/month
A6154	Wound pouch, each	3/week
0101		J, WCCK

\* No published Medicare usual maximum. Monthly quantity subject to medical necessary documentation.

\*\*Elastic Barrier Strips - Y-shape: 60 strips/month

## Medicare Utilization - Wound Care

Usual maximum limits from the Medicare Local Coverage Determinations (LCD) are shown below. Suppliers should examine their applicable Medicare contractor policies for full examination of policy requirements, including medical necessity and documentation.

HCPCS*	Short Description	Usual Maximum
A6196	Alginate dressing ≤ 16 sq. in.	1/day
A6197	Alginate dressing > 16 ≤ 48 sq. in.	1/day
A6199	Alginate dressing wound filler, per 6 inches	1/day
A6203	Sterile composite dressing, $\leq$ 16 sq. in. with border	3/week
A6204	Sterile composite dressing >16 $\leq$ 48 sq. in. with border	3/week
A6206	Sterile contact layer, ≤ 16 sq. in.	1/week
A6207	Sterile contact layer > $16 \le 48$ sq. in.	1/week
A6209	Sterile foam dressing ≤ 16 sq. in. without border	3/week
A6210	Sterile foam dressing > $16 \le 48$ sq. in. without border	3/week
A6211	Sterile foam dressing > 48 sq. in. without border	3/week
A6212	Sterile foam dressing $\leq 16$ sq. in with border	3/week
A6213	Sterile foam dressing > $16 \le 48$ sq. in. with border	3/week
A6234	Sterile hydrocolloid dressing $\leq$ 6 sq. in. without border	3/week
A6235	Sterile hydrocolloid dressing > $16 \le 48$ sq. in. without border	3/week
A6236	Sterile hydrocolloid dressing > 48 sq. in. without border	3/week
A6237	Sterile hydrocolloid dressing $\leq$ 16 sq. in. with border	3/week
A6238	Sterile hydrocolloid dressing > $16 \le 48$ sq. in. with border	3/week
A6240	Sterile hydrocolloid dressing filler, paste, per ounce	3/week
A6248	Hydrogel dressing, wound filler, gel, per fluid ounce	3 fl. oz. per wound/30 days
A6251	Sterile specialty absorptive dressing, wound cover, ≤ 16 sq. in. without border	1/day
A6252	Sterile specialty absorptive dressing, wound cover, > $16 \le 48$ without border	1/day
A6260	Sterile wound cleanser, any type/size	Not covered under the surgical dressing benefit

Reimbursement Disclaimer: Coloplast Corp. provides this information for your general reference and related to the reimbursement of Coloplast products only. Reimbursement, coverage and payment policies can vary from one insurer and region to another, and may change over time. Coloplast does not guarantee coverage or payment of products.

## HCPCS Code Reference

A4322	Irrigation syringe, bulb or piston, each
A4331	Extension drainage tubing, any type, any length, with connector/adapter, for use with urinary leg bag or urostomy pouch, each
A4332	Lubricant, individual sterile packet, each
A4334	Urinary catheter anchoring device, leg strap, each
A4335	Incontinence supply; miscellaneous
A4338	Indwelling catheter; Foley type, two-way latex with coating (Teflon, silicone, silicone elastomer, or hydrophilic, etc.), each
A4340	Indwelling catheter, specialty type, (e.g., Coudé, mushroom, wing, etc.), each
A4344	Indwelling catheter; Foley type, two-way, all silicone, each
A4346	Indwelling catheter; Foley type, three-way for continuous irrigation, each
A4349	Male external catheter, with or without adhesive, disposable, each
A4351	Intermittent urinary catheter; straight tip, with or without coating (Teflon, silicone, silicone elastomer, or hydrophilic, etc.), each
A4352	Intermittent urinary catheter; Coudé (curved) tip, with or without coating (Teflon, silicone, silicone elastomer, or hydrophilic, etc.), each
A4353	Intermittent urinary catheter, with insertion supplies
A4357	Bedside drainage bag, day or night, with or without anti-reflux device, with or without tube, each
A4358	Urinary drainage bag, leg or abdomen, vinyl, with or without tube, with straps, each
A4361	Ostomy faceplate, each
A4362	Solid skin barrier; solid, 4x4 or equivalent; each
A4363	Ostomy clamp, any type, replacement only, each
A4366	Ostomy vent, any type, each
A4367	Ostomy belt, each
A4368	Ostomy filter, any type, each
A4369	Ostomy skin barrier liquid (spray, brush, etc.), per ounce
A4371	Ostomy skin barrier powder, per ounce
A4372	Ostomy skin barrier, solid 4x4 or equivalent, standard wear, with built-in convexity, each
A4373	Ostomy skin barrier, with flange (solid, flexible or accordion), with built-in convexity, any size, each
A4385	Ostomy skin barrier, solid 4 x 4 or equivalent, extended wear, without built-in convexity, each
A4388	Ostomy pouch, drainable, with extended wear barrier attached, (1 piece), each
A4389	Ostomy pouch, drainable, with barrier attached, with built-in convexity (1-piece), each
A4390	Ostomy pouch, drainable, with extended wear barrier attached, with built-in convexity (1-piece), each

## HCPCS Code Reference

A4391	Ostomy pouch, urinary, with extended wear barrier attached (1-piece), each
A4392	Ostomy pouch, urinary with standard wear barrier attached, with built-in convexity (1-piece), each
A4393	Ostomy pouch, urinary, with extended wear barrier attached, with built-in convexity (1 pc), each
A4394	Ostomy deodorant, with or without lubricant, for use in ostomy pouch, per fluid ounce
A4396	Ostomy belt with peristomal hernia support
A4397	Irrigation supply; sleeve, each
A4398	Ostomy irrigation supply; bag, each
A4399	Ostomy Irrigation supply; cone /catheter, with or without brush
A4402	Lubricant, per ounce
A4404	Ostomy ring, each
A4405	Ostomy skin barrier, non-pectin based, paste, per ounce
A4406	Ostomy skin barrier, pectin-based, paste, per ounce
A4407	Ostomy skin barrier, with flange (solid, flexible, or accordion), extended wear, with built-in convexity, 4x4 inches or smaller, each
A4408	Ostomy skin barrier, with flange (solid, flexible or accordion), extended wear, with built-in convexity, larger than 4x4 inches, each
A4409	Ostomy skin barrier, with flange (solid, flexible or accordion), extended wear, without built in convexity, 4x4 inches or smaller, each
A4410	Ostomy skin barrier with flange (solid, flexible or accordion), extended wear, without built-in convexity, larger than 4x4 inches, each
A4411	Ostomy skin barrier, solid 4x4 or equivalent, extended wear, with built-in convexity, each
A4412	Ostomy pouch, drainable, high output, for use on a barrier with flange (2-piece system), without filter, each
A4413	Ostomy pouch, drainable, high output, for use on a barrier with flange (2 piece system) with filter, each
A4414	Ostomy skin barrier, with flange (solid, flexible or accordion), without built-in convexity, 4x4 inches or smaller, each
A4415	Ostomy skin barrier, with flange (solid, flexible or accordion), without built-in convexity, larger than 4x4 inches, each
A4416	Ostomy pouch, closed, with barrier attached, with filter (1-piece,) each
A4417	Ostomy pouch, closed, with barrier attached, with built-in convexity, with filter (1-piece), each
A4418	Ostomy pouch, drainable, high output, for use on a barrier with flange (2-piece system), without filter, each
A4419	Ostomy pouch, closed; for use on barrier with non-locking flange, with filter (2-piece), each $% \left( 2,1,2,2,2,2,2,2,2,2,2,2,2,2,2,2,2,2,2,2$
A4421	Ostomy supply; miscellaneous
A4422	Ostomy pouch absorbent material (sheet, pad, crystal packet) for use in ostomy pouch to thicken liquid stomal output, each
A4423	Ostomy pouch, closed; for use on barrier with locking flange, with filter (2-piece), each
A4424	Ostomy pouch, drainable, with barrier attached, with filter (1-piece), each

HCPCS Code Reference

### Ostomy pouch, drainable; for use on barrier with non-A4425 locking flange, with filter (2-piece system), each A4426 Ostomy pouch, drainable; for use on barrier with locking flange (2-piece system), each Ostomy pouch, drainable; for use on barrier with locking A4427 flange, with filter (2 piece system), each A4428 Ostomy pouch, urinary, with extended wear barrier attached, with faucet-type tap with valve (1-piece), each A4429 Ostomy pouch, urinary, with barrier attached, with built-in convexity, with faucet-type tap with value (1-piece), each A4430 Ostomy pouch, urinary, with extended wear barrier attached, with built-in convexity, with faucet-type tap with valve (1-piece), each A4431 Ostomy pouch, urinary; with barrier attached, with faucet-type tap with valve (1- piece), each A4432 Ostomy pouch, urinary; with barrier attached, with faucet-type tap with value (1- piece), each A4433 Ostomy pouch, urinary; for use on barrier with locking flange (2-piece), each A4434 Ostomy pouch, urinary; for use on barrier with locking flange, with faucet-type tap with valve (2-piece), each Ostomy pouch, drainable, high output, with extended A4435 wear barrier (1-piece system), with or without filter, each A4455 Adhesive remover or solvent (for tape, cement or other adhesive), per ounce A4456 Adhesive remover, wipes, any type, each A4459 Manual pump-operated enema system, includes balloon, catheter and all accessories, reusable, any type A5051 Ostomy pouch, closed; with barrier attached (1-piece), each A5052 Ostomy pouch, closed; without barrier attached (1-piece), each A5053 Ostomy pouch, closed; for use on faceplate, each A5054 Ostomy pouch, closed; for use on barrier with flange (2-piece), each A5055 Stoma cap Ostomy pouch, drainable, with extended wear barrier A5056 attached, with filter (1-piece), each A5057 Ostomy pouch, drainable, with extended wear barrier attached, with built in convexity, with filter, (1-piece), each A5061 Ostomy pouch, drainable; with barrier attached, (1-piece), each A5062 Ostomy pouch, drainable; without barrier attached (1-piece), each A5063 Ostomy pouch, drainable; for use on barrier with flange (2-piece system), each A5071 Ostomy pouch, urinary; with barrier attached (1-piece), each A5072 Ostomy pouch, urinary; without barrier attached (1-piece), each A5073 Ostomy pouch, urinary; for use on barrier with flange (2-piece), each A5081 Stoma plug or seal, any type

A5093 Ostomy accessory; convex insert

## HCPCS Codes

HCPCS (	Codes
A5113	Leg strap; latex, replacement only, per set
A5114	Leg strap; foam or fabric, replacement only, per set
A5120	Skin barrier, wipes or swabs, each
A5121	Skin barrier; solid, 6x6 or equivalent, each
A5122	Skin barrier; solid, 8x8 or equivalent, each
A5131	Appliance cleaner, incontinence and ostomy appliances, per 16 ounces
A6154	Wound pouch, each
A6196	Alginate or other fiber gelling dressing, wound cover, sterile, pad size 16 sq. in. or less, each dressing
A6197	Alginate or other fiber gelling dressing wound cover, sterile, pad size more than 16 sq. in. but less than or equal to 48 sq. in, each dressing
A6199	Alginate or other fiber gelling dressing, wound filler, sterile, per 6 inches
A6203	Composite dressing, sterile, pad size 16 sq. in. or less, with any size adhesive border, each dressing
A6204	Composite dressing, sterile, pad size more than 16 sq. in. but less than or equal to 48 sq. in., with any size adhesive border, each dressing
A6206	Contact layer, sterile, 16 sq. in. or less, each dressing
A6207	Contact layer, sterile, more than 16 sq. in. but less than or equal to 48 sq. in., each dressing
A6209	Foam dressing, wound cover, sterile, pad size 16 sq. in. or less, without adhesive border, each dressing
A6210	Foam dressing, wound cover, sterile, pad size more than 16 sq. in but less than or equal to 48 sq. in. without adhesive border, each dressing
A6211	Foam dressing, wound cover, sterile, pad size more than 48 sq. in. without adhesive border, each dressing
A6212	Foam dressing, wound cover, sterile, pad size 16 sq. in. or less, with any size adhesive border, each dressing
A6213	Foam dressing, wound cover, sterile, pad size more than 16 sq. in but less than or equal to 48 sq. in. with any size adhesive border, each dressing
A6234	Hydrocolloid dressing, wound cover, sterile, pad size 16 sq. in. or less without adhesive border, each dressing
A6235	Hydrocolloid dressing, wound cover, sterile, pad size more than 16 but less than or equal to 48 sq. in. without adhesive border, each dressing
A6236	Hydrocolloid dressing, wound cover, sterile, pad size more than 48 sq. in. without adhesive border, each dressing
A6237	Hydrocolloid dressing, wound cover, sterile, pad size 16 sq. in. or less, with any size adhesive border, each dressing
A6238	Hydrocolloid dressing, wound cover, sterile, pad size more than 16 sq. in. but less than or equal to 48 sq. in., with any size adhesive border, each dressing
A6239	Hydrocolloid dressing, wound cover, sterile, pad size more than 48 sq. in., with any size adhesive border, each dressing
A6240	Hydrocolloid dressing, wound filler, paste, sterile, per ounce
A6248	Hydrogel dressing, wound filler, gel, per fluid ounce
A6251	Specialty absorptive dressing, wound cover, sterile, pad size 16 sq. in. or less, without border adhesive border, each dressing

HCPCS Codes	
A6252	Specialty absorptive dressing, wound cover, sterile, pad size more than 16 sq. in but less than or equal to 48 sq. in without adhesive border, each dressing
A6260	Wound cleansers, any type, any size
A6413	Adhesive bandage, first-aid type, any size, each Not covered under the surgical dressing benefit

## **Product Guide**

Index | Effective June 2020

Index 

# Making life **easier**

Code Pa	ge#	Code	Ρα	je #	Code	Pa	ge #
	<u> </u>						
Continence		4210		12	6430		23
27-4200-В	24	4212	•••••	12	6500	•••••	23
110	13	4214		12	6530		23
112	13	4306		12	8000		24
114	13	4308		12	8030		24
116	13	4310		12	8100		24
	13	4408		12	8200		24
10	13	4410		12	8205		24
212	13	4412		12	8230		24
14	13	4414		12	8235		24
40	13	4416		12			24
					8300		
05	13	4418		12	8305		24
06	13	4608		15	8400		24
08	13	4610	•••••	15	8430		24
10	13	4612		15	8500		24
08 80	13	4614		15	20010		14
10	13	4808		15	20012		14
12	13	4810		15	20014		14
14	13	4812		15	20016		14
16	13	4814		15	20021		14
18	13	4816		15	20022		14
150	13	5062	,	26	20022		14
160	13	5070		30	21026		25
75	30	5100		23	21027		25
08	15	5110		23	21029		25
10	15	5170		25	21034		25
12	15	5171		25	21053		25
14	15	5174		25	21054		25
06	15	5190		23	21346		26
08	15	5200		23	21356		26
10	15	5210		23	22025		22
L2	15	5221		22	22028		22
14	15	5225		22	22030		22
	15			22			22
			1		22035		
18	15	5230		22	22040		22
.008	19	5230H		22	22121		22
.010	19	5235		22	22125		22
.012	19	5235H	I	22	22130		22
.014	19	5240		22	22135		22
	19	5240⊦	I	22	25502		26
108	19	5290		23	27010		11
112	19	5300		23	27012		11
.114	19	5310					11
116	19	5390		23			11
	19						
2114						S	18
2214	19	5410		23		S	18
614	19	5490		23		S	18
2814	19	5500		23		S	18
2816	19	5590		23	28009		18
3114	19	6100		23	28009	S	18
3214	19	6130		23	28022	s	18
3614	19	6200		23	28023	s	18
4110	12	6230		23		s	18
4112	12	6300		23		S	18
4114	12	6330				· · · · · · · · · · · · · · · · · · ·	18
4116	12						
	12	6400		23	2002/	S	18

Code	Ρας	je #
28028		18
280285		18
280295		18
28032		18
280325		18
280335		18
28036		18
28036S		18
28037		18
28037S		18
28408		10
28410		10
28412		10
28414		10
28416		10
28418		10
28422		16
28481		17
28482		17
28483		17
28484		17
28485		17
28486		17
28490		14
28492		14
28494		14
28496		14
28499		17
28500		17
28501		17
28502		17
28503		17
28504		17
28506		10
28508		10
28510		10
28512		10
28514		10
28516		10
28520		16
28522		16
28524		16
28578		11
28580		11
28582		11
28584		11
28606		10
28608		10
28610		10
28612		10 16
28702 28706		
		10 10
28708 28710		10 10
28710		10
28812		11
20012		т <i>Т</i>

Code Page#	Code Page#
<b>28814</b> 11	<b>1757</b> 101
<b>29122</b> 21 <b>29127</b> 21	<b>1758</b> 101 <b>1759</b> 101
<b>29127</b> 21 <b>29129</b> 21	<b>1759</b> 101 <b>1789</b> 180
	<b>2041</b> 71
<b>50501</b> 30 <b>68001</b> 27	<b>2041</b>
<b>68004</b> 27	<b>2115</b> 182 <b>2125</b> 182
68006 R or L 27	<b>2125</b> 182 <b>2156</b> 180
68007 27	<b>2161</b> 180
<b>90000</b> 30	<b>2182</b> 180
AA 6C18 28	<b>2310</b> 71
AA 6C20 28	<b>2315</b> 71
AA 6C22 28	<b>2320</b> 71
AA 6C24 28	<b>2325</b> 71
AA 6106 28	<b>2330</b> 71
AA 6108 28	<b>2340</b> 71
AA 6110 28	<b>2350</b> 71
AA 6112 28	<b>2501</b> 156
AA 6114 28	<b>2650</b> 69, 171
AA 6116 28	<b>2801</b> 100, 180
AA 6118 28	<b>2802</b> 100
AA 6120 28	<b>2803</b> 100
AA 6122 28	<b>2807</b> 100, 181
AA 6124 28	<b>2808</b> 100
AA 6308 28	<b>2831</b> 93
AA 6310 28	<b>2832</b> 93
AA 6312 29	<b>2833</b> 93
AA 6314 29	<b>2836</b> 99
AA 6316 29	<b>2837</b> 99
AA 6318 29	<b>2846</b> 98
AA 6320 29 AA 6322 29	<b>2847</b> 98 <b>2881</b> 92
AA 6324 29	<b>2881</b> 92 <b>2882</b> 92
AA 6408 28	<b>2883</b> 92
AA 6410 28	<b>4215</b> 73
AA 6412 29	<b>4220</b> 73
AA 6414 29	<b>4237</b> 72
AA 6416 29	<b>4247</b> 72
AA 6418 29	<b>5070</b> 73
AA 6422 29	<b>5570</b> 157
AA 6424 29	<b>5575</b> 157
<b>C1110</b> 19	<b>5580</b> 157
C3814 19	<b>5585</b> 157
C3816 19	<b>7583</b> 71
Ostomy	8009 182
<b>0509</b> 72	<b>8114</b> 98
<b>0925</b> 71	<b>8115</b> 98 <b>8116</b> 98
<b>1003</b> 75 <b>1110</b> 76	<b>8116</b> 98 <b>9500</b> 72, 168
<b>1110</b>	<b>10011</b> 86
1500H 74	<b>10011</b> 86
<b>1751</b> 102	<b>10012</b> 86
<b>1752</b> 102	<b>10015</b> 86
<b>1754</b> 101	<b>10016</b> 86
<b>1755</b> 101	<b>10017</b> 86
<b>1756</b> 101	<b>10021</b> 86

Code Page #	Code Page #	Code Page #
10022 96	10505 70	
<b>10022</b> 86 <b>10023</b> 86	<b>10505</b> 79	<b>11009</b> 103
	<b>10506</b> 79	<b>11010</b> 103
<b>10025</b> 86	<b>10512</b> 79	<b>11011</b> 87
<b>10026</b>	<b>10513</b> 79	<b>11011K</b> 119
<b>10027</b> 86	<b>10514</b> 79	<b>11012</b> 87
<b>10031</b> 86	<b>10515</b> 79	<b>11015</b> 87
<b>10032</b> 86	<b>10516</b> 79	<b>11016</b> 87
<b>10035</b> 86 <b>10036</b> 86	<b>10522</b> 79 <b>10523</b> 79	<b>11021</b> 87
		<b>11022</b> 87
<b>10041</b> 86	<b>10524</b> 79	<b>11023</b> 87
<b>10045</b> 86	<b>10525</b> 79 <b>10531</b> 79	<b>11024</b> 87
<b>10101</b> 111		<b>11025</b> 87 <b>11026</b> 87
<b>10102</b> 111		
<b>10103</b> 111	<b>10561</b> 107	<b>11027</b> 87
<b>10104</b> 111	<b>10562</b> 107	<b>11028</b> 87
<b>10105</b> 111	<b>10563</b> 107	<b>11031</b> 87
<b>10106</b> 111 <b>10108</b> 111	<b>10564</b> 107	<b>11032</b> 87 <b>11035</b> 87
	<b>10565</b> 107 <b>10566</b> 107	
<b>10154</b> 90 <b>10155</b> 90		<b>11036</b> 87
<b>10155</b> 90 <b>10156</b> 90	<b>10567</b> 107 <b>10568</b> 107	<b>11041</b> 87 <b>11045</b> 87
	<b>10568</b> 107 <b>10571</b> 107	<b>11045</b> 87 <b>11114</b> 88
<b>10164</b> 90 <b>10165</b> 90	<b>10571</b> 107 <b>10572</b> 107	<b>11114</b> 88 <b>11115</b> 88
<b>10165</b>	<b>10573</b> 107	<b>11113</b> 88 <b>11124</b> 88
<b>10168</b>	<b>10575</b> 107 <b>10586</b> 137	
<b>10187</b>	<b>10300</b> 137 <b>10811</b> 137	<b>11125</b> 88 <b>11126</b> 88
<b>10167</b> 90 <b>10461</b> 126	<b>10811</b> 137 <b>10821</b> 137	<b>11120</b> 88 <b>11127</b> 88
<b>10471</b> 125	<b>10822</b> 137	<b>11127</b> 88 <b>11135</b> 89
<b>10471</b> 125 <b>10472</b> 125	<b>10822</b> 137 <b>10823</b> 137	<b>11135</b> 89 <b>11136</b> 89
<b>10472</b> 125 <b>10473</b> 125	<b>10824</b> 137	<b>11130</b> 89
<b>10474</b> 125	<b>10824</b> 137 <b>10831</b> 136	<b>11186</b> 89
<b>10475</b> 125	<b>10833</b> 136	<b>11194</b> 88
<b>10476</b> 125	<b>10834</b> 136	<b>11195</b> 88
<b>10477</b> 125	<b>10835</b> 136	<b>11196</b> 88
<b>10478</b> 125	<b>10870</b> 133	<b>11197</b> 88
<b>10479</b> 125	<b>10871</b> 132	<b>11301</b> 111
<b>10480</b> 125	<b>10872</b> 132	<b>11302</b> 111
<b>10481</b> 125	<b>10873</b> 132	<b>11303</b> 111
<b>10482</b> 125	<b>10874</b> 132	<b>11304</b> 112
<b>10483</b> 125	<b>10881</b> 132	<b>11305</b> 112
<b>10484</b> 125	<b>10883</b> 132	<b>11306</b> 112
<b>10485</b> 125	<b>10884</b> 132	<b>11307</b> 111
<b>10486</b> 125	<b>10885</b> 132	<b>11308</b> 112
<b>10487</b> 125	<b>10911</b> 113	<b>11312</b> 111
<b>10488</b> 125	<b>10912</b> 113	<b>11313</b> 111
<b>10489</b> 126	<b>10914</b> 113	<b>11314</b> 111
<b>10490</b> 125	<b>10915</b> 113	<b>11315</b> 111
<b>10491</b> 125	<b>10917</b> 113	<b>11401</b> 84
<b>10492</b> 126	<b>10918</b> 113	<b>11402</b> 84
<b>10493</b> 126	<b>11001</b> 104	<b>11412</b> 84
<b>10494</b> 126	<b>11002</b> 104	<b>11413</b> 84
<b>10495</b> 126	<b>11003</b> 121	<b>11421</b> 84
<b>10496</b> 126	<b>11004</b> 104	<b>11422</b> 84
<b>10502</b> 79	<b>11005</b> 159	<b>11425</b> 84
<b>10503</b> 79	<b>11006</b> 121	<b>11426</b> 84
<b>10504</b> 79	<b>11008</b> 159	<b>11461</b> 82

Code	Page #	Code	Page #	Code	Page #
11462	82	12033	68, 170	12356	
11463		12035	68, 170	12374	
11471	82	12035	68, 170	12375	
11472	82	12030	68, 170	12375	
11473	82	12038	68, 170		
11474		12039	68, 170	12385	
11482	82	12045	68, 170	12386	
11483	82	12046	68, 170		153
11484	82	12047	68, 170	12419	
11486	82	12048	68, 170	12432	
11487	82	12049	68, 170		158
11488	82	12050	69, 170	12471	
11491	85	12060	69	12472	158
11492	85	12061	69	12474	157
11497	85	12076	67, 169	12527	
11498	85	12080	69	12528	
11499	85	12110	154	12529	
11500	85	12130	154	12533	152
11512	112	12134	154	12577	
11514	112	12135	154	12578	
11515	112	12136	154	12579	
11517	112	12140	154	12595	158
11518	112	12144	154	12596	158
11590	112	12145	154	12620	152
11591	112	12146	154	12630	152
11612	112	12147	154	12660	151
11690	112	12149	154	12670	151
11803	144	12170	153	12701	95
11804	144	12174	153	12702	95
11805	145	12175	153	12703	95
11806	145	12176	153	12704	95
11807	145	12177	153	12705	95
11808	145	12202	110	12706	95
11809	145	12203	110	12707	95
11814	145	12213	110	12708	95
11816	146	12216	110	12709	95
11817	146	12222	110	12716	
11818	146	12223	110	12719	
11819	146	12262	108	12800	165
11841	92	12272	108	12802	165
11844	91	12273	108	12805	166
11845	91	12274	108	12808	165
11846	91	12281	108	12810	165
11851	92	12282	108	12818	166
11854	91	12283	108	12820	
11855	91	12284	108	12830	
11856	91	12291	110	12834	
12003	72	12292	110	12835	
12004	72	12296	110	12836	
12005	72	12297	110	12841	
12006	72	12298	110	12842	
12007	72	12344		12844	
12008	72	12345		12845	
12009	72	12346		12847	
12014	175	12354		12848	
12032	68, 170	12355		12991	158

Code	Page #	Code Page	# Code	Page #
12992	158	<b>14229</b> 10	1 <b>14329</b>	117
		<b>14232</b> 10		
	158 158	<b>14232</b>		117 117
		<b>14234</b>		
	158 158			
	158	<b>14235</b>		117 116
	158	<b>14230</b>		116
	152	<b>14238</b>		116
	148	<b>14239</b>		116
	148	<b>14243</b> 9		116
	148	<b>14246</b>		116
	148	<b>14249</b>		115
	148	<b>14250</b>		115
	148	<b>14251</b> 9		115
	148	<b>14252</b>		115
	148	<b>14254</b>		115
	148	<b>14255</b>		115
		<b>14257</b>		114
		<b>14258</b>		114
		<b>14261</b>		114
		<b>14262</b>		149
13966		<b>14263</b>	3 <b>14405</b>	149
13974		<b>14271</b> 9	3 <b>14406</b>	114
13975		<b>14272</b> 9	3 <b>14408</b>	114
13976		<b>14273</b> 9	3 <b>14411</b>	149
13984		<b>14274</b>	3 <b>14412</b>	149
13985		<b>14275</b> 9	3 <b>14413</b>	151
13986		<b>14276</b> 9	3 <b>14414</b>	149
14010	73, 168	<b>14277</b> 9	3 <b>14415</b>	149
14015	73, 168	<b>14278</b>	3 <b>14416</b>	151
	. 74, 168	<b>14281</b> 9		151
	74, 168	<b>14282</b> 9		149
	74, 168	<b>14283</b> 9		149
	167	<b>14291</b> 9		149
	167	<b>14292</b> 9		149
	167	<b>14293</b> 9		155
	150 150	<b>14294</b>		155 155
	150 150	<b>14295</b> 9 <b>14296</b> 9		155 155
	150	<b>14297</b> 9		155
	150	<b>14297</b>		155
	150	<b>14301</b> 11		155
	150	<b>14302</b> 11		
	150	<b>14303</b> 11		
	150	<b>14304</b> 11		
	150	<b>14305</b> 11		155
	151	<b>14306</b> 11		
	100, 181	<b>14307</b> 17		
14205	100	<b>14308</b> 17		148
14206	100	<b>14309</b> 11		116
14222	156	<b>14316</b> 11	7 <b>14525</b>	116
14224	101	<b>14317</b> 11	7 <b>14527</b>	115
14225	101	<b>14319</b> 11	7 <b>14552</b>	118
14226	101	<b>14326</b> 11	7 <b>14553</b>	118
14227	101	<b>14327</b> 11		
14228	101	<b>14328</b> 11	7 <b>14555</b>	118

Code	Page #	Code	P	age #	Code	P	age #
14601		15993			16735		127
14602		15994			16737		127
14603		16301			16738		127
14642		16305			16739		127
14643	114	16306			16740		127
14644		16310			16741		127
14645 14646	114	16311 16312			16742 16748		127 127
14681		16313			16755		127
14682					16756		128
14691	179	16321			16757		128
14692		16322		134	16760		128
14717					16761		128
15480		16325			16762		128
15484					16763		128
15521		16327			16764		128
15522	141	16330		134	16765		128
15523	141	16331		134	16766		128
15524	141	16333		134	16767		128
15531	141	16334		134	16768		128
15532	141	16335		134	16770		128
15533		16341		135	16771		128
15534	141	16342		135	16772		128
15606	142				16773		128
15608		16345			16778		128
15623					16801		138
15624		16347			16805		138
15633	142	16350			16810		138 139
15634 15695	142	16351 16353			16821 16822		139
15696	143	16354			16823		139
15697		16355		135	16825		139
15698	143	16471			16826		139
15821	141	16472			16827		139
15860	147	16481		107	16830		139
15869	147	16483		107	16831		139
15870	147	16484		107	16832		139
15871	147	16485		107	16833		139
15874	147	16486		107	16835		139
15875	147	16487		107	16836		139
15876	147	16488		107	16837		139
	147	16491		107			
				107	16841		139
15955		16494		107	16842		139
15956		16705		126	16855		140
		16706			16856		140
15975		16715			16857		140
15976 15981				126 127	16860 16861		140 140
15982		16726		127	16862		140
15983		16727			16863		140
15984		16730			16865		140
15985		16731			16866		140
15986					16867		140
15991	142				16870		140
15992	142	16734		127	16871		140

110         110         111         111         112           16872         140         18403         140         18710         178           16903         80         18405         140         18711         178           16913         80         18501         81         18713         178           16914         80         18510         81         18720         176           16923         80         18513         81         18721         177           16924         80         18520         81         18725         177           16924         80         18520         81         18731         177           16934         80         18523         81         18731         177           16934         80         18524         81         18747         176           16941         80         18523         81         19001         164           16953         80         18601         83         19010         164           16954         80         18603         31         19010         164           16964         80         18603         31         19041 <td< th=""><th>Code Page#</th><th>Code Page #</th><th>Code Page#</th></td<>	Code Page#	Code Page #	Code Page#
16901       .79       18404       .140       18711       .178         16903       .80       18405       .140       18712       .178         16911       .79       18500       .81       18712       .176         16914       .80       18510       .81       18720       .176         16914       .80       18510       .81       18722       .177         16921       .79       18512       .81       18724       .177         16924       .80       18520       .81       18725       .177         16933       .80       18520       .81       18731       .177         16934       .80       18523       .81       18731       .177         16934       .80       18524       .81       18731       .177         16941       .80       18523       .81       19001       .164         16951       .80       18530       .81       19010       .164         16954       .80       18600       .83       19011       .164         16955       .80       18601       .83       19021       .163         16964       .80       18612			
16903       80       18405       140       18712       178         16911       79       18500       81       18713       178         16913       80       18510       81       18721       176         16914       80       18511       81       18722       177         16921       79       18512       81       18723       177         16923       80       18513       81       18725       177         16931       79       18521       81       18735       177         16933       80       18522       81       18735       176         16941       80       18523       81       19701       164         16954       80       18530       81       19001       164         16955       80       18601       83       19020       163         16954       80       18602       83       19021       163         16955       80       18602       83       19021       163         16961       80       18602       83       19041       89         16973       80       18612       83       19042 <t< th=""><th></th><th></th><th></th></t<>			
16911       .79       18500       .81       18713       .178         16913       .80       18501       .81       18720       .176         16914       .80       18510       .81       18721       .177         16923       .80       18513       .81       18722       .177         16924       .80       18520       .81       18730       .177         16934       .80       18523       .81       18730       .177         16934       .80       18523       .81       18730       .177         16934       .80       18523       .81       18730       .177         16934       .80       18523       .81       18730       .177         16934       .80       18530       .81       19001       .164         16951       .80       18600       .83       19010       .164         16954       .80       18601       .83       19020       .163         16964       .80       18610       .83       19040       .89         16971       .80       18613       .83       19042       .89         18170       .36       18620			
16913       80       18501       81       18720       176         16914       80       18510       81       18721       176         16915       80       18513       81       18723       177         16924       80       18520       81       18723       177         16933       79       18512       81       18730       177         16934       80       18523       81       18731       177         16934       80       18524       81       18735       176         16941       80       18530       81       19001       164         16953       80       18531       81       19004       143         16954       80       18600       83       19011       164         16955       80       18601       83       19020       163         16961       80       18602       83       19020       163         16971       80       18611       83       19040       89         16973       80       18612       83       19042       89         16974       80       18613       83       19042			
16914       80       18510       81       18721       176         16915       80       18511       81       18723       177         16921       79       18512       81       18723       177         16923       80       18513       81       18725       177         16931       79       18521       81       18735       176         16934       80       18523       81       18735       176         16941       80       18524       81       18737       176         16942       80       18530       81       19001       164         16951       80       18531       81       19004       143         16953       80       18601       83       19010       164         16954       80       18602       83       19014       69         16954       80       18601       83       19040       89         16954       80       18613       83       19042       89         16974       80       18613       83       19042       89         16974       80       18613       83       10042       12			
16915       80       18511       81       18722       177         16921       79       18512       81       18723       177         16924       80       18513       81       18724       177         16933       80       18520       81       18730       177         16933       80       18523       81       18731       177         16934       80       18523       81       18735       176         16941       80       18530       81       19004       143         16953       80       18531       81       19004       164         16954       80       18600       83       19010       164         16955       80       18602       83       19021       163         16956       80       18602       83       19040       89         16971       80       18613       83       19041       89         16973       80       18613       83       19041       89         16974       80       18613       83       19051       89         18170       136       18620       84       19051       8			
16921       .79       18512       .81       18723       .177         16923       .80       18513       .81       18724       .177         16931       .79       18521       .81       18730       .177         16933       .80       18522       .81       18730       .177         16934       .80       18523       .81       18735       .176         16941       .80       18524       .81       18747       .176         16953       .80       18530       .81       19001       .164         16954       .80       18503       .81       19004       .143         16955       .80       18601       .83       19011       .164         16954       .80       18602       .83       19020       .163         16961       .80       18601       .83       19020       .163         16964       .80       18612       .83       19040       .89         16973       .80       18612       .83       19042       .89         16973       .80       18613       .83       19050       .89         18170       .136       18622			
16923       80       18513       81       18724       177         16934       80       18520       81       18730       177         16933       80       18522       81       18731       177         16934       80       18522       81       18735       176         16941       80       18524       81       18747       176         16942       80       18530       81       19001       164         16955       80       18532       81       19010       164         16955       80       18601       83       19011       164         16955       80       18602       83       19013       164         16956       80       18603       83       19024       89         169673       80       18613       83       19040       89         16974       80       18613       83       19042       89         16974       80       18613       83       19042       89         18170       136       18620       84       19057       70.171         18172       136       18621       83       21004			
16924       80       18520       81       18725       177         16931       79       18521       81       18730       177         16933       80       18522       81       18731       177         16934       80       18523       81       18735       176         16941       80       18523       81       18747       176         16942       80       18530       81       19001       164         16951       80       18532       81       19010       164         16953       80       18601       83       19020       163         16955       80       18601       83       19021       163         169563       80       18602       83       19030       164         16964       80       18610       83       19040       89         16971       80       18612       83       19042       89         16974       80       18613       83       19042       89         1870       136       18622       84       19075       70, 171         18173       136       18623       83       21002			
16931			
16933       80       18522       81       18731       177         16934       80       18523       81       18735       176         16941       80       18523       81       1973       176         16942       80       18530       81       19001       164         16951       80       18531       81       19010       164         16954       80       18600       83       19011       164         16955       80       18601       83       19020       163         16964       80       18602       83       19040       89         16973       80       18613       83       19040       89         16974       80       18613       83       19042       89         16974       80       18613       83       19052       89         18170       136       18620       84       19052       89         18171       136       18621       83       21004       120         18172       136       18631       83       21004       121         18172       136       18622       84       19075			
16934       80       18523       81       18735       176         16941       80       18524       81       18747       176         16942       80       18530       81       19001       164         16951       80       18531       81       19004       143         16953       80       18532       81       19010       164         16955       80       18601       83       19020       163         16961       80       18602       83       19020       163         16964       80       18610       83       19040       89         16974       80       18612       83       19041       89         16974       80       18612       83       19042       89         16974       80       18613       83       19050       89         18152       136       18620       84       19052       89         18170       136       18621       84       19052       89         18171       136       18631       83       21004       120         18180       136       18632       83       21001			
16941       80       18524       81       18747       176         16942       80       18530       81       19001       164         16951       80       18532       81       19004       143         16953       80       18532       81       19010       164         16954       80       18600       83       19020       163         16961       80       18602       83       19020       163         16963       80       18602       83       19020       163         16964       80       18613       83       19040       89         16971       80       18612       83       19040       89         16974       80       18612       83       19042       89         16974       80       18613       83       19051       89         18152       136       18622       84       19052       89         18170       136       18621       83       21002       120         18173       136       18631       83       21002       120         18173       136       18642       83       21005 <td< th=""><th></th><th></th><th></th></td<>			
16942       80       18530       81       19001       164         16951       80       18531       81       19004       143         16953       80       18532       81       19010       164         16954       80       18600       83       19011       164         16955       80       18601       83       19020       163         16961       80       18602       83       19030       164         16963       80       18603       83       19040       89         16973       80       18611       83       19042       89         16974       80       18612       83       19050       89         18170       136       18620       84       19051       89         18171       136       18621       84       19075       70, 171         18172       136       18631       83       21004       120         18181       136       18632       83       21004       121         18182       136       18641       83       21005       159         18302       131       18652       109       21010			
16951       80       18531       81       19004       143         16953       80       18532       81       19010       164         16954       80       18600       83       19011       164         16955       80       18601       83       19020       163         16961       80       18602       83       19020       163         16964       80       18610       83       19020       89         16971       80       18611       83       19040       89         16973       80       18612       83       19042       89         16974       80       18613       83       19050       89         18170       136       18620       84       19057       70,171         18172       136       18630       83       21001       120         18180       136       18631       83       21002       120         18181       136       18632       83       21004       121         18182       136       18632       83       21004       121         18183       136       18652       109       21007			
16954       80       18600       83       19011       164         16955       80       18601       83       19020       163         16961       80       18603       83       19020       163         16963       80       18603       83       19020       163         16964       80       18610       83       19040       89         16971       80       18612       83       19042       89         16974       80       18613       83       19050       89         18152       136       18620       84       19051       89         18170       136       18621       84       19075       70,171         18172       136       18622       84       19075       70,171         18173       136       18630       83       21004       120         18180       136       18630       83       21004       121         18182       136       18640       83       21005       159         18305       131       18650       109       21007       159         18305       131       18651       109       21010 </th <th></th> <th></th> <th></th>			
16955       80       18601       83       19020       163         16961       80       18602       83       19021       163         16963       80       18603       83       19030       164         16964       80       18610       83       19040       89         16971       80       18612       83       19042       89         16973       80       18613       83       19051       89         18152       136       18620       84       19052       89         18170       136       18621       84       19057       70,171         18172       136       18630       83       21001       120         18173       136       18631       83       21002       120         18180       136       18632       83       21004       121         18182       136       18641       83       21005       159         18301       131       18650       109       21007       159         18302       131       18651       109       21004       121         18303       131       18652       109       21007 <th><b>16953</b> 80</th> <th><b>18532</b> 81</th> <th><b>19010</b> 164</th>	<b>16953</b> 80	<b>18532</b> 81	<b>19010</b> 164
16961       80       18602       83       19021       163         16963       80       18603       83       19030       164         16964       80       18610       83       19040       89         16971       80       18611       83       19040       89         16973       80       18612       83       19042       89         16974       80       18613       83       19050       89         18172       136       18620       84       19051       89         18171       136       18621       84       19052       89         18172       136       18631       83       21001       120         18180       136       18632       83       21002       120         18181       136       18640       83       21004       121         18182       136       18641       83       21005       159         18301       131       18652       109       21007       159         18302       131       18651       109       21008       103         18312       131       18655       109       21011	<b>16954</b> 80	<b>18600</b> 83	<b>19011</b> 164
16963       80       18603       83       19030       164         16964       80       18610       83       19040       89         16971       80       18611       83       19041       89         16973       80       18612       83       19042       89         16974       80       18613       83       19050       89         18152       136       18620       84       19051       89         18170       136       18622       84       19052       89         18171       136       18622       84       19075       70,171         18172       136       18630       83       21001       120         18173       136       18631       83       21002       120         18180       136       18641       83       21004       121         18182       136       18641       83       21005       159         18302       131       18651       109       21007       159         18302       131       18651       109       21010       119         18312       131       18652       109       21011 <th><b>16955</b> 80</th> <th>18601 83</th> <th><b>19020</b> 163</th>	<b>16955</b> 80	18601 83	<b>19020</b> 163
16964       80       18610       83       19040       89         16971       80       18611       83       19041       89         16973       80       18612       83       19042       89         16974       80       18613       83       19050       89         18152       136       18620       84       19051       89         18170       136       18622       84       19075       70,171         18172       136       18630       83       21001       120         18173       136       18631       83       21002       120         18180       136       18632       83       21004       121         18182       136       18641       83       21005       159         18301       131       18642       83       21006       159         18302       131       18650       109       21007       159         18302       131       18651       109       21010       119         18312       131       18652       109       21011       120         18314       131       18655       109       21012<	<b>16961</b> 80	<b>18602</b> 83	<b>19021</b> 163
16971       80       18611       83       19041       89         16973       80       18612       83       19042       89         16974       80       18613       83       19050       89         18152       136       18620       84       19051       89         18170       136       18621       84       19055       70, 171         18172       136       18622       84       19075       70, 171         18173       136       18630       83       21001       120         18180       136       18632       83       21002       120         18180       136       18640       83       21004       121         18182       136       18641       83       21005       159         18301       131       18650       109       21007       159         18302       131       18651       109       21007       159         18305       131       18652       109       21010       119         18312       131       18653       109       21011       120         18314       131       18656       109 <t< th=""><th><b>16963</b> 80</th><th><b>18603</b> 83</th><th><b>19030</b> 164</th></t<>	<b>16963</b> 80	<b>18603</b> 83	<b>19030</b> 164
16973       80       18612       83       19042       89         16974       80       18613       83       19050       89         18152       136       18620       84       19051       89         18170       136       18621       84       19052       89         18171       136       18622       84       19075       70,171         18172       136       18630       83       21001       120         18173       136       18631       83       21002       120         18180       136       18632       83       21003       120         18181       136       18640       83       21004       121         18182       136       18641       83       21005       159         18301       131       18650       109       21007       159         18302       131       18651       109       21007       159         18303       131       18652       109       21010       119         18314       131       18655       109       21011       120         18315       131       18656       129       32	<b>16964</b> 80	<b>18610</b> 83	<b>19040</b> 89
16974       80       18613       83       19050       89         18152       136       18620       84       19051       89         18170       136       18621       84       19052       89         18171       136       18622       84       19075       70,171         18172       136       18630       83       21001       120         18173       136       18631       83       21002       120         18180       136       18632       83       21003       120         18181       136       18640       83       21004       121         18182       136       18641       83       21005       159         18301       131       18650       109       21007       159         18302       131       18651       109       21007       159         18303       131       18652       109       21007       159         18304       131       18655       109       21010       119         18313       131       18656       109       21011       120         18314       131       18656       129 <td< th=""><th><b>16971</b> 80</th><th><b>18611</b> 83</th><th><b>19041</b> 89</th></td<>	<b>16971</b> 80	<b>18611</b> 83	<b>19041</b> 89
18152       136       18620       84       19051       89         18170       136       18621       84       19052       89         18171       136       18622       84       19075       70,171         18172       136       18630       83       21001       120         18173       136       18631       83       21002       120         18180       136       18632       83       21003       120         18181       136       18640       83       21004       121         18182       136       18641       83       21005       159         18301       131       18650       109       21007       159         18302       131       18651       109       21007       159         18303       131       18652       109       21008       103         18312       131       18655       109       21010       119         18313       131       18656       109       21011       120         18314       131       18656       109       21013       119         18320       132       18657       109	<b>16973</b> 80	<b>18612</b> 83	<b>19042</b> 89
18170       136       18621       84       19052       89         18171       136       18622       84       19075       70, 171         18172       136       18630       83       21001       120         18173       136       18631       83       21002       120         18180       136       18632       83       21003       120         18181       136       18640       83       21004       121         18182       136       18641       83       21005       159         18301       131       18642       83       21007       159         18302       131       18650       109       21007       159         18303       131       18651       109       21007       159         18304       131       18652       109       21008       103         18312       131       18653       109       21010       119         18314       131       18656       109       21011       120         18315       131       18656       109       21012       119         18320       132       18657       109	<b>16974</b> 80	18613 83	<b>19050</b> 89
18171       136       18622       84       19075       70, 171         18172       136       18630       83       21001       120         18173       136       18631       83       21002       120         18180       136       18632       83       21003       120         18181       136       18640       83       21004       121         18182       136       18641       83       21005       159         18301       131       18642       83       21006       159         18302       131       18650       109       21007       159         18303       131       18651       109       21008       103         18312       131       18652       109       21010       119         18313       131       18655       109       21011       120         18314       131       18656       109       21012       119         18320       132       18657       109       21013       119         18321       131       18665       129       32105       70         18323       131       18665       129	<b>18152</b> 136	<b>18620</b> 84	<b>19051</b> 89
18172       136       18630       83       21001       120         18173       136       18631       83       21002       120         18180       136       18632       83       21003       120         18181       136       18640       83       21004       121         18182       136       18641       83       21005       159         18301       131       18642       83       21006       159         18302       131       18650       109       21007       159         18305       131       18651       109       21007       159         18304       131       18652       109       21008       103         18312       131       18651       109       21010       119         18313       131       18655       109       21011       120         18314       131       18656       109       21012       119         18320       132       18657       109       21013       119         18321       131       18656       129       32155       70         18322       131       18665       129	<b>18170</b> 136	<b>18621</b> 84	<b>19052</b> 89
18173       136       18631       83       21002       120         18180       136       18632       83       21003       120         18181       136       18640       83       21004       121         18182       136       18640       83       21005       159         18301       131       18642       83       21006       159         18302       131       18650       109       21007       159         18305       131       18651       109       21008       103         18312       131       18652       109       21009       119         18313       131       18655       109       21010       119         18314       131       18655       109       21011       120         18315       131       18656       109       21012       119         18320       132       18657       109       21013       119         18321       131       18656       129       32155       70         18323       131       18661       129       32155       70         18330       131       18665       129	<b>18171</b> 136		
18180       136       18632       83       21003       120         18181       136       18640       83       21004       121         18182       136       18641       83       21005       159         18301       131       18642       83       21006       159         18302       131       18650       109       21007       159         18305       131       18651       109       21008       103         18312       131       18652       109       21009       119         18313       131       18653       109       21010       119         18314       131       18655       109       21011       120         18315       131       18656       109       21012       119         18320       132       18657       109       21013       119         18321       131       18658       109       21365       73         18322       131       18661       129       32105       70         18333       131       18665       129       32105       70         18334       131       18666       129			
18181       136       18640       83       21004       121         18182       136       18641       83       21005       159         18301       131       18642       83       21006       159         18302       131       18650       109       21007       159         18305       131       18651       109       21007       159         18312       131       18652       109       21008       103         18313       131       18653       109       21010       119         18314       131       18655       109       21011       120         18315       131       18656       109       21012       119         18320       132       18657       109       21013       119         18321       131       18658       109       21365       73         18322       131       18661       129       26555       68, 169         18323       131       18665       129       32105       70         18330       131       18665       129       32205       70         18331       131       18665       129			
18182       136       18641       83       21005       159         18301       131       18642       83       21006       159         18302       131       18650       109       21007       159         18305       131       18651       109       21008       103         18312       131       18652       109       21009       119         18313       131       18653       109       21010       119         18314       131       18655       109       21011       120         18315       131       18656       109       21012       119         18320       132       18657       109       21013       119         18321       131       18658       109       21365       73         18322       131       18661       129       26555       68, 169         18323       131       18665       129       32105       70         18330       131       18665       129       32205       70         18331       131       18665       129       3205       71         18332       131       18675       130			
18301       131       18642       83       21006       159         18302       131       18650       109       21007       159         18305       131       18651       109       21008       103         18312       131       18652       109       21009       119         18313       131       18653       109       21010       119         18314       131       18655       109       21011       120         18315       131       18656       109       21012       119         18320       132       18657       109       21013       119         18321       131       18656       109       21365       73         18322       131       18661       129       26555       68, 169         18323       131       18665       129       32105       70         18324       131       18665       129       32205       70         18331       131       18665       129       32205       70         18332       131       18675       130       110072       122         18340       131       18675       130			
18302       131       18650       109       21007       159         18305       131       18651       109       21008       103         18312       131       18652       109       21009       119         18313       131       18653       109       21010       119         18314       131       18655       109       21011       120         18315       131       18656       109       21012       119         18320       132       18657       109       21013       119         18321       131       18658       109       21365       73         18322       131       18660       129       26555       68, 169         18323       131       18661       129       32105       70         18324       131       18665       129       32205       70         18330       131       18666       129       32205       70         18331       131       18675       130       110072       122         18333       131       18675       130       120105       69, 171         18333       131       18676       13			
18305       131       18651       109       21008       103         18312       131       18652       109       21009       119         18313       131       18653       109       21010       119         18314       131       18655       109       21010       119         18315       131       18656       109       21011       120         18320       132       18657       109       21012       119         18321       131       18658       109       21365       73         18322       131       18660       129       26555       68, 169         18323       131       18661       129       32105       70         18324       131       18665       129       32205       70         18330       131       18666       129       32205       70         18331       131       18671       129       67589       71         18332       131       18675       130       110072       122         18340       131       18676       130       120105       69, 171         18341       131       18680       163			
18312       131       18652       109       21009       119         18313       131       18653       109       21010       119         18314       131       18655       109       21011       120         18315       131       18656       109       21012       119         18320       132       18657       109       21013       119         18321       131       18658       109       21365       73         18322       131       18660       129       26555       68, 169         18323       131       18661       129       32105       70         18324       131       18665       129       32205       70         18330       131       18666       129       32205       70         18331       131       18670       129       67589       71         18332       131       18671       129       95050       71         18333       131       18675       130       110072       122         18340       132       18676       130       120105       69, 171         18341       131       18680       163<			
18313       131       18653       109       21010       119         18314       131       18655       109       21011       120         18315       131       18656       109       21012       119         18320       132       18657       109       21013       119         18321       131       18658       109       21365       73         18322       131       18660       129       26555       68, 169         18323       131       18661       129       32105       70         18324       131       18665       129       32105       70         18330       131       18666       129       32205       70         18331       131       18667       129       67589       71         18332       131       18671       129       95050       71         18333       131       18675       130       110072       122         18340       132       18676       130       120105       69, 171         18341       131       18680       163       12015       70, 171         18342       131       18681			
18314       131       18655       109       21011       120         18315       131       18656       109       21012       119         18320       132       18657       109       21013       119         18321       131       18658       109       21365       73         18322       131       18660       129       26555       68, 169         18323       131       18661       129       32105       70         18324       131       18665       129       32155       70         18330       131       18666       129       32205       70         18331       131       18670       129       67589       71         18332       131       18671       129       95050       71         18333       131       18675       130       110072       122         18340       132       18676       130       120105       69, 171         18341       131       18680       163       120115       70, 171         18342       131       18681       163       120205       70, 171         18343       131       18690			
18315       131       18656       109       21012       119         18320       132       18657       109       21013       119         18321       131       18658       109       21365       73         18322       131       18660       129       26555       68, 169         18323       131       18661       129       32105       70         18324       131       18665       129       32155       70         18330       131       18666       129       32205       70         18331       131       18670       129       67589       71         18332       131       18671       129       95050       71         18333       131       18675       130       110072       122         18340       132       18676       130       120105       69, 171         18341       131       18680       163       120115       70, 171         18342       131       18681       163       120205       70, 171         18343       131       18690       163       120215       70, 171         18344       131       18691			
18320			
18321       131       18658       109       21365       73         18322       131       18660       129       26555       68, 169         18323       131       18661       129       32105       70         18324       131       18665       129       32105       70         18324       131       18665       129       32205       70         18330       131       18666       129       32205       70         18331       131       18670       129       67589       71         18332       131       18671       129       95050       71         18333       131       18675       130       110072       122         18340       132       18676       130       120105       69,171         18341       131       18680       163       12015       70,171         18342       131       18681       163       120205       70,171         18343       131       18690       163       120215       70,171         18344       131       18691       163       120307       68,170         18401       140       18700			
18323       131       18661       129       32105       70         18324       131       18665       129       32155       70         18330       131       18666       129       32205       70         18331       131       18670       129       67589       71         18332       131       18671       129       95050       71         18333       131       18675       130       110072       122         18340       132       18676       130       120105       69,171         18341       131       18680       163       12015       70,171         18342       131       18681       163       120205       70,171         18343       131       18690       163       120205       70,171         18344       131       18691       163       120205       70,171         18344       131       18691       163       120205       70,171         18344       131       18691       163       120307       68,170         18401       140       18700       175       120427       68,170	<b>18321</b> 131	<b>18658</b> 109	<b>21365</b> 73
18323        131       18661        129       32105        70         18324        131       18665        129       32155        70         18330        131       18666        129       32205        70         18331        131       18670	<b>18322</b> 131	<b>18660</b> 129	<b>26555</b> 68, 169
18330       131       18666       129       32205       70         18331       131       18670       129       67589       71         18332       131       18671       129       95050       71         18333       131       18675       130       110072       122         18340       132       18676       130       120105       69,171         18342       131       18680       163       120115       70,171         18342       131       18681       163       120205       70,171         18343       131       18690       163       120215       70,171         18344       131       18691       163       120307       68,170         18401       140       18700       175       120427       68,170	<b>18323</b> 131		
18331       131       18670       129       67589       71         18332       131       18671       129       95050       71         18333       131       18675       130       110072       122         18340       132       18676       130       120105       69,171         18342       131       18680       163       120115       70,171         18342       131       18681       163       120205       70,171         18343       131       18690       163       120215       70,171         18344       131       18691       163       120307       68,170         18401       140       18700       175       120427       68,170	<b>18324</b> 131	<b>18665</b> 129	<b>32155</b> 70
18332        131       18671        129       95050        71         18333        131       18675        130       110072        122         18340        132       18676        130       120105       69, 171         18341        131       18680        163       120115       70, 171         18342        131       18681        163       120205       70, 171         18343        131       18690        163       120215       70, 171         18344        131       18691        163       120307       68, 170         18401        140       18700	<b>18330</b> 131	<b>18666</b> 129	<b>32205</b> 70
18333        131       18675        130       110072        122         18340        132       18676        130       120105       69, 171         18341        131       18680        163       120115       70, 171         18343        131       18681        163       120205       70, 171         18343        131       18690        163       120215       70, 171         18344        131       18690        163       120215       70, 171         18344        131       18691        163       120307       68, 170         18401        140       18700        175       120427       68, 170	<b>18331</b> 131	<b>18670</b> 129	<b>67589</b> 71
18340          132         18676          130         120105         69, 171           18341          131         18680          163         120115         70, 171           18342          131         18681          163         120205         70, 171           18343          131         18690          163         120215         70, 171           18344          131         18691          163         120307         68, 170           18401          140         18700          175         120427         68, 170	<b>18332</b> 131	<b>18671</b> 129	<b>95050</b> 71
18341        131       18680        163       120115       70, 171         18342        131       18681        163       120205       70, 171         18343        131       18690        163       120215       70, 171         18344        131       18691        163       120307       68, 170         18401        140       18700        175       120427       68, 170	<b>18333</b> 131	<b>18675</b> 130	<b>110072</b> 122
18342        131       18681        163       120205       70, 171         18343        131       18690        163       120215       70, 171         18344        131       18691        163       120307       68, 170         18401        140       18700        175       120427       68, 170	<b>18340</b> 132	<b>18676</b> 130	
18343          131         18690          163         120215         70, 171           18344          131         18691          163         120307         68, 170           18401          140         18700          175         120427         68, 170	<b>18341</b> 131	<b>18680</b> 163	
18344          131         18691          163         120307         68, 170           18401          140         18700          175         120427         68, 170			<b>120205</b> 70, 171
<b>18401</b> 140 <b>18700</b> 175 <b>120427</b> 68, 170			
<b>18402</b> 140 <b>18705</b> 175 <b>120700</b> . 67, 169			
	<b>18402</b> 140	<b>18705</b> 175	<b>120700</b> . 67, 169

Index

Code		P	age #
12072	1 .	67.	169
12074			
		,	
Skin ( 0089	Ca		59
0283			
0203			
0402			
0400			
0407			
0408			
• • • •			
0505			
0506			
1005			
1006			
1300			
1337			
1607			
1611			
1622			
1644			
1645			
1802			
1814			
1843	•••••		
1873	•••••		
1877	••••		
1880	••••		
1944	••••		
1947	••••		
7015	••••		
7041	••••		
7055	••••		
7056	••••		
7057	••••		
7058	••••		
7067			53
7068			53
7069			53
7090			
7091			53
7092			53
7095			53
7229			52
7230			52
7231			52
7233			52
7234			52
7235			52
7251			59
7262			59
7300			50
7301			50
7565			54
7566			54
7567			54

Code Po	age #	Code Page#
<b>120721</b> . 67,	169	<b>7570</b> 56
<b>120740</b> . 67,		<b>7571</b> 56
		<b>7572</b> 56
Skin Care	. 59	<b>7583</b> 59
0089 0283		<b>7710</b> 55
0203		<b>7712</b> 55
0406		<b>7725</b> 55
0407		<b>7910</b> 57
0408		<b>7912</b> 57
0424		<b>7915</b> 57
0505		<b>61760</b> 50
0506		<b>61761</b> 50
1005		<b>61762</b> 50
1006		<b>61768</b> 50
1300		<b>67143</b> 50
1337		<b>67145</b> 50
1607		<b>67146</b> 50
1611		<b>67147</b> 50
1622		<b>67589</b> 59
1644		Wound Care
1645	. 58	<b>1061</b> 42
1802	. 53	<b>1063</b> 42
1814	. 53	<b>1166</b> 42
1843	. 53	<b>1964</b> 42
1873	. 55	<b>1967</b> 42
1877	. 55	<b>3410</b> 35
1880	. 55	<b>3413</b> 35
1944	. 54	<b>3416</b> 35
1947	. 54	<b>3420</b> 36
7015	. 54	<b>3423</b> 36
7041	. 58	<b>3430</b> 36
7055	. 51	<b>3444</b> 35
7056	. 51	<b>3470</b> 36
7057		<b>3473</b> 36
7058		<b>3475</b> 36
7067		<b>3485</b> 36
7068		<b>3488</b> 36
7069		<b>3705</b> 38
7090		<b>3710</b> 38
7091		<b>3715</b> 38
7092		<b>3740</b> 38
7095	. 53	<b>3755</b> 39
7229		<b>3760</b> 39
7230 7231		<b>3765</b> 39
7231 7233		<b>3780</b> 39 <b>3900</b> 42
7233		<b>3900</b> 42 <b>3903</b> 42
7235		<b>3906</b> 42
7251		<b>6105</b> 35
7262		<b>7690</b> 42
7300		<b>9622</b> 37
7301		<b>9625</b> 37
7565		<b>9632</b> 37
7566		<b>9635</b> 37
7567	. 54	<b>9641</b> 37

Code	Ρα	ge #
9643 .		37
33110		40
33115		40
33120		40
33146		40
33280		41
33283		41
33285		41
33400		35
33401		35
33406		35
33434		35
33435		35
33436		35
33437		35
33438		35
33444		35
33445		35
33446		35
33452		35
33453		35
33530		40
33533		40
33536		40
33539		40
33560		43
33561		43
33562		43
33563		43
33564		43
39000		35
39001		35
39010		35
46100		44
46120 46150		44 44
46150		44
46250		44
46250		44
46320		44
46350		44
46390		44
46450		44
.0 /50		. –

## Making life **easier**

From getting to the grocery store to biking up a mountain, users' hopes and aspirations fuel our passion of finding new ways to help make their lives easier.

We constantly question the status quo across every aspect of intimate healthcare and seek to understand our users, to better anticipate and advocate for their needs. Ever since the founding of Coloplast in 1957, we have worked closely with healthcare professionals and the users of our products to discover new solutions as well as new ways to create solutions that are sensitive to their special needs. Together, we move care forward.

Our business includes Ostomy Care, Continence Care, Wound & Skin Care and Interventional Urology. We operate globally and employ about 13,000 employees.

For additional information, contact us at:

1-800-533-0464 info-us@coloplast.com www.coloplast.us



Coloplast Corp. Minneapolis, MN 55411 1-800-533-0464 | www.coloplast.us