



Comfort, Health and Style!



NEW!

JOBST® Compriz2

The 1st layer is for patient comfort
The 2nd layer is for consistent reliable compression

- Moisture permeable design provides high level of patient comfort
- Indicator promotes correct and sustained compression 100% of the time
- Effective edema reduction and recurrence prevention
- Works well with patients on previous compression therapies
- Easily integrated into daily clinical and home routine

THERAPIES. HAND IN HAND.

JOBST® Compri 2

Effective Edema Reduction. Easy-to-Apply.

The new JOBST® Compri2 is a 2-Layer Bandage System specifically designed for **easy application**. It provides **graduated and reliable compression** for up to seven days to help achieve **effective edema reduction** and heal the wound while skin stays dry and comfortable.

Easy Application

- 1st layer of padding is customizable and easy to apply
- 2nd compression layer can be re-applied or adjusted as needed
- Indicator ensures correct compression level is achieved

High Patient Comfort

- Padding can be added in sensitive areas
- Excellent mobility increases compliance and satisfaction
- Patients can wear with most footwear



1st Layer - Non Woven Padding

- **Stretchable** - simplifies application
- **Active Moisture Management** - keeps skin dry and comfortable
- **Customized Padding** - allows for anatomical adjustments and additional protection to sensitive areas



For application instructions of 1st layer



2nd Layer - Short Stretch Bandage

- **Cohesive** - delivers consistent and reliable compression
- **Adjustable** - allows modification during application
- **Visual** - indicator offers simple and effective application of correct compression



For application instructions of 2nd layer

Scannable QR codes and easy to read IFU offer precise instructions for achieving correct application, to deliver high working pressure and low resting pressure, 24 hours a day.

ORDERING INFORMATION:

Available in two lengths and compression levels: regular and lite; ensures precise compression to suit most patients.

	Ankle size	REF No.	Kits per box	HCPCS	
Regular compression (ABPI >0.8)					
JOBST® Compri2	7 1/8 - 9 3/4 inch	7627100	1	A6441 / A6452	
JOBST® Compri2	9 3/4 - 12 1/2 inch	7627101	1	A6441 / A6452	
Lite compression (ABPI 0.5 - 0.8)					
JOBST® Compri2 lite	7 1/8 - 9 3/4 inch	7627102	1	A6441 / A6452	
JOBST® Compri2 lite	9 3/4 - 12 1/2 inch	7627103	1	A6441 / A6452	

REIMBURSEMENT DISCLAIMER: This document includes codes that might be used to bill for BSN medical products and related services. Each provider will have to verify the appropriate codes for each patient. It is the provider's sole responsibility to determine and submit appropriate codes, charges, and modifiers for services rendered. Providers should contact insurers to verify correct coding procedures prior to submitting claims related to any product or service. BSN medical Inc. cannot guarantee coverage or reimbursement with the codes listed in this billing guide. In all cases, providers will need to follow local payer policies for billing and reimbursement.

Please refer to the product label and / or package insert for full instructions on the safe use of these products.



BSN medical Inc. 5825 Carnegie Blvd. Charlotte, NC 28209-4633
 Tel. (+1) 704 554 9933 Fax (+1) 800 835 4325 www.bsnmedical.com
 To order toll-free: BSN medical (+1) 800 552 1157



Comfort, Health and Style!