# **Conforming Stretch Gauze Bandages**

McKesson Conforming Stretch Gauze Bandages are designed as a secondary dressing to help secure placement of primary dressings. They are made of a stretchable, soft cushioning which enables them to provide light compression while conforming to most of the body's natural contours. Ideal as a secondary dressing, they allow for movement while accommodating any increase or decrease in swelling. Their high absorbency helps ensure patient comfort during the healing process. McKesson Conforming Stretch Gauze Bandages feature non-frayed finished edges to avoid fraying and linting. They are available in different lengths and widths for caregiver convenience.

### **Features and Benefits**

- Breathable rayon and polyester blend to help provide added comfort, reduced slippage and shape retention
- Self-adhesive to minimize need for tape
- · Available in sterile or non-sterile
- Not made with natural rubber latex



#### ORDERING INFORMATION

DESCRIPTION	MFR #	QTY
Conforming Stretch Gauze Bandages, Non-Sterile, 1" x 1.7yds, Stretched	16-010	24/bg, 4bg/cs
Conforming Stretch Gauze Bandages, Non-Sterile, 2" x 4.1yds, Stretched	16-011	12/bg, 8bg/cs
Conforming Stretch Gauze Bandages, Non-Sterile, 3" x 4.1yds, Stretched	16-012	12/bg, 8bg/cs
Conforming Stretch Gauze Bandages, Non-Sterile, 4" x 4.1yds, Stretched	16-013	12/bg, 8bg/cs
Conforming Stretch Gauze Bandages, Non-Sterile, 6" x 4.1yds, Stretched	16-014	6/bg, 8bg/cs
Conforming Stretch Gauze Bandages, Sterile, 1" x 1.7yds, Stretched	16-016	1ea, 24pk/bg, 4bg/cs
Conforming Stretch Gauze Bandages, Sterile, 2" x 4.1yds, Stretched	16-017	1ea, 12pk/bg, 8bg/cs
Conforming Stretch Gauze Bandages, Sterile, 3" x 4.1yds, Stretched	16-018	1ea, 12pk/bg, 8bg/cs
Conforming Stretch Gauze Bandages, Sterile, 4" x 4.1yds, Stretched	16-019	1ea, 12pk/bg, 8bg/cs
Conforming Stretch Gauze Bandages, Sterile, 6" x 4.1yds, Stretched	16-020	1ea, 6pk/bg, 8bg/cs

## Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

#### **Never Compromise on Care SM**

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them — making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

© 2016 McKesson Medical-Surgical Inc.