

NEW! aXs2 Pool Lift



The aXs2 is a low profile pool lift that is an ideal access solution for community and hospitality swimming pools.

- Third-party tested & verified ADA compliant
- 360° rotation for transfer safety and flexibility
- Compact design requires minimal deck space
- Integrated armrests
- Sturdy and generous rotomolded seat
- 300 lb/136kg lifting capacity
- LiftOperator® Intelligent Controller with 24V battery
- Full range of operation from handset and control box
- Powder-coated stainless steel and aluminum construction
- Conveniently ships in two boxes



Lift Color



GRAY MIST

Due to printing technology actual color may differ

LiftLock

The LiftLock Locking Square Deck Anchor ensures a smooth ride.




Model No.	Description	Shipping (2 boxes)				
		Weight	Length	Width	Height	
310-0000	aXs2 Pool Lift	Box 1	102 lbs	55"	19"	14"
			46kg	140cm	48cm	36cm
310-0000N	aXs2 Pool Lift without anchor	Box 2	36 lbs*	24"	24"	27"
			16kg	61cm	61cm	69cm

Notes:

- Model 310-0000 includes battery, charger, battery console cover, water resistant hand control, LiftLock stainless steel anchor socket with cover, footrest, armrests and seat belt assembly
- Model 310-0000N comes without LiftLock Anchor, and weight of Box 2 is 34 lbs/15kg.

Parts & Accessories:

- 1001495 LiftOperator Battery
- 300-6200L LiftLock Anchor & Cap
- 970-5000 Seat Saver Cover
- 170-3000A Armrest Replacement (pair)
- 900-6000 Stability Strap
- 900-1000 Seat Belt
- 300-6200-2L LiftLock Anchor Replacement Cap



A completed Deck Profile form is required with your pool lift order.

