



Gravity

Gravity set Roller Clamp(s) 2 SmartSite needle-free valve(s) 9", 86" from 2-piece Male Luer Lock. 60 Drop L: 104 in L: 264 cm PV: 17 mL Fluid Path Sterile

Summary

SKU	42133E
Product Name	Gravity
Product Type	



Weight (lbs)	12.0000	Country of Manufacture	Mexico	Anti-Siphon Valve	No
Approximate Length (cm)	264	Approximate Length (in)	104	Approximate Priming Volume (mL)	17
Bag Access	No	Ball Valve in Drip Chamber	No	Blood Set	No
Burette	No	Case Dimensions	21.75in x 15in x 10in 12lb	Clamp	Roller
Drop	60	Flow Regulator	No	Natural Rubber Latex is not part of the material formulation	Yes
Low Sorbing	No	Luer Lock	2-piece Male Luer Lock	Manifold	None
Number of Back Check Valves	0	Number of Connectors	0	Number of Fuses	0
Pressure Rated	No	Pressure Sensing Disc	No	Sterility	Fluid Path Sterile
Stopcock	None	Texium	No	Color of Tubing	None
T-Connector Type	None	Number of y-sites	2	Y-Site Port Type	SmartSite
Type of Connector	None	Filter	None	Tubing Size	Standardbore
CE Marking	No	PVC Free	No	Quantity	100
Unit Measure	Per case				