



DermaFilm™ Hydrocolloid Dressing

DermaFilm is a hydrocolloid wound dressing for non-infected wounds and shallow pressure ulcers. DermaFilm maintains a moist wound environment. Available in a variety of convenient forms.

Benefits:

- Mild adhesive surface area adheres to periwound skin, but won't stick to wound bed
- Supports moist wound healing & autolytic debridement
- Provides thermal insulation and protection
- Comfortable and conformable
- Can be safely worn for several days, depending on wound condition and exudate levels
- Latex free

Indications:

May be used as a primary or secondary dressing. For the management of low exudating wounds, including:

- non-infected wounds
- partial thickness ulcers
- full-thickness ulcers

Product	Size	Reorder #	Units
DermaFilm Thin	2x2	00229E	20/Box
DermaFilm Thin	4x4	00259E	10/Box
DermaFilm Thin	6x6	00306E	5/Box
DermaFilm HD (reg.)	4x4	00239E	5/Box

Product (cont.)	Size	Reorder #	Units
DermaFilm Sacral	6x7	00279E	5/Box
DermaFilm Heel/Elbow	3.5x5	00299E	10/Box
DermaFilm X–Thin Clear	4x4	00219E	10/Box
DermaFilm X–Thin Clear	6x6	00312E	5/Box

