



Ostomy & Wound Care Product Guide



Introduction...

Safe n Simple™ is a medical care innovator, specializing in the promotion of healthy skin, complete security and confident ostomy and wound care management. We view our product development through the eyes of our customers, who want active, productive lives unimpeded by their conditions. Our goal is to make our customer's ostomy and wound care safe and simple.

We follow two overriding principles: create innovative, high-performance products, and deliver the greatest possible value. We continually push our technicians to develop better ideas, more advanced formulations and superior products, while at the same time making those products convenient and affordable.

We've spent years in **Ostomy care** learning from our customers what they really want. The result is a complete line of ostomate care products designed to *Cleanse, Protect, Secure, and Assure.*

And, a complete line of wound care products created to promote healing and promote cost savings.

Above all, our customers want safety and simplicity. With Safe n Simple, that's our promise.

Table of Contents

Ostomy Care

Page 1	Introduction
Page 2-4	Ostomy & Wound Care - Table of Contents
Page 5	Peri-Stoma Cleanser & Adhesive Remover
.....	Adhesive Remover
Page 6	No-Sting Skin Barrier Film
.....	No-Sting Skin Barrier Wand
.....	No-Sting Skin Barrier Spray
Page 7	Skin Barrier Film
.....	Ostomy Skin Barrier Powder
Page 8	Skin Barrier Ring Conforming Seal
.....	Integrity Skin Barrier Seal
Page 9	Skin Barrier Arcs
.....	1/2 Arc Skin Barrier Strips
.....	Skin Barrier Sheet
Page 10.....	No-Sting Skin Barrier Paste
.....	Skin Barrier Paste
.....	Skin Barrier Paste-Fast Drying
Page 11.....	Assure C Odor Eliminator
.....	Ostomy Pouch Deodorant
Page 12.....	Pouches
.....	Stoma Caps
Page 13.....	Medicare Allowables
Page 14-16.....	Security Support Belts

Table of Contents

Wound Care

Page17 & 18.....Wound Care Intro Page

Page 19.....Alginate
.....Silver Alginate

Page 20.....Fibergel
.....Silver Fibergel

Page 21.....Collagen
.....Silver Collagen

Page 22.....Borderless Foam
.....Bordered Foam

Page 23.....Ionic Silver Foam
.....Borderless Silicone Foam

Page 24.....Bordered Silicone Foam
.....Borderless Silver Silicone Foam

Page 25.....Gauze Sponges
.....Bordered Gauze Sponges

Page 26.....Hydrocolloids
.....Polymer Pads

Page 27.....Hydrogel Impregnated Gauze
.....Hydrogel Sheets

Page 28.....DermaPro Waterproof Plastic Tape
.....DermaPro Waterproof Silicone Tape

Page 29 & 30.....Wound Care Kits

Page 31.....Wound Care Allowables

Special Information

Page 13..... Ostomy Accessory Allowables

Page 31..... Wound Care Allowables

Page 32-35 Stay Connected with SNS



Safe n Simple provides focused care that you can trust in. Our customer service staff are dedicated specialists who are trained to support you with product recommendations to help you improve your quality of life.

We are proud to offer our SimpleStart Ostomy Care Kits to give you an opportunity to trial our ostomy accessory product line.

Call our toll free customer service line to request a kit:

844-767-6334

Email Us: simplestart@sns-medical.com

or fax us at: 248-575-4172

Cleanse

Peri-Stoma Cleanser & Adhesive Remover

**ALCOHOL
FREE!**



Use these soft wipes to breakdown and remove adhesives, dirt and oils from peri-stomal skin. Does not leave an oily residue. No-sting, alcohol-free, and no-rinse formulation. Gentle, convenient and scent-free. Convenient re-sealable pouch.



Item	Product#	Size	HCPC
Large Wipes, 50 Pack	SNS00525	5"x7"	A4456
Large Individual Wipes 5 Box	SNS00505	5"x7"	A4456
Large Individual Wipes 50 Box	SNS00550	5"x7"	A4456
Large Individual Wipes 75 Box	SNS00575	5"x7"	A4456

Customer Testimonial:

"I had some red and irritated skin using other wipes and when I received the Safe n Simple wipes from my medical supply co., they gave me instant relief and the skin healed. I like also, that they are big and so soft. I have nothing but praise for them and I will always use them, nothing compares to them." Sincerely, Hildegard.

Adhesive Remover



If a more aggressive adhesive remover is needed, then this is the wipe to use. The alcohol based wipe has a very gentle formula to clean and remove debris from the skin. Rinse is required.

Item	Product#	Size	HCPC
Individual Wipes 50 Box	SNS00650	2"x2"	A4456

No-Sting Skin Barrier Film

**ALCOHOL
FREE!**



This wipe forms a barrier to protect the skin from adhesives, friction and body fluids. Increases adhesion of tapes and wafers. Gentle, convenient and scent-free.

Item	Product#	Size	HCPC
Large Wipes, 25 Box	SNS00807	5"x7"	A5120
Individual Wipes, 25 Box	SNS80725	2"x2"	A5120
Individual Wipes, 75 Box	SNS80775	2"x2"	A5120
Individual Wipes, 100 Box	SNS80744	2"x2"	A5120

No-Sting Skin Barrier Wand

**ALCOHOL
FREE!**



Applicator wand provides a no-sting, alcohol-free transparent protective coating to the area applied. Increases adhesion of tapes and wafers. Conveniently packaged applicator wands are gentle on the skin, sterile and scent-free.

Item	Product#	Size	HCPC
Wands, 25 Box	SNS80711	1mL	A4369

No-Sting Skin Barrier Spray

**ALCOHOL
FREE!**



No-sting, alcohol-free spray that provides a transparent protective coating onto the area sprayed. Increase adhesion of tapes and wafers. Gentle convenient and scent-free. Ideal for use over powder to create a crusting technique to protect traumatized skin.

Item	Product#	Size	HCPC
Spray Bottle	SNS80792	2 oz.	A4369

Skin Barrier Film



Alcohol based wipe has a very gentle formula that dries quickly. Helps to protect the skin against irritation and adhesive buildup.

Item	Product#	Size	HCPC
Individual Wipes, 50 Box	SNS81850	2"x2"	A5120

Ostomy Skin Barrier Powder



Protects and promotes healing of irritated or broken skin. Dusted onto the skin it absorbs moisture and enhances comfort and wear time of the pouching system. Perfect for crusting.

Item	Product#	Size	HCPC
Applicator Bottle	SNS92305	5 oz.	A4371
Applicator Bottle	SNS92301	1 oz.	A4371



Skin Barrier Ring Conforming Seal



The conforming seals are ideal for any type of ostomy. It will create a secure seal at the base of the stoma and is highly absorbent. These rings “Hold Up Better”, promoting longer pouch wear time and do not leave a mess when removed from the skin. They are very moldable and stretchable with a strong tack and require no cutting. Great for replacing paste.

Item	Product#	Size	HCPC
Individual Seals, 10 Box	SNS684U2	2"	A4385
Individual Seals, 10 Box	SNS684U4	4"	A4385

Integrity Skin Barrier Ring



These durable skin barrier rings maintain their integrity and do not melt down from excessive moisture or body heat. Integrity skin barrier rings are ideal for those with a Colostomy that do not need absorption in a ring. These rings provide a comfortable protective seal to prevent leakage. Sticky, thin, low profile seals stretch and mold easily to fit around the stoma.

Item	Product#	Size	HCPC
Individual Seals, 20 Box	SNS68002	2"	A4385



Customer Testimonial:

"I am using three of your products now, Integrity Seals for fill-in areas, Conforming Seals for around my stoma and also your deodorant and have been getting three days longer wear time now over one day wear time from the other products. You've kinda saved my life."
David B.

Skin Barrier Arcs



Under the pouching system: The skin barrier arcs can be placed directly on the irritated skin around the stoma to create a smooth platform to adhere the pouch or wafer. The arcs protect the skin from trauma caused by repetitive removal of tapes and adhesives.

Over the pouching system: The skin barrier arcs can be placed around the wafer to secure the edges from lifting up, creasing or loosening. Use 1-3 arcs around the edges. Skin barrier arcs provide an excellent alternative to waterproof tapes and give you a waterproof seal around the wafer when bathing, showering or swimming.

Item	Product#	Size	HCPC
Tray - Regular, 30 Arcs	SNS20630	1" Width	A4362
Tray - X-tra Wide, 30 Arcs	SNS21130	2" Width	A4362

1/2 ARC Skin Barrier Strips

**NEW SIZE
NOW AVAILABLE!**



A hydrocolloid 1/2 circle arc that is water resistant and perfect for picture framing around the ostomy pouch. Thin profile for comfortable wear. Non-sterile.

Item	Product#	Size	HCPC
Tray - 1/2 Arc, 20 Strips	SNS22222	1.25" Width	A4362

Skin Barrier Sheet



A thin hydrocolloid sheet that protects the skin from stoma output. Highly absorbent. Cut to fit - no starter hole. Non-sterile.

Item	Product#	Size	HCPC
Individual Packaged 5 Box	SNS21605	4"x4"	A4362

No-Sting Skin Barrier Paste



An alcohol-free formulation to reduce the sting and pain caused by conventional alcohol based pastes. Easy to squeeze tube helps you get the paste out smooth and quick. Creates an optimum seal for a better fitting pouch and longer wear time.

Item	Product#	Size	HCPC
Screw Cap Tube	SNS92802	2 oz.	A4406

Skin Barrier Paste



A protective pectin-based paste that helps to prevent leakage and skin irritation around the stoma providing an optimum seal for a better fit and longer wear time. Easy to squeeze tube helps you get the paste out smooth and quick.

Item	Product#	Size	HCPC
Screw cap Tube	SNS90502	2 oz.	A4406

Skin Barrier Paste - Fast Drying



A protective pectin-based paste that hardens upon application. Contains alcohol for faster drying. Seals and prevents skin irritation around the stoma and is great as a filler for those problem areas.

Item	Product#	Size	HCPC
Screw cap Tube	SNS90516	2 oz.	A4406

Ostomy Pouch Deodorant



Fantastic odor control with no perfumes or artificial scents added for discreet protection. Blue formula helps to see the liquid inside the pouch.

Item	Product#	Size	HCPC
Push Top Dispenser Bottle	SNS40202	2 oz.	A4394
Push Top Dispenser Bottle	SNS40208	8 oz.	A4394

Assure C Odor Eliminator



Assure C is an odor eliminator that also contains a lubricant to help empty pouch contents. This clear formulation eliminates the fear or concern of staining associated with the blue colored deodorants. No perfumes or artificial scents added.


Item	Product#	Size	HCPC
Single Serve Packet	SNS41425	1/4 oz.	A4394
Tilt Tip Dispenser Bottle	SNS41404	4 oz.	A4394
Tilt Tip Dispenser Bottle	SNS41408	8 oz.	A4394



**Be Confident
and Secure!**

Pouches

One-Piece Drainable Pouch



Standard wear drainable skin barrier with a comfortable fit panel and dual adhesive tapes. Cut-to-fit, clamp closure. Latex free.

Item	Product#	Size	HCPC
Transparent Pouch 10 Box	SNS12T100	12"	A5061
Opaque Pouch 10 Box	SNS12B100	12"	A5061

Two (2) clips per box!

One-Piece Drainable Pouch with Filter



Standard wear drainable skin barrier with a comfortable fit panel and dual adhesive tapes. Odor controlling filter to reduce gas buildup and ballooning. Cut-to-fit, clamp closure. Latex free.

Item	Product#	Size	HCPC
Transparent Pouch 10 Box	SNS12T10F	12"	A4424
Opaque Pouch 10 Box	SNS12B10F	12"	A4424

Stoma Cap



A one-piece, lightweight flexible system. Easy to manage. A temporary option to assist with showering, swimming or moments of intimacy. Latex free.

Item	Product#	Size	HCPC
Stoma Cap w/Hydrocolloid Collar 30 Box	SNS14502	3/4-2"	A5055
Stoma Cap w/Tape Collar 30 Box	SNS14506	3/4-2"	A5055

Ostomy Allowables

Medicare Allowables...30 day supply

Peri-Stoma Cleanser & Adhesive Remover	50
Adhesive Remover	50
No-Sting Skin Barrier Film	25
No-Sting Skin Barrier Spray.....	2 oz.
Skin Barrier Film	25
Ostomy Skin Barrier Powder.....	5 oz.
Skin Barrier Ring Conforming Seal	20
Skin Barrier Arcs	20
Skin Barrier Sheet	20
Skin Barrier Paste	4 oz.
No-Sting Skin Barrier Paste	4 oz.
Ostomy Appliance Deodorant.....	as necessary
Assure C Odor Eliminator.....	as necessary
Pouches - Drainable.....	20
Stoma Cap.....	31
Tape 1" x 10 yds.....	2
Tape 2" x 10 yds.....	1

Safe n Simple has provided the above information as a guideline of Medicare policies and assumes no responsibility for inconsistencies. All medical and claims determinations will be made by Medicare and your insurance plan.





Hernia/Ostomy
Support Belt

- ★ Revolutionary 2-way stretch fabric is supportive yet comfortable
- ★ The fabric conforms to fit an individual's shape during daily activity
- ★ A unique glove design allows easy control of the belt

Hernia/Ostomy
Support Belt

"As an ostomate and ET nurse, I've been a long time user and advocate of support belts. In my personal experience, both versions of the Safe n Simple Security support belts are amazingly comfortable while providing the support I need in my physically active lifestyle. The conforming fit glides with body movement and isn't visible when wearing snug fitting clothing. Whether the goal is prevention or hernia support, I'd highly recommend this product."

- K. Wickham, MI

Hernia/Ostomy Support Belts

Made from a unique fabric that gives equal and even stretch. This ensures pressure is evenly distributed throughout the garment. The fabric provides complete support to peristomal hernias while still having the flexibility to ensure the normal function of the stoma.

- Security belts were found to be the most comfortable while providing excellent support
- Security belts can be worn for longer periods of time without feeling restricted
- Security belts are fully adjustable so can be used at all times
- Security belt has a soft Velcro fastener which is ideal for easy fitting and can be adjusted as required

Sizes	Measurement Width
Extra Small	up to 29.5 inches up to 75cms
Small	29.5 to 36.5 inches 75cms to 90cms
Medium	36.5 to 41.25 inches 90cms to 105cms
Large	41.25 to 47.25 inches 105cms to 120cms
Extra Large	47.25 to 53.25 inches 120cms to 135cms
XX Large	53.25 to 59 inches 135cms to 150 cms

Visit our website to find out where to order:

www.sns-medical.com

Questions? Please call: **844-767-6334**

Safe n Simple SECUREWEAR
Clean. Comfortable. Confident.

Security Support Belts


Security
Hernia/Ostomy
Support Belt



Support Belt without Hole

Size	Depth	Color	Product#	Size	Depth	Color	Product#
XS	6"/16cm	Beige	SNS36BG1	XS	6"/16cm	Black	SNS36BK1
S	6"/16cm	Beige	SNS36BG2	S	6"/16cm	Black	SNS36BK2
M	6"/16cm	Beige	SNS36BG3	M	6"/16cm	Black	SNS36BK3
L	6"/16cm	Beige	SNS36BG4	L	6"/16cm	Black	SNS36BK4
XL	6"/16cm	Beige	SNS36BG5	XL	6"/16cm	Black	SNS36BK5
XXL	6"/16cm	Beige	SNS36BG6	XXL	6"/16cm	Black	SNS36BK6
XS	8"/20cm	Beige	SNS3BG21	XS	8"/20cm	Black	SNS3BK21
S	8"/20cm	Beige	SNS3BG22	S	8"/20cm	Black	SNS3BK22
M	8"/20cm	Beige	SNS3BG23	M	8"/20cm	Black	SNS3BK23
L	8"/20cm	Beige	SNS3BG24	L	8"/20cm	Black	SNS3BK24
XL	8"/20cm	Beige	SNS3BG25	XL	8"/20cm	Black	SNS3BK25
XXL	8"/20cm	Beige	SNS3BG26	XXL	8"/20cm	Black	SNS3BK26
XS	10.25"/26cm	Beige	SNS3BG61	XS	10.25"/26cm	Black	SNS3BK61
S	10.25"/26cm	Beige	SNS3BG62	S	10.25"/26cm	Black	SNS3BK62
M	10.25"/26cm	Beige	SNS3BG63	M	10.25"/26cm	Black	SNS3BK63
L	10.25"/26cm	Beige	SNS3BG64	L	10.25"/26cm	Black	SNS3BK64
XL	10.25"/26cm	Beige	SNS3BG65	XL	10.25"/26cm	Black	SNS3BK65
XXL	10.25"/26cm	Beige	SNS3BG66	XXL	10.25"/26cm	Black	SNS3BK66

Need a larger or smaller size? Two holes? Special measurements? We can make it custom for you!
HCPC: A4396

Also Available on Amazon

	Size	Depth	Color	Hole	Product#	Size	Depth	Color	Hole	Product#
Support Belt with Hole - 6"/16cm	XS	6"/16cm	Beige	Right	SNS36GR1	XS	6"/16cm	Beige	Left	SNS36GL1
	S	6"/16cm	Beige	Right	SNS36GR2	S	6"/16cm	Beige	Left	SNS36GL2
	M	6"/16cm	Beige	Right	SNS36GR3	M	6"/16cm	Beige	Left	SNS36GL3
	L	6"/16cm	Beige	Right	SNS36GR4	L	6"/16cm	Beige	Left	SNS36GL4
	XL	6"/16cm	Beige	Right	SNS36GR5	XL	6"/16cm	Beige	Left	SNS36GL5
	XXL	6"/16cm	Beige	Right	SNS36GR6	XXL	6"/16cm	Beige	Left	SNS36GL6
	XS	6"/16cm	Black	Right	SNS36KR1	XS	6"/16cm	Black	Left	SNS36KL1
	S	6"/16cm	Black	Right	SNS36KR2	S	6"/16cm	Black	Left	SNS36KL2
	M	6"/16cm	Black	Right	SNS36KR3	M	6"/16cm	Black	Left	SNS36KL3
	L	6"/16cm	Black	Right	SNS36KR4	L	6"/16cm	Black	Left	SNS36KL4
	XL	6"/16cm	Black	Right	SNS36KR5	XL	6"/16cm	Black	Left	SNS36KL5
	XXL	6"/16cm	Black	Right	SNS36KR6	XXL	6"/16cm	Black	Left	SNS36KL6
Support Belt with Hole - 8"/20cm	XS	8"/20cm	Beige	Right	SNS3BGR1	XS	8"/20cm	Beige	Left	SNS3BGL1
	S	8"/20cm	Beige	Right	SNS3BGR2	S	8"/20cm	Beige	Left	SNS3BGL2
	M	8"/20cm	Beige	Right	SNS3BGR3	M	8"/20cm	Beige	Left	SNS3BGL3
	L	8"/20cm	Beige	Right	SNS3BGR4	L	8"/20cm	Beige	Left	SNS3BGL4
	XL	8"/20cm	Beige	Right	SNS3BGR5	XL	8"/20cm	Beige	Left	SNS3BGL5
	XXL	8"/20cm	Beige	Right	SNS3BGR6	XXL	8"/20cm	Beige	Left	SNS3BGL6
	XS	8"/20cm	Black	Right	SNS3BKR1	XS	8"/20cm	Black	Left	SNS3BKL1
	S	8"/20cm	Black	Right	SNS3BKR2	S	8"/20cm	Black	Left	SNS3BKL2
	M	8"/20cm	Black	Right	SNS3BKR3	M	8"/20cm	Black	Left	SNS3BKL3
	L	8"/20cm	Black	Right	SNS3BKR4	L	8"/20cm	Black	Left	SNS3BKL4
	XL	8"/20cm	Black	Right	SNS3BKR5	XL	8"/20cm	Black	Left	SNS3BKL5
	XXL	8"/20cm	Black	Right	SNS3BKR6	XXL	8"/20cm	Black	Left	SNS3BKL6
Support Belt with Hole - 10.25"/26cm	XS	10.25"/26cm	Beige	Right	SNS3GR61	XS	10.25"/26cm	Beige	Left	SNS3GL61
	S	10.25"/26cm	Beige	Right	SNS3GR62	S	10.25"/26cm	Beige	Left	SNS3GL62
	M	10.25"/26cm	Beige	Right	SNS3GR63	M	10.25"/26cm	Beige	Left	SNS3GL63
	L	10.25"/26cm	Beige	Right	SNS3GR64	L	10.25"/26cm	Beige	Left	SNS3GL64
	XL	10.25"/26cm	Beige	Right	SNS3GR65	XL	10.25"/26cm	Beige	Left	SNS3GL65
	XXL	10.25"/26cm	Beige	Right	SNS3GR66	XXL	10.25"/26cm	Beige	Left	SNS3GL66
	XS	10.25"/26cm	Black	Right	SNS3KR61	XS	10.25"/26cm	Black	Left	SNS3KL61
	S	10.25"/26cm	Black	Right	SNS3KR62	S	10.25"/26cm	Black	Left	SNS3KL62
	M	10.25"/26cm	Black	Right	SNS3KR63	M	10.25"/26cm	Black	Left	SNS3KL63
	L	10.25"/26cm	Black	Right	SNS3KR64	L	10.25"/26cm	Black	Left	SNS3KL64
	XL	10.25"/26cm	Black	Right	SNS3KR65	XL	10.25"/26cm	Black	Left	SNS3KL65
	XXL	10.25"/26cm	Black	Right	SNS3KR66	XXL	10.25"/26cm	Black	Left	SNS3KL66



Safe n Simple focuses on product quality, patient safety, and cost containment. With the introduction of the Simpurity Wound Care product line we are proud to offer wound care products that follow evidence based medicine, facilitate moist wound healing and improve patient outcomes while lowering costs.

Your patients deserve the best care available, with the proven quality of Simpurity wound care product solutions, you can help them achieve greater outcomes and save money doing it.

Simpurity Wound Care Solutions:

- Alginate - Dressing Pads & Rope
- Alginate Silver - Dressing Pads & Rope
- Collagen - Dressing Pads & Powder
- Collagen Silver - Dressing Pads
- DermaPro Waterproof Plastic Tape
- DermaPro Waterproof Silicone Tape
- Fibergel - Dressing Pads
- Fibergel Silver - Dressing Pads
- Foam - Dressing Pads
- Foam - Adhesive Bordered Dressing Pads
- Foam - Silicone Adhesive Dressing Pads
- Foam - Silicone Bordered Dressing Pads
- Foam - Silicone Silver Foam Dressing Pads
- Foam - Silver Dressing Pads
- Hydrocolloid Dressing Pads
- Hydrogel Impregnated Gauze Dressings
- Hydrogel Dressing Sheets
- Polymer - Dressing Pads
- Wound Care Kits

For more information or to request samples please call our toll-free number or visit our website.

844-767-6334 | www.SnS-Medical.com

Enhancing quality of life with exceptional products at an affordable price



Testimonial.....

"I personally use Safe n Simple wound care products in our wound center for several reasons. The Simpurity line of products offers high quality alginate, foam and collagen dressings with and without silver at a very cost effective price. When everyone is looking to use evidenced based treatments and work within our current reimbursement environment this is the PERFECT product line. The products are user friendly. The packaging is appealing and instructions clear. Most importantly our healing rates and results have been excellent."

J. Christopher Brandys, MD FRCS FACS
General and Vascular Surgeon
Director of The Wound Center
Medical Consultant - Safe n Simple

Alginates

Alginate Wound Dressing



Simpurity Alginate Dressings are derived from seaweed. The products are composed of soft, calcium alginate fibers shaped as pads (fibrous mats) or ropes (twisted fibers). When packed in a wound, the alginate interacts with wound exudate to form a soft gel that maintains a moist healing environment. The dressings release of calcium ions can activate prothrombin and improve hemostasis. Has a strong tensile strength.

Item	Product#	Size	HCPC
Dressing Pad, 10 Box	SNS50702	2"x2"	A6196
Dressing Pad, 10 Box	SNS50704	4"x4"	A6196
Dressing Pad, 10 Box	SNS50720	4"x5"	A6197
Dressing Pad, 5 Box	SNS50732	4"x8"	A6197
Dressing Pad, 10 Box	SNS50764	8"x8"	A6198
Rope, 30 Bag	SNS50712	12"	A6199

Silver Alginate Wound Dressing




Simpurity Silver Alginate Dressings are made of natural fibers derived from seaweed through a series of special processes. Simpurity Silver Alginate Dressings can inhibit bacteria growing in the wound bed reducing risk of infection while maintaining a moist environment for wound healing. The dressing has a strong tensile strength, a high absorption capacity and is easy to remove from the wound bed.

Item	Product#	Size	HCPC
Dressing Pad, 10 Box	SNS51702	2"x2"	A6196
Dressing Pad, 10 Box	SNS51720	4"x5"	A6197
Dressing Pad, 10 Box	SNS51764	8"x8"	A6198
Rope, 10 Box	SNS51712	12"	A6199

Fibergels

Fibergel Wound Dressing



Simpurity Fibergel Wound Dressing is a soft, absorbent material that forms into a gel when in contact with exudate. The gelling action creates an optimal environment for moist wound healing. It conforms easily to various wound shapes, supports debridement and can easily be removed in one piece with minimal damage to newly forming granulating tissues. Use as a primary dressing.

Item	Product#	Size	HCPC
Dressing Pad, 10 Box	SNS56702	2"x2"	A6196
Dressing Pad, 10 Box	SNS56704	4"x4"	A6196
Dressing Pad, 10 Box	SNS56745	4"x5"	A6197

Silver Fibergel Wound Dressing



Simpurity Fibergel Ag Wound Dressing has all the benefits of a hydrofiber dressing plus ionic silver to to inhibit bacteria from growing in the wound bed, reducing the risk of infection. The gelling action of Fibergel AG creates an optimal environment for moist wound healing. It conforms easily to various wound shapes, supports debridement and can easily be removed in one piece with minimal damage to newly forming granulating tissues. Use as a primary dressing.


Item	Product#	Size	HCPC
Dressing Pad, 10 Box	SNS56712	2"x2"	A6196
Dressing Pad, 10 Box	SNS56716	4"x4.75"	A6197

Collagens

Collagen Wound Dressing



Smpurity Collagen provides ingredients for tissue regeneration in phase 3 wounds and the 3D scaffold structure of Smpurity Collagen allows fibroblasts to migrate across the wound bed facilitating granulation, angiogenesis and tissue repair. Smpurity Collagen dressings are 100% pure collagen, non-bleached, native undigested bovine collagen. The dressings provide high exudate absorption and retention providing a moist wound healing environment. Easy application and superior conformity to the wound.



Item	Product#	Size	HCPC
Dressing Pad, 10 Box	SNS50002	2"x2"	A6021
Dressing Pad, 5 Box	SNS52222	2"x2"	A6021
Dressing Pad, 5 Box	SNS52244	4"x4"	A6021
Vial Particles, 5 Box	SNS5001G	1 Gram	A6010
Packet Particles, 10 Box	SNS5221G	1 Gram	A6010

Note: SNS50002 is a premium brand collagen wound dressing.

Silver Collagen Wound Dressing



Coming in 2018!

Smpurity Collagen dressing are 100% pure collagen, non-bleached, native undigested bovine collagen, now with Ionic Silver added to the formulation. The dressings provide high exudate absorption and retention providing a moist wound healing environment. Easy application and superior conformity to the wound.

Item	Product#	Size	HCPC
Dressing Pad, 5 Box	SNS51222	2"x2"	A6021
Dressing Pad, 5 Box	SNS51244	4"x4"	A6021

Foams

Foam Wound Dressing



Smpurity Foam Wound Dressings are designed to help manage moderate-to-highly draining wounds. The dressing absorbs wound exudate, promotes moist wound healing, protects the wound from contaminants and permits moisture vapor exchange. These foam dressings have a protective colored film backing to protect against strike through.

Item	Product#	Size	HCPC
Dressing Pad, 12 Box	SNS51W02	2"x2"	A6209
Dressing Pad, 12 Box	SNS51W04	4"x4"	A6209
Dressing Pad, 12 Box	SNS51W24	4"x6"	A6210

Foam Wound Dressing - Adhesive Border



Smpurity Foam wound dressings are composed of a soft and highly absorbent polyurethane foam pad to absorb wound drainage and a transparent adhesive border to hold the dressing in place. This island foam dressing helps maintain a moist environment, which accelerates wound healing. It is especially suitable for moderate to highly exuding wounds and has a protective backing to protect against strike through. The dressings are comfortable, easy to apply, conformable and easy to remove from the wound.

Item	Product#	Size	HCPC
Dressing Pad, 12 Box	SNS72302	2"x2"	A6212
Dressing Pad, 12 Box	SNS72320	4"x5"	A6213

Foams

Foam Wound Dressing - Ionic Silver



Smpurity Foam Dressing is a silver ion-containing polyurethane foam dressing that protects wounds from bacterial penetration while reducing the bioburden in a wound bed. It mainly consists of hydrophilic polyurethane foam substrate and silver compound. It is indicated for the management of moderate-to-high exudate for chronic and acute wounds and helps maintain a moist healing environment. Controlled release of ionic silver provides sustained antibacterial effect. Inhibits a broad-spectrum of infection up to 99.99% for a wide variety of bacteria. Provides longer wear times and pain-free removal.

Item	Product#	Size	HCPC
Dressing Pad, 12 Box	SNS72502	2"x2"	A6209
Dressing Pad, 12 Box	SNS72520	4"x5"	A6210

Foam Wound Dressing - Silicone



Smpurity Silicone Foam Dressing is a superior absorbent dressing for moderate to high wound exudate management. It contains a highly absorbent foam pad and gentle adhesive silicone contact layer. These dressings help maintains a moist wound healing environment with no irritation upon removal, minimizes trauma to the wound and reduces pain during dressing changes.

Item	Product#	Size	HCPC
Dressing Pad, 12 Box	SNS74502	2"x2"	A6209
Dressing Pad, 12 Box	SNS74520	4"x5"	A6210

Foams

Foam Wound Dressing - Bordered Silicone



Simpurity Silicone Foam Bordered Dressing is a superior absorbent dressing for moderate to high wound exudate management. It contains a highly absorbent foam pad and gentle adhesive silicone contact layer and an adhesive border to hold the dressing in place. These dressings help maintain a moist wound healing environment with no irritation upon removal, minimizes trauma to the wound and reduces pain during dressing changes.

Item	Product#	Size	HCPC
Dressing Pad, 12 Box	SNS77733	3"x3"	A6212
Dressing Pad, 12 Box	SNS77744	4"x4"	A6212

Foam Wound Dressing - Silver Silicone



Simpurity Silicone Foam Dressing is a silver ion containing polyurethane foam pad that has a gentle adhesive silicone contact layer that protects wounds from bacterial penetration while reducing the bioburden in the wound bed. It is indicated for the management of moderate-to-high exudate for chronic and acute wounds and helps maintain a moist healing environment. Controlled release of ionic silver provides sustained antibacterial effect. Provides longer wear times and pain-free removal.

Item	Product#	Size	HCPC
Dressing Pad, 12 Box	SNS74422	2"x2"	A6209
Dressing Pad, 12 Box	SNS74445	4"x5"	A6210
Dressing Pad, 12 Box	SNS74466	6"x6"	A6210
Dressing Pad, 12 Box	SNS74477	7"x7"	A6210

Gauze Sterile Sponge



Simpurity Gauze Sterile Sponges are ideal for general wound cleansing, debriding, minor prepping, covering and packing wounds. It is a woven sponge made of 12 ply 100% USP cotton. Contains 2 sponges per packet.

Item	Product#	Size	HCPC
Dressing Pad, 2 Pcs/pkt, 25 pkts	SNS54444	4"x4"	A6402

Bordered Gauze Sterile Sponge



Simpurity Bordered Gauze Sterile Sponges are a non-adherent gauze sponge island dressing with an adhesive border. These woven sponges are made of 12ply 100% USP cotton for use as a secondary dressing. The sponge absorbs exudate for a cleaner wound bed and the adhesive border protects the wound from outside contaminants. Note: 4"x4" has gauze pad size of 2"x2", and 6"x6" has gauze pad size of 4"x4".

Item	Product#	Size	HCPC
Dressing Pad, 30 Pads	SNS52316	4"x4"	A6219
Dressing Pads, 30 Pads	SNS52324	6"x6"	A6220

Hydrocolloids

Hydrocolloid Wound Dressing



Simpurity Hydrocolloid Wound Dressing is a thin hydrocolloid sheet. Simpurity Hydrocolloids are absorbent, and delivers a gentle and secure adhesion that provides a barrier to external debris and bacteria when securely adhered to the skin. Individually packaged, sterile and latex free. For use on light-to-moderately exuding wounds.

Item	Product#	Size	HCPC
Dressing Pad, 12 Box	SNS55402	2"x2"	A6234
Dressing Pad, 12 Box	SNS55404	4"x4"	A6234
Dressing Pad, 10 Box	SNS21149	7"x7"	A6236

Polymer Pads

Polymer Wound Dressing - Highly Absorbent



Simpurity Polymer Foam Wound Dressing is a thin, sterile, superabsorbent single-use dressing, for use on highly exuding wounds. Simpurity Polymer Foam absorbs the exudate and locks it into the super-absorbent core, allowing for significantly more fluid to be handled by the dressing while keeping it away from the wound and surrounding tissue. This can help to reduce maceration, as well as reduce exposure to harmful components of wound fluid.

Item	Product#	Size	HCPC
Dressing Pad, 10 Box	SNS59020	4"x5"	A6251
Dressing Pad, 10 Box	SNS59080	8"x10"	A6251

Hydrogels

Hydrogel - Impregnated Gauze



Smpurity Hydrogel Impregnated Gauze Wound Dressing is designed to provide a moist wound healing environment for wounds with minimal or no exudate. Contains a high viscosity gel impregnated into the gauze pad. Easy to apply.

Item	Product#	Size	HCPC
Dressing Pad, 10 Box	SNS58802	2"x2"	A6231
Dressing Pad, 10 Box	SNS58820	4"x5"	A6232

Hydrogel Dressing Sheet



Smpurity Hydrogel Absorbent Dressing Sheets are designed to provide a moist wound healing environment for wounds with minimal to no exudate. It is permeable to air and moisture vapor for faster wound healing. Easy to remove from the skin without damaging tissues. Conforms to the wound bed preventing bacterial infection and contamination. The dressing can be cut to fit the wound site. Latex free.

Item	Product#	Size	HCPC
Dressing Pad, 10 Box	SNS58302	2"x2"	A6242
Bordered Dressing Pad, 12 Box	SNS58312	2"x2"	A6245

DermaPro Tapes

DermaPro Waterproof Plastic Tape



Smpurity DermaPro Waterproof Plastic Tape is superior for securing wound dressings and repelling moisture from outside contaminants. DermaPro has a latex free adhesive and is gentle to remove from sensitive skin areas. The roll contains a start strip for easy dispensing.

Item	Product#	Size	HCPC
Plastic Tape, 1 Roll	SNS57220	2"x15'	A4452

DermaPro Waterproof Silicone Tape



Smpurity DermaPro Waterproof Silicone Tape is latex free and has a silicone adhesive that is gentle to remove from sensitive skin areas. Superior for securing wound dressings and repelling moisture from outside contaminants. Contains 12-1" rolls per box, and 6-2" rolls per box.

Item	Product#	Size	HCPC
Silicone Tape, 12 Rolls	SNS57230	1"x15'	A4452
Silicone Tape, 6 Rolls	SNS57232	2"x15'	A4452



Wound Care Kits



Eliminate confusion and ensure your patients receive the products you prescribe

Simpurity wound care kits are the perfect answer to ensure patient compliance. Typical compliance rates are 60-70% higher when patients leave the practice with the product in hand vs. traveling to the pharmacy. These pre-packaged kits contain our premium, pre-selected dressings designed to address most wound types. The easy to dispense wound care kits include a 15 or 30-Day supply of primary dressings, secondary dressings, gauze sponges and a saline wound wash.

RECOMMENDED FOR: Diabetic ulcers, pressure ulcers, stasis ulcers, post-surgical wounds and donor sites.

Simpurity wound care kits help you keep control of the dressing choice to achieve optimum wound outcomes, see page 30 for a list of our representatives to call for pricing.

COMPLETE KIT CONTENTS: Each Kit Includes 15/30 Primary Dressings of your Choice along with 15/30 Simpurity Bordered Gauze Sponge Dressings, 15/30 Packets of Simpurity Gauze Sponges and a bottle of Saline Wound Wash.



KIT CONTAINS:

PRIMARY DRESSINGS:

Smpurity Alginate 2"x2" Dressing Pads
 Smpurity Silver Alginate 2"x2" Dressing Pads
 Smpurity Silver Fibergel 2"x2" Dressing Pads
 Smpurity Collagen 2"x2" Dressing Pads
 Smpurity Collagen Particles 1g Packets

Item	Product#	Amount
Alginate 2"x2" Pad Wound Care Kit	SNSWCKA1	15 Day Supply
Alginate 2"x2" Pad Wound Care Kit	SNSWCKA2	30 Day Supply
Silver Alginate 2"x2" Wound Care Kit	SNSWCKS1	15 Day Supply
Silver Alginate 2"x2" Wound Care Kit	SNSWCKS2	30 Day Supply
Silver Fibergel 2"x2" Wound Care Kit	SNSWCKF3	15 Day Supply
Silver Fibergel 2"x2" Wound Care Kit	SNSWCKF4	30 Day Supply
Collagen 2"x2" Wound Care Kit	SNSWCKC1	15 Day Supply
Collagen 2"x2" Wound Care Kit	SNSWCKC2	30 Day Supply
Collagen 1g Pkt Wound Care Kit	SNSWCKP1	15 Day Supply
Collagen 1g Pkt Wound Care Kit	SNSWCKP2	30 Day Supply

Wound Care Allowables

Medicare Allowables...30 day supply

Medicare Part B - Indications and Limitations of Coverage for Surgical Dressings

Alginate Pad	30 each
Alginate Rope (per 6").....	60 inches
Collagen Pad	30 each
Collagen Powder	30 grams
Foam Pads.....	12 each
Gauze Pad - Sterile with Border.....	30 each
Hydrocolloid Pad	12 each
Hydrogel Impregnated Gauze	30 each
Hydrogel Pad	30 each
Hydrogel Pad with Border	12 each
Specialty Absorbent Dressing	30 each

Safe n Simple has provided the above information as a guideline of Medicare policies and assumes no responsibility for inconsistencies. All medical and claims determinations will be made by Medicare and your insurance plan.

Stay Connected....

Contact Our Team...

1. Call Safe n Simple Customer Service
Toll Free 844-767-6334 or
Amy Ham 248-214-4877
Mon-Fri 9am-5pm EST
2. Find a Supplier or Order from our E-Store at:
www.sns-medical.com
3. Mailing Address: 8112 Rosebud Lane
Clarkston, MI 48348
Shipping Address: 5827 Terex
Clarkston, MI 48346

Inside Sales Team...

Holly Loos	248-496-2508
Michelle Pitylak	248-212-2243

Outside Sales Team...

Midwest:	Joseph LeRoy	248-974-3245
Northeast:	Michele Endreszl	248-977-8324
Southeast:	Beth Peter	561-526-6718
Ind. Rep. CA:	Dan Penrose	909-730-0707
Ind. Rep. IN/OH:	Joseph DeHeer	317-771-4319
Ind. Rep. MD/VA:	Jim Spradlin	804-752-2897
Ind. Rep. Canada:	Peter West	519-726-5959

Samples Team...

To request samples call our customer service team
or email us at: samples@sns-medical.com

Clinicians/New Ostomy Patients:
Request a **SimpleStart Ostomy Sample Kit**
email us at: simplestart@sns-medical.com
or fax us at: 248-575-4172

Stay Connected....

Supplier Orders:

1. Call Safe n Simple Customer Service
Toll Free 844-767-6334 or
Amy Ham 248-214-4877
Mon-Fri 9am-5pm EST
2. Email Us: purchaseorder@sns-medical.com
3. FAX Us: 248-698-3454
4. Mailing Address: Shipping Address:
8112 Rosebud Lane 5827 Terex
Clarkston, MI 48348 Clarkston, MI 48346

Marketing Team:

Have a special event or sponsorship opportunity you would like to share with Safe n Simple?

Please feel free to email us: marketing@sns-medical.com

Call Us:
Michelle Bliszack - Marketing Manager
248-408-9501 Direct
844-767-6334 Toll Free

At Safe n Simple we are committed to providing the highest quality products while keeping them affordable to the consumer and providing them the best possible outcomes.

Log on and “like” our Facebook page.

You will have easy access to promotions, new product updates and much more!

VISIT US TODAY!

<https://www.facebook.com/Safe-n-Simple>

Join our Mailing List - Visit: www.sns-medical.com



Our Mission:

Safe n Simple is an innovator, developer and marketer of new medical products. We enhance the patient's quality of life with superior products at an affordable price bringing value to the communities in which we live.

Safe n Simple Ostomy Care Products

Adhesive Remover
Peri-Stoma Cleanser & Adhesive Remover
No-Sting Skin Barrier Film
Skin Barrier Film
Ostomy Skin Barrier Powder
Skin Barrier Ring - Conforming Seal
Skin Barrier Ring - Integrity Seal
Skin Barrier Arc
X-Tra Wide Skin Barrier Arc
1/2 Arc Skin Barrier Strips
Skin Barrier 4"x4" Sheet
Skin Barrier Paste
No-Sting Skin Barrier Paste
Skin Barrier Paste - Fast Drying
Ostomy Pouch Deodorant
Assure C Odor Eliminator
One-Piece, CTF Pouches
One-Piece, CTF Pouches with Filters
Stoma Caps
Security Support Belts

Simpurity Wound Care Products

Alginate - Dressing Pads & Rope
Alginate Silver - Dressing Pads & Rope
Collagen - Dressing Pads & Powder
Collagen Silver - Dressing Pads
DermaPro Waterproof Plastic Tape
DermaPro Waterproof Silicone Tape
Fibergel - Dressing Pads
Fibergel Silver - Dressing Pads
Foam - Dressing Pads
Foam - Adhesive Bordered Dressing Pads
Foam - Silicone Adhesive Dressing Pads
Foam - Silicone Bordered Dressing Pads
Foam - Silicone Silver Foam Dressing Pads
Foam - Silver Dressing Pads
Hydrocolloid Dressing Pads
Hydrogel Dressing Sheets
Hydrogel Gauze Dressing Pads
Polymer Dressing Pads
Wound Care Kits



Safe n Simple™
Clean. Comfortable. Confident.

Clarkston, MI 48346 • P: 844.767.6334 • F: 248.698.3454

www.sns-medical.com

Email: samples@sns-medical.com

Like us on Facebook:

<https://www.facebook.com/Safe-n-Simple>