

D3W

3MHz Waterproof Probe

Newman Medical

Newman Medical ultrasound probes have been designed to provide excellent sensitivity and durability to withstand everyday office use. The **DigiDop 3 MHz Waterproof Obstetrical Probe** has a large ultrasound beam optimized for deep penetration for late term pregnancy or larger patients. The probe has been designed for water applications and for comfort for both the provider and the patient. **Made in the USA in Colorado.** Each Newman Medical Doppler and probe must pass a series of comprehensive tests to assure quality and the durability of each product.



Features:

- ✦ **Best Sound** - Digitally Optimized Processing – DOP™ produces up to 6x signal of other Dopplers
- ✦ **Best Warranty** – 5 year \$0 cost warranty – **ALL** Parts. **ALL** Labor.
- ✦ **Best Support & Service** – Knowledgeable Staff & Guaranteed Service

- ✦ 3 MHz waterproof large beam, obstetrical probe
- ✦ Optimized for early fetal heartbeat detection and use throughout pregnancy
- ✦ For use with DigiDop series of Dopplers
- ✦ User replaceable cord
- ✦ Dimensions: 100 x 35 mm (4 x 1.5 in)
- ✦ Weight: 20 grams (0.75 ounces)



Newman Medical
5350 Vivian Street | Unit C | Arvada | CO | 80002
Toll-free 800-267-5549 | Fax-888-202-3627
www.newman-medical.com

NM0004B