

D5

5MHz Vascular Probe

Newman Medical

Newman Medical ultrasound probes have been designed to provide excellent sensitivity and durability to withstand everyday office use. The **DigiDop 5 MHz Vascular Probe** has a larger ultrasound beam optimized for deeper vessels in vascular applications. The probe has been designed for comfort for both the provider and the patient. **Made in the USA in Colorado.** Each Newman Medical Doppler and probe must pass a series of comprehensive tests to assure quality and the durability of each product.

DigiDop 5MHz Probe



Features:

- ✦ **Best Sound** - Digitally Optimized Processing – DOP™ produces up to 6x signal of other Dopplers
- ✦ **Best Warranty** – 5 year \$0 cost warranty – **ALL** Parts. **ALL** Labor.
- ✦ **Best Support & Service** – Knowledgeable Staff & Guaranteed Service

- ✦ 5 MHz larger beam, vascular probe
- ✦ Optimized for deeper vessels in vascular applications
- ✦ For use with DigiDop series of Dopplers
- ✦ User replaceable cord
- ✦ Dimensions: 115 x 25 mm (4.5 x 1 in)
- ✦ Weight: 20 grams (0.75 ounces)



Newman Medical
5350 Vivian Street | Unit C | Arvada | CO | 80002
Toll-free 800-267-5549 | Fax-888-202-3627
www.newman-medical.com

NM0004B