

# DigiDop Dopplers

"Built to Last"

**TOUGH** and *Sensitive*



## Configure to your application -

Interchangeable vascular & obstetrical probes. Smooth. Comfortable. Easy clean.

## Five Year Warranty- Includes Probes

Zero Cost Warranty includes parts and labor on Doppler and probes.

## Sensitive - More listening, less searching

Unique digitally Optimized Processing (Dop™) provides up to 6x signal of other Dopplers.

## Minimize Service Returns

Flexible, magnetic probe holder won't crack or break. Tough stainless/Polycarbonate case.

## Save \$100's in service costs

User replaceable cords, great warranty

MADE IN  
  
THE USA



## About Newman Medical

The management of Newman Medical has been in the Doppler/vascular business since 1976, starting with founding Imex Medical Systems. They were also the founders of the Vascular Disease Foundation, a non-profit patient education foundation.

Call us at 800.267.5549

[www.newman-medical.com](http://www.newman-medical.com)

Newman Medical Dopplers are available on Premier & MediGroup GPO contracts

# What makes the DigiDop difference!

## ① Best Sound – Digitally Optimized Processing

The DigiDop provides a unique powerful combination of tough durability along with proven, measured sensitivity. We know that the real key to a great Doppler is performance – sensitivity and sound quality. While others may say this; we measure it – and have found that our signal is more than six times greater than some other Dopplers. Our Digitally Optimized Processing - DOP™ - eliminates most background noise while letting the Doppler sounds through bright and clear. With DOP each interchangeable probe can be optimized with different filtering in the main unit.

## ② Best Warranty -5 year \$0 cost warranty -ALL parts. ALL labor.

**True Lower Costs of Ownership – Five Year \$0 Warranty – Including Probes** Our 5 year warranty covers parts and labor, including our probes and cords. This can easily be a savings of hundreds of dollars compared to other Dopplers on a single service issue. We are able to offer this warranty because our design and tough materials significantly reduce the number and cost of warranty issues.

## ③ Best Support & Service – Knowledgeable Staff & Guaranteed Service

Description	Model	5 Year \$0 Warranty	Unbreakable Probe Holder	Interchangeable Probes?	Display?	Rechargeable?
DigiDop 300	DD-300-**	✓	✓	✓		
DigiDop 301	DD-301-**	✓	✓	✓		✓
DigiDop 700	DD-700-**	✓	✓	✓	✓	
DigiDop 701	DD-701-**	✓	✓	✓	✓	✓

\*\* When ordering, place the desired probe number from the probe section below in this location.

### Probe Options

D2	2MHz obstetrical probe - for deeper penetration
D2W	2MHz waterproof obstetrical probe for deeper penetration/water applications
D3	3MHz obstetrical probe - early heartbeat detection & general use
D3W	3MHz waterproof obstetrical probe for general use water applications
D5	5MHz vascular probe - pen style vascular for deeper vessels
D8	8MHz vascular probe - pen style for superficial vessels & general use

### Technical Specifications

Weight	340 grams (12 ounces)
Dimensions	150 x 65 x 35 mm (6.0 x 2.5 x 1.25 in)
Battery Life	900 minutes Alkaline/720 minutes NiMH
Warranty	5 years <b>ALL</b> parts <b>ALL</b> labor

### Obstetrical Probes – Sensitive, Durable, Versatile

Our interchangeable obstetrical probes are easy to hold and clean. Each Doppler probe is individually tuned for maximum sensitivity. The probe cords are user replaceable.



#### D2

2MHz – Deeper penetration



#### D2W

Waterproof 2MHz for water bath/birth



#### D3

3MHz – Very sensitive to early fetal heart, best all around OB probe



#### D3W

Waterproof 3MHz probe for water bath/birth

### Vascular probes are also available:

D5 – 5MHz, D8 – 8MHz, PPG – audio ppg

### DigiDop Dock Options

Secure your DigiDop by mounting on the wall, an IV pole, a roll-stand, or a table dock. Theft resistant fastener keeps it there.



Roll Stand Option  
includes basket



Table Dock Option



IV Pole Option



Wall Option

### How to order!

Simply add probe choice to end of main unit choice  
For example, DD-700-D3 is a Display Doppler with 3MHz probe