

DPPG

Photoplethysmography Probe

Newman Medical

Newman Medical PPG probes, like our Doppler probes, have been designed to provide excellent sensitivity and durability for everyday use. The unique **DigiDop PPG Vascular Probe** can be used on any DigiDop and provides an audible signal/beep for each pulse. It is an excellent choice for taking a toe pressure for the Toe/Brachial Index exam and also in any application for monitoring of blood flow at the extremities. The PPG probe's sensor is independently replaceable. **Made in the USA in Colorado.** Each Newman Medical Doppler and probe must pass a series of comprehensive tests to assure quality and the durability of each product.



Features:

- ✦ **Best Sound** - Digitally Optimized Processing – DOP™ produces up to 6x signal of other Dopplers
- ✦ **Best Warranty** – 5 year \$0 cost warranty – **ALL** Parts. **ALL** Labor.
- ✦ **Best Support & Service** – Knowledgeable Staff & Guaranteed Service
- ✦ Replaceable PPG sensor
- ✦ For use with patients with incompressible arteries
- ✦ For use with DigiDop series of Dopplers
- ✦ User replaceable cord
- ✦ Dimensions: 100 x 25 mm (4 x 1 in)
- ✦ PPG Sensor with 3.5 foot cord
- ✦ Weight: 20 grams (0.75 ounces)



Newman Medical
5350 Vivian Street | Unit C | Arvada | CO | 80002
Toll-free 800-267-5549 | Fax-888-202-3627
www.newman-medical.com

NM0004B